# 3.0T 16ch Shoulder Coil

## Service Manual

M7100SH / NC143002

Distributed by GE Healthcare

**NeoCoil**  
Innovations to **sharpen** your image

NC143SCM Rev B  
June 2020  
©2020 NeoCoil, LLC  
All rights reserved.

---

## Medical Device Directive

This product conforms to the requirements of the council directive 93/42/EEC concerning medical devices when it bears the following CE Marking of Conformity; first year of CE Mark is 2020:

The CE Mark logo consists of the letters "CE" in a bold, sans-serif font, enclosed within a black square.

## European Representative

**EC REP**

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

**CAUTION:** The waste of electrical and electronic equipment must not be disposed as unsorted municipal waste and must be collected separately. Contact an authorized representative of the manufacturer for information concerning the decommissioning of your equipment.

## US Federal Law 21 CFR Part 801

**CAUTION:** Federal law restricts this device to sale by or on the order of a physician.

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
---

# ATTENTION!

Installation of this coil **MUST** be performed by a GE field engineer. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration can result in system errors.

---

THIS PAGE LEFT  
INTENTIONALLY BLANK.

# 1 Overview

## Device Description

The 16ch Shoulder Coil is a receive-only phased array RF coil designed for optimum signal-to-noise ratio (SNR) and uniform coverage of the shoulder anatomy for use with GE Healthcare Magnetic Resonance Imaging (MRI) scanners. The coil receives magnetic resonance signals generated in hydrogen nuclei (protons) in the shoulder while blocking the high-frequency magnetic field applied by the MRI scanner at specified timings.


The received signals are amplified before being transferred to the MRI scanner through the coil's system cable. The amplified signals are processed into tomographic images of the shoulder anatomy by the MRI scanner.

Images are typically generated as axial, sagittal, coronal oblique slices and include coverage of the Humerus, Humeral Head, Labrum, Labral Tear, Glenoid, Scapula, Clavicle, and Rotator Cuff regions of the shoulder anatomy.

During a shoulder MRI exam, the patient lays in the supine position on the MRI scanner table. The rigid Posterior Array is positioned beneath the patient's affected shoulder. The flexible Anterior Array is positioned on top of the patient's affected shoulder. The 16ch Shoulder Posterior Housing Pad and 16ch Shoulder Pad provided with the coil promote patient comfort and facilitate patient positioning. Table Straps provided with the MRI scanner may be used to minimize unintended patient movement during the exam.

## Shipping Contents

**NOTE:** Items listed below are included with your original coil purchase. For replacement parts, refer to the *FRUs and Replacement Components* section of this manual.

No.	Item	Dimensions (cm) L x W x H	Weight (kg)
1	<b>16ch Shoulder Coil</b> 	112.2 cm x 28.2 cm x 18.2 cm	3.0 kg
2	<b>16ch Shoulder Posterior Housing Pad</b> 	37.4 cm x 28.0 cm x 0.6 cm	0.2 kg
3	<b>16ch Shoulder Pad</b> 	52.0 cm x 52.0 cm x 5.7 cm	0.6 kg

---

## Related Documentation

- 16ch Shoulder Coil Instructions for Use (English), NeoCoil Part # NC143IFU-EN

---

## Damage In Transportation

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or signed for by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier immediately upon discovery or, in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within a 14-day period.

To file a report with GE:

- Call 1-800-548-3366 and use option 6.
- Contact your local Service Coordinator for more information on this process.

## Product Locator

Each serialized component is supplied with a product locator card (PLC). The number of PLCs for any given site will depend upon the MCAT that's been ordered. Make sure to submit the complete number of installation cards for each catalog.

The Global Install Base Database (also known as the Product Locator System) tracks shipment, transshipment, and field location of the serialized models. There are two methods for submitting Product Locator information.

**Method 1:** For **U.S. ONLY** at this time, the preferred method of submitting information is the FE. Site Verification website at: <http://gib.gehealthcare.com/gib/>

The FE Site Verification has two components available on the web from the main menu:

- Installation/deinstallation product locator model and serial numbers.
- Ability to add/modify the ship to address information.

**Method 2:** Once the "Shipping Card" is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the "Installation Card" and extra shipment cards are attached:

- Verify that serial and model number on each rating plate matches the installation card numbers before removing the installation card.
- Process just the installation card and discard any extra shipment cards and labels. (There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination.)

## Environmental Requirements

The Shoulder Coil shall be stored in a room that meets the temperature and humidity parameters listed below.

	Transport/Storage	Operating Conditions
Humidity	5% to 95% Non-condensing	10% to 90% Non-condensing
Temperature	-40° C to 70° C -40° F to 158° F	15° C to 21° C 59° F to 69.8° F
Atmospheric Pressure	50 kPa - 106 kPa	

The pads for the Shoulder Coil shall be stored in a room that meets the temperature and humidity parameters listed below.

	Transport/Storage	Operating Conditions
Humidity	5% to 95% Non-condensing	10% to 90% Non-condensing
Temperature	-40° C to 55° C -40° F to 131° F	15° C to 21° C 59° F to 69.8° F
Atmospheric Pressure	50 kPa - 106 kPa	

## Coil Specifications


Item	Dimensions (cm) L x W x H	Weight (kg)
16ch Shoulder Coil	112.2 cm x 28.2 cm x 18.2 cm	3.0 kg
16ch Shoulder Posterior Housing Pad	37.4 cm x 28.0 cm x 0.6 cm	0.2 kg
16ch Shoulder Pad	52.0 cm x 52.0 cm x 5.7 cm	0.6 kg

## Use of Symbols

Coil users should understand the following symbols used for caution and warning explanations:

Symbol	Description
	<b>CAUTION:</b> Indicates a hazard that may cause a minor injury or equipment damage.
	<b>WARNING:</b> Indicates a hazard that may cause serious injury.
	Accessory for an MRI scanner to specifically and directly assist the medical functionality of the MRI scanner in terms of its intended purpose.
	Do not loop coil cable.
	Dipose of RF coil properly.
	Manufacturer
	Date of manufacture
	Batch code
	Catalogue number

## Use of Symbols (continued)

Symbol	Meaning
	Serial number
	Consult Instructions for Use
	Instructions for Use are available in an electronic format.
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Receive only RF coil
	CE Marking of Conformity
	Authorized representative in the European Community
	ETL Recognized Component (Canada and U.S.)

## Coil Classification

Symbol	Classification & Ratings
	<b>Class II equipment</b> Coil is Class II equipment in accordance with the international safety standard IEC 60601-1.
	<b>Type BF applied</b> Coil is type BF applied parts, complying with the specified requirements of the international safety standard IEC 60601-1.

This equipment is acceptable for Continuous Operation, having ordinary protection against ingress of liquids or fluids (IPX0).

All components of this system are considered non-sterile and are not for use in the presence of flammable anesthetics.

---

## Coil Care

**WARNING:** Electric shock may occur if the coil is attached to the system during cleaning or when it is still wet. Detach cable from the scanner before attempting to clean. Do not touch connectors with bare fingers. Never press sharp objects against connector surface. Do not reattach connector after cleaning until the coil has dried completely.

**CAUTION:** Do not spray or pour cleaning solution directly onto the coil. Do not submerge any of the parts in the solution. Parts of the equipment contain sensitive electronic components that may be damaged by the solution. Under no circumstances should the coil be placed into any type of sterilizer. All parts should be cleaned only according to the procedure described below.

Thoroughly clean the coil and accessories using the following procedure:

- Wipe with a cloth that has been dampened in a solution of 10% bleach and 90% tap water or 30% isopropyl alcohol and 70% tap water.
- Let the surfaces dry before use.

**WARNING:** Do not clean coil while it is in use with a patient or connected to the MRI scanner.

After cleaning, visually inspect the following to ensure that the coil and its accessories remain safe and effective:

- surfaces of the coil and accessories are dry.
- soil and contaminants have been removed.

Should any part of the coil need to be returned for service, wipe it down with a 10% bleach solution (as described above) to minimize risk of exposure to potentially infectious agents.

Store the Shoulder Coil in an air conditioned scan or equipment room that meets the temperature and humidity parameters listed under “Environmental Requirements”. Do not store other coils on top of the Shoulder Coil.

Perform maintenance and inspection according to your site’s protocols.

# 2 Installation

Installation and configuration must be performed by your service representative. No calibrations are required upon installation or replacement of the 16ch Shoulder Coil by a field service engineer.

## Multi-Coil Quality Assurance (MCQA)

### MCQA requirements

Required Persons	Preliminary Requirements	Procedure	Finalization
1	0 minutes	20 minutes	0 minutes

### MCQA setup

**NOTE:** The device system cable must be routed away from the bore wall.


**NOTE:** Any other removable surface coils and pads must be removed from the table.

**CAUTION:** All RF coil-related tests must be run on a well-calibrated MRI scanner that has passed all system tests, particularly white pixel and correlated noise.

### MCQA test procedure

1. Place the 16ch Shoulder Coil on the MRI scanner table as shown in Figure 1.

**Figure 1:** 16ch Shoulder Coil setup

2. Plug the 16ch Shoulder Coil into port P2 of the MRI scanner.
3. Place the Head TLT SiOil Phantom (GEHC PN# 2359877) on the divot of the patient contacting surface of the Posterior Array and place the flexible Anterior Array on top of the phantom as shown in Figure 2.


**Figure 2:** Phantom setup

**CAUTION:** To prevent inaccurate results, the 16ch Shoulder Coil's system cable must be routed away from the bore wall.

4. Position the 16ch Shoulder Coil on the MRI scanner table such that MRI scanner landmark aligns with the coil's landmark lines as shown in Figure 3.

**Figure 3:** Landmark alignment

5. Route the 16ch Shoulder Coil system cable down the edge of the MRI scanner table; do not loop the coil's system cable.

**NOTE:** Improper routing of the coil's system cable could result in inaccurate results.

6. From the MR Service desktop window, Image quality tab, launch the Multi-Coil QA Tool.
7. In the Multi-Coil QA Tool:
  - a. Populate the required fields, such as Coil Serial #.
  - b. Press start.
  - c. Wait for the test to complete.
8. Confirm MCQA test status; refer to Figure 4.

**NOTE:** Values shown in Figure 4 are for reference only and may not match the actual values observed.


Figure 4: Example MCQA user interface

9. If PASS is displayed, all coil elements are functioning properly.
10. If FAIL is displayed, or if any issues were encountered with the 16ch Shoulder Coil while performing MCQA, discontinue use of the coil and contact your service representative.

---

THIS PAGE LEFT  
INTENTIONALLY BLANK.

# 3 Replacement Parts

## FRUs and Replacement Components

For U.S. and Canada, replacement parts marked with a star (\*) are available through NeoCoil at Tel: 1.888.MRI.COIL (674-2645), Email: [contactus@neocoil.com](mailto:contactus@neocoil.com) or the Web: [www.neocoil.com](http://www.neocoil.com).

Outside of the U.S. and Canada and for all other parts contact your GEHC.

Description	NeoCoil Part #	GE Part #
3.0T 16ch Shoulder Coil	NC143002	5849000-2
Pad Kit Includes: <ul style="list-style-type: none"><li>• 16ch Shoulder Posterior Housing Pad</li><li>• 16ch Shoulder Pad</li></ul>	NC143005	5849000-7
GE Replacement Straps ( <i>available only from GE</i> ) <ul style="list-style-type: none"><li>• 14" W x 46" Long (Gray) - Accessory Order # E8802MC</li><li>• 6" W x 46" Long (Gray) - Accessory Order # E8802MD</li></ul>	N/A	46-306594P1 46-306593P1
SiOil Spherical Phantom	N/A	2359877

---

## Disposal of Components

When any part of this product is to be disposed of, it is necessary to consider the possible impact on the environment and to observe all applicable laws. Consult your GE service representative for instructions concerning proper disposal.

---

# 4 Troubleshooting

## General Notes

When encountering a problem with the Shoulder Coil, attempt these troubleshooting tips.

Problem	Possible Cause	Solution
Unable to Scan with Coil	Other coils are plugged in	Unplug all other surface coils from the scanner.
Bias Fault	Bent pin	Visually inspect the system cable for a bent pin. If there is a bent pin, return the coil for service.
Coil Not Recognized by Scanner	Bent pin	Visually inspect the system cable for a bent pin. If there is a bent pin, return the coil for service.

---

# 5 Contact Information

## **GE Healthcare Americas**

USA 800-558-5102

Canada 800-668-0732

**GE Healthcare Europe** (33) 1-41-19-76-76

## **GE Healthcare Asia/Australia**

China 86-21-62192228

Taiwan 886-2-2505-7900

Singapore 65-291-8528

Australia 61-2-9975-5501

Japan 81-120-48-2630

Korea 82-31-740-6119

India 91-80-845-2923

**GE Healthcare Website** [www.gehealthcare.com](http://www.gehealthcare.com)

**NeoCoil, LLC** 1-888-674-2645

**NeoCoil Website** [www.neocoil.com](http://www.neocoil.com)

## **Manufacturer's name and address**

NeoCoil, LLC

N27 W23910A Paul Rd.

Pewaukee, WI 53072 USA