



GE Medical Systems

Technical Publications

Vivid q N

C € 0123

Reference Manual

GEMS #: FQ092028

Rev. 02

Copyright© 2015 By General Electric Co.
Reference Documentation



Regulatory Requirements

This product complies with regulatory requirements of the following European Directive 93/42/EEC concerning medical devices.



This manual is a reference for Vivid q N. It applies to all versions of the 12.x software for the Vivid q N ultrasound system.



GE Medical Systems

MANUAL STATUS
FQ092028 Rev. 02
12 August 2012

COMPANY DATA

GE VINGMED ULTRASOUND A/S
Strandpromenaden 45
N-3191 Horten, Norway
Tel.: (+47) 3302 1100 Fax: (+47) 3302 1350

Revision History

| Revision | Date | Reason for Change |
|----------|----------------|---------------------------------|
| 1 | 12 August 2012 | Initial release |
| 2 | 10 July 2015 | Change CE Notified Body to 0123 |

List of affected pages

| Revision | Affected Pages |
|----------|----------------|
| 1 | All |
| 2 | All |

Please verify that you are using the latest revision of this document. Information pertaining to this document is maintained on ePDM (GE Medical Systems electronic Product Data Management). If you need to know the latest revision, contact your distributor, local GE Sales Representative or in the USA call the GE Ultrasound Clinical Answer Center at 1 800 682 5327 or 1 262 524 5698.

Revision History

Regulatory information



CAUTION *Any changes to accessories, peripheral units or any other part of the system must be approved by the manufacturer. Ignoring this advice may compromise the regulatory approvals obtained for the product.*

Directives

The GE Healthcare Ultrasound product families are tested to meet all applicable requirements in relevant EU Directives and European/International standards.

- Council Directive 93/42/EEC concerning MDD (Medical Devices Directive): the CE label affixed to the product testifies compliance to this Directive.

The location of the CE marking is shown in the User Manual, Direction 5400907-1000, Chapter 1 - "Safety" as specified on the "Device Labels".

- European registered place of business:

GE Medical Systems Information Technologies GmbH,
Munzinger Strasse 5 D-79111 Freiburg, Germany
Tel: (+49) 761 45 43 0 Fax: (+49) 76145 43 233

Product Classifications

The Vivid q N ultrasound unit confirms to the following classifications:

- According to 93/42/EEC Medical Device Directive, this is Class IIa Medical Device.
- According to IEC/EN 60601-1, Equipment is Class I, Type B with BF or CF Applied Parts.
- According to CISPR 11, this is Group 1, Class B ISM Equipment.
- Classification according to the degree of protection against ingress of water as detailed in the current edition of IEC 60529 (section 6.1.1).
- The system is non-protected (IPX0).
- The footswitch rate IPX8 is suitable for use in surgical rooms.
- IPX1 for external use transducers
- IPX7 for internal use transducers
- IPX8 for ICE transducers

Conformity to Standards

To fulfill the requirements of relevant EC directives and/or European Harmonized/ International standards, the following documents/standards have been used:

- International Electrotechnical Commission (IEC).
 - IEC/EN 60601-1: Medical Electrical Equipment, Part 1 General Requirements for Safety
 - IEC/EN 60601-1-2: Medical electrical equipment - Part 1-2: General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests
 - IEC/EN 60601-1-4: Medical electrical equipment - Part 1-4: General requirements for safety - Collateral standard: Programmable electrical medical systems
 - IEC/EN 60601-1-6: Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability
 - IEC/EN 60601-2-37: Medical Electrical Equipment Part 2-37: Particular Requirements for the Safety of Ultrasonic Medical Diagnostic and Monitoring Equipment
 - IEC/EN 60601-1-1: Medical Electrical Equipment Part 1-1: General Requirements for Safety Collateral Standard: Safety Requirements for Medical Electrical Systems
- IEC/EN 62304: Medical device software - Software life-cycle processes
- Canadian Standards Association (CSA).
 - CSA 22.2, 601.1 Medical Electrical Equipment, Part 1 General Requirements for Safety.
- NEMA Standards Publication UD 2: "Acoustic Output Measurement Standard For Diagnostic Ultrasound Equipment".
- NEMA Standards Publication UD 3: "Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices On Diagnostic Ultrasound Equipment".
- Medical Device Good Manufacturing Practice Manual issued by the FDA (Food and Drug Administration, Department of Health, USA).

Certifications

- Quality management standards for medical devices: General Electric Medical Systems is ISO 9001 and ISO13485 certified.

Software License Acknowledgements

- WindowBlinds™ OCX © Stardock®

Table of Contents

Revision History

| | |
|------------------------------------|----------|
| List of affected pages..... | i |
| | ii |

Regulatory information

| | |
|--|------------|
| Directives | iii |
| Product Classifications..... | iii |
| Conformity to Standards | iv |
| Certifications..... | v |
| Software License Acknowledgements | v |

Introduction

| | |
|----------------------------------|----------|
| Contact information | 2 |
|----------------------------------|----------|

Chapter 1

Measurements

| | |
|---------------------------------------|-----------|
| Measurement overview | 6 |
| Cardiac measurements..... | 6 |
| AFI Measurements | 19 |
| Auto-EF measurements..... | 22 |
| Measurement formulas | 24 |
| Formulas–Cardiac | 24 |
| Formulas–Generic | 55 |
| Formulas–Vascular..... | 57 |
| Formulas–OB..... | 58 |
| Formulas–GYN..... | 62 |
| Measurement accuracy | 64 |
| General | 64 |
| Sources of error | 64 |
| Optimizing Measurement Accuracy | 66 |
| Measurement Uncertainties..... | 66 |
| DICOM SR Measurements | 69 |

Chapter 2

OB Tables

| | |
|-----------------------------|------------|
| ASUM..... | 72 |
| Berkowitz | 74 |
| Brenner | 75 |
| Campbell | 75 |
| Eriksen | 76 |
| Goldstein..... | 77 |
| Hadlock | 78 |
| Hansmann..... | 85 |
| Hellman | 94 |
| Hill..... | 94 |
| Hohler..... | 95 |
| Jeanty..... | 95 |
| JSUM | 108 |
| Kurtz..... | 112 |
| Mayden..... | 113 |
| Mercer | 114 |
| Merz..... | 115 |
| Moore | 125 |
| Nelson | 125 |
| Osaka | 126 |
| Paris | 130 |
| Rempen | 133 |
| Robinson..... | 138 |
| Tokyo..... | 138 |
| Tokyo Shinozuka..... | 142 |
| Williams..... | 149 |
| Yarkoni..... | 149 |

Chapter 3

Acoustic information

| | |
|--|------------|
| The real-time display of acoustic output indices..... | 153 |
|--|------------|

| | |
|--|------------|
| Thermal Index..... | 153 |
| Mechanical Index:..... | 154 |
| Concerns Surrounding the Use of Diagnostic Ultrasound.... | 155 |
| Default Settings and Output Levels..... | 155 |
| Controls Affecting Acoustic Output | 156 |
| Track 3 Summary Table | 158 |
| Acoustic Parameters as Measured in Water..... | 161 |
| Definitions, symbols and abbreviations | 161 |
| Translations of definitions, symbols, and abbreviations ... | 162 |
| Acoustic Output Reporting Tables | |
| for Track 3/IEC 60601-2-37..... | 186 |
| Explanation of Footnotes | 186 |
| Multiple focal-zones | 186 |
| Operating Conditions | 187 |
| Transducer Model: 3Sc-RS | 188 |
| Transducer Model: M4S-RS | 194 |
| Transducer Model: 5S-RS | 200 |
| Transducer Model: 6S-RS | 206 |
| Transducer Model: 7S-RS | 212 |
| Transducer Model: 10S-RS | 218 |
| Transducer Model: 12S-RS | 224 |
| Transducer Model: e8C-RS | 230 |
| Transducer Model: 3C-RS | 235 |
| Transducer Model: 4C-RS | 240 |
| Transducer Model: 8C-RS | 245 |
| Transducer Model: 8L-RS..... | 250 |
| Transducer Model: 9L-RS..... | 255 |
| Transducer Model: 12L-RS..... | 260 |
| Transducer Model: i12L-RS..... | 265 |
| Transducer Model: 6T-RS | 270 |
| Transducer Model: 6Tc-RS..... | 276 |
| Transducer Model: 9T-RS | 282 |
| Transducer Model: P2D-RS..... | 288 |
| Transducer Model: P6D-RS..... | 289 |
| Transducer Model: AcuNav8 | 290 |
| Transducer Model: AcuNav™ 10/SoundStar™ 3D 10FG / eco 10FG..... | 296 |

Chapter 4

Electromagnetic Compatibility

Guidance and manufacturer's declaration303

Introduction

The Vivid q N ultrasound unit is a compact, high performance portable digital ultrasound imaging system.

The system provides image generation in 2D (B Mode), Color Flow Mapping (CFM), Power Doppler (Angio), M-Mode, Color M-Mode (CMM), PW and CW Doppler spectra, Tissue Velocity Imaging / Tissue Tracking (TVI/TT), Tissue Synchronization Imaging (TSI), B Flow Imaging (BFI/Bflow) and LVO Contrast option applications.

The fully digital architecture of the Vivid q N unit allows optimal usage of all scanning modes and probe types, throughout the full spectrum of operating frequencies.

Contact information

If additional information or assistance is needed, please contact the local distributor or the appropriate support resource listed below:

| | |
|---|--|
| Europe GE Ultraschall KG Deutschland GmbH & Co. Beethovenstraße 239 Postfach 11 05 60 D-42655 Solingen | Tel: 0130 81 6370 Tel: (49)(0) 212-28-02-208 |
| USA GE Medical Systems Ultrasound Service Engineering 4855 W. Electric Avenue Milwaukee, WI 53219 On-line Applications Support | Tel: (1) 800-437-1171 Fax: (1) 414-647-4090 Tel: (1) 800-682-5327 or (262) 524-5698 |
| Canada GE Medical Systems Ultrasound Service Engineering 4855 W. Electric Avenue Milwaukee, WI 53219 On-line Applications Support | Tel: (1) 800-664-0732 Tel: (1) 800-682-5327 or (262) 524-5698 |
| Asia GE Ultrasound Asia Service Department Ultrasound 298 Tiong Gahru Road # 15-01/06 Central Plaza Singapore 168730 | Tel: (65) 291-8528 Fax: (65) 272-3997 |

| | |
|--|---|
| Latin and South America GE Medical Systems Ultrasound Service Engineering 4855 W. Electric Avenue Milwaukee, WI 53219 On-line Applications Support | Tel: (1) 305-735-2304 Tel: (1) 800-682-5327 or (262) 524-5698 |
| Brazil GE Ultrasound Rua Tomas Carvalhal, 711 Paraiso Cep: 04006-002 - São Paulo, SP | Tel: (55.11) 887-8099 Fax: (55.11) 887-9948 |

Chapter 1

Measurements

This chapter includes the following information:

- **Measurement overview 6**
 - Cardiac measurements 6
 - AFI Measurements 19
 - Auto-EF measurements 22
- **Measurement formulas..... 24**
 - Formulas–Cardiac 24
 - Formulas–Generic 55
 - Formulas–Vascular 57
 - Formulas–OB 58
 - Formulas–GYN 62
- **Measurement accuracy 64**
 - General 64
 - Sources of error 64
 - Optimizing Measurement Accuracy 66
 - Measurement Uncertainties 66
- **DICOM SR Measurements 69**

Measurements

Measurement overview

The following table shows the cardiac measurements available on the Vivid q N ultrasound unit.

Cardiac measurements

| Abbreviation | Definition | Unit |
|--------------|---|-----------------|
| %FS | LV Fractional Shortening, 2D | % |
| %FS | LV Fractional Shortening, M-mode | % |
| %IVS Thck | IVS Fractional Shortening, 2D | % |
| %IVS Thck | IVS Fractional Shortening, M-mode | % |
| %LVPW Thck | LV Posterior Wall Fractional Shortening, 2D | % |
| %LVPW Thck | LV Posterior Wall Fractional Shortening, M-mode | % |
| Ao Abd AP | Abdominal Aorta AP diameter, 2D | cm |
| Ao Arch Diam | Aortic Arch Diameter | cm |
| Ao asc | Ascending Aortic Diameter | cm |
| Ao Desc Diam | Descending Aortic Diameter | cm |
| Ao Isthmus | Aortic Isthmus | cm |
| Ao Root Diam | Aortic Root Diameter | cm |
| Ao Root Diam | Aortic Root Diameter, M-mode | cm |
| Ao Sinus | Aortic Root at the Sinus, 2D | cm |
| AR ERO | PISA: Regurgitant Orifice Area | cm ² |
| AR Flow | PISA: Regurgitant Flow | ml/s |
| AR PHT | AV Insuf. Pressure Half Time | ms |
| AR Rad | PISA: Radius of Aliased Point | cm |
| AR RV | PISA: Regurgitant Volume Flow | ml |
| AR Vel | PISA: Aliased Velocity | m/s |
| AR Vmax | Aortic Insuf. Peak Velocity | m/s |
| AR VTI | Aortic Insuf. Velocity Time Integral | cm |

| Abbreviation | Definition | Unit |
|-----------------|--|---------------------------------|
| ARed max PG | Aortic Insuf. End-Diastole Pressure Gradient | mm Hg |
| ARed Vmax | Aortic Insuf. End-Diastolic Velocity | m/s |
| AV Acc Slope | Aortic Valve Flow Acceleration | m/s ² |
| AV Acc Time | Aortic Valve Acceleration Time | ms |
| AV AccT/ET | AV Acceleration to Ejection Time Ratio | |
| AV CO | Cardiac Output by Aortic Flow | l/min |
| AV Cusp | Aortic Valve Cusp Separation, 2D | cm |
| AV Cusp | Aortic Valve Cusp Separation, M-mode | cm |
| AV Dec Time | Aortic Valve Deceleration Time | ms |
| AV Diam | Aortic Diameter, 2D | cm |
| AV max PG | Aortic Valve Peak Pressure Gradient | mm Hg |
| AV mean PG | Aortic Valve Mean Pressure Gradient | mm Hg |
| AV SV | Stroke Volume by Aortic Flow | ml |
| AV Vmax | Aortic Valve Peak Velocity | m/s |
| AV Vmean | AV Mean Velocity | m/s |
| AV VTI | Aortic Valve Velocity Time Integral | cm |
| AVA (Vmax) | AV Area by Continuity Equation by Peak V | cm ² |
| AVA (VTI) | AV Area by Continuity Equation VTI | cm ² |
| AVA Planimetry | Aortic Valve Area | cm ² |
| AVAI Planimetry | Aortic Valve Area Index | cm ² /m ² |
| AVET | Aortic Valve Ejection Time | ms |
| AVET | Aortic Valve Ejection Time, M-mode | ms |
| CO (A-L A2C) | CO 2CH, Single Plane, Area-Length | l/min |
| CO (A-L A4C) | CO 4CH, Single Plane, Area-Length | l/min |
| CO (Biplane) | CO, Bi-Plane, MOD | l/min |
| CO (bullet) | CO, Bi-Plane, Bullet | l/min |

Measurements

| Abbreviation | Definition | Unit |
|---------------|--|-------|
| CO (MOD A2C) | CO 2CH, Single Plane, MOD(Simpson) | l/min |
| CO (MOD A4C) | CO 4CH, Single Plane, 4CH, MOD(Simpson) | l/min |
| CO(Cube) | Cardiac Output, 2D, Cubic | l/min |
| CO(Cube) | Cardiac Output, M-mode, Cubic | l/min |
| CO(Teich) | Cardiac Output, 2D, Teicholtz | l/min |
| CO(Teich) | Cardiac Output, M-mode, Teicholtz | l/min |
| D-E Excursion | MV Anterior Leaflet Excursion | cm |
| D-E Excursion | Mitral Valve D-E Slope | cm |
| EDV (bullet) | LV Volume, Diastolic, Bi-Plane, Bullet | ml |
| EDV(Cube) | Left Ventricle Volume, Diastolic, 2D, Cubic | ml |
| EDV(Cube) | Left Ventricle Volume, Diastolic, M-mode, Cubic | ml |
| EDV(Teich) | Left Ventricle Volume, Diastolic, 2D, Teicholz | ml |
| EDV(Teich) | Left Ventricle Volume, Diastolic, M-mode, Teicholz | ml |
| EF (A-L A2C) | Ejection Fraction 2CH, Single Plane, Area-Length | % |
| EF (A-L A4C) | Ejection Fraction 4CH, Single Plane, Area-Length | % |
| EF (Biplane) | Ejection Fraction, Bi-Plane, MOD | % |
| EF (bullet) | Ejection Fraction 2CH, Bi-Plane, Bullet | % |
| EF (MOD A2C) | Ejection Fraction 2CH, Single Plane, MOD(Simpson) | % |
| EF (MOD A4C) | Ejection Fraction 4CH, Single Plane, 4CH, MOD(Simpson) | % |
| E-F Slope | Mitral Valve E-F Slope | m/s |
| EF(Cube) | Ejection Fraction, 2D, Cubic | % |
| EF(Cube) | Ejection Fraction, M-mode, Cubic | % |
| EF(Teich) | Ejection Fraction, 2D, Teicholtz | % |
| EF(Teich) | Ejection Fraction, M-mode, Teicholtz | % |
| EPSS | E-Point-to-Septum Separation, M-mode | cm |

| Abbreviation | Definition | Unit |
|--------------|---|-----------------|
| EPSS 2D | E-Point-to-Septum Separation, 2D | cm |
| ERO | Effective Regurgitant Orifice | cm ² |
| ESV (bullet) | LV Volume, Systolic, Bi-Plane, Bullet | ml |
| ESV(Cube) | Left Ventricle Volume, Systolic, 2D, Cubic | ml |
| ESV(Cube) | Left Ventricle Volume, Systolic, M-mode, Cubic | ml |
| ESV(Teich) | Left Ventricle Volume, Systolic, 2D, Teicholz | ml |
| ESV(Teich) | Left Ventricle Volume, Systolic, M-mode, Teicholz | ml |
| HR | AV Heart Rate, Dop | BPM |
| HR | Heart Rate, 2D, Teicholtz | bpm |
| HR | Heart Rate for 2CH study | bpm |
| HR | Heart Rate for 4CH study | bpm |
| HR | Heart Rate for 2CH AL study | bpm |
| HR | Heart Rate for 2CH MOD study | bpm |
| HR | Heart Rate for 4CH AL study | bpm |
| HR | Heart Rate for 4CH MOD study | bpm |
| HR | Heart Rate for Bullet study | bpm |
| HR | Heart Rate for BiPlane MOD study | bpm |
| HR | LV Heart Rate, Dop | bpm |
| HR | Heart Rate, M-mode, Teicholtz | bpm |
| HR | Heart Rate | bpm |
| IVC | Inferior Vena Cava | cm |
| IVCT | Isovolumic Contraction Time | ms |
| IVRT | Isovolumic Relaxation Time | ms |
| IVSd | Interventricular Septum Thickness, Diastolic, 2D | cm |
| IVSd | IVS Thickness, Diastolic, M-mode | cm |
| IVSs | Interventricular Septum Thickness, Systolic, 2D | cm |

Measurements

| Abbreviation | Definition | Unit |
|-----------------|--|------------------|
| IVSs | IVS Thickness, Systolic, M-mode | cm |
| LA Diam | Left Atrium Diameter, 2D | cm |
| LA Diam | Left Atrium Diameter, M-mode | cm |
| LA Diam | Right Atrium Diameter, 2D | cm |
| LA Major | Left Atrium Major | cm |
| LA Minor | Left Atrium Minor | cm |
| LA/Ao | LA Diameter to AoRoot Diameter Ratio, 2D | |
| LA/Ao | LA Diameter to AoRoot Diameter Ratio, M-mode | |
| LAEDV (MOD A4C) | LA Volume, Diastolic Single Plane, MOD | ml |
| LAEDV (MOD BP) | LA Volume, Diastolic, Bi-Plane, MOD | ml |
| LAESV (MOD A4C) | LA Volume, Systolic, Single Plane, MOD | ml |
| LAESV (MOD BP) | LA Volume, Systolic, Bi-Plane, MOD | ml |
| LIMP | Left Index of Mysocardial Performance | |
| LVA (s) | Left Ventricular Area, Systolic, 2CH | cm ² |
| LVAd (A2C) | Left Ventricular Area, Diastolic, 2CH | cm ² |
| LVAd (A4C) | Left Ventricular Area, Diastolic, 4CH | cm ² |
| LVAd(sax) | LV area, SAX, Diastolic | cm ² |
| LVAend (d) | LV Endocardial Area, SAX | cm ² |
| LVAepi (d) | LV Epicardial Area, SAX | cm ² |
| LVAs (A4C) | Left Ventricular Area, Systolic, 4CH | cm ² |
| LVAs(sax) | LV area, SAX, Systolic | cm ² |
| LVd Mass | LV Mass, Diastolic, 2D | g |
| LVd Mass | LV Mass, Diastolic, M-mode | g |
| LVd Mass Index | LV Mass Index, Diastolic, 2D | g/m ² |
| LVd Mass Index | LV Mass Index, Diastolic, M-mode | g/m ² |

| Abbreviation | Definition | Unit |
|-----------------|--|-------------------|
| LVEDV (A-L A2C) | LV Volume, Diastolic, 2CH, Area-Length | ml |
| LVEDV (A-L A4C) | LV Volume, Diastolic, 4CH, Area-Length | ml |
| LVEDV (MOD A2C) | LV Volume, Diastolic, Single Plane, 2CH, MOD | ml |
| LVEDV (MOD A4C) | LV Volume, Diastolic, Single Plane, 4CH, MOD | ml |
| LVEDV (MOD BP) | LV Volume, Diastolic, Bi-Plane, MOD | ml |
| LVESV (A-L A2C) | LV Volume, Systolic, 2CH, Area-Length | ml |
| LVESV (A-L A4C) | LV Volume, Systolic, 4CH, Area-Length | ml |
| LVESV (MOD A2C) | LV Volume, Systolic, Single Plane, 2CH, MOD | ml |
| LVESV (MOD A4C) | LV Volume, Systolic, Single Plane, 4CH, MOD | ml |
| LVESV (MOD BP) | LV Volume, Systolic, Bi-Plane, MOD | ml |
| LVESV (MOD LAX) | LV Volume, Diastolic, Apical View, LAX, MOD | ml |
| LVESV (MOD LAX) | LV Volume, Systolic, Apical View, LAX, MOD | ml |
| LVET | Left Ventricle Ejection Time | ms |
| LVIDd | LV Internal Dimension, Diastolic, 2D | cm |
| LVIDd | LV Internal Dimension, Diastolic, M-mode | cm |
| LVIDs | LV Internal Dimension, Systolic, 2D | cm |
| LVIDs | LV Internal Dimension, Systolic, M-mode | cm |
| LVLd (apical) | Left Ventricular Length, Diastolic, 2D | cm |
| LVLs (apical) | Left Ventricular Length, Systolic, 2D | cm |
| LVOT Area | Left Ventricle Outflow Tract Area | cm ² |
| LVOT CO | Cardiac Output by Aortic Flow | l/min |
| LVOT Diam | Left Ventricular Outflow Tract Diameter | cm |
| LVOT max PG | LVOT Peak Pressure Gradient | mm Hg |
| LVOT mean PG | LVOT Mean Pressure Gradient | mm Hg |
| LVOT SI | Stroke Volume Index by Aortic Flow | ml/m ² |
| LVOT SV | Stroke Volume by Aortic Flow | ml |

Measurements

| Abbreviation | Definition | Unit |
|------------------------|--|------------------|
| LVOT Vmax | LVOT Peak Velocity | m/s |
| LVOT Vmean | LVOT Mean Velocity | m/s |
| LVOT VTI | LVOT Velocity Time Integral | cm |
| LVPWd | Left Ventricular Posterior Wall Thickness, Diastolic, 2D | cm |
| LVPWd | Left Ventricular Posterior Wall Thickness, Diastolic, M-mode | cm |
| LVPWs | Left Ventricular Posterior Wall Thickness, Systolic, 2D | cm |
| LVPWs | Left Ventricular Posterior Wall Thickness, Systolic, M-mode | cm |
| LVs Mass | LV Mass, Systolic, 2D | g |
| LVs Mass | LV Mass, Systolic, M-mode | g |
| LVs Mass Index | LV Mass Index, Systolic, 2D | g/m ² |
| LVs Mass Index | LV Mass Index, Systolic, M-mode | g/m ² |
| LAA _d (A2C) | Left Atrium Area, Apical 2C | cm ² |
| LAA _d (A4C) | Left Atrium Area, Apical 4C | cm ² |
| MCO | Mitral Valve closure to Opening | ms |
| MP Area | Mitral Valve Prosthesis | cm ² |
| MR Acc Time | MV Regurg. Flow Acceleration | s |
| MR ERO | PISA: Regurgitant Orifice Area | cm ² |
| MR Flow | PISA: Regurgitant Flow | ml/s |
| MR max PG | Mitral Regurg. Peak Pressure Gradient | mm Hg |
| MR Rad | PISA: Radius of Aliased Point | cm |
| MR RV | PISA: Regurgitant Volume Flow | ml |
| MR Vel | PISA: Aliased Velocity | m/s |
| MR Vmax | Mitral Regurg. Peak Velocity | m/s |

| Abbreviation | Definition | Unit |
|-----------------|---------------------------------------|-------------------|
| MR Vmax | PISA: CW Peak Velocity | m/s |
| MR Vmean | Mitral Regurg. Mean Velocity | m/s |
| MR VTI | Mitral Regurg. Velocity Time Integral | cm |
| MR VTI | PISA: CW Velocity Time Integral | cm |
| MV A Dur | Mitral Valve A-Wave Duration | ms |
| MV A Velocity | MV Velocity Peak A | m/s |
| MV Acc Slope | Mitral Valve Flow Acceleration | m/s ² |
| MV Acc Time | Mitral Valve Acceleration Time | ms |
| MV Acc/Dec Time | MV: Acc.Time/Decel.Time Ratio | |
| MV an diam | Mitral Valve Annulus Diameter, 2D | cm |
| MV CO | Cardiac Output by Mitral Flow | l/min |
| MV Dec Slope | Mitral Valve Flow Deceleration | m/s ² |
| MV Dec Time | Mitral Valve Deceleration Time | ms |
| MV E Velocity | MV Velocity Peak E | m/s |
| MV E/A Ratio | Mitral Valve E-Peak to A-Peak Ratio | |
| MV max PG | Mitral Valve Peak Pressure Gradient | mm Hg |
| MV mean PG | Mitral Valve Mean Pressure Gradient | mm Hg |
| MV PHT | Mitral Valve Pressure Half Time | ms |
| MV SI | Stroke Volume Index by Mitral Flow | ml/m ² |
| MV SV | Stroke Volume by Mitral Flow | ml |
| MV Time to Peak | Mitral Valve Time to Peak | ms |
| MV Vmax | Mitral Valve Peak Velocity | m/s |
| MV Vmean | MV Mean Velocity | m/s |
| MV VTI | Mitral Valve Velocity Time Integral | cm |
| MVA | Mitral Valve Area | cm ² |
| MVA By PHT | Mitral Valve Area according to PHT | cm ² |

Measurements

| Abbreviation | Definition | Unit |
|----------------------|--|------------------|
| MVA by plan | Mitral Valve Area, 2D | cm ² |
| MVET | Mitral Valve Ejection Time | ms |
| P Vein A | Pulmonary Vein Velocity Peak A (reverse) | m/s |
| P Vein A Dur | Pulmonary Vein A-Wave Duration | ms |
| P Vein D | Pulmonary Vein End-Diastolic Peak Velocity | m/s |
| P Vein S | Pulmonary Vein Systolic Peak Velocity | m/s |
| PAEDP | Pulmonary Artery Diastolic Pressure | mm Hg |
| PE(d) | Pericard Effusion, M-mode | cm |
| PEs | Pericard Effusion, 2D | cm |
| PR max PG | Pulmonic Insuf. Peak Pressure Gradient | mm Hg |
| PR mean PG | Pulmonic Insuf. Mean Pressure Gradient | mm Hg |
| PR PHT | Pulmonic Insuf. Pressure Half Time | ms |
| PR Vmax | Pulmonic Insuf. Peak Velocity | m/s |
| PR VTI | Pulmonic Insuf. Velocity Time Integral | cm |
| PRend max PG | Pulmonic Insuf. End-Diastole Pressure Gradient | mm Hg |
| PRend Vmax | Pulmonic Insuf. End-Diastolic Velocity | m/s |
| Pulmonic Diam | Pulmonary Artery Diameter, 2D | cm |
| PV Acc Slope | Pulmonic Valve Flow Acceleration | m/s ² |
| PV Acc Time | Pulmonic Valve Acceleration Time | ms |
| PV Acc Time/ET Ratio | PV Acceleration to Ejection Time Ratio | |
| PV an diam | Pulmonic Valve Annulus Diameter, 2D | cm |
| PV Ann Area | Pulmonic Valve Area | cm ² |
| PV CO | Cardiac Output by Pulmonic Flow | l/min |
| PV CO | Cardiac Output by Pulmonic Flow | l/min |
| PV max PG | Pulmonic Valve Peak Pressure Gradient | mm Hg |
| PV mean PG | Pulmonic Valve Mean Pressure Gradient | mm Hg |

| Abbreviation | Definition | Unit |
|-----------------|---|-----------------|
| PV SV | Stroke Volume by Pulmonic Flow | ml |
| PV Vmax | Pulmonary Artery Peak Velocity | m/s |
| PV Vmax | Pulmonic Valve Peak Velocity | m/s |
| PV Vmean | PV Mean Velocity | m/s |
| PV VTI | Pulmonic Valve Velocity Time Integral | cm |
| PVA (VTI) | Pulmonary Artery Velocity Time Integral | cm ² |
| PVein S/D Ratio | Pulmonary Vein SD Ratio | |
| PVET | Pulmonic Valve Ejection Time | ms |
| PVPEP | Pulmonic Valve Pre-Ejection Period | ms |
| PVPEP/ET Ratio | PV Pre-Ejection to Ejection Time Ratio | |
| Qp/Qs | Pulmonic-to-Systemic Flow Ratio | |
| RA Major | Right Atrium Major, 2D | cm |
| RA Minor | Right Atrium Minor, 2D | cm |
| RAEDV A2C | Right Atrium End Diastolic Volume, Apical 2 chamber | cm ³ |
| RAEDV A-L | RA End Diastolic Volume (A-L) | ml |
| RAEDV MOD | RA Volume Diastolic, Single Plan, MOD | ml |
| RAEDV MOD | RA End Diastolic Volume (MOD) | ml |
| RAESV A-L | RA End Systole Volume (A-L) | ml |
| RAESV MOD | RA Volume, Systolic, Single Plane, MOD | ml |
| RAESV MOD | RA End Systole Volume (MOD) | ml |
| RALd | Right Atrium Length, Diastole | cm |
| RALs | RA Length, systole | cm |
| RIMP | Right Index of Myocardial Performance | |
| RJA (A4C) | Regurgitant jet area | cm ² |
| RJA/LAA | Regurgitant jet area ratio RJA/LAA | |
| RV Major | Right Ventricle Major | cm |

Measurements

| Abbreviation | Definition | Unit |
|--------------|---|-------------------|
| RV Minor | Right Ventricle Minor | cm |
| RVAWd | Right Ventricle Wall Thickness, Diastolic, 2D | cm |
| RVAWs | Right Ventricle Wall Thickness, Systolic, 2D | cm |
| RVET | Right Ventricle Ejection Time | s |
| RVIDd | Right Ventricle Diameter, Diastolic, 2D | cm |
| RVIDd | Right Ventricle Diameter, Diastolic, M-mode | cm |
| RVIDs | Right Ventricle Diameter, Systolic, 2D | cm |
| RVIDs | Right Ventricle Diameter, Systolic, M-mode | cm |
| RVOT Area | Right Ventricle Outflow Tract Area | cm ² |
| RVOT Diam | RV Output Tract Diameter, 2D | cm |
| RVOT Diam | RV Output Tract Diameter, M-Mode | cm |
| RVOT max PG | RVOT Peak Pressure Gradient | mm Hg |
| RVOT meanPG | RVOT Mean Pressure Gradient | mm Hg |
| RVOT SI | LV Stroke Volume Index by Pulmonic Flow | ml/m ² |
| RVOT SV | Stroke Volume by Pulmonic Flow | ml |
| RVOT Vmax | RVOT Peak Velocity | m/s |
| RVOT Vmean | RVOT Mean Velocity | m/s |
| RVOT VTI | RVOT Velocity Time Integral | cm |
| RVSP | Right Ventricle Systolic Pressure | mm Hg |
| RVWd | Right Ventricle Wall Thickness, Diastolic, M-mode | cm |
| RVWs | Right Ventricle Wall Thickness, Systolic, M-mode | cm |
| RAA (d) | Right Atrium Area, 2D, Diastole | cm ² |
| RAA (s) | Right Atrium Area, 2D, Systole | cm ² |
| SI (A-L A2C) | LV Stroke Index, Single Plane, 2CH, Area-Length | ml/m ² |
| SI (A-L A4C) | LV Stroke Index, Single Plane, 4CH, Area-Length | ml/m ² |
| SI (Biplane) | LV Stroke Index, Bi-Plane, MOD | ml/m ² |

| Abbreviation | Definition | Unit |
|---------------|---|-------------------|
| SI (bullet) | LV Stroke Index, Bi-Plane, Bullet | ml/m ² |
| SI (MOD A2C) | LV Stroke Index, Single Plane, 2CH, MOD | ml/m ² |
| SI (MOD A4C) | LV Stroke Index, Single Plane, 4CH, MOD | ml/m ² |
| SI (Teich) | LV Stroke Index, Teicholz, 2D | ml/m ² |
| SI (Teich) | LV Stroke Index, Teicholz, M-mode | ml/m ² |
| SV (A-L A2C) | LV Stroke Volume, Single Plane, 2CH, Area-Length | ml |
| SV (A-L A4C) | LV Stroke Volume, Single Plane, 4CH, Area-Length | ml |
| SV (Biplane) | LV Stroke Volume, Bi-Plane, MOD | ml |
| SV (bullet) | LV Stroke Volume, Bi-Plane, Bullet | ml |
| SV (MOD A2C) | LV Stroke Volume, Single Plane, 2CH, MOD(Simpson) | ml |
| SV (MOD A4C) | LV Stroke Volume, Single Plane, 4CH, MOD(Simpson) | ml |
| SV(Cube) | LV Stroke Volume, 2D, Cubic | ml |
| SV(Cube) | LV Stroke Volume, M-mode, Cubic | ml |
| SV(Teich) | LV Stroke Volume, 2D, Teicholtz | ml |
| SV(Teich) | LV Stroke Volume, M-mode, Teicholtz | ml |
| Systemic Diam | Systemic Vein Diameter, 2D | cm |
| Systemic Vmax | Systemic Vein Peak Velocity | m/s |
| Systemic VTI | Systemic Vein Velocity Time Integral | cm |
| TCO | Tricuspid Valve Closure to Opening | ms |
| TR max PG | Tricuspid Regurg. Peak Pressure Gradient | mm Hg |
| TR mean PG | Tricuspid Regurg. Mean Pressure Gradient | mm Hg |
| TR Vmax | Tricuspid Regurg. Peak Velocity | m/s |
| TR Vmean | Tricuspid Regurg. Mean Velocity | m/s |
| TR VTI | Tricuspid Regurgitation Velocity Time Integral | cm |

Measurements

| Abbreviation | Definition | Unit |
|---------------|--|------------------|
| TV A dur | Tricuspid Valve A-Wave Duration | ms |
| TV A Velocity | Tricuspid Valve A Velocity | m/s |
| TV Acc Time | Tricuspid Valve Time to Peak | ms |
| TV Ann Area | Tricuspid Valve Area | cm ² |
| TV ann diam | Tricuspid Valve Annulus Diameter, 2D | cm |
| TV Area | Tricuspid Valve Area, 2D | cm ² |
| TV CO | Cardiac Output by Tricuspid Flow | l/min |
| TV Dec Slope | Tricuspid Valve Flow Deceleration | m/s ² |
| TV E Velocity | Tricuspid Valve E Velocity | m/s |
| TV E/A Ratio | Tricuspid Valve E-Peak to A-Peak Ratio | |
| TV max PG | Tricuspid Valve Peak Pressure Gradient | mm Hg |
| TV mean PG | Tricuspid Valve Mean Pressure Gradient | mm Hg |
| TV mean PG | Tricuspid Valve Mean Pressure Gradient | mm Hg |
| TV PHT | Tricuspid Valve Pressure Half Time | ms |
| TV SV | Stroke Volume by Tricuspid Flow | ml |
| TV Vmean | TV Mean Velocity | m/s |
| TV VTI | Tricuspid Valve Velocity Time Integral | cm |
| VSD max PG | VSD Peak Pressure Gradient | mm Hg |
| VSD Vmax | VSD Peak Velocity | m/s |

AFI Measurements

| AFI shown parameters | The AFI shown parameters are shown only on screen | |
|----------------------|---|------|
| Abbreviation | Definition | Unit |
| PSS | Peak Systolic Longitudinal Strain Segmental most negative peak longitudinal strain during systole. Replaced with a positive systolic peak if positive peak exceeds $3/4 \times \text{abs}(\text{negative})$. Formula: $S = (L - L_0) / L_0$, where L is the instantaneous length of the segment and L_0 is the initial length at start systole (QRS). The length is measured along the ROI center line. | % |
| PSI | Post Systolic Index (based on longitudinal strain) Segmental, percentage of post-systolic contraction Formula: $PSI = 100 \times (PS - ESS) / PS$ where ESS = strain at AVC and PS = peak strain after AVC. | % |
| HR | HR Heart Rate Heart rate based on the time between left and right QRS markers in 2D Strain. Formula: $HR = 60 / \text{Cycle time}$ | BPM |

AFI parameters exported to the Worksheet

| Abbreviation | Definition | Unit |
|---|---|------|
| AVC | Aortic Valve Closure Time interval between QRS (start systole) and Aortic valve closure. | ms |
| A4C: | Apical 4 chamber view | |
| BS peak sys SL | Peak Systolic Longitudinal Strain - Basal Septum ¹ | % |
| MS peak sys SL | Peak Systolic Longitudinal Strain - Mid Septum ¹ | % |
| ¹ Segmental most negative peak longitudinal strain during systole. Replaced with a positive systolic peak if positive peak exceeds $3/4 \times \text{abs}(\text{negative})$. | | |

Measurements

| Abbreviation | Definition | Unit |
|---|---|------|
| AS peak sys SL | Peak Systolic Longitudinal Strain - Apical Septum+A19:C55 ¹ | % |
| BL peak sys SL | Peak Systolic Longitudinal Strain - Basal Lateral segment ¹ | % |
| ML peak sys SL | Peak Systolic Longitudinal Strain - Mid Lateral segment ¹ | % |
| AL peak sys SL | Peak Systolic Longitudinal Strain - Apical Lateral segment ¹ | % |
| A2C: | Apical 2 chamber view | |
| BI peak sys SL | Peak Systolic Longitudinal Strain - Basal Inferior segment ¹ | % |
| MI peak sys SL | Peak Systolic Longitudinal Strain - Mid Inferior segment ¹ | % |
| AI peak sys SL | Peak Systolic Longitudinal Strain - Apical Inferior segment ¹ | % |
| BA peak sys SL | Peak Systolic Longitudinal Strain - Basal Anterior segment ¹ | % |
| MA peak sys SL | Peak Systolic Longitudinal Strain - Mid Anterior segment ¹ | % |
| AA peak sys SL | Peak Systolic Longitudinal Strain - Apical Anterior segment ¹ | % |
| APLAX: | Apical long axis view | |
| BP peak sys SL | Peak Systolic Longitudinal Strain - Basal Posterior segment ¹ | % |
| MP peak sys SL | Peak Systolic Longitudinal Strain - Mid Posterior segment ¹ | % |
| AP peak sys SL | Peak Systolic Longitudinal Strain - Apical Posterior segment ¹ | % |
| ¹ Segmental most negative peak longitudinal strain during systole. Replaced with a positive systolic peak if positive peak exceeds 3/4*abs(negative). | | |

| Abbreviation | Definition | Unit |
|---|--|------|
| BAS peak sys SL | Peak Systolic Longitudinal Strain - Basal Antero-septal segment ¹ | % |
| MAS peak sys SL | Peak Systolic Longitudinal Strain - Mid Antero-septal segment ¹ | % |
| GLPS_LAX | Global longitudinal strain in apical long axis view Peak contraction of the entire myocardial view wall. Only calculated if at least 5 segments have acceptable TQ. Formula: $(\text{MinWallLength} - \text{MaxWallLength}) / \text{MaxWallLength}$ | % |
| GLPS_A4C | Global longitudinal strain in apical 4 chamber view Peak contraction of the entire myocardial view wall. Only calculated if at least 5 segments have acceptable TQ. Formula: $(\text{MinWallLength} - \text{MaxWallLength}) / \text{MaxWallLength}$ | % |
| GLPS_A2C | Global longitudinal strain in apical 2 chamber view Peak contraction of the entire myocardial view wall. Only calculated if at least 5 segments have acceptable TQ. Formula: $(\text{MinWallLength} - \text{MaxWallLength}) / \text{MaxWallLength}$ | % |
| GLPS_Avg | Average global longitudinal strain. Averaged between 3 views. Formula: $(\text{GLPS_LAX} + \text{GLPS_A4C} + \text{GLPS_A2C}) / 3$ | % |
| ¹ Segmental most negative peak longitudinal strain during systole. Replaced with a positive systolic peak if positive peak exceeds $3/4 * \text{abs}(\text{negative})$. | | |

Measurements

Auto-EF measurements

| AutoEF parameters exported to the Worksheet | | |
|---|--|-------|
| Abbreviation | Definition | Unit |
| Apical 4 chamber view | | |
| HR 4Ch Q | Heart rate, calculated from 4CH clip | BPM |
| LVVED 4Ch Q | LV volume, calculated using the modified Simpson method (Method of discs) from 4CH clip | ml |
| LVVES 4Ch Q | LV volume, calculated using the modified Simpson method (Method of discs) from 4CH clip | ml |
| LVEF 4Ch Q | Ejection fraction, calculated from 4CH clip Formula: $EF = 100 * (LVVED\ 4Ch\ Q - LVVES\ 4Ch\ Q) / LVVED\ 4Ch\ Q$ | % |
| LVSV 4Ch Q | LV Stroke Volume, calculated using the modified Simpson method (Method of discs) from 4CH clip | ml |
| LVCO 4Ch Q | LV Cardiac Output, calculated using the modified Simpson method (Method of discs) from 4CH clip | l/min |
| LVLs 4Ch Q | LV main axis size, calculated from 4CH clip | cm |
| LVLd 4Ch Q | LV main axis size, calculated from 4CH clip | cm |
| Apical 2 chamber view | | |
| HR 2Ch Q | Heart rate, calculated from 2Ch clip | BPM |
| LVVED 2Ch Q | LV volume, calculated using the modified Simpson method (Method of discs) from 2Ch clip | ml |
| LVVES 2Ch Q | LV volume, calculated using the modified Simpson method (Method of discs) from 2Ch clip | ml |
| LVEF 2Ch Q | Ejection fraction, calculated from 2Ch clip Formula: $EF = 100 * (LVVED\ 2Ch\ Q - LVVES\ 2Ch\ Q) / LVVED\ 2Ch\ Q$ | % |
| LVSV 2Ch Q | LV Stroke Volume, calculated using the modified Simpson method (Method of discs) from 2Ch clip | ml |
| LVCO 2Ch Q | LV Cardiac Output, calculated using the modified Simpson method (Method of discs) from 2Ch clip | l/min |
| LVLs 2Ch Q | LV main axis size, calculated from 2Ch clip | cm |

| AutoEF parameters exported to the Worksheet | | |
|---|---|-------|
| Abbreviation | Definition | Unit |
| LVLd 2Ch Q | LV main axis size, calculated from 2Ch clip | cm |
| Bi-plane | | |
| LVVED BiP Q | LV volume, calculated using the modified Simpson method (Method of discs) from 4Ch and 2Ch clips | ml |
| LVVES BiP Q | LV volume, calculated using the modified Simpson method (Method of discs) from 4Ch and 2Ch clips | ml |
| LVEF BiP Q | Ejection fraction, calculated from 4Ch and 2Ch clips Formula: $EF = 100 * (LVVED \text{ BiP Q} - LVVES \text{ BiP Q}) / LVVED \text{ BiP Q}$ | % |
| LVSV BiP Q | LV Stroke Volume, calculated using the modified Simpson method (Method of discs) from 4Ch and 2Ch clips | ml |
| LVCO BiP Q | LV Cardiac Output, calculated using the modified Simpson method (Method of discs) from 4Ch and 2Ch clips | l/min |

Measurement formulas

Formulas-Cardiac

The following table lists the cardiac calculations. The folders where to find the calculations and related measurements are indicated in brackets "[]".

| |
|--|
| %FS [Dimension, Cube/Teicholz] Mode: 2D:CF Formula: $((\{LVIDd\}-\{LVIDs\})/\{LVIDd\})$ Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz] Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC] |
| %FS [Generic, Dimension] Mode: MM:CM:AMM Formula: $((\{LVIDd\}-\{LVIDs\})/\{LVIDd\})$ Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension] Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER] |
| %IVS Thck [Dimension] Mode: 2D:CF Formula: $((\{IVSs\}-\{IVSd\})/\{IVSd\})$ Needs measurement: IVSs [Dimension], IVSd [Dimension] Measured by: LVs [2DLV], IVSs [2DCALIPER] |
| %IVS Thck [Dimension] Mode: MM:CM:AMM Formula: $((\{IVSs\}-\{IVSd\})/\{IVSd\})$ Needs measurement: IVSs [Dimension], IVSd [Dimension] Measured by: IVSs [MMDISCALIPER] |
| %LVPW Thck [Dimension] Mode: 2D:CF Formula: $((\{LVPWs\}-\{LVPWd\})/\{LVPWd\})$ Needs measurement: LVPWs [Dimension], LVPWd [Dimension] Measured by: LVs [2DLV], LVPWs [2DCALIPER] |
| %LVPW Thck [Dimension] Mode: MM:CM:AMM Formula: $((\{LVPWs\}-\{LVPWd\})/\{LVPWd\})$ Needs measurement: LVPWs [Dimension], LVPWd [Dimension] Measured by: LVPWs [MMDISCALIPER] |

| |
|--|
| Ao st junct/Ao [Dimension] Mode: 2D:CF Formula: {Ao st junct}/{Ao Diam} Needs measurement: Ao st junct [Dimension], Ao Diam [Dimension] Measured by: Ao st junct [2DCALIPER] |
| Ao/LA [Generic, Dimension] Mode: MM:CM:AMM Formula: {Ao Diam}/{LA Diam} Needs measurement: Ao Diam [Generic, Dimension], LA Diam [Generic, Dimension] Measured by: LA/Ao [MMLAAO] |
| AP Area [Aortic] Mode: CW:PW Formula: {LVOT Diam} ² *0.785*({LVOT VTI}/{AP VTI}) Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], AP VTI [Aortic] Measured by: AP Area [SDMANTRACE] |
| AR ERO [PISA] Mode: CF:CW:PW Formula: {AR Flow}/{AR Vmax} Needs measurement: AR Flow [PISA], AR Vmax [PISA] Measured by: AR Trace [AUTOCALC] |
| AR RV [PISA] Mode: CF:CW:PW Formula: {AR Flow}/{AR Vmax}*{AR VTI} Needs measurement: AR Flow [PISA], AR Vmax [PISA], AR VTI [PISA] Measured by: AR Trace [AUTOCALC] |
| AV Acc Time/ET Ratio [Aortic] Mode: CW:PW Formula: {AV AccT}/{AVET} Needs measurement: AV AccT [Aortic], AVET [Aortic] Measured by: AVET [SDTIMECALIPER] |
| AV Area [Dimension] Mode: 2D:CF Formula: 3.14/4*{AV Diam} ² Needs measurement: AV Diam [Dimension] Measured by: AV Diam [2DCALIPER] |
| AV CI [Aortic] Mode: CW:PW Formula: (({AV Diam} ² *0.785*{AV VTI})*{HR}/60)/{BSA} Needs measurement: AV Diam [Aortic], AV VTI [Aortic], HR [Aortic] Measured by: AV Trace [SDMANTRACE] |

Measurements

| |
|---|
| AV CO [Aortic] Mode: CW:PW Formula: $(\{AV\ Diam\}^2 * 0.785 * \{AV\ VTI\}) * \{HR\} / 60$ Needs measurement: AV Diam [Aortic], AV VTI [Aortic], HR [Aortic] Measured by: AV Trace [SDMANTRACE] |
| AV SI [Aortic] Mode: CW:PW Formula: $(\{AV\ Diam\}^2 * 0.785 * \{AV\ VTI\}) / \{BSA\}$ Needs measurement: AV Diam [Aortic], AV VTI [Aortic] Measured by: AV Trace [SDMANTRACE] |
| AV SV [Aortic] Mode: CW:PW Formula: $\{AV\ Diam\}^2 * 0.785 * \{AV\ VTI\}$ Needs measurement: AV Diam [Aortic], AV VTI [Aortic] Measured by: AV Trace [SDMANTRACE] |
| AVA (VTI) [Aortic] Mode: 2D:CW:PW Formula: $3.14 / 4 * \{LVOT\ Diam\}^2 * \{LVOT\ VTI\} / \{AV\ VTI\}$ Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], AV VTI [Aortic] Measured by: AV Trace [AUTOCALC] |
| AVA Vmax [Aortic] Mode: 2D:CW:PW Formula: $3.14 / 4 * \{LVOT\ Diam\}^2 * \text{abs}(\{LVOT\ Vmax\} / \{AV\ Vmax\})$ Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic] Measured by: AV Vmax [AUTOCALC] |
| AVA Vmax [Aortic] Mode: 2D:CW:PW Formula: $3.14 / 4 * \{LVOT\ Diam\}^2 * \text{abs}(\{LVOT\ Vmax\} / \{AV\ Vmax\})$ Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic] Measured by: AV Trace [AUTOCALC] |
| AVA Vmax, Pt [Aortic] Mode: 2D:CW:PW Formula: $3.14 / 4 * \{LVOT\ Diam\}^2 * \text{abs}(\{LVOT\ Vmax\} / \{AV\ Vmax\})$ Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic] Measured by: AV Vmax [AUTOCALC] |
| AVA Vmax, Pt [Aortic] Mode: 2D:CW:PW Formula: $3.14 / 4 * \{LVOT\ Diam\}^2 * \text{abs}(\{LVOT\ Vmax\} / \{AV\ Vmax\})$ Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic] Measured by: AV Trace [AUTOCALC] |

AVAI (VTI) [Aortic]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{LVOT \text{ Diam}\}^2 \times \{LVOT \text{ VTI}\} / \{AV \text{ VTI}\} \times \{BSA\}$

Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], AV VTI [Aortic]

Measured by: AV Trace [AUTOCALC]

AVAI Vmax [Aortic]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{LVOT \text{ Diam}\}^2 \times \text{abs}(\{LVOT \text{ Vmax}\} / \{AV \text{ Vmax}\} \times \{BSA\})$

Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic]

Measured by: AV Vmax [AUTOCALC]

AVAI Vmax [Aortic]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{LVOT \text{ Diam}\}^2 \times \text{abs}(\{LVOT \text{ Vmax}\} / \{AV \text{ Vmax}\} \times \{BSA\})$

Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic]

Measured by: AV Trace [AUTOCALC]

AVAI Vmax, Pt [Aortic]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{LVOT \text{ Diam}\}^2 \times \text{abs}(\{LVOT \text{ Vmax}\} / \{AV \text{ Vmax}\} \times \{BSA\})$

Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic]

Measured by: AV Vmax [AUTOCALC]

AVAI Vmax, Pt [Aortic]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{LVOT \text{ Diam}\}^2 \times \text{abs}(\{LVOT \text{ Vmax}\} / \{AV \text{ Vmax}\} \times \{BSA\})$

Needs measurement: LVOT Diam [Aortic], LVOT Vmax [Aortic], AV Vmax [Aortic]

Measured by: AV Trace [AUTOCALC]

CI A-L A2C [Single Plane A2C, AutoBiplane]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ A-L A2C}\} - \{LVESV \text{ A-L A2C}\}) \times \{HR\} / 60) / \{BSA\}$

Needs measurement: LVEDV A-L A2C [Single Plane A2C, AutoBiplane], LVESV A-L A2C [Single Plane A2C, AutoBiplane], HR [Single Plane A2C, AutoBiplane]

Measured by: R-R [2DCALIPER], A2C [2DAUTOVOLUME]

CI A-L A2C [Single Plane A2C]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ A-L A2C}\} - \{LVESV \text{ A-L A2C}\}) \times \{HR\} / 60 / \text{Auto}) / \{BSA\}$

Needs measurement: LVEDV A-L A2C [Single Plane A2C], LVESV A-L A2C [Single Plane A2C], HR [Single Plane A2C]

Measured by: LVESV A2C [2DVOLUMETRACE]

Measurements

| |
|--|
| <p>CI A-L A4C [Single Plane A4C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $((\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}) * \{HR\} / 60) / \{BSA\}$</p> <p>Needs measurement: LVEDV A-L A4C [Single Plane A4C, AutoBiplane], LVESV A-L A4C [Single Plane A4C, AutoBiplane], HR [Single Plane A4C, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], A4C [2DAUTOVOLUME]</p> |
| <p>CI A-L A4C [Single Plane A4C]</p> <p>Mode: 2D:CF</p> <p>Formula: $((\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}) * \{HR\} / 60) / \{BSA\}$</p> <p>Needs measurement: LVEDV A-L A4C [Single Plane A4C], LVESV A-L A4C [Single Plane A4C], HR [Single Plane A4C]</p> <p>Measured by: LVESV A4C [2DVOLUMETRACE]</p> |
| <p>CI A-L LAX [Single Plane LAX, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $((\{LVEDV\ A-L\ LAX\} - \{LVESV\ A-L\ LAX\}) * \{HR\} / 60) / \{BSA\}$</p> <p>Needs measurement: LVEDV A-L LAX [Single Plane LAX, AutoBiplane], LVESV A-L LAX [Single Plane LAX, AutoBiplane], HR [Single Plane LAX, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], AutoVolume [2DAUTOVOLUME]</p> |
| <p>CI Biplane [Biplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $d = \text{biplane}(\{LVLd\ A4C\}, \{LVD\text{Disks}\}, \{LVLd\ A2C\}, \{LVD\text{Disks}\})$</p> <p>Needs measurement: LVLd A4C [Biplane], LVLd A2C [Biplane], LVLs A4C [Biplane], LVLs A2C [Biplane], HR [Biplane]</p> <p>Measured by: R-R [2DCALIPER]</p> |
| <p>CI bp el [Biplane Ellipse]</p> <p>Mode: 2D:CF</p> <p>Formula: $((d-s) * \{ECG/HeartRate\} / 60) / \{BSA\}$ where: $s = (8 / (3 * 3.14159)) * \{LVAs(A4C)\} * \{LVAs(sax\ MV)\} / \{2D/LVIDs\}$ $d = (8 / (3 * 3.14159)) * \{LVAd\ A4C\} * \{LVAd\ (sax\ MV)\} / \{LVIDd\}$</p> <p>Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse], LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse], HR [Biplane Ellipse]</p> <p>Measured by: R-R [2DCALIPER]</p> |
| <p>CI bullet [Bullet]</p> <p>Mode: 2D:CF</p> <p>Formula: $((d-s) * \{ECG/HeartRate\} / 60) / \{BSA\}$ where: $s = 5/6 * \{LVAs(sax)\} * \{LVLs(apical)\}$ $d = 5/6 * \{LVAd\ sax\} * \{LVLd\ apical\}$</p> <p>Needs measurement: LVAd sax [Bullet], LVLd apical [Bullet], LVAs sax [Bullet], LVLs apical [Bullet], HR [Bullet]</p> <p>Measured by: R-R [2DCALIPER]</p> |
| <p>CI MOD A2C [Single Plane A2C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $((\{LVEDV\ MOD\ A2C\} - \{LVESV\ MOD\ A2C\}) * \{HR\} / 60) / \{BSA\}$</p> <p>Needs measurement: LVEDV MOD A2C [Single Plane A2C, AutoBiplane], LVESV MOD A2C [Single Plane A2C, AutoBiplane], HR [Single Plane A2C, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], A2C [2DAUTOVOLUME]</p> |

CI MOD A2C [Single Plane A2C]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ MOD A2C}\} - \{LVESV \text{ MOD A2C}\}) * \{HR\} / 60) / \{BSA\}$

Needs measurement: LVEDV MOD A2C [Single Plane A2C], LVESV MOD A2C [Single Plane A2C], HR [Single Plane A2C]

Measured by: LVESV A2C [2DVOLUMETRACE]

CI MOD A4C [Single Plane A4C, AutoBiplane]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ MOD A4C}\} - \{LVESV \text{ MOD A4C}\}) * \{HR\} / 60) / \{BSA\}$

Needs measurement: LVEDV MOD A4C [Single Plane A4C, AutoBiplane], LVESV MOD A4C [Single Plane A4C, AutoBiplane], HR [Single Plane A4C, AutoBiplane]

Measured by: R-R [2DCALIPER], A4C [2DAUTOVOLUME]

CI MOD A4C [Single Plane A4C]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ MOD A4C}\} - \{LVESV \text{ MOD A4C}\}) * \{HR\} / 60) / \{BSA\}$

Needs measurement: LVEDV MOD A4C [Single Plane A4C], LVESV MOD A4C [Single Plane A4C], HR [Single Plane A4C]

Measured by: LVESV A4C [2DVOLUMETRACE]

CI MOD LAX [Single Plane LAX, AutoBiplane]

Mode: 2D:CF

Formula: $((\{LVEDV \text{ MOD LAX}\} - \{LVESV \text{ MOD LAX}\}) * \{HR\} / 60) / \{BSA\}$

Needs measurement: LVEDV MOD LAX [Single Plane LAX, AutoBiplane], LVESV MOD LAX [Single Plane LAX, AutoBiplane], HR [Single Plane LAX, AutoBiplane]

Measured by: R-R [2DCALIPER], AutoVolume [2DAUTOVOLUME]

CI mod sim [Modified Simpson]

Mode: 2D:CF

Formula: $((d-s) * \{ECG/HeartRate\} / 60) / \{BSA\}$ where: $s = (\{LVLs(apical)\} / 9) * ((4 * \{LVAs(sax \text{ MV})\}) + (2 * \{LVAs(sax \text{ PM})\}) + \sqrt{\{LVAs(sax \text{ MV})\} * \{LVAs(sax \text{ PM})\}})$ $d = (\{LVLD \text{ apical}\} / 9) * ((4 * \{LVAd(sax \text{ MV})\}) + (2 * \{LVAd(sax \text{ PM})\}) + \sqrt{\{LVAd(sax \text{ MV})\} * \{LVAd(sax \text{ PM})\}})$

Needs measurement: LVLD apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson], LVLs apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson], HR [Modified Simpson]

Measured by: R-R [2DCALIPER]

CI(Cube) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: $((d-s) * \{ECG/HeartRate\} / 60) / \{BSA\}$ where: $s = \{2D/LVIDs\}^3$ $d = \{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz], HR [Dimension, Cube/Teicholz]

Measured by: R-R [2DCALIPER]

CI(Cube) [Generic, Dimension]

Mode: MM:CM:AMM

Formula: $((dv-sv) * \{MM/HeartRate\} / 60) / \{BSA\}$ where: $sv = \{MM/LVIDs\}^3$ $dv = \{LVIDd\}^3$

Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension], HR [Generic, Dimension]

Measured by: Heartrate [MMTIMECALIPER]

Measurements

| |
|--|
| <p>CI(Teich) [Dimension, Cube/Teicholz]</p> <p>Mode: 2D:CF</p> <p>Formula: $((d-s) \cdot \{ECG/HeartRate\}/60) / \{BSA\}$ where: $s = 7 / (2.4 + \{2D/LVIDs\}) \cdot \{2D/LVIDs\}^3$ $d = 7 / (2.4 + \{LVIDd\}) \cdot \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz], HR [Dimension, Cube/Teicholz]</p> <p>Measured by: R-R [2DCALIPER]</p> |
| <p>CI(Teich) [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $((dv-sv) \cdot \{MM/HeartRate\}/60) / \{BSA\}$ where: $sv = 7 / (2.4 + \{MM/LVIDs\}) \cdot \{MM/LVIDs\}^3$ $dv = 7 / (2.4 + \{LVIDd\}) \cdot \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension], HR [Generic, Dimension]</p> <p>Measured by: Heartrate [MMTIMECALIPER]</p> |
| <p>CO A-L A2C [Single Plane A2C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $(\{LVEDV\ A-L\ A2C\} - \{LVESV\ A-L\ A2C\}) \cdot \{HR\} / 60$</p> <p>Needs measurement: LVEDV A-L A2C [Single Plane A2C, AutoBiplane], LVESV A-L A2C [Single Plane A2C, AutoBiplane], HR [Single Plane A2C, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], A2C [2DAUTOVOLUME]</p> |
| <p>CO A-L A2C [Single Plane A2C]</p> <p>Mode: 2D:CF</p> <p>Formula: $(\{LVEDV\ A-L\ A2C\} - \{LVESV\ A-L\ A2C\}) \cdot \{HR\} / 60$</p> <p>Needs measurement: LVEDV A-L A2C [Single Plane A2C], LVESV A-L A2C [Single Plane A2C], HR [Single Plane A2C]</p> <p>Measured by: LVESV A2C [2DVOLUMETRACE]</p> |
| <p>CO A-L A4C [Single Plane A4C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $(\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}) \cdot \{HR\} / 60$</p> <p>Needs measurement: LVEDV A-L A4C [Single Plane A4C, AutoBiplane], LVESV A-L A4C [Single Plane A4C, AutoBiplane], HR [Single Plane A4C, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], A4C [2DAUTOVOLUME]</p> |
| <p>CO A-L A4C [Single Plane A4C]</p> <p>Mode: 2D:CF</p> <p>Formula: $(\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}) \cdot \{HR\} / 60$</p> <p>Needs measurement: LVEDV A-L A4C [Single Plane A4C], LVESV A-L A4C [Single Plane A4C], HR [Single Plane A4C]</p> <p>Measured by: LVESV A4C [2DVOLUMETRACE]</p> |
| <p>CO A-L LAX [Single Plane LAX, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $(\{LVEDV\ A-L\ LAX\} - \{LVESV\ A-L\ LAX\}) \cdot \{HR\} / 60$</p> <p>Needs measurement: LVEDV A-L LAX [Single Plane LAX, AutoBiplane], LVESV A-L LAX [Single Plane LAX, AutoBiplane], HR [Single Plane LAX, AutoBiplane]</p> <p>Measured by: R-R [2DCALIPER], AutoVolume [2DAUTOVOLUME]</p> |

CO A-L LAX [Single Plane LAX]

Mode: 2D:CF

Formula: $(\{LVEDV \text{ A-L LAX}\} - \{LVESV \text{ A-L LAX}\}) * \{HR\} / 60$

Needs measurement: LVEDV A-L LAX [Single Plane LAX], LVESV A-L LAX [Single Plane LAX], HR [Single Plane LAX]

Measured by: LVESV LAX [2DVOLUMETRACE]

CO Biplane [Biplane]

Mode: 2D:CF

Formula: $d = \text{biplane}(\{LVLd \text{ A4C}\}, \{LVDIsks\}, \{LVLd \text{ A2C}\}, \{LVDIsks\})$

Needs measurement: LVLd A4C [Biplane], LVLd A2C [Biplane], LVLs A4C [Biplane], LVLs A2C [Biplane], HR [Biplane]

Measured by: R-R [2DCALIPER]

CO bp el [Biplane Ellipse]

Mode: 2D:CF

Formula: $(d-s) * \{ECG/HeartRate\} / 60$ where: $s = (8 / (3 * 3.14159)) * \{LVAs(A4C)\} * \{LVAs(sax \text{ MV})\} / \{2D/LVIDs\}$ $d = (8 / (3 * 3.14159)) * \{LVAd \text{ A4C}\} * \{LVAd (sax \text{ MV})\} / \{LVIDd\}$

Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse], LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse], HR [Biplane Ellipse]

Measured by: R-R [2DCALIPER]

CO bullet [Bullet]

Mode: 2D:CF

Formula: $(d-s) * \{ECG/HeartRate\} / 60$ where: $s = 5/6 * \{LVAs(sax)\} * \{LVLs(apical)\}$ $d = 5/6 * \{LVAd \text{ sax}\} * \{LVLd \text{ apical}\}$

Needs measurement: LVAd sax [Bullet], LVLd apical [Bullet], LVLs apical [Bullet], HR [Bullet]

Measured by: R-R [2DCALIPER]

CO MOD A2C [Single Plane A2C, AutoBiplane]

Mode: 2D:CF

Formula: $(\{LVEDV \text{ MOD A2C}\} - \{LVESV \text{ MOD A2C}\}) * \{HR\} / 60$

Needs measurement: LVEDV MOD A2C [Single Plane A2C, AutoBiplane], LVESV MOD A2C [Single Plane A2C, AutoBiplane], HR [Single Plane A2C, AutoBiplane]

Measured by: R-R [2DCALIPER], A2C [2DAUTOVOLUME]

CO MOD A2C [Single Plane A2C]

Mode: 2D:CF

Formula: $(\{LVEDV \text{ MOD A2C}\} - \{LVESV \text{ MOD A2C}\}) * \{HR\} / 60$

Needs measurement: LVEDV MOD A2C [Single Plane A2C], LVESV MOD A2C [Single Plane A2C], HR [Single Plane A2C]

Measured by: LVESV A2C [2DVOLUMETRACE]

CO MOD A4C [Single Plane A4C, AutoBiplane]

Mode: 2D:CF

Formula: $(\{LVEDV \text{ MOD A4C}\} - \{LVESV \text{ MOD A4C}\}) * \{HR\} / 60$

Needs measurement: LVEDV MOD A4C [Single Plane A4C, AutoBiplane], LVESV MOD A4C [Single Plane A4C, AutoBiplane], HR [Single Plane A4C, AutoBiplane]

Measured by: R-R [2DCALIPER], A4C [2DAUTOVOLUME]

Measurements

| |
|---|
| CO MOD A4C [Single Plane A4C] Mode: 2D:CF Formula: $\{(\text{LVEDV MOD A4C}) - (\text{LVESV MOD A4C})\} \times \{\text{HR}\} / 60$ Needs measurement: LVEDV MOD A4C [Single Plane A4C], LVESV MOD A4C [Single Plane A4C], HR [Single Plane A4C] Measured by: LVESV A4C [2DVOLUMETRACE] |
| CO MOD LAX [Single Plane LAX, AutoBiplane] Mode: 2D:CF Formula: $\{(\text{LVEDV MOD LAX}) - (\text{LVESV MOD LAX})\} \times \{\text{HR}\} / 60$ Needs measurement: LVEDV MOD LAX [Single Plane LAX, AutoBiplane], LVESV MOD LAX [Single Plane LAX, AutoBiplane], HR [Single Plane LAX, AutoBiplane] Measured by: R-R [2DCALIPER], AutoVolume [2DAUTOVOLUME] |
| CO MOD LAX [Single Plane LAX] Mode: 2D:CF Formula: $\{(\text{LVEDV MOD LAX}) - (\text{LVESV MOD LAX})\} \times \{\text{HR}\} / 60$ Needs measurement: LVEDV MOD LAX [Single Plane LAX], LVESV MOD LAX [Single Plane LAX], HR [Single Plane LAX] Measured by: LVESV LAX [2DVOLUMETRACE] |
| CO mod sim [Modified Simpson] Mode: 2D:CF Formula: $(d-s) \times \{\text{ECG/HeartRate}\} / 60$ where: $s = \{(\text{LVLs(apical)}) / 9\} \times \{((4 \times \{\text{LVAs(sax MV)})\}) + (2 \times \{\text{LVAs(sax d = } (\text{LVLd apical}) / 9\} \times ((4 \times \{\text{LVAd (sax MV)})\}) + (2 \times \{\text{LVAd sax PM}\}) + \text{sqrt}(\{\text{LVAd (sax MV)}\} \times \{\text{LVAd sax PM}\})\})$ Needs measurement: LVLd apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson], LVLs apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson], HR [Modified Simpson] Measured by: R-R [2DCALIPER] |
| CO(A-L) [Generic] Mode: 2D:CF Formula: $\{(\text{EDV(A-L)}) - (\text{ESV(A-L)})\} \times \{\text{HR}\} / 60$ Needs measurement: ESV(A-L) [Generic], HR [Generic] Measured by: R-R [2DCALIPER] |
| CO(Cube) [Dimension, Cube/Teicholz] Mode: 2D:CF Formula: $(d-s) \times \{\text{ECG/HeartRate}\} / 60$ where: $s = \{2D/\text{LVIDs}\}^3$ $d = \{\text{LVIDd}\}^3$ Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz], HR [Dimension, Cube/Teicholz] Measured by: R-R [2DCALIPER] |
| CO(Cube) [Generic, Dimension] Mode: MM:CM:AMM Formula: $(dv-sv) \times \{\text{MM/HeartRate}\} / 60$ where: $sv = \{\text{MM/LVIDs}\}^3$ $dv = \{\text{LVIDd}\}^3$ Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension], HR [Generic, Dimension] Measured by: Heartrate [MMTIMECALIPER] |

CO(Teich) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: $(d-s) \cdot \{ECG/HeartRate\}/60$ where: $s = 7/(2.4 + \{2D/LVIDs\}) \cdot \{2D/LVIDs\}^3$ $d = 7/(2.4 + \{LVIDd\}) \cdot \{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz], HR [Dimension, Cube/Teicholz]

Measured by: R-R [2DCALIPER]

CO(Teich) [Generic, Dimension]

Mode: MM:CM:AMM

Formula: $(dv-sv) \cdot \{MM/HeartRate\}/60$ where: $sv = 7/(2.4 + \{MM/LVIDs\}) \cdot \{MM/LVIDs\}^3$ $dv = 7/(2.4 + \{LVIDd\}) \cdot \{LVIDd\}^3$

Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension], HR [Generic, Dimension]

Measured by: Heartrate [MMTIMECALIPER]

EDV bp el [Biplane Ellipse]

Mode: 2D:CF

Formula: $(8/(3 \cdot 3.14159)) \cdot \{LVAd A4C\} \cdot \{LVAd (sax MV)\} / \{LVIDd\}$

Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse]

Measured by: LVEF BP-EL [AUTOCALC]

EDV bullet [Bullet]

Mode: 2D:CF

Formula: $5/6 \cdot \{LVAd sax\} \cdot \{LVLd apical\}$

Needs measurement: LVAd sax) [Bullet], LVLd apical [Bullet]

Measured by: LVEF Bullet [AUTOCALC]

EDV mod sim [Modified Simpson]

Mode: 2D:CF

Formula: $(\{LVLd apical\}/9) \cdot ((4 \cdot \{LVAd (sax MV)\}) + (2 \cdot \{LVAd sax PM\}) + \sqrt{\{LVAd (sax MV)\} \cdot \{LVAd sax PM\}})$

Needs measurement: LVLd apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson]

Measured by: EF mod sim [AUTOCALC]

EDV(Cube) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: $\{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz]

Measured by: LVd [2DLV], LVIDd [2DCALIPER], EF(Cube) [AUTOCALC]

EDV(Cube) [Generic, Dimension]

Mode: MM:CM:AMM

Formula: $\{LVIDd\}^3$

Needs measurement: LVIDd [Generic, Dimension]

Measured by: LV Study [MMLV], LVIDd [MMDISCALIPER]

Measurements

| |
|--|
| EDV(Teich) [Dimension, Cube/Teicholz] Mode: 2D:CF Formula: $7/(2.4 + \{LVIDd\}) * \{LVIDd\}^3$ Needs measurement: LVIDd [Dimension, Cube/Teicholz] Measured by: LVd [2DLV], LVIDd [2DCALIPER], EF(Cube) [AUTOCALC] |
| EDV(Teich) [Generic, Dimension] Mode: MM:CM:AMM Formula: $7/(2.4 + \{LVIDd\}) * \{LVIDd\}^3$ Needs measurement: LVIDd [Generic, Dimension] Measured by: LV Study [MMLV], LVIDd [MMDISCALIPER] |
| EF A-L A2C [Biplane, Single Plane A2C, AutoBiplane] Mode: 2D:CF Formula: $(\{LVEDV\ A-L\ A2C\} - \{LVESV\ A-L\ A2C\}) / \{LVEDV\ A-L\ A2C\}$ Needs measurement: LVEDV A-L A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV A-L A2C [Biplane, Single Plane A2C, AutoBiplane] Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME] |
| EF A-L A4C [Biplane, Single Plane A4C, AutoBiplane] Mode: 2D:CF: Formula: $(\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}) / \{LVEDV\ A-L\ A4C\}$ Needs measurement: LVEDV A-L A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV A-L A4C [Biplane, Single Plane A4C, AutoBiplane] Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME] |
| EF A-L LAX [Single Plane LAX, AutoBiplane] Mode: 2D:CF Formula: $(\{LVEDV\ A-L\ LAX\} - \{LVESV\ A-L\ LAX\}) / \{LVEDV\ A-L\ LAX\}$ Needs measurement: LVEDV A-L LAX [Single Plane LAX, AutoBiplane], LVESV A-L LAX [Single Plane LAX, AutoBiplane] Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME] |
| EF Biplane [Biplane, AutoBiplane] Mode: 2D:CF Formula: $d = \text{biplane}(\{LVLd\ A4C\}, \{LVDisks\}, \{LVLd\ A2C\}, \{LVDisks\})$ Needs measurement: LVLd A4C [Biplane, AutoBiplane], LVLd A2C [Biplane, AutoBiplane], LVLs A4C [Biplane, AutoBiplane], LVLs A2C [Biplane, AutoBiplane] Measured by: EF Biplane [AUTOCALC] |
| EF mod sim [Modified Simpson] Mode: 2D:CF Formula: $(\{LVLd\ apical\}/9) * ((4 * \{LVAd\ (sax\ MV)\}) + (2 * \{LVAd\ sax\ PM\}) + \text{sqrt}(\{LVAd\ (sax\ MV)\} * \{LVAd\ sax\ PM\}))$ Needs measurement: LVLd apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson], LVLs apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson] Measured by: EF mod sim [AUTOCALC] |

| |
|---|
| <p>EF(A-L) [Generic]</p> <p>Mode: 2D:CF</p> <p>Formula: $\frac{({EDV(A-L)} - {ESV(A-L)})}{{EDV(A-L)}}$</p> <p>Needs measurement: ESV(A-L) [Generic], EDV(A-L) [Generic]</p> <p>Measured by: EF Volume [AUTOCALC]</p> |
| <p>EF(Cube) [Dimension, Cube/Teicholz]</p> <p>Mode: 2D:CF</p> <p>Formula: $(d-s)/d$ where: $s = \{2D/LVIDs\}^3$ $d = \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz]</p> <p>Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]</p> |
| <p>EF(Cube) [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $(dv-sv)/dv$ where: $sv = \{MM/LVIDs\}^3$ $dv = \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension]</p> <p>Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER]</p> |
| <p>EF(MOD) [Generic]</p> <p>Mode: 2D:CF</p> <p>Formula: $\frac{({EDV(MOD)} - {ESV(MOD)})}{{EDV(MOD)}}$</p> <p>Needs measurement: EDV(MOD) [Generic], ESV(MOD) [Generic]</p> <p>Measured by: EF Volume [AUTOCALC]</p> |
| <p>EF(Teich) [Dimension, Cube/Teicholz]</p> <p>Mode: 2D:CF</p> <p>Formula: $(d-s)/d$ where: $s = 7/(2.4 + \{2D/LVIDs\}) * \{2D/LVIDs\}^3$ $d = 7/(2.4 + \{LVIDd\}) * \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz]</p> <p>Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]</p> |
| <p>EF(Teich) [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $(dv-sv)/dv$ where: $sv = 7/(2.4 + \{MM/LVIDs\}) * \{MM/LVIDs\}^3$ $dv = 7/(2.4 + \{LVIDd\}) * \{LVIDd\}^3$</p> <p>Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension]</p> <p>Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER]</p> |
| <p>ESV bp el [Biplane Ellipse]</p> <p>Mode: 2D:CF</p> <p>Formula: $(8/(3 * 3.14159)) * \{LVAs A4C\} * \{LVAs sax MV\} / \{LVIDs\}$</p> <p>Needs measurement: LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse]</p> <p>Measured by: LVEF BP-EL [AUTOCALC]</p> |
| <p>ESV bullet [Bullet]</p> <p>Mode: 2D:CF</p> <p>Formula: $5/6 * \{LVAs sax\} * \{LVLs apical\}$</p> <p>Needs measurement: LVAs sax [Bullet], LVLs apical [Bullet]</p> <p>Measured by: LVEF Bullet [AUTOCALC]</p> |

Measurements

| |
|---|
| <p>ESV mod sim [Modified Simpson]</p> <p>Mode: 2D:CF</p> <p>Formula: $((\{LVls\ apical\}/9)*((4*\{LVAs\ sax\ MV\})+(2*\{LVAs\ sax\ PM\}))+sqrt(\{LVAs\ sax\ MV\}*\{LVAs\ sax\ PM\}))$</p> <p>Needs measurement: LVls apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson]</p> <p>Measured by: EF mod sim [AUTOCALC]</p> |
| <p>ESV(Cube) [Dimension, Cube/Teicholz]</p> <p>Mode: 2D:CF</p> <p>Formula: $\{LVIDs\}^3$</p> <p>Needs measurement: LVIDs [Dimension, Cube/Teicholz]</p> <p>Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]</p> |
| <p>ESV(Cube) [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $\{LVIDs\}^3$</p> <p>Needs measurement: LVIDs [Generic, Dimension]</p> <p>Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER]</p> |
| <p>ESV(Teich) [Dimension, Cube/Teicholz]</p> <p>Mode: 2D:CF</p> <p>Formula: $7/(2.4+\{LVIDs\})*\{LVIDs\}^3$</p> <p>Needs measurement: LVIDs [Dimension, Cube/Teicholz]</p> <p>Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]</p> |
| <p>ESV(Teich) [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $7/(2.4+\{LVIDs\})*\{LVIDs\}^3$</p> <p>Needs measurement: LVIDs [Generic, Dimension]</p> <p>Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER]</p> |
| <p>HR (Generic, Dimension, Biplane, Modified Simpson, Cube/Teicholz, Single Plane A4C, Single Plane A2C, Single Plane LAX, Bullet, Biplane Ellipse)</p> <p>Mode: 2D:CF</p> <p>Formula: $60/\{R-R\}$</p> <p>Needs measurement: R-R [Generic, Dimension, Biplane, Modified Simpson, Cube/Teicholz, Single Plane A4C, Single Plane A2C, Single Plane LAX, Bullet, Biplane Ellipse]</p> <p>Measured by: R-R [2DCALIPER]</p> <p>Used to calculate: CO(A-L),CO(Teich),CI(Teich),CO(Cube),CI(Cube),CO Biplane,CI Biplane,CO mod sim,CI mod sim,CI A-L A4C,CO MOD A4C,CI MOD A4C,CI A-L A2C,CO A-L A2C,CI A-L A2C,CO MOD A2C,CI MOD A2C,CO A-L LAX,CI A-L LAX,CO MOD LAX,CI MOD LAX,CO bullet,CI bullet,CO bp el,CI bp el</p> |
| <p>HR [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $60/\{Time\}$</p> <p>Needs measurement: Time [Generic, Dimension]</p> <p>Measured by: Heartrate [MMTIMECALIPER]</p> <p>Used to calculate: CO(Cube),CO(Teich),CI(Teich),CI(Cube)</p> |

| |
|---|
| HR [Generic] Mode: CW:PW Formula: $60/\{\text{Time}\}$ Needs measurement: Time [Generic] Measured by: Heartrate [SDTIMECALIPER] |
| IVSd/LVPWd [Dimension] Mode: MM:CM:AMM Formula: $\{\text{IVSd}\}/\{\text{LVPWd}\}$ Needs measurement: IVSd [Dimension], LVPWd [Dimension] Measured by: LVPWd [MMDISCALIPER] |
| LA/Ao [Generic, Dimension] Mode: MM:CM:AMM Formula: $\{\text{LA Diam}\}/\{\text{Ao Diam}\}$ Needs measurement: LA Diam [Generic, Dimension], Ao Diam [Generic, Dimension] Measured by: LA/Ao [MMLAAO] |
| LIMP [Mitral Valve, Aortic] Mode: CW:PW Formula: $(\{\text{MCO}\}-\{\text{AVET}\})/\{\text{AVET}\}$ Needs measurement: MCO [Mitral Valve, Aortic], AVET [Mitral Valve, Aortic] Measured by: LIMP [AUTOCALC] |
| LVCI Dopp [Aortic] Mode: PW Formula: $(\{\text{LVOT Diam}\}^2 \times 0.785 \times \{\text{LVOT VTI}\} \times \{\text{HR}\} / 60) / \{\text{BSA}\}$ Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], HR [Aortic], Measured by: LVOT Trace [SDMANTRACE] |
| LVCO Dopp [Aortic] Mode: PW Formula: $\{\text{LVOT Diam}\}^2 \times 0.785 \times \{\text{LVOT VTI}\} \times \{\text{HR}\} / 60$ Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], HR [Aortic] Measured by: LVOT Trace [SDMANTRACE] |
| LVd Mass (ASE) [Generic] Mode: MM:CM:AMM Formula: $((1.04 \times (\{\text{IVSd}\} + \{\text{LVIDd}\} + \{\text{LVPWd}\})^3 - (\{\text{LVIDd}\})^3) \times 0.8 + 0.6) / 1000$ Needs measurement: IVSd [Generic], LVIDd [Generic], LVPWd [Generic] Measured by: LV Study [MMLV] |
| LVd Mass [Dimension] Mode: 2D:CF Formula: $((1.04 \times (\{\text{IVSd}\} + \{\text{LVIDd}\} + \{\text{LVPWd}\})^3 - (\{\text{LVIDd}\})^3) - 13.6) / 1000$ Needs measurement: IVSd [Dimension], LVIDd [Dimension], LVPWd [Dimension], LVIDd [Dimension] Measured by: LVPWd [2DCALIPER] |

Measurements

| |
|--|
| LVd Mass [Generic, Dimension] Mode: MM:CM:AMM Formula: $((1.04 * ((IVSd) + \{LVIDd\} + \{LVPWd\})^3 - \{LVIDd\}^3) - 13.6) / 1000$ Needs measurement: IVSd [Generic, Dimension], LVPWd [Generic, Dimension], LVIDd [Generic, Dimension] Measured by: LV Study [MMLV], LVPWs [MMDISCALIPER] |
| LVd Mass A-L [Mass] Mode: 2D:CF Formula: $1.05 * 5/6 * (\{LVAd(sax\ epi)\} * (\{LVLd(apical)\} + t) - \{LVAd(sax\ PM)\} * \{LVLd(apical)\}) / 1000$ where: $t = \sqrt{(\{LVAd\ sax\ EPI\} / 3.14159) - \sqrt{(\{LVAd\ sax\ PM\} / 3.14159)}}$ Needs measurement: LVAd sax EPI [Mass], LVAd sax PM [Mass], LVLd apical [Mass] Measured by: LVMass(d) [AUTOCALC] |
| LVd Mass I A-L [Mass] Mode: 2D:CF Formula: $m / \{BSA\}$ where: $m = 1.05 * 5/6 * (\{LVAd(sax\ epi)\} * (\{LVLd(apical)\} + t) - \{LVAd(sax\ PM)\} * \{LVLd(apical)\}) / 1000$ $t = \sqrt{(\{LVAd\ sax\ EPI\} / 3.14159) - \sqrt{(\{LVAd\ sax\ PM\} / 3.14159)}}$ Needs measurement: LVAd sax EPI [Mass], LVAd sax PM [Mass], LVLd apical [Mass] Measured by: LVMass(d) [AUTOCALC] |
| LVd Mass Ind (ASE) [Generic] Mode: MM:CM:AMM Formula: $((1.04 * ((IVSd) + \{LVIDd\} + \{LVPWd\})^3 - \{LVIDd\}^3) * 0.8 + 0.6) / 1000 / \{BSA\}$ Needs measurement: IVSd [Generic], LVIDd [Generic], LVPWd [Generic] Measured by: LV Study [MMLV] |
| LVd Mass Index [Dimension] Mode: 2D:CF Formula: $m / \{BSA\}$ where $m = ((1.04 * ((IVSd) + \{LVIDd\} + \{LVPWd\})^3 - \{LVIDd\}^3) - 13.6) / 1000$ Needs measurement: IVSd [Dimension], LVIDd [Dimension], LVPWd [Dimension], LVIDd [Dimension] Measured by: LVPWd [2DCALIPER] |
| LVd Mass Index [Generic, Dimension] Mode: MM:CM:AMM Formula: $((1.04 * ((IVSd) + \{LVIDd\} + \{LVPWd\})^3 - \{LVIDd\}^3) - 13.6) / 1000 / \{BSA\}$ Needs measurement: IVSd [Generic, Dimension], LVIDd [Generic, Dimension], LVPWd [Generic, Dimension] Measured by: LV Study [MMLV], LVPWs [MMDISCALIPER] |
| LVEDV MOD BP [Biplane, AutoBiplane] Mode: 2D:CF Formula: $biplane(\{LVLd\ A4C\}, \{LVDIsks\}, \{LVLd\ A2C\}, \{LVDIsks\})$ Needs measurement: LVLd A4C [Biplane, AutoBiplane], LVLd A2C [Biplane, AutoBiplane] Measured by: EF Biplane [AUTOCALC] |

LVEF BP-EL [Biplane Ellipse]

Mode: 2D:CF

Formula: $(d-s)/d$ where: $s = (8/(3*3.14159))*\{LVAs(A4C)\}*\{LVAs(sax\ MV)\}/\{2D/LVIDs\}$ $d = (8/(3*3.14159))*\{LVAd\ A4C\}*\{LVAd\ (sax\ MV)\}/\{LVIDd\}$

Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse], LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse]

Measured by: LVEF BP-EL [AUTOCALC]

LVEF Bullet [Bullet]

Mode: 2D:CF

Formula: $(d-s)/d$ where: $s = 5/6*\{LVAs(sax)\}*\{LVls(apical)\}$ $d = 5/6*\{LVAd\ sax\}*\{LVld\ apical\}$

Needs measurement: LVAd sax [Bullet], LVld apical [Bullet], LVls apical [Bullet]

Measured by: LVEF Bullet [AUTOCALC]

LVEF MOD A2C [Biplane, Single Plane A2C, AutoBiplane]

Mode: 2D:CF

Formula: $(\{LVEDV\ MOD\ A2C\} - \{LVESV\ MOD\ A2C\})/\{LVEDV\ MOD\ A2C\}$

Needs measurement: LVEDV MOD A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV MOD A2C [Biplane, Single Plane A2C, AutoBiplane]

Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME]

LVEF MOD A4C [Biplane, Single Plane A4C, AutoBiplane]

Mode: 2D:CF

Formula: $(\{LVEDV\ MOD\ A4C\} - \{LVESV\ MOD\ A4C\})/\{LVEDV\ MOD\ A4C\}$

Needs measurement: LVEDV MOD A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV MOD A4C [Biplane, Single Plane A4C, AutoBiplane]

Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME]

LVEF MOD LAX [Single Plane LAX, AutoBiplane]

Mode: 2D:CF

Formula: $(\{LVEDV\ MOD\ LAX\} - \{LVESV\ MOD\ LAX\})/\{LVEDV\ MOD\ LAX\}$

Needs measurement: LVEDV MOD LAX [Single Plane LAX, AutoBiplane], LVESV MOD LAX [Single Plane LAX, AutoBiplane]

Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME]

LVESV MOD BP [Biplane, AutoBiplane]

Mode: 2D:CF

Formula: $biplane(\{LVls\ A4C\}, \{LVDisks\}, \{LVls\ A2C\}, \{LVDisks\})$

Needs measurement: LVls A4C [Biplane, AutoBiplane], LVls A2C [Biplane, AutoBiplane]

Measured by: EF Biplane [AUTOCALC]

LVIDd Index [Dimension]

Mode: 2D:CF

Formula: $\{LVIDd\}/\{BSA\}$

Needs measurement: LVIDd [Dimension],

Measured by: LVIDd [2DCALIPER]

Measurements

| |
|---|
| LVIDd Index [Dimension] Mode: MM:CM:AMM Formula: $\{\text{LVIDd}\}/\{\text{BSA}\}$ Needs measurement: LVIDd [Dimension] Measured by: LVIDd [MMDISCALIPER] |
| LVIDs Index [Dimension] Mode: 2D:CF Formula: $\{\text{LVIDs}\}/\{\text{BSA}\}$ Needs measurement: LVIDs [Dimension] Measured by: LVIDs [2DCALIPER] |
| LVIDs Index [Dimension] Mode: MM:CM:AMM Formula: $\{\text{LVIDs}\}/\{\text{BSA}\}$ Needs measurement: LVIDs [Dimension] Measured by: LVIDs [MMDISCALIPER] |
| LVOT Area [Dimension] Mode: 2D:CF Formula: $3.14/4 * \{\text{LVOT Diam}\}^2$ Needs measurement: LVOT Diam [Dimension] Measured by: LVOT Diam [2DCALIPER] |
| LVOT Diam [Aortic] Mode: CW:PW Formula: $\{\text{LVOT Diam}\}$ Needs measurement: LVOT Diam [Aortic] Measured by: AP Area [SDMANTRACE] Used to calculate: AP Area |
| LVOT Diam [Mitral Valve] Mode: CW:PW Formula: $\{\text{LVOT Diam}\}$ Needs measurement: LVOT Diam [Mitral Valve] Measured by: MP Area [SDMANTRACE] Used to calculate: MP Area |
| LVOT VTI [Aortic] Mode: CW:PW Formula: $\{\text{LVOT VTI}\}$ Needs measurement: LVOT VTI [Aortic] Measured by: AP Area [SDMANTRACE] Used to calculate: AP Area |

| |
|--|
| <p>LVOT VTI [Mitral Valve]</p> <p>Mode: CW:PW</p> <p>Formula: {LVOT VTI}</p> <p>Needs measurement: LVOT VTI [Mitral Valve]</p> <p>Measured by: MP Area [SDMANTRACE]</p> <p>Used to calculate: MP Area</p> |
| <p>LVPEP/ET [Aortic]</p> <p>Mode: CW:PW</p> <p>Formula: {LVPEP}/{LVET}</p> <p>Needs measurement: LVPEP [Aortic], LVET [Aortic]</p> <p>Measured by: LVET [SDTIMECALIPER]</p> |
| <p>LVPEP/ET [Time]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: {LVPEP}/{LVET}</p> <p>Needs measurement: LVPEP [Time], LVET [Time]</p> <p>Measured by: LVET [MMTIMECALIPER]</p> |
| <p>LVs Mass (ASE) [Generic]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $((1.04 * ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - (\{LVIDs\})^3)) * 0.8 + 0.6) / 1000$</p> <p>Needs measurement: IVSs [Generic], LVIDs [Generic], LVPWs [Generic]</p> <p>Measured by: LV Study [MMLV]</p> |
| <p>LVs Mass [Dimension]</p> <p>Mode: 2D:CF</p> <p>Formula: $((1.04 * ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - (\{LVIDs\})^3)) - 13.6) / 1000$</p> <p>Needs measurement: IVSs [Dimension], LVIDs [Dimension], LVPWs [Dimension]</p> <p>Measured by: LVPWs [2DCALIPER]</p> |
| <p>LVs Mass [Generic, Dimension]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $((1.04 * ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - (\{LVIDs\})^3)) - 13.6) / 1000$</p> <p>Needs measurement: IVSs [Generic, Dimension], LVIDs [Generic, Dimension], LVPWs [Generic, Dimension]</p> <p>Measured by: LV Study [MMLV], LVPWs [MMDISCALIPER]</p> |
| <p>LVs Mass A-L [Mass]</p> <p>Mode: 2D:CF</p> <p>Formula: $1.05 * 5/6 * (\{LVAs(sax\ epi)\} * (\{LVLs(apical)\} + t) - \{LVAs(sax\ PM)\} * \{LVLs(apical)\}) / 1000$ where: $t = \sqrt{(\{LVAs\ sax\ EPI\} / 3.14159) - \sqrt{(\{LVAs\ sax\ PM\} / 3.14159)}}$</p> <p>Needs measurement: LVAs sax EPI [Mass], LVAs sax PM [Mass], LVLs apical [Mass]</p> <p>Measured by: LVMass(s) [AUTOCALC]</p> |
| <p>LVs Mass Ind (ASE) [Generic]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $((1.04 * ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - (\{LVIDs\})^3)) * 0.8 + 0.6) / 1000 / \{BSA\}$</p> <p>Needs measurement: IVSs [Generic], LVIDs [Generic], LVPWs [Generic]</p> <p>Measured by: LV Study [MMLV]</p> |

Measurements

| |
|--|
| LVs Mass Ind A-L [Mass] Mode: 2D:CF Formula: $m/\{BSA\}$ where: $m = 1.05 \times 5/6 \times (\{LVAs(sax\ epi)\} \times (\{LVLs(apical)\} + t) - \{LVAs(sax\ PM)\} \times \{LVLs(apical)\}) / 1000$ $t = \sqrt{(\{LVAs\ sax\ EPI\} / 3.14159) - \sqrt{(\{LVAs\ sax\ PM\} / 3.14159)}}$ Needs measurement: LVAs sax EPI [Mass], LVAs sax PM [Mass], LVLs apical [Mass] Measured by: LVMass(s) [AUTOCALC] |
| LVs Mass Index [Dimension] Mode: 2D:CF Formula: $m/\{BSA\}$ where: $m = ((1.04 \times ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - \{LVIDs\}^3)) - 13.6) / 1000$ Needs measurement: IVSs [Dimension], LVIDs [Dimension], LVPWs [Dimension] Measured by: LVPWs [2DCALIPER] |
| LVs Mass Index [Generic, Dimension] Mode: MM:CM:AMM Formula: $((1.04 \times ((\{IVSs\} + \{LVIDs\} + \{LVPWs\})^3 - \{LVIDs\}^3)) - 13.6) / 1000 / \{BSA\}$ Needs measurement: IVSs [Generic, Dimension], LVIDs [Generic, Dimension], LVPWs [Generic, Dimension] Measured by: LV Study [MMLV], LVPWs [MMDISCALIPER] |
| LVSI Dopp [Aortic] Mode: PW Formula: $\{LVOT\ Diam\}^2 \times 0.785 \times \{LVOT\ VTI\} / \{BSA\}$ Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic], Measured by: LVOT Trace [SDMANTRACE] |
| LVSV Dopp [Aortic] Mode: PW Formula: $\{LVOT\ Diam\}^2 \times 0.785 \times \{LVOT\ VTI\}$ Needs measurement: LVOT Diam [Aortic], LVOT VTI [Aortic] Measured by: LVOT Trace [SDMANTRACE] |
| MP Area [Mitral Valve] Mode: CW:PW Formula: $\{LVOT\ Diam\}^2 \times 0.785 \times (\{LVOT\ VTI\} / \{MP\ VTI\})$ Needs measurement: LVOT Diam [Mitral Valve], LVOT VTI [Mitral Valve], MP VTI [Mitral Valve] Measured by: MP Area [SDMANTRACE] |
| MR ERO [PISA] Mode: CF:CW:PW Formula: $\{MR\ Flow\} / \{MR\ Vmax\}$ Needs measurement: MR Flow [PISA], MR Vmax [PISA] Measured by: MR Trace [AUTOCALC] |
| MR RV [PISA] Mode: CF:CW:PW Formula: $\{MR\ Flow\} / \{MR\ Vmax\} \times \{MR\ VTI\}$ Needs measurement: MR Flow [PISA], MR Vmax [PISA], MR VTI [PISA] Measured by: MR Trace [AUTOCALC] |

| |
|--|
| MV AccT/DecT [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV AccT}\} / \{\text{MV DecT}\}$ Needs measurement: MV AccT [Mitral Valve], MV DecT [Mitral Valve] Measured by: MV AccT [SDCALIPER] |
| MV Area [Dimension] Mode: 2D:CF Formula: $3.14/4 * \{\text{MV Ann Diam}\}^2$ Needs measurement: MV Ann Diam [Dimension] Measured by: MV Ann Diam [2DCALIPER] |
| MV CI [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV Ann Diam}\}^2 * 0.785 * \{\text{MV VTI}\} * \{\text{HR}\} / 60 / \{\text{BSA}\}$ Needs measurement: MV Ann Diam [Mitral Valve], MV VTI [Mitral Valve], HR [Mitral Valve] Measured by: MV Trace [SDMANTRACE] |
| MV CO [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV Ann Diam}\}^2 * 0.785 * \{\text{MV VTI}\} * \{\text{HR}\} / 60$ Needs measurement: MV Ann Diam [Mitral Valve], MV VTI [Mitral Valve], HR [Mitral Valve] Measured by: MV Trace [SDMANTRACE] |
| MV E/A Ratio [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV E Vel}\} / \{\text{MV A Vel}\}$ Needs measurement: MV E Vel [Mitral Valve], MV A Vel [Mitral Valve] Measured by: MV A Vel [SDPTCALIPER], MV A Vel [AUTOCALC] |
| MV SI [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV Ann Diam}\}^2 * 0.785 * \{\text{MV VTI}\} / \{\text{BSA}\}$ Needs measurement: MV Ann Diam [Mitral Valve], MV VTI [Mitral Valve], Measured by: MV Trace [SDMANTRACE] |
| MV SV [Mitral Valve] Mode: CW:PW Formula: $\{\text{MV Ann Diam}\}^2 * 0.785 * \{\text{MV VTI}\}$ Needs measurement: MV Ann Diam [Mitral Valve], MV VTI [Mitral Valve] Measured by: MV Trace [SDMANTRACE] |
| MVA (VTI) [Mitral Valve] Mode: 2D:CW:PW Formula: $3.14/4 * \{\text{LVOT Diam}\}^2 * \{\text{LVOT VTI}\} / \{\text{MV VTI}\}$ Needs measurement: LVOT Diam [Mitral Valve], LVOT VTI [Mitral Valve], MV VTI [Mitral Valve] Measured by: MV Trace [AUTOCALC] |

Measurements

| |
|--|
| MVA By PHT [Mitral Valve] Mode: CW:PW Formula: $22/(\{MV\ PHT\})$ Needs measurement: MV PHT [Mitral Valve] Measured by: MV E/A Velocity [SDEA3], MV PHT [SDCALIPER] |
| P Vein S/D Ratio [Pulmonary Vein] Mode: PW Formula: $\{P\ Vein\ S\}/\{P\ Vein\ D\}$ Needs measurement: P Vein S [Pulmonary Vein], P Vein D [Pulmonary Vein] Measured by: P Vein D [SDPTCALIPER] |
| PAEDP [Pulmonic] Mode: CW:PW Formula: $\{PRend\ PG\} + \{RAP\}$ Needs measurement: PRend PG [Pulmonic], RAP [Pulmonic] Measured by: PRend Vmax [AUTOCALC] |
| PR ERO [PISA] Mode: CF:CW:PW Formula: $\{PR\ Flow\}/\{PR\ Vmax\}$ Needs measurement: PR Flow [PISA], PR Vmax [PISA] Measured by: PR Trace [AUTOCALC] |
| PR RV [PISA] Mode: CF:CW:PW Formula: $\{PR\ Flow\}/\{PR\ Vmax\} * \{PR\ VTI\}$ Needs measurement: PR Flow [PISA], PR Vmax [PISA], PR VTI [PISA] Measured by: PR Trace [AUTOCALC] |
| Pulmonic CO [Shunts, Congenital Heart] Mode: CW:PW Formula: $\{Pulmonic\ SV\} * \{Pulmonic\ HR\} / 60$ Needs measurement: Pulmonic SV [Shunts, Congenital Heart], Pulmonic HR [Shunts, Congenital Heart] Measured by: Pulmonic VTI [SDMANTRACE] |
| Pulmonic SV [Shunts, Congenital Heart] Mode: CW:PW Formula: $3.14159/4 * \{Pulmonic\ Diam\}^2 * \{Pulmonic\ VTI\}$ Needs measurement: Pulmonic Diam [Shunts, Congenital Heart], Pulmonic VTI [Shunts, Congenital Heart] Measured by: Pulmonic VTI [SDMANTRACE], Pulmonic VTI [SDMANTRACE] Used to calculate: Pulmonic CO |
| PV A/MV A Dur [Pulmonary Vein] Mode: PW Formula: $\{P\ Vein\ A\ Dur\}/\{MV\ A\ Dur\}$ Needs measurement: P Vein A Dur [Pulmonary Vein], MV A Dur [Pulmonary Vein] Measured by: P Vein A Dur [SDTIMECALIPER] |

| |
|--|
| PV A/MV VTI [Pulmonary Vein] Mode: PW Formula: $\{P \text{ Vein A Dur}\} / \{MV \text{ VTI}\}$ Needs measurement: P Vein A Dur [Pulmonary Vein], MV VTI [Pulmonary Vein] Measured by: P Vein A Dur [SDTIMECALIPER] |
| PV AccT/ET [Pulmonic] Mode: CW:PW Formula: $\{PV \text{ AccT}\} / \{PVET\}$ Needs measurement: PV AccT [Pulmonic], PVET [Pulmonic] Measured by: PVET [SDTIMECALIPER] |
| PV A-MV A Dur [Pulmonary Vein] Mode: PW Formula: $\{P \text{ Vein A Dur}\} - \{MV \text{ A Dur}\}$ Needs measurement: P Vein A Dur [Pulmonary Vein], MV A Dur [Pulmonary Vein] Measured by: P Vein A Dur [SDTIMECALIPER] |
| PV Area [Dimension] Mode: 2D:CF Formula: $3.14/4 * \{PV \text{ Ann Diam}\}^2$ Needs measurement: PV Ann Diam [Dimension] Measured by: PV Ann Diam [2DCALIPER] |
| PV CI [Pulmonic, Valvular PS] Mode: CW:PW Formula: $\{(\{PV \text{ Ann Diam}\}^2 * 0.785 * \{PV \text{ VTI}\}) * \{HR\} / 60\} / \{BSA\}$ Needs measurement: PV Ann Diam [Pulmonic, Valvular PS], PV VTI [Pulmonic, Valvular PS], HR [Pulmonic, Valvular PS] Measured by: PV Trace [SDMANTRACE] |
| PV CO [Pulmonic, Valvular PS] Mode: CW:PW Formula: $\{(\{PV \text{ Ann Diam}\}^2 * 0.785 * \{PV \text{ VTI}\}) * \{HR\} / 60\}$ Needs measurement: PV Ann Diam [Pulmonic, Valvular PS], PV VTI [Pulmonic, Valvular PS], HR [Pulmonic, Valvular PS] Measured by: PV Trace [SDMANTRACE] |
| PV SI [Pulmonic, Valvular PS] Mode: CW:PW Formula: $\{(\{PV \text{ Ann Diam}\}^2 * 0.785 * \{PV \text{ VTI}\}) / \{BSA\}\}$ Needs measurement: PV Ann Diam [Pulmonic, Valvular PS], PV VTI [Pulmonic, Valvular PS] Measured by: PV Trace [SDMANTRACE] |
| PV SV [Pulmonic, Valvular PS] Mode: CW:PW Formula: $\{PV \text{ Ann Diam}\}^2 * 0.785 * \{PV \text{ VTI}\}$ Needs measurement: PV Ann Diam [Pulmonic, Valvular PS], PV VTI [Pulmonic, Valvular PS] Measured by: PV Trace [SDMANTRACE] |

Measurements

| |
|--|
| <p>PVA (Vmax) [Pulmonic]</p> <p>Mode: 2D:CW:PW</p> <p>Formula: $3.14/4 \times \{RVOT\ Diam\}^2 \times \{RVOT\ Vmax\} / \{PV\ Vmax\}$</p> <p>Needs measurement: RVOT Diam [Pulmonic], RVOT Vmax [Pulmonic], PV Vmax [Pulmonic]</p> <p>Measured by: PV Vmax [AUTOCALC]</p> |
| <p>PVA (Vmax) [Pulmonic]</p> <p>Mode: 2D:CW:PW</p> <p>Formula: $3.14/4 \times \{RVOT\ Diam\}^2 \times \{RVOT\ Vmax\} / \{PV\ Vmax\}$</p> <p>Needs measurement: RVOT Diam [Pulmonic], RVOT Vmax [Pulmonic], PV Vmax [Pulmonic]</p> <p>Measured by: PV Trace [AUTOCALC]</p> |
| <p>PVA (VTI) [Pulmonic]</p> <p>Mode: 2D:CW:PW</p> <p>Formula: $3.14/4 \times \{RVOT\ Diam\}^2 \times \{RVOT\ VTI\} / \{PV\ VTI\}$</p> <p>Needs measurement: RVOT Diam [Pulmonic], RVOT VTI [Pulmonic], PV VTI [Pulmonic]</p> <p>Measured by: PV Trace [AUTOCALC]</p> |
| <p>Qp/Qs [Shunts, Congenital Heart]</p> <p>Mode: CW:PW</p> <p>Formula: $3.14159/4 \times \{Pulmonic\ Diam\}^2 \times \{Pulmonic\ VTI\} / (3.14159/4 \times \{Systemic\ Diam\}^2 \times \{Systemic\ VTI\})$</p> <p>Needs measurement: Pulmonic Diam [Shunts, Congenital Heart], Pulmonic VTI [Shunts, Congenital Heart], Systemic Diam [Shunts, Congenital Heart], Systemic VTI [Shunts, Congenital Heart]</p> <p>Measured by: Qp/Qs [AUTOCALC]</p> |
| <p>RIMP [Pulmonic, Tricuspid Valve]</p> <p>Mode: CW:PW</p> <p>Formula: $(\{TCO\} - \{PVET\}) / \{PVET\}$</p> <p>Needs measurement: TCO [Pulmonic, Tricuspid Valve], PVET [Pulmonic, Tricuspid Valve], PVET [Pulmonic, Tricuspid Valve]</p> <p>Measured by: RIMP [AUTOCALC]</p> |
| <p>RVOT Area [Dimension]</p> <p>Mode: 2D:CF</p> <p>Formula: $3.14/4 \times \{RVOT\ Diam\}^2$</p> <p>Needs measurement: RVOT Diam [Dimension]</p> <p>Measured by: RVOT Diam [2DCALIPER]</p> |
| <p>RVOT CI [Pulmonic, Valvular PS]</p> <p>Mode: PW</p> <p>Formula: $((\{RVOT\ Diam\}^2 \times 0.785 \times \{RVOT\ VTI\}) \times \{HR\} / 60) / \{BSA\}$</p> <p>Needs measurement: RVOT Diam [Pulmonic, Valvular PS], RVOT VTI [Pulmonic, Valvular PS], HR [Pulmonic, Valvular PS],</p> <p>Measured by: RVOT Trace [SDMANTRACE]</p> |

RVOT CO [Pulmonic, Valvular PS]

Mode: PW

Formula: $\{(\text{RVOT Diam})^2 \times 0.785 \times \{\text{RVOT VTI}\} \times \{\text{HR}\} / 60$

Needs measurement: RVOT Diam [Pulmonic, Valvular PS], RVOT VTI [Pulmonic, Valvular PS], HR [Pulmonic, Valvular PS]

Measured by: RVOT Trace [SDMANTRACE]

RVOT SI [Pulmonic, Valvular PS]

Mode: PW

Formula: $\{(\text{RVOT Diam})^2 \times 0.785 \times \{\text{RVOT VTI}\} / \{\text{BSA}\}$

Needs measurement: RVOT Diam [Pulmonic, Valvular PS], RVOT VTI [Pulmonic, Valvular PS],

Measured by: RVOT Trace [SDMANTRACE]

RVOT SV [Pulmonic, Valvular PS]

Mode: PW

Formula: $\{\text{RVOT Diam}\}^2 \times 0.785 \times \{\text{RVOT VTI}\}$

Needs measurement: RVOT Diam [Pulmonic, Valvular PS], RVOT VTI [Pulmonic, Valvular PS]

Measured by: RVOT Trace [SDMANTRACE]

RVPEP/ET [Pulmonic]

Mode: CW:PW

Formula: $\{\text{RVPEP}\} / \{\text{RVET}\}$

Needs measurement: RVPEP [Pulmonic], RVET [Pulmonic]

Measured by: RVET [SDTIMECALIPER]

RVPEP/ET [Time]

Mode: MM:CM:AMM

Formula: $\{\text{RVPEP}\} / \{\text{RVET}\}$

Needs measurement: RVPEP [Time], RVET [Time]

Measured by: RVET [MMTIMECALIPER]

RVSP [Tricuspid Valve]

Mode: CW:PW

Formula: $\{\text{TR maxPG}\} + \{\text{RAP}\}$

Needs measurement: TR maxPG [Tricuspid Valve], RAP [Tricuspid Valve]

Measured by: TR Vmax [AUTOCALC]

SI A-L A2C [Biplane, Single Plane A2C, AutoBiplane]

Mode: 2D:CF

Formula: $\{(\text{LVEDV A-L A2C}) - \{\text{LVESV A-L A2C}\} / \{\text{BSA}\}$

Needs measurement: LVEDV A-L A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV A-L A2C [Biplane, Single Plane A2C, AutoBiplane]

Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME]

Measurements

| |
|---|
| SI A-L A4C [Biplane, Single Plane A4C, AutoBiplane] Mode: 2D:CF Formula: $\frac{(\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\})}{\{BSA\}}$ Needs measurement: LVEDV A-L A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV A-L A4C [Biplane, Single Plane A4C, AutoBiplane] Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME] |
| SI A-L LAX [Single Plane LAX, AutoBiplane] Mode: 2D:CF Formula: $\frac{(\{LVEDV\ A-L\ LAX\} - \{LVESV\ A-L\ LAX\})}{\{BSA\}}$ Needs measurement: LVEDV A-L LAX [Single Plane LAX, AutoBiplane], LVESV A-L LAX [Single Plane LAX, AutoBiplane] Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME] |
| SI Biplane [Biplane, AutoBiplane] Mode: 2D:CF Formula: $d = \text{biplane}(\{LVLd\ A4C\}, \{LVDisks\}, \{LVLd\ A2C\}, \{LVDisks\})$ Needs measurement: LVLd A4C [Biplane, AutoBiplane], LVLd A2C [Biplane, AutoBiplane], LVLs A4C [Biplane, AutoBiplane], LVLs A2C [Biplane, AutoBiplane] Measured by: EF Biplane [AUTOCALC] |
| SI bp el [Biplane Ellipse] Mode: 2D:CF Formula: $(d-s)/\{BSA\}$ where: $s = (8/(3*3.14159)) * \{LVAs(A4C)\} * \{LVAs(sax\ MV)\} / \{2D/LVIDs\}$ $d = (8/(3*3.14159)) * \{LVAd\ A4C\} * \{LVAd\ (sax\ MV)\} / \{LVIDd\}$ Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse], LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse] Measured by: LVEF BP-EL [AUTOCALC] |
| SI bullet [Bullet] Mode: 2D:CF Formula: $(d-s)/\{BSA\}$ where: $s = 5/6 * \{LVAs(sax)\} * \{LVLs(apical)\}$ $d = 5/6 * \{LVAd\ sax\} * \{LVLd\ apical\}$ Needs measurement: LVAd sax [Bullet], LVLd apical [Bullet], LVLs apical [Bullet] Measured by: LVEF Bullet [AUTOCALC] |
| SI MOD A2C [Biplane, Single Plane A2C, AutoBiplane] Mode: 2D:CF Formula: $\frac{(\{LVEDV\ MOD\ A2C\} - \{LVESV\ MOD\ A2C\})}{\{BSA\}}$ Needs measurement: LVEDV MOD A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV MOD A2C [Biplane, Single Plane A2C, AutoBiplane] Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME] |
| SI MOD A4C [Biplane, Single Plane A4C, AutoBiplane] Mode: 2D:CF Formula: $\frac{(\{LVEDV\ MOD\ A4C\} - \{LVESV\ MOD\ A4C\})}{\{BSA\}}$ Needs measurement: LVEDV MOD A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV MOD A4C [Biplane, Single Plane A4C, AutoBiplane] Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME] |

SI MOD LAX [Single Plane LAX, AutoBiplane]

Mode: 2D:CF

Formula: $\frac{\{LVEDV \text{ MOD LAX}\} - \{LVESV \text{ MOD LAX}\}}{\{BSA\}}$

Needs measurement: LVEDV MOD LAX [Single Plane LAX, AutoBiplane], LVESV MOD LAX [Single Plane LAX, AutoBiplane]

Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME]

SI mod sim [Modified Simpson]

Mode: 2D:CF

Formula: $d - s / \{BSA\}$ where: $s = \frac{\{LVLs(apical)\}}{9} * ((4 * \{LVAs(sax \text{ MV})\}) + (2 * \{LVAs(sax \text{ PM})\}) + \sqrt{\{LVAs(sax \text{ MV})\} * \{LVAs(sax \text{ PM})\}})$ $d = \frac{\{LVLd \text{ apical}\}}{9} * ((4 * \{LVAd(sax \text{ MV})\}) + (2 * \{LVAd(sax \text{ PM})\}) + \sqrt{\{LVAd(sax \text{ MV})\} * \{LVAd(sax \text{ PM})\}})$

Needs measurement: LVLd apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson], LVLs apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson]

Measured by: EF mod sim [AUTOCALC]

SI(A-L) [Generic]

Mode: 2D:CF

Formula: $\frac{\{EDV(A-L)\} - \{ESV(A-L)\}}{\{BSA\}}$

Needs measurement: ESV(A-L) [Generic]

Measured by: EF Volume [AUTOCALC]

SI(Cube) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: $(d - s) / \{BSA\}$ where: $s = \{2D/LVIDs\}^3$ $d = \{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz]

Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]

SI(Cube) [Generic]

Mode: MM:CM:AMM

Formula: $(dv - sv) / \{BSA\}$ where: $sv = \{MM/LVIDs\}^3$ $dv = \{LVIDd\}^3$

Needs measurement: LVIDd [Generic], LVIDs [Generic],

Measured by: LV Study [MMLV]

SI(MOD) [Generic]

Mode: 2D:CF

Formula: $\frac{\{EDV(MOD)\} - \{ESV(MOD)\}}{\{BSA\}}$

Needs measurement: EDV(MOD) [Generic], ESV(MOD) [Generic]

Measured by: EF Volume [AUTOCALC]

SI(Teich) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: $(d - s) / \{BSA\}$ $s = 7 / (2.4 + \{2D/LVIDs\}) * \{2D/LVIDs\}^3$ $d = 7 / (2.4 + \{LVIDd\}) * \{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz]

Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]

Measurements

| |
|--|
| <p>SI(Teich) [Generic]</p> <p>Mode: MM:CM:AMM</p> <p>Formula: $(dv-sv)/\{BSA\}$ where: $sv = 7/(2.4+\{MM/LVIDs\})^3$ $dv = 7/(2.4+\{LVIDd\})^3$</p> <p>Needs measurement: LVIDd [Generic], LVIDd [Generic], LVIDs [Generic]</p> <p>Measured by: LV Study [MMLV]</p> |
| <p>SV A-L A2C [Biplane, Single Plane A2C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $\{LVEDV\ A-L\ A2C\} - \{LVESV\ A-L\ A2C\}$</p> <p>Needs measurement: LVEDV A-L A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV A-L A2C [Biplane, Single Plane A2C, AutoBiplane]</p> <p>Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME]</p> |
| <p>SV A-L A4C [Biplane, Single Plane A4C, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $\{LVEDV\ A-L\ A4C\} - \{LVESV\ A-L\ A4C\}$</p> <p>Needs measurement: LVEDV A-L A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV A-L A4C [Biplane, Single Plane A4C, AutoBiplane]</p> <p>Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME]</p> |
| <p>SV A-L LAX [Single Plane LAX, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $\{LVEDV\ A-L\ LAX\} - \{LVESV\ A-L\ LAX\}$</p> <p>Needs measurement: LVEDV A-L LAX [Single Plane LAX, AutoBiplane], LVESV A-L LAX [Single Plane LAX, AutoBiplane]</p> <p>Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME]</p> |
| <p>SV Biplane [Biplane, AutoBiplane]</p> <p>Mode: 2D:CF</p> <p>Formula: $d = \text{biplane}(\{LVLd\ A4C\}, \{LVDisk\}, \{LVLd\ A2C\}, \{LVDisk\})$</p> <p>Needs measurement: LVLd A4C [Biplane, AutoBiplane], LVLd A2C [Biplane, AutoBiplane], LVLs A4C [Biplane, AutoBiplane], LVLs A2C [Biplane, AutoBiplane]</p> <p>Measured by: EF Biplane [AUTOCALC]</p> |
| <p>SV bp el [Biplane Ellipse]</p> <p>Mode: 2D:CF</p> <p>Formula: $d = (8/(3 \times 3.14159)) \times \{LVAd\ A4C\} \times \{LVAd\ (sax\ MV)\} / \{LVIDd\}$</p> <p>Needs measurement: LVAd A4C [Biplane Ellipse], LVAd (sax MV) [Biplane Ellipse], LVIDd [Biplane Ellipse], LVAs A4C [Biplane Ellipse], LVAs sax MV [Biplane Ellipse], LVIDs [Biplane Ellipse]</p> <p>Measured by: LVEF BP-EL [AUTOCALC]</p> |
| <p>SV bullet [Bullet]</p> <p>Mode: 2D:CF</p> <p>Formula: $d-s$ where: $s = 5/6 \times \{LVAs(sax)\} \times \{LVLs(apical)\}$ $d = 5/6 \times \{LVAd\ sax\} \times \{LVLd\ apical\}$</p> <p>Needs measurement: LVAd sax [Bullet], LVLd apical [Bullet], LVLs apical [Bullet]</p> <p>Measured by: LVEF Bullet [AUTOCALC]</p> |

SV MOD A2C [Biplane, Single Plane A2C, AutoBiplane]

Mode: 2D:CF

Formula: {LVEDV MOD A2C}-{LVESV MOD A2C}

Needs measurement: LVEDV MOD A2C [Biplane, Single Plane A2C, AutoBiplane], LVESV MOD A2C [Biplane, Single Plane A2C, AutoBiplane]

Measured by: EF SP A2C [AUTOCALC], LVESV A2C [2DVOLUMETRACE], A2C [2DAUTOVOLUME]

SV MOD A4C [Biplane, Single Plane A4C, AutoBiplane]

Mode: 2D:CF

Formula: {LVEDV MOD A4C}-{LVESV MOD A4C}

Needs measurement: LVEDV MOD A4C [Biplane, Single Plane A4C, AutoBiplane], LVESV MOD A4C [Biplane, Single Plane A4C, AutoBiplane]

Measured by: EF SP A4C [AUTOCALC], LVESV A4C [2DVOLUMETRACE], A4C [2DAUTOVOLUME]

SV MOD LAX [Single Plane LAX, AutoBiplane]

Mode: 2D:CF

Formula: {LVEDV MOD LAX}-{LVESV MOD LAX}

Needs measurement: LVEDV MOD LAX [Single Plane LAX, AutoBiplane], LVESV MOD LAX [Single Plane LAX, AutoBiplane]

Measured by: LVESV LAX [2DVOLUMETRACE], EF SP LAX [AUTOCALC], AutoVolume [2DAUTOVOLUME]

SV mod sim [Modified Simpson]

Mode: 2D:CF

Formula: $((\{LVLd\ apical\}/9)*((4*\{LVAd\ (sax\ MV)\})+(2*\{LVAd\ sax\ PM\}))+sqrt(\{LVAd\ (sax\ MV)\}*\{LVAd\ sax\ PM\}))$

Needs measurement: LVLd apical [Modified Simpson], LVAd (sax MV) [Modified Simpson], LVAd sax PM [Modified Simpson], LVLs apical [Modified Simpson], LVAs sax MV [Modified Simpson], LVAs sax PM [Modified Simpson]

Measured by: EF mod sim [AUTOCALC]

SV(A-L) [Generic]

Mode: 2D:CF

Formula: {EDV(A-L)}-{ESV(A-L)}

Needs measurement: ESV(A-L) [Generic]

Measured by: EF Volume [AUTOCALC]

SV(Cube) [Dimension, Cube/Teicholz]

Mode: 2D:CF

Formula: d-s where: $s = \{2D/LVIDs\}^3$ $d = \{LVIDd\}^3$

Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz]

Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC]

SV(Cube) [Generic, Dimension]

Mode: MM:CM:AMM

Formula: dv-sv where: $sv = \{MM/LVIDs\}^3$ $dv = \{LVIDd\}^3$

Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension]

Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER]

Measurements

| |
|--|
| SV(MOD) [Generic] Mode: 2D:CF Formula: {EDV(MOD)}-{ESV(MOD)} Needs measurement: EDV(MOD) [Generic], ESV(MOD) [Generic] Measured by: EF Volume [AUTOCALC] |
| SV(Teich) [Dimension, Cube/Teicholz] Mode: 2D:CF Formula: d-s where: $s = 7 / (2.4 + \{2D/LVIDs\} * \{2D/LVIDs\}^3)$ $d = 7 / (2.4 + \{LVIDd\} * \{LVIDd\}^3)$ Needs measurement: LVIDd [Dimension, Cube/Teicholz], LVIDs [Dimension, Cube/Teicholz] Measured by: LVs [2DLV], LVIDs [2DCALIPER], EF(Cube) [AUTOCALC] |
| SV(Teich) [Generic, Dimension] Mode: MM:CM:AMM Formula: dv-sv where: $sv = 7 / (2.4 + \{MM/LVIDs\} * \{MM/LVIDs\}^3)$ $dv = 7 / (2.4 + \{LVIDd\} * \{LVIDd\}^3)$ Needs measurement: LVIDd [Generic, Dimension], LVIDs [Generic, Dimension] Measured by: LV Study [MMLV], LVIDs [MMDISCALIPER] |
| Systemic CO [Shunts, Congenital Heart] Mode: CW:PW Formula: {Systemic SV}*{Systemic HR}/60 Needs measurement: Systemic SV [Shunts, Congenital Heart], Systemic HR [Shunts, Congenital Heart] Measured by: Systemic VTI [SDMANTRACE] |
| Systemic SV [Shunts, Congenital Heart] Mode: CW:PW Formula: $3.14159/4 * \{Systemic Diam\}^2 * \{Systemic VTI\}$ Needs measurement: Systemic Diam [Shunts, Congenital Heart], Systemic VTI [Shunts, Congenital Heart] Measured by: Systemic VTI [SDMANTRACE], Systemic VTI [SDMANTRACE] Used to calculate: Systemic CO |
| TR ERO [PISA] Mode: CF:CW:PW Formula: {TR Flow}/{TR Vmax} Needs measurement: TR Flow [PISA], TR Vmax [PISA] Measured by: TR Trace [AUTOCALC] |
| TR RV [PISA] Mode: CF:CW:PW Formula: {TR Flow}/{TR Vmax}*{TR VTI} Needs measurement: TR Flow [PISA], TR Vmax [PISA], TR VTI [PISA] Measured by: TR Trace [AUTOCALC] |
| TV AccT/DecT [Tricuspid Valve] Mode: CW:PW Formula: {TV AccT}/{TV Dec Time} Needs measurement: TV AccT [Tricuspid Valve], TV Dec Time [Tricuspid Valve] Measured by: TV AccT [SDCALIPER] |

TV Area [Dimension]

Mode: 2D:CF

Formula: $3.14/4 \times \{TV \text{ Ann Diam}\}^2$

Needs measurement: TV Ann Diam [Dimension]

Measured by: TV Ann Diam [2DCALIPER]

TV CI [Tricuspid Valve]

Mode: CW:PW

Formula: $((\{TV \text{ Ann Diam}\}^2 \times 0.785 \times \{TV \text{ VTI}\}) \times \{HR\} / 60) / \{BSA\}$

Needs measurement: TV Ann Diam [Tricuspid Valve], TV VTI [Tricuspid Valve], HR [Tricuspid Valve]

Measured by: TV Trace [SDMANTRACE]

TV CO [Tricuspid Valve]

Mode: CW:PW

Formula: $(\{TV \text{ Ann Diam}\}^2 \times 0.785 \times \{TV \text{ VTI}\}) \times \{HR\} / 60$

Needs measurement: TV Ann Diam [Tricuspid Valve], TV VTI [Tricuspid Valve], HR [Tricuspid Valve]

Measured by: TV Trace [SDMANTRACE]

TV E/A Ratio [Tricuspid Valve]

Mode: CW:PW

Formula: $\{TV \text{ E Vel}\} / \{TV \text{ A Vel}\}$

Needs measurement: TV E Vel [Tricuspid Valve], TV A Vel [Tricuspid Valve]

Measured by: TV A Vel [SDPTCALIPER]

TV SI [Tricuspid Valve]

Mode: CW:PW

Formula: $(\{TV \text{ Ann Diam}\}^2 \times 0.785 \times \{TV \text{ VTI}\}) / \{BSA\}$

Needs measurement: TV Ann Diam [Tricuspid Valve], TV VTI [Tricuspid Valve]

Measured by: TV Trace [SDMANTRACE]

TV SV [Tricuspid Valve]

Mode: CW:PW

Formula: $\{TV \text{ Ann Diam}\}^2 \times 0.785 \times \{TV \text{ VTI}\}$

Needs measurement: TV Ann Diam [Tricuspid Valve], TV VTI [Tricuspid Valve]

Measured by: TV Trace [SDMANTRACE]

TVA (Vmax) [Tricuspid Valve]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{RVOT \text{ Diam}\}^2 \times \{RVOT \text{ Vmax}\} / \{TV \text{ Vmax}\}$

Needs measurement: RVOT Diam [Tricuspid Valve], RVOT Vmax [Tricuspid Valve], TV Vmax [Tricuspid Valve]

Measured by: TV Vmax [AUTOCALC]

TVA (Vmax) [Tricuspid Valve]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{RVOT \text{ Diam}\}^2 \times \{RVOT \text{ Vmax}\} / \{TV \text{ Vmax}\}$

Needs measurement: RVOT Diam [Tricuspid Valve], RVOT Vmax [Tricuspid Valve], TV Vmax [Tricuspid Valve]

Measured by: TV Trace [AUTOCALC]

Measurements

TVA (VTI) [Tricuspid Valve]

Mode: 2D:CW:PW

Formula: $3.14/4 \times \{RVOT\ Diam\}^2 \times \{RVOT\ VTI\} / \{TV\ VTI\}$

Needs measurement: RVOT Diam [Tricuspid Valve], RVOT VTI [Tricuspid Valve], TV VTI [Tricuspid Valve]

Measured by: TV [AUTOCALC]

Vcf mean [Dimension]

Mode: MM:CM:AMM

Formula: $(\{LVIDd\} - \{LVIDs\}) / (\{LVIDd\} \times \{LVET\})$

Needs measurement: LVIDd [Dimension], LVIDs [Dimension], LVET [Dimension]

Measured by: Vcf [MMTIMECALIPER]

Vcf mn (corr) [Dimension]

Mode: MM:CM:AMM

Formula: $(\{LVIDd\} - \{LVIDs\}) / (\{LVIDd\} \times \{LVET\} / \sqrt{\{Time\}})$

Needs measurement: LVIDd [Dimension], LVIDs [Dimension], LVET [Dimension], Time [Dimension]

Measured by: Vcf [MMTIMECALIPER]

Formulas-Generic

| Calc Mnemonic | Calc Name | Input Measurements | Formula |
|---------------|---|--|---|
| BSA | Body Surface Area | Patient weight (kg) and height (m) | $0.007184 \times \text{Weight}^{0.425} \times \text{Height}^{0.725}$ |
| BSA | Body Surface Area | Patient weight (kg) | $0.1 \times \text{Weight}^{0.667}$ |
| MaxPG | Maximum Pressure Gradient | two Doppler blood flow peak velocities | $\text{MaxPG}[\text{mmHg}] = 4 \times (v_1^2 - v_2^2)$ |
| MeanPG | Mean Pressure Gradient | flow velocities from one time marker to another time marker in a Doppler display | $\text{MeanPG}[\text{mmHg}] = \frac{4 \times \sum_{i=1}^n (V_i^2)}{n}$ |
| % Stenosis | Stenosis Ratio | two areas (by ellipse, trace, circle or distance) | $\% \text{ Stenosis} = [1 - (A_{\text{residual}} / A_{\text{lumen}})] \times 100$ |
| PI | Pulsatility Index | two Doppler blood flow peak velocities and TAMAX | $PI = (V_{\text{max}} - V_{\text{diastole}}) / \text{TAMAX}^a$ |
| RI | Resistivity Index | two Doppler blood flow peak velocities | $RI = (V_{\text{max}} - V_{\text{diastole}}) / V_{\text{max}}^a$ |
| HR | Heart Rate (beats/minute) | one 2 beat time interval | $HR[\text{BPM}] = 120[\text{sec}] / 2\text{beat time}[\text{sec}]$ |
| A/B Ratio | Velocities Ratio | two Doppler blood flow peak velocities | $A/B = V_1 / V_2$ |
| TAMAX | Time Averaged Maximum Velocity (Trace Method is Peak or manual) | two time marks in a Doppler display | $\text{TAMAX} = \sum \{V_t\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$ |
| TAMIN | Time Averaged Minimum Velocity (Trace method is Floor) | two time marks in a Doppler display | $\text{TAMIN} = \sum \{V_b\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$ |
| TAMEAN | Time Averaged Mean Velocity (Trace method is Mean) | two time marks in a Doppler display | $\text{TAMEAN} = \sum \{V_b\} \text{ from } t_1 \text{ to } t_2 / (t_2 - t_1) [\text{cm/s or m/s}]$ |

a. Vdiastole = Vmin or Vmid-diastole (depends on preset selection)

| Calc Name | Input Measurements | Formula |
|-----------------------------|---|---|
| Volume (spherical) | one distance | $\text{Vol}[\text{ml}] = (\pi/6) \times d^3$ |
| Volume (prolate spheroidal) | two distances, $d_1 > d_2$ | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2^2$ |
| Volume (prolate spheroidal) | one ellipse, d_1 major axis, d_2 minor axis | $\text{Vol}[\text{ml}] = (\pi/6) \times d_1 \times d_2^2$ |

Measurements

| Calc Name | Input Measurements | Formula |
|---------------------|--|---|
| Volume (spheroidal) | three distances | $Vol[ml] = (\pi/6) \times d1 \times d2 \times d3$ |
| Volume (spheroidal) | one distance d1, one ellipse, d2 major axis, d3 minor axis | $Vol[ml] = (\pi/6) \times d1 \times d2 \times d3$ |
| Volume Flow | two time marks in a Doppler display, and one distance | $Volume\ Flow = 60 \times T_{Amean} \times V_{Fdiam} \times V_{Fdiam} \times (\pi/4)$ |

Formulas-Vascular

Vascular Calculation Formulas

| Calc Mnemonic | Calc Name | Input Measurements | Formula |
|---|---|---|---|
| RT ECA | Right External Carotid Artery Velocity | one Doppler blood flow peak velocity | RT ECA=v1 [cm/s or m/s] |
| RT CCA | Right Common Carotid Artery Velocity | one Doppler blood flow peak velocity | RT CCA=v1 [cm/s or m/s] |
| RT BIFURC | Right Carotid Bifurcation Velocity | one Doppler blood flow peak velocity | RT BIFURC=v1 [cm/s or m/s] |
| RT ICA | Right Internal Carotid Artery Velocity | one Doppler blood flow peak velocity | RT ICA=v1 [cm/s or m/s] |
| RT ICA/CCA | Right Internal Carotid Artery Velocity/Common Carotid Artery Velocity Ratio | two Doppler blood flow peak velocities | RT ICA/CCA= V_{ICA}/V_{CCA} |
| LT ECA, LT CCA, LT BIFURC, LT ICA, LT ICA/CCA | Same as above, for Left Carotid Artery | Same as above | Same as above |
| A/B Ratio | Velocities Ratio | two Doppler blood flow peak velocities | $A/B = V_1/V_2$ |
| % Stenosis | Stenosis Ratio | two areas (by ellipse, trace, circle or distance) | % Stenosis= $[1 - (A_{residual}/A_{lumen})] \times 100$ |
| S/D Ratio | Systolic Velocity/Diastolic Velocities Ratio | two Doppler blood flow peak velocities | $S/D = V_{systole}/V_{diastole}^a$ |
| PI | Pulsatility Index | two Doppler blood flow peak velocities and TAMAX | $PI = (V_{max} - V_{diastole})/TAMAX^a$ |
| RI | Resistivity Index | two Doppler blood flow peak velocities | $RI = (V_{max} - V_{diastole})/V_{max}^a$ |
| HR | Heart Rate (beats/minute) | one 2 beat time interval (measured manually or automatically) | HR[BPM]=120[sec]/2 beat time[sec] |

a. $V_{diastole} = V_{min}$ or $V_{end-diastole}$ (depends on preset selection)

Formulas-OB

OB Calculation Formulas

| Calc Mnemonic | Calc Name | Input Measurements | Formula | Author Reference |
|---------------|-------------------------|--|---|---|
| AC | Abdominal Circumference | circumference by trace, ellipse, circle or two distances | $AC = 13.3 + 1.61 (GA) - 0.00998 (GA)^2$ | Hadlock et al, Radiology, 152:497-501, 1984 |
| BPD | Biparietal Diameter | one distance | $BPD = -3.08 + 0.41 (GA) - 0.000061 (GA)^3$ | |
| CRL | Crown Rump Length | one distance | $CRL = 1.684969 + 0.315646 \times d1 + 0.049306 \times d1^2 + 0.004057 \times d1^3 + 0.000120456 \times d1^4$ | |
| FL | Femur Length | one distance | $FL = -3.91 + 0.427 (GA) - 0.0034 (GA)^2$ | |
| GS | Gestational Sac | three distances | $GS [wk] = 1.42450142 * (d1 + d2 + d3) / 3 + 3.6225$ | |
| HC | Head Circumference | circumference by trace, ellipse, circle or two distances | $HC = -11.48 + 1.56 (GA) - 0.0002548 (GA)^3$ | Hansmann, Ultrasound Diagnosis in Obstetrics and Gynecology 438-9, 1985 |
| HC | Head Circumference | one ellipse | $HC [mm] = 2.325 * (BPD [mm]^2 + OFD [mm]^2)^{0.5}$ | |
| EF | Ejection Fraction | two distances on M-Mode (End-diastolic dimension and End-systolic dimension on M-Mode) | $EF = (1 - Ds^3 / Dd^3)$ | n/a |

CUA Hadlock Formulas

| Calc Mnemonic | Calc Name | Formula |
|--|--------------------------|--|
| CUA ^a | Composite Ultrasound Age | <ol style="list-style-type: none"> CUA (BPD) = $9.54 + 1.482 * \text{BPD} + 0.1676 * \text{BPD}^2$ CUA (HC) = $8.96 + 0.540 * \text{HC} + 0.0003 * \text{HC}^3$ CUA (AC) = $8.14 + 0.753 * \text{AC} + 0.0036 * \text{AC}^2$ CUA (FL) = $10.35 + 2.460 * \text{FL} + 0.170 * \text{FL}^2$ CUA (BPD, HC) = $10.32 + 0.009 * \text{HC}^2 + 1.3200 * \text{BPD} + 0.00012 * \text{HC}^3$ CUA (BPD, AC) = $9.57 + 0.524 * \text{AC} + 0.1220 * \text{BPD}^2$ CUA (BPD, FL) = $10.50 + 0.197 * \text{BPD} * \text{FL} + 0.9500 * \text{FL} + 0.7300 * \text{BPD}$ CUA (HC, AC) = $10.31 + 0.012 * \text{HC}^2 + 0.3850 * \text{AC}$ CUA (HC, FL) = $11.19 + 0.070 * \text{HC} * \text{FL} + 0.2630 * \text{HC}$ CUA (AC, FL) = $10.47 + 0.442 * \text{AC} + 0.3140 * \text{FL}^2 - 0.0121 * \text{FL}^3$ CUA (BPD, HC, AC) = $10.58 + 0.005 * \text{HC}^2 + 0.3635 * \text{AC} + 0.02864 * \text{BPD} * \text{AC}$ CUA (BPD, HC, FL) = $11.38 + 0.070 * \text{HC} * \text{FL} + 0.9800 * \text{BPD}$ CUA (BPD, AC, FL) = $10.61 + 0.175 * \text{BPD} * \text{FL} + 0.2970 * \text{AC} + 0.7100 * \text{FL}$ CUA (HC, AC, FL) = $10.33 + 0.031 * \text{HC} * \text{FL} + 0.3610 * \text{HC} + 0.0298 * \text{AC} * \text{FL}$ CUA (BPD, HC, AC, FL) = $10.85 + 0.060 * \text{HC} * \text{FL} + 0.6700 * \text{BPD} + 0.1680 * \text{AC}$ |
| Author Reference: Hadlock, Radiology, 1984 152:497-501 | | |
| ^a) Formulas are used only if Hadlock HC, FL, AC and BPD are used and CUA is selected as the preset in the CUA/AUA for Hadlock preset in the System M&A Preset Menu. If other authors are used, CUA automatically changes to AUA and an average value is displayed. | | |

Measurements

EFW Calculation Formulas


| Calc Mnemonic | Calc Name | Input Measurements | Formula | Author Reference |
|---------------|------------------------|----------------------|--|-----------------------------------|
| EFW | Estimated Fetal Weight | AC and HC | $EFW [kg] = 10^{(1.182 + (0.0273 * HC [cm]) + (0.07057 * AC [cm] - (0.00063 * AC [cm]^2) - (0.0002184 * AC [cm] * HC [cm]))}$ | Hadlock, Radiology, 150:535, 1984 |
| EFW | Estimated Fetal Weight | BPD, AC, and FL | $EFW [g] = 10^{(1.335 - (0.0034 * AC [cm] * FL [cm]) + (0.0316 * BPD [cm]) + (0.0457 * AC [cm]) + (0.1623 * FL [cm])}$ | Hadlock, AJOG, 151:333, 1985 |
| EFW | Estimated Fetal Weight | AC, HC, and FL | $EFW [g] = 10^{(1.326 - (0.00326 * AC [cm] * FL [cm]) + (0.0107 * HC [cm]) + (0.0438 * AC [cm]) + (0.158 * FL [cm]))}$ | Hadlock, AJOG, 151:333, 1985 |
| EFW | Estimated Fetal Weight | AC, HC, BPD, FL | $EFW [g] = 10^{(1.3596 - (0.00386 * AC [cm] * FL [cm]) + (0.0064 * HC [cm]) + (0.00061 * BPD [cm] * AC [cm]) + (0.0424 * AC [cm]) + (0.174 * FL [cm]))}$ | Hadlock, AJOG, 151:333, 1985 |
| EFW | Estimated Fetal Weight | AC and FL | $EFW [g] = 10^{(1.304 + 0.05281 * AC [cm] + 0.1938 * FL [cm] - 0.004 * AC [cm] * FL [cm])}$ | Hadlock, Radiology, 150:535, 1984 |
| EFW | Estimated Fetal Weight | | $EFW [g] = -3200.40479 + 157.07186 * AC [cm] + 15.90391 * BPD [cm]^2$ | Merz |
| EFW | Estimated Fetal Weight | | $EFW [g] = 0.515263 - 0.105775 * BPD [mm] + (0.000930707 * BPD [mm]^2 + 0.0649145 * TAD [mm] - 0.00020562 * TAD [mm]^2$ | German |
| EFW | Estimated Fetal Weight | AC and BPD | $EFW [kg] = 10^{(-1.7492 + 0.166 * BPD [cm] + 0.046 * AC [cm] - 2.646 * AC [cm] * BPD [cm]/1000)}$ | Shepard, AJOG, 142:47, 1982 |
| EFW | Estimated Fetal Weight | BPD [cm] and AC [cm] | $EFW [g] = 10^{(1.7288 + 0.09184 * BPD [cm] + 0.02581 * AC [cm] + 0.00011 * BPD [cm] * AC [cm])}$ | Shepard/Warsoff |

| Calc Mnemonic | Calc Name | Input Measurements | Formula | Author Reference |
|---------------|------------------------|----------------------|--|---------------------|
| EFW | Estimated Fetal Weight | BPD [cm] and AC [cm] | $EFW [g] = 10^{(3-1.7492 + (0.166 * BPD [cm]) + (0.04 * A [cm]) - (0.002646 * AC [cm] * BPD [cm]))}$ | Richards/ Berkowitz |
| EFW | Estimated Fetal Weight | | $EFW [g] = 1.07 * BPD [cm]^3 + 3.42 * APTD [cm] * TTD [cm] * FL [cm]$ | Tokyo University |
| EFW | Estimated Fetal Weight | BPD, AxT, FL [cm] | $EFW1 [g] = 1.07 * BPD [cm]^3 + 3.42 * AxT [cm]^2 * FL [cm]$ | Tokyo Shinozuka |
| EFW | Estimated Fetal Weight | BPD, AC, FL [cm] | $EFW2 [g] = 1.07 * BPD [cm]^3 + 0.30 * AC [cm]^2 * FL [cm]$ | Tokyo Shinozuka |
| EFW | Estimated Fetal Weight | BPD, AxT, LV [cm] | $EFW3 [g] = 1.07 * BPD [cm]^3 + 2.91 * AxT [cm]^2 * LV [cm]$ | Tokyo Shinozuka |

Amniotic Fluid Index (AFI)

The normal values are considered to be:


- 36-40 weeks
- 0-5 cm = very low
- 5.1-8.0 cm = low
- 8.1-18.0 cm = normal
- >18.0 = high



Sagittal

Dr. Rutherford/Dr. Phelan, *Obstetrics and Gynecology*, Volume 70, No. 3, Part 1, p.353-6, Sept. 1987.

- 28-40 weeks
- 15.0 cm = average
- >20.0 - 24.0 = hydramnios
- <5.0-6.0 = Oligohydramnios



Transverse

Dr. C.C. Smith, *The Female Patient*, Volume 15, p.85-97, March 1990.

Formulas-GYN

GYN Calculation Formulas

| Calc Mnemonic | Calc Name | Input Measurements | Formula |
|---------------|--------------------------------------|--|--|
| UT-L | Uterine Length | one distance | Ut-L[cm or mm]=d1 |
| UT-H | Uterine Height | one distance | Ut-H[cm or mm]=d1 |
| UT-W | Uterine Width | one distance | Ut-W[cm or mm]=d1 |
| UT-Volume | Uterine Volume | | |
| UtPFD | Uterus Portio-Fundus Distance | | |
| UtAP | Anterior-Posterior Uterus Diameter | | |
| UtQ | Transverse Uterus Diameter | | |
| Endo | Endometrium Thickness | one distance | Endo[cm or mm]=d1 |
| Lt. Ov-L | Left Ovarian Length | one distance | Lt. Ov-L[cm or mm]=d1 |
| Lt. Ov-H | Left Ovarian Height | one distance | Lt. Ov-H[cm or mm]=d1 |
| Lt. Ov-W | Left Ovarian Width | one distance | Lt. Ov-W[cm or mm]=d1 |
| Lt. Ov-Volume | Left Ovarian Volume | | |
| Rt. Ov-L | Right Ovarian Length | one distance | Rt. Ov-L[cm or mm]=d1 |
| Rt. Ov-H | Right Ovarian Height | one distance | Rt. Ov-H[cm or mm]=d1 |
| Rt. Ov-W | Right Ovarian Width | one distance | Rt. Ov-W[cm or mm]=d1 |
| Rt. Ov-Volume | Right Ovarian Volume | | |
| Lt. Ov-RI | Left Ovarian Vessel Resistive Index | two Doppler blood flow peak velocities | Lt. Ov-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}^a$ |
| Ut-RI | Uterine Vessel Resistive Index | two Doppler blood flow peak velocities | Ut-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}^a$ |
| Rt. Ov-RI | Right Ovarian Vessel Resistive Index | two Doppler blood flow peak velocities | Rt. Ov-RI= $(V_{\max}-V_{\text{diastole}})/V_{\max}^a$ |

| Calc Mnemonic | Calc Name | Input Measurements | Formula |
|---|--|----------------------------|---|
| LtOvFo[ml] | Left Ovary Follicles | One distance ^b | $D1[\text{cm}]^3 \times \pi/6$ |
| | | Two distances ^b | $D1[\text{cm}]^2 \times D2[\text{cm}] \times \pi/6$: (D1 < D2) |
| | | | $D1[\text{cm}] \times D2[\text{cm}] \times 2 \times \pi/6$: (D2 < D1) |
| | | Three distances | $D1[\text{cm}] \times D2[\text{cm}] \times D3[\text{cm}] \times \pi/6$ |
| RtOvFo[ml] | Right Ovary Follicles | One distance ^b | $D1[\text{cm}]^3 \times \pi/6$ |
| | | Two distances ^b | $D1[\text{cm}]^2 \times D2[\text{cm}] \times \pi/6$: (D1 < D2) |
| | | | $D1[\text{cm}] \times D2[\text{cm}] \times 2 \times \pi/6$: (D2 < D1) |
| | | Three distances | $D1[\text{cm}] \times D2[\text{cm}] \times D3[\text{cm}] \times \pi/6$ |
| Lt. Ov-PI | Left Ovarian Vessel Pulsatility Index | | |
| Rt. Ov-PI | Right Ovarian Vessel Pulsatility Index | | |
| ^a) $V_{\text{diastole}} = V_{\text{min}}$ or $V_{\text{end-diastole}}$ (depends on preset selection) ^b) To calculate LtOvFo or RtOvFo with one (or two distances), press the CLEAR key after the first (or second distance) measurement(s). | | | |

Measurement accuracy

General

When using the Measurement and Analysis (M&A) package, it is important to keep in mind the different aspects that affect the accuracy of the measurements. These include acoustical properties, patient echogenicity, measurement tools and algorithms, scanner setup (especially Field-of-view or Range settings), probe type used, and operator inputs.

Sources of error

Image Quality

The accuracy of each measurement is highly dependent on image quality. Image quality is highly dependent on system design, operator variability, and patient echogenicity. The operator variability and patient echogenicity are independent of the ultrasound system

Operator variability

One of the largest potential sources of error is operator variability. A skilled operator can reduce this by optimizing the image quality for each type of measurement. Clear identification of structures, good probe alignment and correct cursor placement is important. Because of pixel resolution, the accuracy of a measurement decreases with decreasing distance on screen. Therefore it is important when scaling the object on the screen to avoid measuring objects that are too small.

See also “Optimizing Measurement Accuracy” below for recommended techniques.

Image measurement

The accuracy in lateral direction is limited by the beam width and the beam positioning. The radial accuracy is mainly limited by the acoustic pulse length.

Doppler alignment

If alignment is not possible, you may use the Angle Correction control to compensate if the flow direction is known.

Errors in velocity measurements increase with the cosine of the angle between the measured flow and the ultrasound beam. For example, an alignment error of 20 degrees, will give a 6% under-estimation of the velocities, while an error of 40 degrees will cause the under-estimation to be 24%. Optimize transducer position to align the beam with the flow direction.

Screen pixel resolution

The display screen is composed of an array of square picture elements (pixels). The smallest resolvable unit is +/- 1 pixel. This pixel error is only significant when measuring short distances on the screen. By observing good scanning practices, the settings of the field of view should be such that the measured distance covers a relatively large portion of the screen. When such scaling is impossible, the pixel error may come into play. The pixel error is +/- 0.2% of the full ultrasound area in the User Screen.

Algorithms

Some formulae used in clinical calculations are based on assumptions or approximations. For example the volume calculations from 2D or M mode assume a certain, 'ideal' shape of the heart chamber, while the actual shape can vary quite much between individuals. Also, formulae taking several "raw" measurements as inputs are prone to increased errors, depending on the combination of input variable accuracies. For example, the Cardiac Output formula from Doppler is sensitive to errors in the entered Diameter, since this will be squared in the formula.

Speed of Sound in Tissue

The average value 1540 meters / second is used for all calculations. Depending on the tissue structures, this generalization may give errors from 2% (typical) to 5% (much fatty tissue layers present).

Optimizing Measurement Accuracy

Probe selection

Select a transducer appropriate for the application, and optimize the transducer frequencies used. Higher imaging frequencies give better resolution, but less penetration than lower frequencies. Lower Doppler frequencies can measure higher max velocities, and at greater depths, but with less velocity resolution than higher Doppler frequencies.

Field of View

All display modes should be adjusted so that the area of interest covers as large portion of the display as possible. Use **Depth, Angle, Zoom, Horizontal Sweep** and **Velocity** controls to optimize the different modes.

Cursor Placement

Avoid placement of the cursor near the array edges when using convex or linear probes.

All measurements are dependent on the accuracy of their “input” data. Consistency and precision in placing cursors and drawing traces correctly on the images are important.

Measurement Uncertainties

The accuracy percentages reported below are based on data taken with optimum control settings, using calibrated phantoms and test equipment. The table below only includes errors related to the system with probes.

The calibration was done for the basic measurable parameters: Distance, Time and Velocity.

Independent sources of uncertainty contribute to a total uncertainty by a RMS (Root Mean Square) combination of the sources. Refer to the discussions above regarding measurement accuracy and sources of error when reading the table below.

| Measurement | Range | Accuracy | Comments |
|---------------------------------|--------------------------|------------------------|----------------------------------|
| 2D Calipers | | | |
| Distance | 1 - 10 cm | 7% | |
| | > 10 cm | 5% | |
| Area | >1 cm ² | 10% | |
| Volume (area + distance) | 20 - 150 cm ³ | 10 ml | |
| M-mode Calipers | | | |
| Distance | 1 - 10 cm | 7% | |
| dt | 0.5 - 1.5 s | 0.5% | With optimal sweep speed setting |
| ECG alignment with M-mode data | All ranges | +/- 10 ms | |
| Spectrum Calipers | | | |
| Velocity | 0.2 - 1.5 m/s | 6% | |
| dt | 0.5 - 1.5 s | 0.5% | With optimal sweep speed setting |
| ECG alignment with Doppler data | All ranges | +/- 10 ms | |
| Q Analysis | | | |
| Velocity from TVI data | | 10% | |
| Strain | -5 to -25% | +/-3 percentage points | |
| Strain rate | -2.0 to + 2.0 1/sec | 20% | |

| Measurement | Range | Accuracy | Comments |
|-------------------|-------------|----------|----------|
| 2D Auto EF | | | |
| ES Volume | 20 - 120 ml | +/-35 ml | |

Measurements

| Measurement | Range | Accuracy | Comments |
|------------------------------|-------------|-------------------------|----------|
| ED Volume | 20 - 120 ml | +/-45 ml | |
| Ejection Fraction | 30 - 80% | +/-15 percentage points | |
| ESV_BiP or EDV_BiP | 20 - 120 ml | +/-20 ml | |
| AFI | | | |
| Global Longitudinal Strain | -5 to -25% | +/-3 percentage points | |
| Regional Longitudinal Strain | 0 to -25% | +/-9 percentage points | |

DICOM SR Measurements

DICOM Structured Reporting (SR) is a standardized format for medical results. Vivid q N and EchoPAC PC supports the specialized form for Adult Echo and Vascular Ultrasound ("Supplement 72" and "Supplement 71" respectively, issued by the DICOM Standards Committee) for M&A results.

Note: For a full comprehensive table of DICOM SR measurements, please refer to the relevant Conformance Statement document located on the following GE Healthcare WEB site:

http://www.gehealthcare.com/us/en/interoperability/dicom/products/ultrasound_dicom.html

"Supplements 71 and 72" do not support all M&A results from Vivid q N and EchoPAC PC. "Supplement 72" limits the information that is possible to send to the following:

- Publicly coded parameters, no pediatric or fetal cardiac and no unassigned measurements.
- Basic modes: 2D, M-mode, Color Flow, PW and CW Doppler.
- Publicly coded methods, no Modified Simpson or Bullet methods.
- Basic derivations (Average and Last); no references between derived measurements and the ones they were made from.
- Wall Motion Scoring: individual segment scores only according to 16-segment model; no graded Hypokinesis (only Hypokinesis is used).

Chapter 2

OB Tables

- ASUM 72
- Berkowitz 74
- Brenner 75
- Campbell 75
- Eriksen 76
- Goldstein 77
- Hadlock 78
- Hansmann..... 85
- Hellman 94
- Hill 94
- Hohler..... 95
- Jeanty..... 95
- JSUM 108
- Kurtz..... 112
- Mayden..... 113
- Mercer 114
- Merz..... 115
- Moore 125
- Nelson 125
- Osaka 126
- Paris 130
- Rempen 133
- Robinson 138
- Tokyo 138
- Tokyo Shinozuka 142
- Williams 149
- Yarkoni..... 149

ASUM

ASUM

Table 2-1: AC: ASUM, Deler (Fetal Age)
Unit: AC (mm); Age (Days); 2SD (Days)

| AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|-----|
| <35 | n/a | — | 126 | 126 | 10 | 228 | 189 | 14 | 331 | 252 | 18 |
| 35 | 70 | 8 | 137 | 133 | 10 | 240 | 196 | 14 | 342 | 259 | 18 |
| 46 | 77 | 8 | 149 | 140 | 10 | 251 | 203 | 14 | 354 | 266 | 20 |
| 57 | 84 | 8 | 160 | 147 | 10 | 263 | 210 | 14 | 365 | 273 | 20 |
| 69 | 91 | 8 | 171 | 154 | 10 | 274 | 217 | 14 | 377 | 280 | 20 |
| 80 | 98 | 9 | 183 | 161 | 10 | 285 | 224 | 16 | >377 | n/a | — |
| 92 | 105 | 9 | 194 | 168 | 12 | 297 | 231 | 16 | | | |
| 103 | 112 | 9 | 206 | 175 | 12 | 308 | 238 | 18 | | | |
| 114 | 119 | 9 | 217 | 182 | 12 | 320 | 245 | 18 | | | |

Table 2-2: BPD: ASUM, Aust NZ, Obstet Gynaecol 1989: 29:26 (Fetal Age)
Unit: BPD (mm); Age (Week); 2SD (Week - * signifies No Data)

| BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <20 | n/a | — | 40 | 123 | 8 | 61 | 171 | 13 | 82 | 225 | 18 |
| 20 | 84 | 4 | 41 | 126 | 9 | 62 | 173 | 13 | 83 | 228 | 18 |
| 21 | 86 | 4 | 42 | 128 | 9 | 63 | 176 | 14 | 84 | 231 | 19 |
| 22 | 88 | 4 | 43 | 130 | 9 | 64 | 178 | 14 | 85 | 234 | 0 |
| 23 | 90 | 4 | 44 | 132 | 9 | 65 | 181 | 14 | 86 | 237 | 0 |
| 24 | 92 | 5 | 45 | 134 | 9 | 66 | 183 | 14 | 87 | 240 | 0 |
| 25 | 94 | 5 | 46 | 136 | 10 | 67 | 186 | 15 | 88 | 244 | 0 |
| 26 | 95 | 5 | 47 | 139 | 10 | 68 | 188 | 15 | 89 | 247 | 0 |
| 27 | 97 | 5 | 48 | 141 | 10 | 69 | 191 | 15 | 90 | 251 | 0 |
| 28 | 99 | 5 | 49 | 143 | 10 | 70 | 193 | 15 | 91 | 255 | 0 |
| 29 | 101 | 6 | 50 | 145 | 11 | 71 | 196 | 16 | 92 | 259 | 0 |
| 30 | 103 | 6 | 51 | 147 | 11 | 72 | 199 | 16 | 93 | 264 | 0 |
| 31 | 105 | 6 | 52 | 149 | 11 | 73 | 201 | 16 | 94 | 270 | 0 |
| 32 | 107 | 6 | 53 | 152 | 11 | 74 | 204 | 16 | 95 | 276 | 0 |
| 33 | 109 | 7 | 54 | 154 | 12 | 75 | 206 | 17 | 96 | 284 | 0 |
| 34 | 111 | 7 | 55 | 157 | 12 | 76 | 209 | 17 | 97 | 292 | 0 |
| 35 | 113 | 7 | 56 | 159 | 12 | 77 | 212 | 17 | 98 | 301 | 0 |
| 36 | 115 | 7 | 57 | 161 | 12 | 78 | 214 | 17 | >98 | n/a | — |
| 37 | 117 | 8 | 58 | 164 | 13 | 79 | 217 | 17 | | | |
| 38 | 119 | 8 | 59 | 166 | 13 | 80 | 220 | 18 | | | |
| 39 | 121 | 8 | 60 | 169 | 13 | 81 | 222 | 18 | | | |

Table 2-3: CRL: ASUM, Silva et al 1991.6 (Fetal Age)
Unit: CRL (mm); Age (Days); 2SD (* No Data available)

| CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <2 | n/a | — | 16 | 56 | * | 34 | 71 | * | 58 | 86 | * |
| 2 | 42 | * | 17 | 57 | * | 36 | 72 | * | 60 | 87 | * |
| 3 | 43 | * | 18 | 58 | * | 37 | 73 | * | 62 | 88 | * |
| 4 | 44 | * | 19 | 59 | * | 38 | 74 | * | 64 | 89 | * |
| 5 | 45 | * | 20 | 60 | * | 40 | 75 | * | 66 | 90 | * |
| 6 | 46 | * | 22 | 61 | * | 41 | 76 | * | 68 | 91 | * |
| 7 | 47 | * | 23 | 62 | * | 43 | 77 | * | 70 | 92 | * |
| 8 | 48 | * | 24 | 63 | * | 45 | 78 | * | 72 | 93 | * |
| 9 | 49 | * | 25 | 64 | * | 46 | 79 | * | 74 | 94 | * |
| 10 | 50 | * | 26 | 65 | * | 48 | 80 | * | 76 | 95 | * |
| 11 | 51 | * | 27 | 66 | * | 50 | 81 | * | 78 | 96 | * |
| 12 | 52 | * | 29 | 67 | * | 51 | 82 | * | 80 | 97 | * |
| 13 | 53 | * | 30 | 68 | * | 53 | 83 | * | 82 | 98 | * |
| 14 | 54 | * | 31 | 69 | * | 55 | 84 | * | >82 | n/a | — |
| 15 | 55 | * | 33 | 70 | * | 57 | 85 | * | | | |

Berkowitz

Table 2-4: BD: Berkowitz (Fetal Age)
Unit: BD (mm); Age (Day); SD (mm)

| BD | Age | SD | BD | Age | SD | BD | Age | SD | BD | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <13 | n/a | — | 25 | 112 | 0 | 38 | 155 | 0 | 51 | 217 | 0 |
| 13 | 81 | 0 | 26 | 116 | 0 | 39 | 159 | 0 | 52 | 223 | 0 |
| 14 | 82 | 0 | 27 | 120 | 0 | 40 | 162 | 0 | 53 | 230 | 0 |
| 15 | 84 | 0 | 28 | 124 | 0 | 41 | 166 | 0 | 54 | 237 | 0 |
| 16 | 86 | 0 | 29 | 128 | 0 | 42 | 169 | 0 | 55 | 244 | 0 |
| 17 | 88 | 0 | 30 | 130 | 0 | 43 | 173 | 0 | 56 | 251 | 0 |
| 18 | 91 | 0 | 31 | 132 | 0 | 44 | 179 | 0 | 57 | 258 | 0 |
| 19 | 95 | 0 | 32 | 135 | 0 | 45 | 185 | 0 | 58 | 266 | 0 |
| 20 | 98 | 0 | 33 | 138 | 0 | 46 | 191 | 0 | 59 | 275 | 0 |
| 21 | 102 | 0 | 34 | 142 | 0 | 47 | 197 | 0 | >59 | n/a | — |
| 22 | 105 | 0 | 35 | 145 | 0 | 48 | 202 | 0 | | | |
| 23 | 109 | 0 | 36 | 149 | 0 | 49 | 207 | 0 | | | |
| 24 | 110 | 0 | 37 | 152 | 0 | 50 | 212 | 0 | | | |

Brenner

Table 2-5: EFW: Brenner (Fetal Growth)

GP, Table/Graph Range = 10%: 90%

Age (Weeks); Mini/Mean/Max (grams)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|------|------|------|------|------|------|------|
| 21.0 | 280 | 410 | 860 | 33.0 | 1480 | 2010 | 2690 |
| 22.0 | 320 | 480 | 920 | 34.0 | 1670 | 2220 | 2880 |
| 23.0 | 370 | 550 | 990 | 35.0 | 1870 | 2430 | 3090 |
| 24.0 | 420 | 640 | 1080 | 36.0 | 2190 | 2650 | 3290 |
| 25.0 | 490 | 740 | 1180 | 37.0 | 2310 | 2870 | 3470 |
| 26.0 | 570 | 860 | 1320 | 38.0 | 2510 | 3030 | 3610 |
| 27.0 | 660 | 990 | 1470 | 39.0 | 2680 | 3170 | 3750 |
| 28.0 | 770 | 1150 | 1660 | 40.0 | 2750 | 3280 | 3870 |
| 29.0 | 890 | 1310 | 1890 | 41.0 | 2800 | 3360 | 3980 |
| 30.0 | 1030 | 1460 | 2100 | 42.0 | 2830 | 3410 | 4060 |
| 31.0 | 1180 | 1630 | 2290 | 43.0 | 2840 | 3420 | 4100 |
| 32.0 | 1310 | 1810 | 2500 | 44.0 | 2790 | 3390 | 4110 |

Campbell

Table 2-6: HC/AC Ratio: Campbell, Br J Obstet Gynaecol 1977, 84:165-174

(Fetal Growth)

Unit: GA (Weeks); Min/Max (Index)

| GA | Min | Max | GA | Min | Max | GA | Min | Max |
|-----|------|------|----|------|------|-----|------|------|
| <13 | n/a | — | 23 | 1.05 | 1.21 | 35 | 0.93 | 1.11 |
| 13 | 1.14 | 1.31 | 25 | 1.04 | 1.22 | 37 | 0.92 | 1.05 |
| 15 | 1.05 | 1.39 | 27 | 1.05 | 1.22 | 39 | 0.87 | 1.06 |
| 17 | 1.07 | 1.29 | 29 | 0.99 | 1.21 | 41 | 0.93 | 1.00 |
| 19 | 1.09 | 1.26 | 31 | 0.96 | 1.17 | >42 | n/a | n/a |
| 21 | 1.06 | 1.25 | 33 | 0.96 | 1.11 | | | |

Eriksen

Table 2-7: TAD: Eriksen (Fetal Age)

Unit: TAD (mm); Age (Day); SD (mm)

| TAD | Age | SD | TAD | Age | SD | TAD | Age | SD | TAD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|------|-----|----|
| <23 | n/a | — | 45 | 134 | 0 | 68 | 182 | 0 | 91 | 232 | 0 |
| 23 | 91 | 0 | 46 | 136 | 0 | 69 | 184 | 0 | 92 | 234 | 0 |
| 24 | 93 | 0 | 47 | 138 | 0 | 70 | 186 | 0 | 93 | 236 | 0 |
| 25 | 95 | 0 | 48 | 140 | 0 | 71 | 188 | 0 | 94 | 239 | 0 |
| 26 | 97 | 0 | 49 | 142 | 0 | 72 | 190 | 0 | 95 | 241 | 0 |
| 27 | 99 | 0 | 50 | 144 | 0 | 73 | 192 | 0 | 96 | 243 | 0 |
| 28 | 101 | 0 | 51 | 146 | 0 | 74 | 195 | 0 | 97 | 245 | 0 |
| 29 | 103 | 0 | 52 | 148 | 0 | 75 | 197 | 0 | 98 | 247 | 0 |
| 30 | 105 | 0 | 53 | 150 | 0 | 76 | 199 | 0 | 99 | 250 | 0 |
| 31 | 107 | 0 | 54 | 152 | 0 | 77 | 201 | 0 | 100 | 252 | 0 |
| 32 | 109 | 0 | 55 | 154 | 0 | 78 | 203 | 0 | 101 | 254 | 0 |
| 33 | 111 | 0 | 56 | 156 | 0 | 79 | 206 | 0 | 102 | 256 | 0 |
| 34 | 113 | 0 | 57 | 158 | 0 | 80 | 208 | 0 | 103 | 259 | 0 |
| 35 | 115 | 0 | 58 | 161 | 0 | 81 | 210 | 0 | 104 | 261 | 0 |
| 36 | 117 | 0 | 59 | 163 | 0 | 82 | 212 | 0 | 105 | 263 | 0 |
| 37 | 119 | 0 | 60 | 165 | 0 | 83 | 214 | 0 | 106 | 266 | 0 |
| 38 | 120 | 0 | 61 | 167 | 0 | 84 | 217 | 0 | 107 | 268 | 0 |
| 39 | 122 | 0 | 62 | 169 | 0 | 85 | 219 | 0 | 108 | 270 | 0 |
| 40 | 124 | 0 | 63 | 171 | 0 | 86 | 221 | 0 | 109 | 273 | 0 |
| 41 | 126 | 0 | 64 | 173 | 0 | 87 | 223 | 0 | 110 | 275 | 0 |
| 42 | 128 | 0 | 65 | 175 | 0 | 88 | 225 | 0 | 111 | 277 | 0 |
| 43 | 130 | 0 | 66 | 177 | 0 | 89 | 228 | 0 | 112 | 280 | 0 |
| 44 | 132 | 0 | 67 | 179 | 0 | 90 | 230 | 0 | >112 | n/a | — |

Goldstein

Table 2-8: TCD: Goldstein et al, Am J OB/GYN, May 1987 (Fetal Growth)
 Unit: TCD (Weeks); Age/Quat1/Mean/Quat3/Max (mm)

| Age | Min | Quat1 | Mean | Quat3 | Max |
|-----|-----|-------|------|-------|-----|
| 15 | 10 | 12 | 14 | 15 | 16 |
| 16 | 14 | 16 | 16 | 16 | 17 |
| 17 | 16 | 17 | 17 | 18 | 18 |
| 18 | 17 | 18 | 18 | 19 | 19 |
| 19 | 18 | 18 | 19 | 19 | 22 |
| 20 | 18 | 19 | 20 | 20 | 22 |
| 21 | 19 | 20 | 22 | 23 | 24 |
| 22 | 21 | 23 | 23 | 24 | 24 |
| 23 | 22 | 23 | 24 | 25 | 26 |
| 24 | 22 | 24 | 25 | 27 | 28 |
| 25 | 23 | 21.5 | 28 | 28 | 29 |
| 26 | 25 | 28 | 29 | 30 | 32 |
| 27 | 26 | 28.5 | 30 | 31 | 32 |
| 28 | 27 | 30 | 31 | 32 | 34 |
| 29 | 29 | 32 | 34 | 36 | 38 |
| 30 | 31 | 32 | 35 | 37 | 40 |
| 31 | 32 | 35 | 38 | 39 | 43 |
| 32 | 33 | 36 | 38 | 40 | 42 |
| 33 | 32 | 36 | 40 | 43 | 44 |
| 34 | 33 | 38 | 40 | 41 | 44 |
| 35 | 31 | 37 | 40.5 | 43 | 47 |
| 36 | 36 | 29 | 43 | 52 | 55 |
| 37 | 37 | 37 | 45 | 52 | 55 |
| 38 | 40 | 40 | 48.5 | 52 | 55 |
| 39 | 52 | 52 | 52 | 55 | 55 |

Hadlock

Table 2-9: AC: Hadlock, Radiology 1984, Vol. 152:497 (Fetal Age)

Unit: AC (mm); Age (Week); 2SD (Week)

| AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD |
|-----|------|-------|-----|------|-------|-----|------|-------|------|------|-------|
| <50 | n/a | — | 135 | 19.0 | ± 2.1 | 225 | 26.9 | ± 2.2 | 315 | 35.4 | ± 3.0 |
| 50 | 12.0 | ± 1.7 | 140 | 19.4 | ± 2.1 | 230 | 27.4 | ± 2.2 | 320 | 35.9 | ± 3.0 |
| 55 | 12.4 | ± 1.7 | 145 | 19.8 | ± 2.1 | 235 | 27.8 | ± 2.2 | 321 | 36.0 | ± 3.1 |
| 60 | 12.8 | ± 1.7 | 150 | 20.2 | ± 2.1 | 240 | 28.3 | ± 2.2 | 325 | 36.4 | ± 3.1 |
| 65 | 13.2 | ± 1.7 | 155 | 20.7 | ± 2.1 | 245 | 28.7 | ± 2.2 | 330 | 36.9 | ± 3.1 |
| 70 | 13.6 | ± 1.7 | 160 | 21.1 | ± 2.1 | 250 | 29.2 | ± 2.2 | 335 | 37.4 | ± 3.1 |
| 75 | 14.0 | ± 1.7 | 165 | 21.5 | ± 2.1 | 255 | 29.7 | ± 2.2 | 340 | 37.9 | ± 3.1 |
| 80 | 14.4 | ± 1.7 | 170 | 22.0 | ± 2.1 | 258 | 30.0 | ± 2.2 | 345 | 38.4 | ± 3.1 |
| 85 | 14.8 | ± 1.7 | 175 | 22.4 | ± 2.1 | 259 | 30.1 | ± 3.0 | 350 | 38.9 | ± 3.1 |
| 90 | 15.2 | ± 1.7 | 180 | 22.9 | ± 2.1 | 260 | 30.2 | ± 3.0 | 355 | 39.4 | ± 3.1 |
| 95 | 15.6 | ± 1.7 | 185 | 23.3 | ± 2.1 | 265 | 30.6 | ± 3.0 | 360 | 39.9 | ± 3.1 |
| 100 | 16.0 | ± 1.7 | 190 | 23.7 | ± 2.1 | 270 | 31.1 | ± 3.0 | 365 | 40.4 | ± 3.1 |
| 105 | 16.4 | ± 1.7 | 192 | 23.9 | ± 2.1 | 275 | 31.6 | ± 3.0 | 370 | 40.9 | ± 3.1 |
| 110 | 16.9 | ± 1.7 | 193 | 24.0 | ± 2.2 | 280 | 32.0 | ± 3.0 | 375 | 41.4 | ± 3.1 |
| 115 | 17.3 | ± 1.7 | 195 | 24.2 | ± 2.2 | 285 | 32.5 | ± 3.0 | 380 | 42.0 | ± 3.1 |
| 120 | 17.7 | ± 1.7 | 200 | 24.6 | ± 2.2 | 290 | 33.0 | ± 3.0 | 385 | 42.5 | ± 3.1 |
| 123 | 17.9 | ± 1.7 | 205 | 25.1 | ± 2.2 | 295 | 33.5 | ± 3.0 | >385 | n/a | — |
| 124 | 18.0 | ± 2.1 | 210 | 25.5 | ± 2.2 | 300 | 34.0 | ± 3.0 | | | |
| 125 | 18.1 | ± 2.1 | 215 | 26.0 | ± 2.2 | 305 | 34.5 | ± 3.0 | | | |
| 130 | 18.5 | ± 2.1 | 220 | 26.4 | ± 2.2 | 310 | 34.9 | ± 3.0 | | | |

Table 2-10: AC: Hadlock, AJR; 139: 367-370; 1982 (Fetal Age)

Unit: AC (mm); Age (Days); SD (Days)

| AC | Age | SD | AC | Age | SD | AC | Age | SD | AC | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|------|-----|----|
| <47 | n/a | — | 138 | 133 | 14 | 230 | 189 | 15 | 305 | 241 | 21 |
| 47 | 84 | 13 | 144 | 136 | 14 | 235 | 192 | 15 | 310 | 245 | 21 |
| 53 | 87 | 13 | 151 | 140 | 14 | 241 | 196 | 15 | 314 | 248 | 21 |
| 60 | 91 | 13 | 157 | 143 | 14 | 246 | 199 | 15 | 319 | 252 | 21 |
| 67 | 94 | 13 | 163 | 147 | 14 | 251 | 203 | 15 | 323 | 255 | 18 |
| 74 | 98 | 13 | 174 | 154 | 14 | 256 | 206 | 15 | 328 | 259 | 18 |
| 80 | 101 | 13 | 180 | 157 | 14 | 261 | 210 | 15 | 332 | 262 | 18 |
| 87 | 105 | 13 | 186 | 161 | 14 | 266 | 213 | 21 | 337 | 266 | 18 |
| 93 | 106 | 13 | 192 | 164 | 14 | 271 | 217 | 21 | 341 | 269 | 18 |
| 100 | 112 | 13 | 197 | 168 | 14 | 276 | 220 | 21 | 344 | 273 | 18 |
| 106 | 115 | 13 | 203 | 171 | 15 | 281 | 224 | 21 | 349 | 276 | 18 |
| 113 | 119 | 13 | 208 | 175 | 15 | 286 | 227 | 21 | 353 | 280 | 18 |
| 119 | 122 | 13 | 214 | 178 | 15 | 291 | 231 | 21 | >353 | n/a | — |
| 126 | 126 | 13 | 219 | 182 | 15 | 296 | 234 | 21 | | | |
| 132 | 129 | 14 | 225 | 185 | 15 | 300 | 238 | 21 | | | |

Table 2-11: BPD: Hadlock, Radiology 1984, Vol. 152:497^a (Fetal Age)
Unit: BPD (mm); Age (Week); 2SD (Week)

| BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD |
|-----|------|-------|-----|------|-------|-----|------|-------|------|------|-------|
| <14 | n/a | — | 36 | 17.0 | ± 1.2 | 59 | 24.1 | ± 2.2 | 82 | 33.0 | ± 3.1 |
| 14 | 11.9 | ± 1.2 | 37 | 17.3 | ± 1.2 | 60 | 24.5 | ± 2.2 | 83 | 33.4 | ± 3.1 |
| 15 | 12.1 | ± 1.2 | 38 | 17.6 | ± 1.2 | 61 | 24.8 | ± 2.2 | 84 | 33.8 | ± 3.1 |
| 16 | 12.3 | ± 1.2 | 39 | 17.9 | ± 1.2 | 62 | 25.2 | ± 2.2 | 85 | 34.2 | ± 3.1 |
| 17 | 12.5 | ± 1.2 | 40 | 18.1 | ± 1.7 | 63 | 25.5 | ± 2.2 | 86 | 34.7 | ± 3.1 |
| 18 | 12.8 | ± 1.2 | 41 | 18.4 | ± 1.7 | 64 | 25.9 | ± 2.2 | 87 | 35.1 | ± 3.1 |
| 19 | 13.0 | ± 1.2 | 42 | 18.7 | ± 1.7 | 65 | 26.3 | ± 2.2 | 88 | 35.6 | ± 3.1 |
| 20 | 13.2 | ± 1.2 | 43 | 19.0 | ± 1.7 | 66 | 26.6 | ± 2.2 | 89 | 36.0 | ± 3.2 |
| 21 | 13.4 | ± 1.2 | 44 | 19.3 | ± 1.7 | 67 | 27.0 | ± 2.2 | 90 | 36.5 | ± 3.2 |
| 22 | 13.6 | ± 1.2 | 45 | 19.6 | ± 1.7 | 68 | 27.4 | ± 2.2 | 91 | 36.9 | ± 3.2 |
| 23 | 13.8 | ± 1.2 | 46 | 19.9 | ± 1.7 | 69 | 27.7 | ± 2.2 | 92 | 37.4 | ± 3.2 |
| 24 | 14.1 | ± 1.2 | 47 | 20.2 | ± 1.7 | 70 | 28.1 | ± 2.2 | 93 | 37.8 | ± 3.2 |
| 25 | 14.3 | ± 1.2 | 48 | 20.5 | ± 1.7 | 71 | 28.5 | ± 2.2 | 94 | 38.3 | ± 3.2 |
| 26 | 14.5 | ± 1.2 | 49 | 20.8 | ± 1.7 | 72 | 28.9 | ± 2.2 | 95 | 38.7 | ± 3.2 |
| 27 | 14.8 | ± 1.2 | 50 | 21.1 | ± 1.7 | 73 | 29.3 | ± 2.2 | 96 | 39.2 | ± 3.2 |
| 28 | 15.0 | ± 1.2 | 51 | 21.5 | ± 1.7 | 74 | 29.7 | ± 2.2 | 97 | 39.7 | ± 3.2 |
| 29 | 15.2 | ± 1.2 | 52 | 21.8 | ± 1.7 | 75 | 30.1 | ± 3.1 | 98 | 40.2 | ± 3.2 |
| 30 | 15.5 | ± 1.2 | 53 | 22.1 | ± 1.7 | 76 | 30.5 | ± 3.1 | 99 | 40.6 | ± 3.2 |
| 31 | 15.7 | ± 1.2 | 54 | 22.4 | ± 1.7 | 77 | 30.9 | ± 3.1 | 100 | 41.1 | ± 3.2 |
| 32 | 16.0 | ± 1.2 | 55 | 22.8 | ± 1.7 | 78 | 31.3 | ± 3.1 | 101 | 41.6 | ± 3.2 |
| 33 | 16.3 | ± 1.2 | 56 | 23.1 | ± 1.7 | 79 | 31.7 | ± 3.1 | 102 | 42.1 | ± 3.2 |
| 34 | 16.5 | ± 1.2 | 57 | 23.4 | ± 1.7 | 80 | 32.1 | ± 3.1 | 103 | 42.6 | ± 3.2 |
| 35 | 16.8 | ± 1.2 | 58 | 23.8 | ± 1.7 | 81 | 32.5 | ± 3.1 | >103 | n/a | — |

a. Variability of GA estimate by BPD at term is ± 2 SD (6 weeks)

Hadlock

Table 2-12: BPD: Hadlock, J Ultrasound Med 1:97-104, April 1982 (Fetal Age)
Unit: BPD (mm); Age (Days); SD (Days)

| BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|-----|
| <20 | n/a | — | 40 | 126 | 10 | 61 | 175 | 9 | 82 | 233 | 14 |
| 20 | 85 | 6 | 41 | 128 | 10 | 62 | 177 | 9 | 83 | 237 | 14 |
| 21 | 88 | 6 | 42 | 130 | 10 | 63 | 180 | 9 | 84 | 239 | 14 |
| 22 | 90 | 6 | 43 | 132 | 10 | 64 | 183 | 9 | 85 | 243 | 14 |
| 23 | 92 | 6 | 44 | 134 | 10 | 65 | 185 | 9 | 86 | 246 | 14 |
| 24 | 93 | 6 | 45 | 137 | 10 | 66 | 188 | 9 | 87 | 249 | 14 |
| 25 | 95 | 6 | 46 | 139 | 10 | 67 | 190 | 9 | 88 | 253 | 25 |
| 26 | 97 | 6 | 47 | 141 | 10 | 68 | 193 | 9 | 89 | 256 | 25 |
| 27 | 99 | 6 | 48 | 144 | 10 | 69 | 196 | 9 | 90 | 259 | 25 |
| 28 | 102 | 6 | 49 | 146 | 10 | 70 | 198 | 9 | 91 | 263 | 25 |
| 29 | 103 | 6 | 50 | 148 | 10 | 71 | 201 | 9 | 92 | 266 | 25 |
| 30 | 105 | 6 | 51 | 151 | 10 | 72 | 204 | 9 | 93 | 270 | 25 |
| 31 | 107 | 6 | 52 | 153 | 10 | 73 | 207 | 9 | 94 | 272 | 25 |
| 32 | 109 | 6 | 53 | 155 | 10 | 74 | 209 | 9 | 95 | 276 | 25 |
| 33 | 111 | 6 | 54 | 158 | 10 | 75 | 213 | 14 | 96 | 279 | 25 |
| 34 | 113 | 6 | 55 | 160 | 10 | 76 | 216 | 14 | 97 | 284 | 25 |
| 35 | 116 | 6 | 56 | 162 | 10 | 77 | 218 | 14 | 98 | 287 | 25 |
| 36 | 118 | 6 | 57 | 163 | 10 | 78 | 221 | 14 | 99 | 291 | 25 |
| 37 | 120 | 6 | 58 | 167 | 10 | 79 | 224 | 14 | 100 | 294 | 25 |
| 38 | 122 | 6 | 59 | 169 | 9 | 80 | 228 | 14 | >100 | n/a | — |
| 39 | 124 | 6 | 60 | 172 | 9 | 81 | 230 | 14 | | | |

Table 2-13: CI: Hadlock, AJR, 137: 83, 1981 (Fetal Growth)

| Min (Index) | Max (Index) |
|-------------|-------------|
| 70 | 86 |

Table 2-14: CRL: Hadlock, Radiology 1992, Vol. 182:501 (Fetal Age)
Unit: CRL (mm); Age (Week); SD (Week)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|------|-------|-----|------|-------|-----|------|-------|------|------|-------|
| <2 | n/a | — | 32 | 10.1 | ± 0.5 | 63 | 12.7 | ± 0.6 | 94 | 15.3 | ± 0.7 |
| 2 | 5.7 | ± 0.3 | 33 | 10.2 | ± 0.5 | 64 | 12.8 | ± 0.6 | 95 | 15.3 | ± 0.7 |
| 3 | 5.9 | ± 0.3 | 34 | 10.3 | ± 0.5 | 65 | 12.8 | ± 0.6 | 96 | 15.4 | ± 0.7 |
| 4 | 6.1 | ± 0.3 | 35 | 10.4 | ± 0.5 | 66 | 12.9 | ± 0.6 | 97 | 15.5 | ± 0.7 |
| 5 | 6.2 | ± 0.3 | 36 | 10.5 | ± 0.5 | 67 | 13.0 | ± 0.6 | 98 | 15.6 | ± 0.7 |
| 6 | 6.4 | ± 0.3 | 37 | 10.6 | ± 0.5 | 68 | 13.1 | ± 0.6 | 99 | 15.7 | ± 0.7 |
| 7 | 6.6 | ± 0.3 | 38 | 10.7 | ± 0.5 | 69 | 13.1 | ± 0.6 | 100 | 15.9 | ± 0.7 |
| 8 | 6.7 | ± 0.3 | 39 | 10.8 | ± 0.5 | 70 | 13.2 | ± 0.6 | 101 | 16.0 | ± 0.7 |
| 9 | 6.9 | ± 0.3 | 40 | 10.9 | ± 0.5 | 71 | 13.3 | ± 0.6 | 102 | 16.1 | ± 0.7 |
| 10 | 7.1 | ± 0.3 | 41 | 11.0 | ± 0.5 | 72 | 13.4 | ± 0.6 | 103 | 16.2 | ± 0.7 |
| 11 | 7.2 | ± 0.3 | 42 | 11.1 | ± 0.5 | 73 | 13.4 | ± 0.6 | 104 | 16.3 | ± 0.7 |
| 12 | 7.4 | ± 0.3 | 43 | 11.2 | ± 0.5 | 74 | 13.5 | ± 0.6 | 105 | 16.4 | ± 0.7 |
| 13 | 7.5 | ± 0.3 | 44 | 11.2 | ± 0.5 | 75 | 13.6 | ± 0.6 | 106 | 16.5 | ± 0.7 |
| 14 | 7.7 | ± 0.3 | 45 | 11.3 | ± 0.5 | 76 | 13.7 | ± 0.6 | 107 | 16.6 | ± 0.7 |
| 15 | 7.9 | ± 0.4 | 46 | 11.4 | ± 0.5 | 77 | 13.8 | ± 0.6 | 108 | 16.7 | ± 0.7 |
| 16 | 8.0 | ± 0.4 | 47 | 11.5 | ± 0.5 | 78 | 13.8 | ± 0.6 | 109 | 16.8 | ± 0.7 |
| 17 | 8.1 | ± 0.4 | 48 | 11.6 | ± 0.5 | 79 | 13.9 | ± 0.6 | 110 | 16.9 | ± 0.8 |
| 18 | 8.3 | ± 0.4 | 49 | 11.7 | ± 0.5 | 80 | 14.0 | ± 0.6 | 111 | 17.0 | ± 0.8 |
| 19 | 8.4 | ± 0.4 | 50 | 11.7 | ± 0.5 | 81 | 14.1 | ± 0.6 | 112 | 17.1 | ± 0.8 |
| 20 | 8.6 | ± 0.4 | 51 | 11.8 | ± 0.5 | 82 | 14.2 | ± 0.6 | 113 | 17.2 | ± 0.8 |
| 21 | 8.7 | ± 0.4 | 52 | 11.9 | ± 0.5 | 83 | 14.2 | ± 0.6 | 114 | 17.3 | ± 0.8 |
| 22 | 8.9 | ± 0.4 | 53 | 12.0 | ± 0.5 | 84 | 14.3 | ± 0.6 | 115 | 17.4 | ± 0.8 |
| 23 | 9.0 | ± 0.4 | 54 | 12.0 | ± 0.5 | 85 | 14.4 | ± 0.6 | 116 | 17.5 | ± 0.8 |
| 24 | 9.1 | ± 0.4 | 55 | 12.1 | ± 0.5 | 86 | 14.5 | ± 0.6 | 117 | 17.6 | ± 0.8 |
| 25 | 9.2 | ± 0.4 | 56 | 12.2 | ± 0.5 | 87 | 14.6 | ± 0.6 | 118 | 17.7 | ± 0.8 |
| 26 | 9.4 | ± 0.4 | 57 | 12.3 | ± 0.5 | 88 | 14.7 | ± 0.7 | 119 | 17.8 | ± 0.8 |
| 27 | 9.5 | ± 0.4 | 58 | 12.3 | ± 0.5 | 89 | 14.8 | ± 0.7 | 120 | 17.9 | ± 0.8 |
| 28 | 9.6 | ± 0.4 | 59 | 12.4 | ± 0.6 | 90 | 14.9 | ± 0.7 | 121 | 18.0 | ± 0.8 |
| 29 | 9.7 | ± 0.4 | 60 | 12.5 | ± 0.6 | 91 | 15.0 | ± 0.7 | >121 | n/a | — |
| 30 | 9.9 | ± 0.4 | 61 | 12.6 | ± 0.6 | 92 | 15.1 | ± 0.7 | | | |
| 31 | 10.0 | ± 0.4 | 62 | 12.6 | ± 0.6 | 93 | 15.2 | ± 0.7 | | | |

Table 2-15: EFW: Hadlock (Fetal Age)
Unit: EFW (grams); Mean (Weeks); SD (grams)

| EFW | Mean | 2SD | EFW | Mean | 2SD | EFW | Mean | 2SD |
|-----|------|-----|------|------|-----|-------|------|-----|
| <35 | n/a | — | 399 | 21 | 51 | 2162 | 33 | 275 |
| 35 | 10 | 4 | 478 | 22 | 61 | 2377 | 34 | 302 |
| 45 | 11 | 6 | 568 | 23 | 72 | 2595 | 35 | 330 |
| 58 | 12 | 7 | 670 | 24 | 85 | 2813 | 36 | 357 |
| 73 | 13 | 9 | 785 | 25 | 101 | 3028 | 37 | 385 |
| 93 | 14 | 12 | 913 | 26 | 116 | 3236 | 38 | 411 |
| 117 | 15 | 15 | 1055 | 27 | 134 | 3435 | 39 | 436 |
| 146 | 16 | 19 | 1210 | 28 | 154 | 3619 | 40 | 460 |
| 181 | 17 | 23 | 1379 | 29 | 175 | 3787 | 41 | 481 |
| 223 | 18 | 28 | 1559 | 30 | 198 | >3787 | n/a | — |
| 273 | 19 | 35 | 1751 | 31 | 222 | | | |
| 331 | 20 | 42 | 1953 | 32 | 248 | | | |

Table 2-16: FL: Hadlock, Radiology 1984, Vol. 152:497 (Fetal Age)
Unit: FL (mm); Age (Week); 2SD (Week)

| FL | Age | 2SD | FL | Age | 2SD | FL | Age | 2SD | FL | Age | 2SD |
|----|------|-------|----|------|-------|----|------|-------|-----|------|-------|
| <6 | n/a | — | 25 | 17.6 | ± 1.4 | 45 | 24.9 | ± 2.1 | 65 | 33.5 | ± 3.0 |
| 6 | 11.9 | ± 1.4 | 26 | 17.9 | ± 1.4 | 46 | 25.3 | ± 2.1 | 66 | 34.0 | ± 3.0 |
| 7 | 12.2 | ± 1.4 | 27 | 18.2 | ± 1.8 | 47 | 25.7 | ± 2.1 | 67 | 34.5 | ± 3.0 |
| 8 | 12.4 | ± 1.4 | 28 | 18.6 | ± 1.8 | 48 | 26.1 | ± 2.1 | 68 | 34.9 | ± 3.0 |
| 9 | 12.7 | ± 1.4 | 29 | 18.9 | ± 1.8 | 49 | 26.5 | ± 2.1 | 69 | 35.4 | ± 3.0 |
| 10 | 13.0 | ± 1.4 | 30 | 19.3 | ± 1.8 | 50 | 26.9 | ± 2.1 | 70 | 35.9 | ± 3.0 |
| 11 | 13.3 | ± 1.4 | 31 | 19.6 | ± 1.8 | 51 | 27.3 | ± 2.1 | 71 | 36.4 | ± 3.1 |
| 12 | 13.5 | ± 1.4 | 32 | 20.0 | ± 1.8 | 52 | 27.7 | ± 2.1 | 72 | 36.9 | ± 3.1 |
| 13 | 13.8 | ± 1.4 | 33 | 20.3 | ± 1.8 | 53 | 28.2 | ± 2.1 | 73 | 37.4 | ± 3.1 |
| 14 | 14.1 | ± 1.4 | 34 | 20.7 | ± 1.8 | 54 | 28.6 | ± 2.1 | 74 | 37.9 | ± 3.1 |
| 15 | 14.4 | ± 1.4 | 35 | 21.0 | ± 1.8 | 55 | 29.0 | ± 2.1 | 75 | 38.4 | ± 3.1 |
| 16 | 14.7 | ± 1.4 | 36 | 21.4 | ± 1.8 | 56 | 29.5 | ± 2.1 | 76 | 38.9 | ± 3.1 |
| 17 | 15.0 | ± 1.4 | 37 | 21.8 | ± 1.8 | 57 | 29.9 | ± 2.1 | 77 | 39.4 | ± 3.1 |
| 18 | 15.3 | ± 1.4 | 38 | 22.2 | ± 1.8 | 58 | 30.3 | ± 3.0 | 78 | 39.9 | ± 3.1 |
| 19 | 15.6 | ± 1.4 | 39 | 22.5 | ± 1.8 | 59 | 30.8 | ± 3.0 | 79 | 40.4 | ± 3.1 |
| 20 | 16.0 | ± 1.4 | 40 | 22.9 | ± 1.8 | 60 | 31.2 | ± 3.0 | 80 | 40.9 | ± 3.1 |
| 21 | 16.3 | ± 1.4 | 41 | 23.3 | ± 1.8 | 61 | 31.7 | ± 3.0 | 81 | 41.4 | ± 3.1 |
| 22 | 16.6 | ± 1.4 | 42 | 23.7 | ± 1.8 | 62 | 32.1 | ± 3.0 | 82 | 42.0 | ± 3.1 |
| 23 | 16.9 | ± 1.4 | 43 | 24.1 | ± 2.1 | 63 | 32.6 | ± 3.0 | 83 | 42.5 | ± 3.1 |
| 24 | 17.2 | ± 1.4 | 44 | 24.5 | ± 2.1 | 64 | 33.1 | ± 3.0 | >83 | n/a | — |

Table 2-17: FL: Hadlock, AJR 138: 875-878, May 1982 (Fetal Age)
Unit: FL (mm); Age (Days); 2SD (Days)

| FL | Age | SD | FL | Age | SD | FL | Age | SD | FL | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <10 | n/a | — | 27 | 125 | 6 | 45 | 171 | 10 | 63 | 226 | 10 |
| 10 | 90 | 6 | 28 | 127 | 6 | 46 | 174 | 10 | 64 | 230 | 10 |
| 11 | 92 | 6 | 29 | 130 | 6 | 47 | 177 | 10 | 65 | 233 | 10 |
| 12 | 94 | 6 | 30 | 132 | 6 | 48 | 180 | 10 | 66 | 237 | 10 |
| 13 | 95 | 6 | 31 | 134 | 6 | 49 | 183 | 10 | 67 | 239 | 10 |
| 14 | 97 | 6 | 32 | 137 | 6 | 50 | 185 | 10 | 68 | 243 | 10 |
| 15 | 99 | 6 | 33 | 139 | 6 | 51 | 189 | 10 | 69 | 246 | 10 |
| 16 | 101 | 6 | 34 | 142 | 6 | 52 | 192 | 10 | 70 | 250 | 10 |
| 17 | 104 | 6 | 35 | 145 | 6 | 53 | 195 | 10 | 71 | 253 | 11 |
| 18 | 106 | 6 | 36 | 147 | 6 | 54 | 197 | 10 | 72 | 257 | 11 |
| 19 | 108 | 6 | 37 | 150 | 6 | 55 | 201 | 10 | 73 | 260 | 11 |
| 20 | 110 | 6 | 38 | 153 | 6 | 56 | 204 | 10 | 74 | 264 | 11 |
| 21 | 112 | 6 | 39 | 155 | 6 | 57 | 207 | 10 | 75 | 268 | 11 |
| 22 | 114 | 6 | 40 | 157 | 6 | 58 | 210 | 10 | 76 | 272 | 11 |
| 23 | 116 | 6 | 41 | 160 | 6 | 59 | 213 | 10 | 77 | 275 | 11 |
| 24 | 118 | 6 | 42 | 163 | 6 | 60 | 216 | 10 | 78 | 279 | 11 |
| 25 | 120 | 6 | 43 | 166 | 6 | 61 | 220 | 10 | 79 | 283 | 11 |
| 26 | 123 | 6 | 44 | 169 | 10 | 62 | 223 | 10 | >79 | n/a | — |

Table 2-18: HC: Hadlock, Radiology 1984, Vol. 152:497 (Fetal Age)
Unit: HC (mm); Age (Week); 2SD (Week)

| HC | Age | 2SD | HC | Age | 2SD | HC | Age | 2SD | HC | Age | 2SD |
|-----|------|-------|-----|------|-------|-----|------|-------|------|------|-------|
| <55 | n/a | — | 135 | 17.0 | ± 1.2 | 215 | 23.6 | ± 1.5 | 290 | 31.9 | ± 3.0 |
| 55 | 12.0 | ± 1.2 | 140 | 17.3 | ± 1.2 | 219 | 23.9 | ± 1.5 | 295 | 32.6 | ± 3.0 |
| 60 | 12.3 | ± 1.2 | 145 | 17.7 | ± 1.2 | 220 | 24.0 | ± 2.1 | 300 | 33.3 | ± 3.0 |
| 65 | 12.6 | ± 1.2 | 149 | 18.0 | ± 1.2 | 225 | 24.5 | ± 2.1 | 305 | 33.9 | ± 3.0 |
| 70 | 12.8 | ± 1.2 | 150 | 18.1 | ± 1.5 | 230 | 25.0 | ± 2.1 | 310 | 34.6 | ± 3.0 |
| 75 | 13.1 | ± 1.2 | 155 | 18.4 | ± 1.5 | 235 | 25.5 | ± 2.1 | 315 | 35.3 | ± 3.0 |
| 80 | 13.4 | ± 1.2 | 160 | 18.8 | ± 1.5 | 240 | 26.1 | ± 2.1 | 319 | 35.9 | ± 3.0 |
| 85 | 13.7 | ± 1.2 | 165 | 19.2 | ± 1.5 | 245 | 26.6 | ± 2.1 | 320 | 36.1 | ± 2.7 |
| 90 | 14.0 | ± 1.2 | 170 | 19.6 | ± 1.5 | 250 | 27.1 | ± 2.1 | 325 | 36.8 | ± 2.7 |
| 95 | 14.3 | ± 1.2 | 175 | 20.0 | ± 1.5 | 255 | 27.7 | ± 2.1 | 330 | 37.6 | ± 2.7 |
| 100 | 14.7 | ± 1.2 | 180 | 20.4 | ± 1.5 | 260 | 28.3 | ± 2.1 | 335 | 38.3 | ± 2.7 |
| 105 | 15.0 | ± 1.2 | 185 | 20.8 | ± 1.5 | 265 | 28.9 | ± 2.1 | 340 | 39.1 | ± 2.7 |
| 110 | 15.3 | ± 1.2 | 190 | 21.3 | ± 1.5 | 270 | 29.4 | ± 2.1 | 345 | 39.9 | ± 2.7 |
| 115 | 15.6 | ± 1.2 | 195 | 21.7 | ± 1.5 | 274 | 29.9 | ± 2.1 | 350 | 40.7 | ± 2.7 |
| 120 | 16.0 | ± 1.2 | 200 | 22.2 | ± 1.5 | 275 | 30.0 | ± 3.0 | 355 | 41.6 | ± 2.7 |
| 125 | 16.3 | ± 1.2 | 205 | 22.6 | ± 1.5 | 280 | 30.7 | ± 3.0 | 360 | 42.4 | ± 2.7 |
| 130 | 16.6 | ± 1.2 | 210 | 23.1 | ± 1.5 | 285 | 31.3 | ± 3.0 | >360 | n/a | — |

Table 2-19: HC: Hadlock, AJR 138: 649-653, 1982 (Fetal Age)
Unit: HC (mm); Age (Days); 2SD (Days)

| HC | Age | SD | HC | Age | SD | HC | Age | SD | HC | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|------|-----|----|
| <69 | n/a | — | 169 | 136 | 11 | 260 | 196 | 16 | 322 | 252 | 19 |
| 69 | 84 | 9 | 175 | 140 | 11 | 264 | 199 | 16 | 325 | 255 | 24 |
| 75 | 87 | 9 | 181 | 143 | 11 | 269 | 203 | 16 | 328 | 259 | 24 |
| 81 | 91 | 9 | 187 | 147 | 11 | 273 | 206 | 16 | 331 | 262 | 24 |
| 88 | 94 | 9 | 193 | 150 | 11 | 278 | 210 | 16 | 334 | 266 | 24 |
| 96 | 98 | 9 | 198 | 154 | 11 | 282 | 213 | 19 | 337 | 269 | 24 |
| 103 | 101 | 9 | 204 | 157 | 11 | 286 | 217 | 19 | 340 | 273 | 24 |
| 110 | 105 | 9 | 209 | 161 | 11 | 291 | 220 | 19 | 343 | 276 | 24 |
| 117 | 108 | 9 | 215 | 164 | 11 | 294 | 224 | 19 | 345 | 280 | 24 |
| 124 | 112 | 9 | 220 | 168 | 11 | 298 | 227 | 19 | 348 | 286 | 24 |
| 131 | 115 | 9 | 225 | 171 | 16 | 302 | 231 | 19 | 351 | 287 | 24 |
| 137 | 119 | 9 | 230 | 175 | 16 | 306 | 234 | 19 | 353 | 290 | 24 |
| 144 | 122 | 9 | 240 | 182 | 16 | 309 | 238 | 19 | 354 | 294 | 24 |
| 150 | 126 | 9 | 245 | 185 | 16 | 312 | 241 | 19 | >354 | n/a | — |
| 157 | 129 | 11 | 250 | 189 | 16 | 316 | 245 | 19 | | | |
| 163 | 133 | 11 | 255 | 192 | 16 | 319 | 248 | 19 | | | |

Hadlock

Table 2-20: FL/HC Ratio: Hadlock, J Ultrasound Med 1984, 3: 439-442 (Fetal Growth)
Unit: GA (Weeks)

| GA | Min | Max | GA | Min | Max | GA | Min | Max |
|-----|------|------|----|------|------|-----|------|------|
| <15 | n/a | — | 24 | 18.7 | 20.9 | 34 | 19.4 | 21.8 |
| 15 | 15.3 | 17.1 | 25 | 18.7 | 20.3 | 35 | 20.1 | 22.3 |
| 16 | 13.3 | 16.5 | 26 | 18.6 | 20.4 | 36 | 20.1 | 22.1 |
| 17 | 14.6 | 17.6 | 27 | 18.6 | 20.4 | 37 | 20.8 | 22.6 |
| 18 | 15.8 | 18.0 | 28 | 18.8 | 20.6 | 38 | 20.9 | 22.7 |
| 19 | 16.1 | 18.3 | 29 | 19.6 | 20.8 | 39 | 20.6 | 23.4 |
| 20 | 16.8 | 19.8 | 30 | 19.2 | 21.4 | 40 | 20.7 | 22.5 |
| 21 | 15.9 | 20.3 | 31 | 19.3 | 21.3 | 41 | 21.6 | 23.2 |
| 22 | 18.4 | 20.2 | 32 | 19.1 | 21.3 | 42 | 20.1 | 23.9 |
| 23 | 19.2 | 20.8 | 33 | 19.9 | 21.5 | >42 | n/a | n/a |

Table 2-21: FL/AC Ratio: Hadlock (Fetal Growth)
Unit: Age (Weeks)

| Age | Min (Index) | Max (Index) |
|-----|-------------|-------------|
| 21 | 20 | 24 |
| 42 | 20 | 24 |

Hansmann

Table 2-22: AC: Hansmann (Fetal Age)
 (Hansmann:M & Al:Geburtsh, u, Frauenheilk 39:656, 1979)
 Unit: AC (mm); Age (Weeks/Days); SD (mm)

| AC | Age | SD | AC | Age | SD | AC | Age | SD | AC | Age | SD |
|-----|-------|----|-----|-------|----|-----|-------|----|-----|-------|----|
| <53 | n/a | — | 99 | 15w2d | 0 | 146 | 20w2d | 0 | 193 | 25w2d | 0 |
| 53 | 11w1d | 0 | 100 | 15w3d | 0 | 147 | 20w2d | 0 | 194 | 25w3d | 0 |
| 54 | 11w2d | 0 | 101 | 15w4d | 0 | 148 | 20w3d | 0 | 195 | 25w4d | 0 |
| 55 | 11w2d | 0 | 102 | 15w5d | 0 | 149 | 20w3d | 0 | 196 | 25w4d | 0 |
| 56 | 11w3d | 0 | 103 | 15w5d | 0 | 150 | 20w4d | 0 | 197 | 25w5d | 0 |
| 57 | 11w3d | 0 | 104 | 15w6d | 0 | 151 | 20w4d | 0 | 198 | 25w5d | 0 |
| 58 | 11w4d | 0 | 105 | 16w0d | 0 | 152 | 20w5d | 0 | 199 | 25w6d | 0 |
| 59 | 11w4d | 0 | 106 | 16w0d | 0 | 153 | 20w6d | 0 | 200 | 26w0d | 0 |
| 60 | 11w5d | 0 | 107 | 16w1d | 0 | 154 | 21w0d | 0 | 201 | 26w0d | 0 |
| 61 | 11w6d | 0 | 108 | 16w2d | 0 | 155 | 21w1d | 0 | 202 | 26w1d | 0 |
| 62 | 12w0d | 0 | 109 | 16w3d | 0 | 156 | 21w2d | 0 | 203 | 26w2d | 0 |
| 63 | 12w1d | 0 | 110 | 16w3d | 0 | 157 | 21w2d | 0 | 204 | 26w3d | 0 |
| 64 | 12w2d | 0 | 111 | 16w4d | 0 | 158 | 21w3d | 0 | 205 | 26w3d | 0 |
| 65 | 12w2d | 0 | 112 | 16w5d | 0 | 159 | 21w3d | 0 | 206 | 26w4d | 0 |
| 66 | 12w3d | 0 | 113 | 16w6d | 0 | 160 | 21w4d | 0 | 207 | 26w5d | 0 |
| 67 | 12w3d | 0 | 114 | 16w6d | 0 | 161 | 21w4d | 0 | 208 | 26w6d | 0 |
| 68 | 12w4d | 0 | 115 | 17w0d | 0 | 162 | 21w5d | 0 | 209 | 26w6d | 0 |
| 69 | 12w5d | 0 | 116 | 17w1d | 0 | 163 | 21w6d | 0 | 210 | 27w0d | 0 |
| 70 | 12w5d | 0 | 117 | 17w2d | 0 | 164 | 22w0d | 0 | 211 | 27w1d | 0 |
| 71 | 12w6d | 0 | 118 | 17w2d | 0 | 165 | 22w1d | 0 | 212 | 27w2d | 0 |
| 72 | 12w6d | 0 | 119 | 17w3d | 0 | 166 | 22w2d | 0 | 213 | 27w2d | 0 |
| 73 | 13w0d | 0 | 120 | 17w3d | 0 | 167 | 22w3d | 0 | 214 | 27w3d | 0 |
| 74 | 13w0d | 0 | 121 | 17w4d | 0 | 168 | 22w4d | 0 | 215 | 27w4d | 0 |
| 75 | 13w1d | 0 | 122 | 17w4d | 0 | 169 | 22w5d | 0 | 216 | 27w4d | 0 |
| 76 | 13w2d | 0 | 123 | 17w5d | 0 | 170 | 22w5d | 0 | 217 | 27w5d | 0 |
| 77 | 13w2d | 0 | 124 | 17w6d | 0 | 171 | 22w6d | 0 | 218 | 27w5d | 0 |
| 78 | 13w3d | 0 | 125 | 18w0d | 0 | 172 | 23w0d | 0 | 219 | 27w6d | 0 |
| 79 | 13w3d | 0 | 126 | 18w1d | 0 | 173 | 23w1d | 0 | 220 | 28w0d | 0 |
| 80 | 13w4d | 0 | 127 | 18w2d | 0 | 174 | 23w2d | 0 | 221 | 28w0d | 0 |
| 81 | 13w4d | 0 | 128 | 18w3d | 0 | 175 | 23w2d | 0 | 222 | 28w1d | 0 |
| 82 | 13w5d | 0 | 129 | 18w3d | 0 | 176 | 23w3d | 0 | 223 | 28w2d | 0 |
| 83 | 13w6d | 0 | 130 | 18w4d | 0 | 177 | 23w3d | 0 | 224 | 28w3d | 0 |
| 84 | 14w0d | 0 | 131 | 18w5d | 0 | 178 | 23w4d | 0 | 225 | 28w4d | 0 |
| 85 | 14w1d | 0 | 132 | 18w6d | 0 | 179 | 23w4d | 0 | 226 | 28w5d | 0 |
| 86 | 14w2d | 0 | 133 | 18w6d | 0 | 180 | 23w5d | 0 | 227 | 28w5d | 0 |
| 87 | 14w2d | 0 | 134 | 19w0d | 0 | 181 | 23w6d | 0 | 228 | 28w6d | 0 |
| 88 | 14w3d | 0 | 135 | 19w1d | 0 | 182 | 24w0d | 0 | 229 | 29w0d | 0 |
| 89 | 14w3d | 0 | 136 | 19w2d | 0 | 183 | 24w1d | 0 | 230 | 29w1d | 0 |
| 90 | 14w4d | 0 | 137 | 19w2d | 0 | 184 | 24w2d | 0 | 231 | 29w2d | 0 |
| 91 | 14w5d | 0 | 138 | 19w3d | 0 | 185 | 24w3d | 0 | 232 | 29w2d | 0 |
| 92 | 14w5d | 0 | 139 | 19w3d | 0 | 186 | 24w4d | 0 | 233 | 29w3d | 0 |
| 93 | 14w6d | 0 | 140 | 19w4d | 0 | 187 | 24w5d | 0 | 234 | 29w3d | 0 |
| 94 | 14w6d | 0 | 141 | 19w4d | 0 | 188 | 24w5d | 0 | 235 | 29w4d | 0 |
| 95 | 15w0d | 0 | 142 | 19w5d | 0 | 189 | 24w6d | 0 | 236 | 29w4d | 0 |
| 96 | 15w0d | 0 | 143 | 19w6d | 0 | 190 | 25w0d | 0 | 237 | 29w5d | 0 |
| 97 | 15w1d | 0 | 144 | 20w0d | 0 | 191 | 25w1d | 0 | 238 | 29w6d | 0 |
| 98 | 15w2d | 00 | 145 | 20w1d | 0 | 192 | 25w2d | 0 | 239 | 30w0d | 0 |

Table 2-22: AC: Hansmann (Fetal Age)
(Hansmann:M & Al:Geburtsh, u, Frauenheilk 39:656, 1979)(Continued)
Unit: AC (mm); Age (Weeks/Days); SD (mm)

| AC | Age | SD | AC | Age | SD | AC | Age | SD | AC | Age | SD |
|-----|-------|----|-----|-------|----|-----|-------|----|------|-------|----|
| 240 | 30w1d | 0 | 261 | 32w3d | 0 | 282 | 34w4d | 0 | 303 | 36w5d | 0 |
| 241 | 30w2d | 0 | 262 | 32w3d | 0 | 283 | 34w4d | 0 | 304 | 36w6d | 0 |
| 242 | 30w3d | 0 | 263 | 32w4d | 0 | 284 | 34w5d | 0 | 305 | 37w0d | 0 |
| 243 | 30w3d | 0 | 264 | 32w4d | 0 | 285 | 34w6d | 0 | 306 | 37w1d | 0 |
| 244 | 30w4d | 0 | 265 | 32w5d | 0 | 286 | 35w0d | 0 | 307 | 37w2d | 0 |
| 245 | 30w5d | 0 | 266 | 32w6d | 0 | 287 | 35w1d | 0 | 308 | 37w3d | 0 |
| 246 | 30w6d | 0 | 267 | 33w0d | 0 | 288 | 35w2d | 0 | 309 | 37w3d | 0 |
| 247 | 30w6d | 0 | 268 | 33w1d | 0 | 289 | 35w3d | 0 | 310 | 37w4d | 0 |
| 248 | 31w0d | 0 | 269 | 33w2d | 0 | 290 | 35w3d | 0 | 311 | 37w5d | 0 |
| 249 | 31w1d | 0 | 270 | 33w3d | 0 | 291 | 35w4d | 0 | 312 | 37w6d | 0 |
| 250 | 31w2d | 0 | 271 | 33w3d | 0 | 292 | 35w5d | 0 | 313 | 37w6d | 0 |
| 251 | 31w3d | 0 | 272 | 33w4d | 0 | 293 | 35w6d | 0 | 314 | 38w0d | 0 |
| 252 | 31w3d | 0 | 273 | 33w5d | 0 | 294 | 35w6d | 0 | 315 | 38w1d | 0 |
| 253 | 31w4d | 0 | 274 | 33w6d | 0 | 295 | 36w0d | 0 | 316 | 38w2d | 0 |
| 254 | 31w5d | 0 | 275 | 33w6d | 0 | 296 | 36w1d | 0 | 317 | 38w4d | 0 |
| 255 | 31w6d | 0 | 276 | 34w0d | 0 | 297 | 36w2d | 0 | 318 | 38w5d | 0 |
| 256 | 31w6d | 0 | 277 | 34w1d | 0 | 298 | 36w2d | 0 | 319 | 39w0d | 0 |
| 257 | 32w0d | 0 | 278 | 34w2d | 0 | 299 | 36w3d | 0 | 320 | 39w1d | 0 |
| 258 | 32w1d | 0 | 279 | 34w2d | 0 | 300 | 36w3d | 0 | >320 | n/a | — |
| 259 | 32w2d | 0 | 280 | 34w3d | 0 | 301 | 36w4d | 0 | | | |
| 260 | 32w2d | 0 | 281 | 34w3d | 0 | 302 | 36w4d | 0 | | | |

Table 2-23: BPD: Hansmann (Fetal Age)
 Ultrasound Diagnosis in Obstetrics & Gynecology, 438-439, 1985

Known LMP (left)—Unknown LMP (right)

Unit: BPD (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD |
|-----|-------|-----|------|-------|-----|-----|-------|-----|------|-------|-----|
| <14 | n/a | — | 60 | 22w6d | 5 | <14 | n/a | — | 60 | 23w2d | 10 |
| 14 | 10w0d | 0 | 61 | 23w1d | 5 | 14 | 9w1d | 7 | 61 | 23w4d | 10 |
| 15 | 10w1d | 0 | 62 | 23w4d | 5 | 15 | 9w3d | 7 | 62 | 24w0d | 10 |
| 16 | 10w2d | 0 | 63 | 23w6d | 5 | 16 | 9w5d | 7 | 63 | 24w2d | 10 |
| 17 | 10w5d | 0 | 64 | 24w1d | 6 | 17 | 10w0d | 7 | 64 | 24w4d | 10 |
| 18 | 10w6d | 0 | 65 | 24w4d | 6 | 18 | 10w2d | 7 | 65 | 24w6d | 10 |
| 19 | 11w1d | 0 | 66 | 24w6d | 6 | 19 | 10w4d | 7 | 66 | 25w1d | 11 |
| 20 | 11w3d | 0 | 67 | 25w1d | 6 | 20 | 10w6d | 7 | 67 | 25w3d | 12 |
| 21 | 11w5d | 0 | 68 | 25w3d | 6 | 21 | 11w1d | 7 | 68 | 25w6d | 10 |
| 22 | 12w0d | 0 | 69 | 25w5d | 6 | 22 | 11w3d | 7 | 69 | 26w1d | 10 |
| 23 | 12w2d | 0 | 70 | 26w1d | 6 | 23 | 11w5d | 7 | 70 | 26w3d | 10 |
| 24 | 12w4d | 5 | 71 | 26w3d | 6 | 24 | 12w0d | 7 | 71 | 26w5d | 12 |
| 25 | 12w6d | 5 | 72 | 26w6d | 6 | 25 | 12w2d | 7 | 72 | 27w1d | 11 |
| 26 | 13w1d | 5 | 73 | 27w1d | 6 | 26 | 12w4d | 7 | 73 | 27w3d | 13 |
| 27 | 13w2d | 5 | 74 | 27w3d | 6 | 27 | 12w6d | 7 | 74 | 27w6d | 12 |
| 28 | 13w4d | 4 | 75 | 27w6d | 6 | 28 | 13w1d | 7 | 75 | 28w1d | 12 |
| 29 | 13w6d | 4 | 76 | 28w1d | 6 | 29 | 13w3d | 8 | 76 | 28w4d | 13 |
| 30 | 14w1d | 4 | 77 | 28w4d | 6 | 30 | 13w5d | 7 | 77 | 28w6d | 13 |
| 31 | 14w3d | 4 | 78 | 28w6d | 6 | 31 | 14w0d | 8 | 78 | 29w2d | 15 |
| 32 | 14w4d | 4 | 79 | 29w2d | 6 | 32 | 14w2d | 8 | 79 | 29w5d | 16 |
| 33 | 14w6d | 4 | 80 | 29w5d | 6 | 33 | 14w4d | 9 | 80 | 30w0d | 15 |
| 34 | 15w2d | 4 | 81 | 30w0d | 6 | 34 | 15w0d | 9 | 81 | 30w3d | 15 |
| 35 | 15w4d | 4 | 82 | 30w3d | 6 | 35 | 15w2d | 8 | 82 | 31w0d | 15 |
| 36 | 15w6d | 4 | 83 | 30w5d | 6 | 36 | 15w4d | 9 | 83 | 31w2d | 16 |
| 37 | 16w1d | 4 | 84 | 31w2d | 6 | 37 | 16w0d | 8 | 84 | 31w6d | 17 |
| 38 | 16w3d | 4 | 85 | 31w5d | 6 | 38 | 16w2d | 9 | 85 | 32w2d | 17 |
| 39 | 16w5d | 4 | 86 | 32w1d | 6 | 39 | 16w4d | 9 | 86 | 32w5d | 18 |
| 40 | 17w0d | 4 | 87 | 32w4d | 6 | 40 | 17w0d | 9 | 87 | 33w2d | 20 |
| 41 | 17w2d | 4 | 88 | 33w0d | 7 | 41 | 17w2d | 9 | 88 | 33w5d | 19 |
| 42 | 17w4d | 4 | 89 | 33w3d | 7 | 42 | 17w4d | 9 | 89 | 34w2d | 19 |
| 43 | 17w6d | 4 | 90 | 33w6d | 7 | 43 | 17w6d | 9 | 90 | 34w5d | 19 |
| 44 | 18w1d | 4 | 91 | 34w3d | 7 | 44 | 18w1d | 9 | 91 | 35w1d | 25 |
| 45 | 18w3d | 4 | 92 | 34w6d | 7 | 45 | 18w4d | 9 | 92 | 35w6d | 24 |
| 46 | 18w5d | 4 | 93 | 35w3d | 7 | 46 | 18w6d | 9 | 93 | 36w5d | 21 |
| 47 | 19w0d | 4 | 94 | 36w0d | 7 | 47 | 19w1d | 10 | 94 | 37w3d | 19 |
| 48 | 19w2d | 5 | 95 | 36w3d | 7 | 48 | 19w3d | 10 | 95 | 38w3d | 22 |
| 49 | 19w4d | 5 | 96 | 37w1d | 7 | 49 | 19w5d | 10 | 96 | 38w6d | 25 |
| 50 | 19w6d | 5 | 97 | 37w6d | 7 | 50 | 20w0d | 10 | 97 | 39w0d | 22 |
| 51 | 20w1d | 5 | 98 | 38w4d | 7 | 51 | 20w3d | 10 | 98 | 39w2d | 20 |
| 52 | 20w3d | 5 | 99 | 39w3d | 7 | 52 | 20w5d | 10 | 99 | 39w3d | 22 |
| 53 | 20w6d | 5 | 100 | 40w3d | 7 | 53 | 21w0d | 11 | 100 | 39w4d | 20 |
| 54 | 21w1d | 5 | 101 | 41w3d | 7 | 54 | 21w3d | 10 | 101 | 39w5d | 20 |
| 55 | 21w2d | 5 | >101 | n/a | — | 55 | 21w5d | 10 | 102 | 39w6d | 19 |
| 56 | 21w4d | 5 | | | | 56 | 22w0d | 9 | 103 | 40w0d | 19 |
| 57 | 21w6d | 5 | | | | 57 | 22w2d | 9 | 104 | 40w1d | 19 |
| 58 | 22w2d | 5 | | | | 58 | 22w5d | 9 | 105 | 40w2d | 17 |
| 59 | 22w4d | 5 | | | | 59 | 23w0d | 10 | >105 | n/a | — |

Table 2-24: CRL: Hansmann (Fetal Age)

Ultrasound Diagnosis in Obstetrics & Gynecology, 438-439, 1985

Unit: CRL (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD |
|-------------|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| Known LMP | | | | | | | | | | | |
| <13 | n/a | — | 54 | 12w0d | 15 | 96 | 15w3d | 11 | 138 | 19w2d | 15 |
| 13 | 7w4d | 0 | 55 | 12w1d | 16 | 97 | 15w3d | 11 | 139 | 19w3d | 15 |
| 14 | 7w5d | 0 | 56 | 12w1d | 16 | 98 | 15w4d | 11 | 140 | 19w4d | 15 |
| 15 | 8w0d | 0 | 57 | 12w2d | 16 | 99 | 15w4d | 11 | 141 | 19w4d | 16 |
| 16 | 8w1d | 0 | 58 | 12w2d | 16 | 100 | 15w5d | 11 | 142 | 19w5d | 16 |
| 17 | 8w2d | 0 | 59 | 12w3d | 16 | 101 | 15w5d | 10 | 143 | 19w5d | 16 |
| 18 | 8w3d | 0 | 60 | 12w3d | 16 | 102 | 15w6d | 10 | 144 | 19w6d | 16 |
| 19 | 8w4d | 7 | 61 | 12w4d | 15 | 103 | 15w6d | 10 | 145 | 20w0d | 16 |
| 20 | 8w5d | 7 | 62 | 12w4d | 15 | 104 | 16w0d | 10 | 146 | 20w1d | 17 |
| 21 | 8w6d | 8 | 63 | 12w5d | 15 | 105 | 16w1d | 10 | 147 | 20w2d | 17 |
| 22 | 9w0d | 8 | 64 | 12w5d | 15 | 106 | 16w2d | 10 | 148 | 20w2d | 17 |
| 23 | 9w1d | 10 | 65 | 12w6d | 15 | 107 | 16w2d | 10 | 149 | 20w3d | 17 |
| 24 | 9w2d | 10 | 66 | 12w6d | 15 | 108 | 16w3d | 10 | 150 | 20w4d | 17 |
| 25 | 9w3d | 11 | 67 | 13w0d | 15 | 109 | 16w3d | 10 | 151 | 20w4d | 0 |
| 26 | 9w4d | 11 | 68 | 13w1d | 15 | 110 | 16w4d | 10 | 152 | 20w5d | 0 |
| 27 | 9w4d | 11 | 69 | 13w1d | 15 | 111 | 16w4d | 11 | 153 | 20w5d | 0 |
| 28 | 9w5d | 11 | 70 | 13w2d | 15 | 112 | 16w5d | 11 | 154 | 20w6d | 0 |
| 29 | 9w6d | 11 | 71 | 13w3d | 15 | 113 | 16w5d | 11 | 155 | 21w0d | 0 |
| 30 | 10w0d | 12 | 72 | 13w3d | 15 | 114 | 16w6d | 11 | 156 | 21w0d | 0 |
| 31 | 10w0d | 12 | 73 | 13w4d | 15 | 115 | 17w0d | 11 | 157 | 21w1d | 0 |
| 32 | 10w1d | 12 | 74 | 13w4d | 15 | 116 | 17w1d | 12 | 158 | 21w1d | 0 |
| 33 | 10w2d | 12 | 75 | 13w5d | 15 | 117 | 17w2d | 12 | 159 | 21w2d | 0 |
| 34 | 10w3d | 12 | 76 | 13w5d | 15 | 118 | 17w2d | 12 | 160 | 21w3d | 0 |
| 35 | 10w3d | 13 | 77 | 13w6d | 15 | 119 | 17w3d | 12 | 161 | 21w3d | 0 |
| 36 | 10w4d | 13 | 78 | 13w6d | 15 | 120 | 17w3d | 12 | 162 | 21w4d | 0 |
| 37 | 10w5d | 13 | 79 | 14w0d | 15 | 121 | 17w4d | 13 | 163 | 21w4d | 0 |
| 38 | 10w5d | 13 | 80 | 14w0d | 15 | 122 | 17w5d | 13 | 164 | 21w5d | 0 |
| 39 | 10w6d | 13 | 81 | 14w1d | 13 | 123 | 17w5d | 13 | 165 | 21w6d | 0 |
| 40 | 10w6d | 13 | 82 | 14w1d | 13 | 124 | 17w6d | 13 | 166 | 21w6d | 0 |
| 41 | 11w0d | 14 | 83 | 14w2d | 13 | 125 | 18w0d | 13 | 167 | 22w0d | 0 |
| 42 | 11w1d | 14 | 84 | 14w2d | 13 | 126 | 18w1d | 14 | 168 | 22w0d | 0 |
| 43 | 11w1d | 14 | 85 | 14w3d | 13 | 127 | 18w1d | 14 | 169 | 22w1d | 0 |
| 44 | 11w2d | 14 | 86 | 14w3d | 13 | 128 | 18w2d | 14 | 170 | 22w1d | 0 |
| 45 | 11w2d | 14 | 87 | 14w4d | 13 | 129 | 18w2d | 14 | 171 | 22w2d | 0 |
| 46 | 11w3d | 14 | 88 | 14w4d | 13 | 130 | 18w3d | 15 | 172 | 22w2d | 0 |
| 47 | 11w3d | 15 | 89 | 14w5d | 13 | 131 | 18w4d | 15 | 173 | 22w3d | 0 |
| 48 | 11w4d | 15 | 90 | 14w6d | 13 | 132 | 18w4d | 15 | 174 | 22w3d | 0 |
| 49 | 11w4d | 15 | 91 | 14w6d | 12 | 133 | 18w5d | 15 | 175 | 22w4d | 0 |
| 50 | 11w5d | 15 | 92 | 15w0d | 12 | 134 | 18w6d | 15 | >175 | n/a | — |
| 51 | 11w5d | 15 | 93 | 15w0d | 12 | 135 | 19w0d | 15 | | | |
| 52 | 11w6d | 15 | 94 | 15w1d | 12 | 136 | 19w1d | 15 | | | |
| 53 | 11w6d | 15 | 95 | 15w2d | 12 | 137 | 19w1d | 15 | | | |
| Unknown LMP | | | | | | | | | | | |

Table 2-24: CRL: Hansmann (Fetal Age) (Continued)
 Ultrasound Diagnosis in Obstetrics & Gynecology, 438-439, 1985
 Unit: CRL (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD |
|-----|------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <6 | n/a | — | 22 | 9w1d | 7 | 54 | 12w3d | 9 | 106 | 16w2d | 13 |
| 6 | 6w1d | 7 | 23 | 9w2d | 7 | 56 | 12w4d | 9 | 110 | 16w4d | 14 |
| 7 | 6w2d | 7 | 24 | 9w3d | 7 | 58 | 12w5d | 9 | 113 | 17w0d | 14 |
| 8 | 6w4d | 7 | 26 | 9w5d | 7 | 60 | 12w6d | 9 | 116 | 17w2d | 14 |
| 9 | 6w6d | 7 | 28 | 10w0d | 8 | 63 | 13w0d | 10 | 120 | 17w4d | 14 |
| 10 | 7w0d | 7 | 30 | 10w2d | 8 | 66 | 13w2d | 10 | 123 | 18w0d | 14 |
| 11 | 7w2d | 7 | 32 | 10w3d | 8 | 70 | 13w3d | 10 | 126 | 18w2d | 14 |
| 12 | 7w3d | 7 | 34 | 10w5d | 8 | 73 | 13w5d | 11 | 130 | 18w6d | 14 |
| 13 | 7w4d | 7 | 36 | 10w6d | 8 | 76 | 13w6d | 11 | 133 | 19w1d | 15 |
| 14 | 7w6d | 7 | 38 | 11w1d | 8 | 80 | 14w1d | 11 | 136 | 19w4d | 16 |
| 15 | 8w0d | 7 | 40 | 11w2d | 8 | 83 | 14w2d | 12 | 140 | 20w0d | 16 |
| 16 | 8w2d | 7 | 42 | 11w3d | 8 | 86 | 14w4d | 12 | 143 | 20w3d | 16 |
| 17 | 8w3d | 7 | 44 | 11w4d | 9 | 90 | 14w6d | 12 | 146 | 20w6d | 16 |
| 18 | 8w4d | 7 | 46 | 11w6d | 9 | 93 | 15w1d | 12 | 150 | 21w3d | 16 |
| 19 | 8w5d | 7 | 48 | 12w0d | 9 | 96 | 15w3d | 12 | >150 | n/a | — |
| 20 | 8w6d | 7 | 50 | 12w1d | 9 | 100 | 15w5d | 12 | | | |
| 21 | 9w0d | 7 | 52 | 12w2d | 9 | 103 | 16w0d | 13 | | | |

Table 2-25: FL: Hansmann (Fetal Age)
 Ultrasound Diagnosis in Obstetrics and Gynecology, 438-439, 1985
 Known/Unknown LMP; Unit: FL (mm); Age (Weeks/Days); 2SD (Week)

| FL | Age | 2SD | FL | Age | 2SD | FL | Age | 2SD | FL | Age | 2SD |
|-----|-------|-----|----|-------|-----|----|-------|-----|-----|-------|-----|
| <12 | n/a | — | 28 | 18w4d | 4 | 45 | 24w6d | 5 | 62 | 32w1d | 5 |
| 12 | 13w4d | 0 | 29 | 18w6d | 4 | 46 | 25w2d | 5 | 63 | 32w5d | 5 |
| 13 | 13w6d | 0 | 30 | 19w2d | 4 | 47 | 25w4d | 5 | 64 | 33w1d | 6 |
| 14 | 14w1d | 0 | 31 | 19w4d | 4 | 48 | 26w0d | 5 | 65 | 33w5d | 6 |
| 15 | 14w3d | 0 | 32 | 20w0d | 4 | 49 | 26w3d | 5 | 66 | 34w1d | 6 |
| 16 | 14w5d | 5 | 33 | 20w3d | 4 | 50 | 26w6d | 5 | 67 | 34w5d | 6 |
| 17 | 15w1d | 5 | 34 | 20w5d | 4 | 51 | 27w3d | 5 | 68 | 35w1d | 6 |
| 18 | 15w2d | 4 | 35 | 21w1d | 5 | 52 | 27w5d | 5 | 69 | 35w5d | 6 |
| 19 | 15w5d | 4 | 36 | 21w3d | 5 | 53 | 28w1d | 5 | 70 | 36w1d | 6 |
| 20 | 16w0d | 4 | 37 | 21w6d | 5 | 54 | 28w4d | 5 | 71 | 36w5d | 6 |
| 21 | 16w2d | 4 | 38 | 22w1d | 5 | 55 | 29w0d | 5 | 72 | 37w2d | 6 |
| 22 | 16w4d | 4 | 39 | 22w4d | 5 | 56 | 29w3d | 6 | 73 | 37w6d | 6 |
| 23 | 16w6d | 4 | 40 | 22w6d | 5 | 57 | 29w6d | 6 | 74 | 38w3d | 7 |
| 24 | 17w2d | 4 | 41 | 23w2d | 5 | 58 | 30w2d | 6 | 75 | 39w0d | 7 |
| 25 | 17w4d | 4 | 42 | 23w5d | 5 | 59 | 30w5d | 5 | >75 | n/a | — |
| 26 | 17w6d | 4 | 43 | 24w0d | 5 | 60 | 31w2d | 5 | | | |
| 27 | 18w2d | 4 | 44 | 24w3d | 5 | 61 | 31w5d | 5 | | | |

Table 2-26: GS: Hansmann (Fetal Age)
Hansmann: M and AI: Geburtsh, u, Frauenheilk 39: 656, 1979
Unit: GS (mm); Age (Days); SD (mm)

| GS | Age | SD | GS | Age | SD | GS | Age | SD | GS | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <10 | n/a | — | 24 | 47 | 5 | 39 | 61 | 5 | 54 | 76 | 5 |
| 10 | 33 | 5 | 25 | 48 | 5 | 40 | 62 | 5 | 55 | 77 | 5 |
| 11 | 34 | 5 | 26 | 49 | 5 | 41 | 63 | 5 | 56 | 78 | 5 |
| 12 | 35 | 5 | 27 | 50 | 5 | 42 | 64 | 5 | 57 | 79 | 5 |
| 13 | 36 | 5 | 28 | 51 | 5 | 43 | 65 | 5 | 58 | 80 | 5 |
| 14 | 37 | 5 | 29 | 52 | 5 | 44 | 66 | 5 | 59 | 81 | 5 |
| 15 | 38 | 5 | 30 | 53 | 5 | 45 | 67 | 5 | 60 | 82 | 5 |
| 16 | 39 | 5 | 31 | 54 | 5 | 46 | 68 | 5 | 61 | 83 | 5 |
| 17 | 40 | 5 | 32 | 55 | 5 | 47 | 69 | 5 | 62 | 84 | 5 |
| 18 | 41 | 5 | 33 | 56 | 5 | 48 | 70 | 5 | 63 | 85 | 5 |
| 19 | 42 | 5 | 34 | 57 | 5 | 49 | 71 | 5 | 64 | 86 | 5 |
| 20 | 43 | 5 | 35 | 58 | 5 | 50 | 72 | 5 | 65 | 87 | 5 |
| 21 | 44 | 5 | 36 | 58 | 5 | 51 | 73 | 5 | >65 | n/a | — |
| 22 | 45 | 5 | 37 | 59 | 5 | 52 | 74 | 5 | | | |
| 23 | 46 | 5 | 38 | 60 | 5 | 53 | 75 | 5 | | | |

Table 2-27: HC: Hansmann (Fetal Age)
Ultrasound Diagnosis in Obstetrics and Gynecology, 438-439, 1985
Known/Unknown LMP; Unit: HC (mm); Age (Weeks/Days); 2SD (mm)

| HC | Age | 2SD | HC | Age | 2SD | HC | Age | 2SD | HC | Age | 2SD |
|------|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <105 | n/a | — | 165 | 18w4d | 16 | 230 | 23w5d | 18 | 295 | 29w5d | 19 |
| 105 | 13w3d | 0 | 170 | 19w0d | 16 | 235 | 24w1d | 18 | 300 | 30w2d | 19 |
| 110 | 14w0d | 0 | 175 | 19w3d | 16 | 240 | 24w4d | 18 | 305 | 30w5d | 19 |
| 115 | 14w3d | 14 | 180 | 19w5d | 16 | 245 | 25w0d | 18 | 310 | 31w2d | 19 |
| 120 | 14w6d | 14 | 185 | 20w1d | 17 | 250 | 25w3d | 18 | 315 | 32w1d | 20 |
| 125 | 15w3d | 14 | 190 | 20w4d | 17 | 255 | 25w6d | 18 | 320 | 32w5d | 20 |
| 130 | 15w5d | 14 | 195 | 21w0d | 17 | 260 | 26w2d | 18 | 325 | 33w3d | 20 |
| 135 | 16w1d | 14 | 200 | 21w2d | 17 | 265 | 26w5d | 18 | 330 | 34w2d | 20 |
| 140 | 16w4d | 14 | 205 | 21w5d | 17 | 270 | 27w1d | 18 | 335 | 35w1d | 20 |
| 145 | 17w0d | 15 | 210 | 22w1d | 17 | 275 | 27w4d | 19 | 340 | 36w2d | 20 |
| 150 | 17w3d | 15 | 215 | 22w4d | 17 | 280 | 28w1d | 19 | 345 | 37w6d | 20 |
| 155 | 17w6d | 16 | 220 | 23w0d | 17 | 285 | 28w5d | 19 | >345 | n/a | — |
| 160 | 18w1d | 16 | 225 | 23w3d | 17 | 290 | 29w1d | 19 | | | |

Table 2-28: OFD: Hansmann (Fetal Age)
 Ultrasound Diagnosis in Obstetrics and Gynecology, 438-439, 1985
Known/Unknown LMP; Unit: OFD (mm); Age (Weeks/Days); 2SD (mm)

| OFD | Age | 2SD | OFD | Age | 2SD | OFD | Age | 2SD | OFD | Age | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <34 | n/a | — | 54 | 18w4d | 5 | 75 | 23w2d | 7 | 96 | 29w0d | 8 |
| 34 | 13w3d | 0 | 55 | 18w6d | 5 | 76 | 23w4d | 7 | 97 | 29w3d | 8 |
| 35 | 13w5d | 0 | 56 | 19w0d | 6 | 77 | 23w6d | 7 | 98 | 29w5d | 8 |
| 36 | 14w0d | 0 | 57 | 19w2d | 6 | 78 | 24w1d | 7 | 99 | 30w0d | 8 |
| 37 | 14w2d | 5 | 58 | 19w3d | 6 | 79 | 24w2d | 7 | 100 | 30w3d | 8 |
| 38 | 14w4d | 5 | 59 | 19w5d | 6 | 80 | 24w4d | 7 | 101 | 30w5d | 8 |
| 39 | 14w6d | 5 | 60 | 20w0d | 6 | 81 | 24w6d | 7 | 102 | 31w1d | 8 |
| 40 | 15w1d | 5 | 61 | 20w1d | 6 | 82 | 25w1d | 7 | 103 | 31w4d | 8 |
| 41 | 15w3d | 5 | 62 | 20w2d | 6 | 83 | 25w2d | 7 | 104 | 32w0d | 8 |
| 42 | 15w5d | 5 | 63 | 20w4d | 6 | 84 | 25w4d | 7 | 105 | 32w3d | 8 |
| 43 | 16w0d | 5 | 64 | 20w6d | 6 | 85 | 25w6d | 7 | 106 | 32w6d | 8 |
| 44 | 16w1d | 5 | 65 | 21w0d | 6 | 86 | 26w1d | 7 | 107 | 33w3d | 8 |
| 45 | 16w3d | 5 | 66 | 21w2d | 6 | 87 | 26w3d | 7 | 108 | 33w6d | 8 |
| 46 | 16w4d | 5 | 67 | 21w4d | 6 | 88 | 26w5d | 7 | 109 | 34w3d | 8 |
| 47 | 16w6d | 5 | 68 | 21w5d | 6 | 89 | 27w0d | 7 | 110 | 35w0d | 8 |
| 48 | 17w1d | 5 | 69 | 22w0d | 6 | 90 | 27w2d | 7 | 111 | 35w4d | 8 |
| 49 | 17w3d | 5 | 70 | 22w1d | 7 | 91 | 27w4d | 8 | 112 | 36w2d | 8 |
| 50 | 17w4d | 5 | 71 | 22w3d | 7 | 92 | 27w6d | 8 | 113 | 37w0d | 8 |
| 51 | 17w6d | 5 | 72 | 22w4d | 7 | 93 | 28w1d | 8 | 114 | 38w0d | 8 |
| 52 | 18w1d | 5 | 73 | 22w6d | 7 | 94 | 28w3d | 8 | 115 | 39w0d | 8 |
| 53 | 18w2d | 5 | 74 | 23w1d | 7 | 95 | 28w5d | 8 | >115 | n/a | — |

Table 2-29: TAD: Hansmann (Fetal Age)
Hansmann: M and AI: Geburtsh, u, Frauenheilk 39: 656, 1979
Unit: TAD (mm); Age (Days); SD (mm)

| TAD | Age | SD | TAD | Age | SD | TAD | Age | SD | TAD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|------|-----|----|
| <20 | n/a | — | 41 | 130 | 4 | 63 | 179 | 4 | 85 | 232 | 5 |
| 20 | 87 | 4 | 42 | 132 | 4 | 64 | 182 | 4 | 86 | 235 | 5 |
| 21 | 89 | 4 | 43 | 135 | 4 | 65 | 184 | 4 | 87 | 237 | 5 |
| 22 | 91 | 4 | 44 | 137 | 4 | 66 | 186 | 4 | 88 | 240 | 5 |
| 23 | 93 | 4 | 45 | 139 | 4 | 67 | 188 | 4 | 89 | 242 | 5 |
| 24 | 95 | 4 | 46 | 141 | 4 | 68 | 191 | 5 | 90 | 245 | 5 |
| 25 | 97 | 4 | 47 | 143 | 4 | 69 | 193 | 5 | 91 | 247 | 5 |
| 26 | 99 | 4 | 48 | 146 | 4 | 70 | 195 | 5 | 92 | 250 | 5 |
| 27 | 101 | 4 | 49 | 148 | 4 | 71 | 198 | 5 | 93 | 252 | 5 |
| 28 | 103 | 4 | 50 | 150 | 4 | 72 | 200 | 5 | 94 | 255 | 5 |
| 29 | 105 | 4 | 51 | 152 | 4 | 73 | 203 | 5 | 95 | 258 | 5 |
| 30 | 107 | 4 | 52 | 155 | 4 | 74 | 205 | 5 | 96 | 261 | 5 |
| 31 | 109 | 4 | 53 | 157 | 4 | 75 | 208 | 5 | 97 | 264 | 5 |
| 32 | 111 | 4 | 54 | 159 | 4 | 76 | 210 | 5 | 98 | 267 | 5 |
| 33 | 113 | 4 | 55 | 161 | 4 | 77 | 212 | 5 | 99 | 270 | 5 |
| 34 | 115 | 4 | 56 | 164 | 4 | 78 | 215 | 5 | 100 | 273 | 5 |
| 35 | 117 | 4 | 57 | 166 | 4 | 79 | 217 | 5 | 101 | 276 | 5 |
| 36 | 119 | 4 | 58 | 168 | 4 | 80 | 220 | 5 | 102 | 279 | 5 |
| 37 | 122 | 4 | 59 | 170 | 4 | 81 | 222 | 5 | 103 | 282 | 5 |
| 38 | 124 | 4 | 60 | 173 | 4 | 82 | 225 | 5 | >103 | n/a | — |
| 39 | 126 | 4 | 61 | 175 | 4 | 83 | 227 | 5 | | | |
| 40 | 128 | 4 | 62 | 177 | 4 | 84 | 230 | 5 | | | |

Table 2-30: ThD: Hansmann (Fetal Age)
 Ultrasound Diagnosis in Obstetrics and Gynecology, 438-439, 1985
Known/Unknown LMP; Unit: ThD (mm); Age (Weeks/Days); 2SD (mm)

| ThD | Age | 2SD | ThD | Age | 2SD | ThD | Age | 2SD | ThD | Age | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|------|-------|-----|
| <20 | n/a | — | 41 | 18w5d | 5 | 63 | 25w5d | 7 | 85 | 33w1d | 9 |
| 20 | 12w4d | 0 | 42 | 19w0d | 5 | 64 | 26w1d | 7 | 86 | 33w4d | 9 |
| 21 | 12w6d | 0 | 43 | 19w3d | 5 | 65 | 26w3d | 7 | 87 | 33w6d | 9 |
| 22 | 13w1d | 0 | 44 | 19w5d | 5 | 66 | 26w5d | 7 | 88 | 34w2d | 9 |
| 23 | 13w3d | 0 | 45 | 19w6d | 5 | 67 | 27w0d | 7 | 89 | 34w4d | 9 |
| 24 | 13w4d | 4 | 46 | 20w2d | 5 | 68 | 27w3d | 8 | 90 | 35w0d | 9 |
| 25 | 13w6d | 4 | 47 | 20w4d | 6 | 69 | 27w5d | 8 | 91 | 35w3d | 10 |
| 26 | 14w1d | 4 | 48 | 20w6d | 6 | 70 | 28w0d | 8 | 92 | 35w5d | 10 |
| 27 | 14w3d | 4 | 49 | 21w2d | 6 | 71 | 28w3d | 8 | 93 | 36w1d | 10 |
| 28 | 14w6d | 4 | 50 | 21w4d | 6 | 72 | 28w5d | 8 | 94 | 36w3d | 10 |
| 29 | 15w1d | 4 | 51 | 21w6d | 6 | 73 | 29w1d | 8 | 95 | 36w6d | 10 |
| 30 | 15w2d | 4 | 52 | 22w1d | 6 | 74 | 29w3d | 8 | 96 | 37w1d | 10 |
| 31 | 15w4d | 4 | 53 | 22w4d | 6 | 75 | 29w5d | 8 | 97 | 37w4d | 10 |
| 32 | 15w6d | 4 | 54 | 22w6d | 6 | 76 | 30w1d | 8 | 98 | 38w1d | 11 |
| 33 | 16w2d | 4 | 55 | 23w1d | 6 | 77 | 30w3d | 8 | 99 | 38w4d | 11 |
| 34 | 16w4d | 4 | 56 | 23w3d | 6 | 78 | 30w5d | 8 | 100 | 38w6d | 11 |
| 35 | 16w6d | 4 | 57 | 23w6d | 7 | 79 | 31w1d | 8 | 101 | 39w3d | 12 |
| 36 | 17w1d | 5 | 58 | 24w1d | 7 | 80 | 31w3d | 8 | 102 | 39w6d | 14 |
| 37 | 17w3d | 5 | 59 | 24w3d | 7 | 81 | 31w5d | 8 | 103 | 40w2d | 14 |
| 38 | 17w5d | 5 | 60 | 24w6d | 7 | 82 | 32w1d | 9 | 104 | 40w5d | 14 |
| 39 | 18w1d | 5 | 61 | 25w1d | 7 | 83 | 32w4d | 9 | 105 | 41w2d | 14 |
| 40 | 18w3d | 5 | 62 | 25w3d | 7 | 84 | 32w6d | 9 | >105 | n/a | — |

Hellman

Table 2-31: GS: Hellman (Fetal Age)
A/OG 103: 789, 1969
Unit: GS (mm); Age (Week); SD (Week)

| GS | Age | SD | GS | Age | SD | GS | Age | SD | GS | Age | SD |
|-----|-----|-------|----|-----|-------|----|------|-------|-----|------|-------|
| <10 | n/a | — | 23 | 6.9 | ± 1.0 | 37 | 8.9 | ± 1.0 | 51 | 10.9 | ± 1.0 |
| 10 | 5.0 | ± 1.0 | 24 | 7.0 | ± 1.0 | 38 | 9.0 | ± 1.0 | 52 | 11.0 | ± 1.0 |
| 11 | 5.2 | ± 1.0 | 25 | 7.2 | ± 1.0 | 39 | 9.2 | ± 1.0 | 53 | 11.2 | ± 1.0 |
| 12 | 5.3 | ± 1.0 | 26 | 7.3 | ± 1.0 | 40 | 9.3 | ± 1.0 | 54 | 11.3 | ± 1.0 |
| 13 | 5.5 | ± 1.0 | 27 | 7.5 | ± 1.0 | 41 | 9.5 | ± 1.0 | 55 | 11.5 | ± 1.0 |
| 14 | 5.6 | ± 1.0 | 28 | 7.6 | ± 1.0 | 42 | 9.6 | ± 1.0 | 56 | 11.6 | ± 1.0 |
| 15 | 5.8 | ± 1.0 | 29 | 7.8 | ± 1.0 | 43 | 9.7 | ± 1.0 | 57 | 11.7 | ± 1.0 |
| 16 | 5.9 | ± 1.0 | 30 | 7.9 | ± 1.0 | 44 | 9.9 | ± 1.0 | 58 | 11.9 | ± 1.0 |
| 17 | 6.0 | ± 1.0 | 31 | 8.0 | ± 1.0 | 45 | 10.0 | ± 1.0 | 59 | 12.0 | ± 1.0 |
| 18 | 6.2 | ± 1.0 | 32 | 8.2 | ± 1.0 | 46 | 10.2 | ± 1.0 | 60 | 12.2 | ± 1.0 |
| 19 | 6.3 | ± 1.0 | 33 | 8.3 | ± 1.0 | 47 | 10.3 | ± 1.0 | >60 | n/a | — |
| 20 | 6.5 | ± 1.0 | 34 | 8.5 | ± 1.0 | 48 | 10.5 | ± 1.0 | | | |
| 21 | 6.6 | ± 1.0 | 35 | 8.6 | ± 1.0 | 49 | 10.6 | ± 1.0 | | | |
| 22 | 6.8 | ± 1.0 | 36 | 8.8 | ± 1.0 | 50 | 10.7 | ± 1.0 | | | |

Hill

Table 2-32: TCD: Hill (Fetal Age)
Obstet Gyn, 75: 981-984, 1990
Unit: TCD (mm); Age (Weeks); SD (Week)

| TCD | Age | SD | TCD | Age | SD | TCD | Age | SD |
|-----|------|--------|-----|------|--------|-----|------|-------|
| <14 | n/a | — | 28 | 24.9 | ± 1.01 | 43 | 33.9 | ± 1.2 |
| 14 | 15.2 | ± 0.5 | 29 | 25.5 | ± 1.01 | 44 | 34.4 | ± 1.2 |
| 15 | 15.8 | ± 0.5 | 30 | 26.2 | ± 1.01 | 45 | 34.8 | ± 1.2 |
| 16 | 16.5 | ± 0.5 | 31 | 26.9 | ± 1.01 | 46 | 35.3 | ± 1.2 |
| 17 | 17.2 | ± 0.5 | 32 | 27.5 | ± 1.01 | 47 | 35.7 | ± 1.2 |
| 18 | 17.9 | ± 0.5 | 33 | 28.1 | ± 1.01 | 48 | 36.1 | ± 1.6 |
| 19 | 18.6 | ± 0.9 | 34 | 28.8 | ± 1.01 | 49 | 36.5 | ± 1.6 |
| 20 | 19.3 | ± 0.9 | 35 | 29.4 | ± 1.01 | 50 | 36.8 | ± 1.6 |
| 21 | 20.0 | ± 0.9 | 36 | 30.0 | ± 1.2 | 51 | 37.2 | ± 1.6 |
| 22 | 20.7 | ± 0.9 | 37 | 30.6 | ± 1.2 | 52 | 37.5 | ± 1.6 |
| 23 | 21.4 | ± 0.9 | 38 | 31.2 | ± 1.2 | 54 | 38.0 | ± 1.6 |
| 24 | 22.1 | ± 0.9 | 39 | 31.8 | ± 1.2 | 55 | 38.3 | ± 1.6 |
| 25 | 22.8 | ± 0.9 | 40 | 32.3 | ± 1.2 | 56 | 38.5 | ± 1.6 |
| 26 | 23.5 | ± 0.9 | 41 | 32.8 | ± 1.2 | >56 | n/a | — |
| 27 | 24.2 | ± 1.01 | 42 | 33.4 | ± 1.2 | | | |

Hohler

Table 2-33: FL: Hohler (Fetal Growth)
Communications in Brief, 143: 479-481, 1982

| Age (Weeks) | Min (Index) | Max (Index) |
|-------------|-------------|-------------|
| 23 | 71 | 87 |
| 40 | 71 | 87 |

Jeanty

Table 2-34: AC: Jeanty (Fetal Age)
Jeanty, Radiology 143: 513, 1982
Unit: AC (mm); Age (Day); SD (mm)

| AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD | AC | Age | 2SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|-----|
| <50 | n/a | — | 115 | 122 | 22 | 185 | 169 | 22 | 255 | 218 | 22 |
| 50 | 79 | 22 | 120 | 125 | 22 | 190 | 172 | 22 | 260 | 222 | 22 |
| 55 | 82 | 22 | 125 | 129 | 22 | 195 | 176 | 22 | 265 | 226 | 22 |
| 60 | 85 | 22 | 130 | 132 | 22 | 200 | 179 | 22 | 270 | 230 | 22 |
| 65 | 89 | 22 | 135 | 135 | 22 | 205 | 182 | 22 | 275 | 234 | 22 |
| 70 | 92 | 22 | 140 | 139 | 22 | 210 | 186 | 22 | 280 | 239 | 22 |
| 75 | 95 | 22 | 145 | 142 | 22 | 215 | 189 | 22 | 285 | 244 | 22 |
| 80 | 99 | 22 | 150 | 146 | 22 | 220 | 192 | 22 | 290 | 249 | 22 |
| 85 | 102 | 22 | 155 | 149 | 22 | 225 | 196 | 22 | 295 | 254 | 22 |
| 90 | 105 | 22 | 160 | 152 | 22 | 230 | 199 | 22 | 300 | 259 | 22 |
| 95 | 109 | 22 | 165 | 156 | 22 | 235 | 203 | 22 | 305 | 265 | 22 |
| 100 | 112 | 22 | 170 | 159 | 22 | 240 | 206 | 22 | 310 | 272 | 22 |
| 105 | 115 | 22 | 175 | 162 | 22 | 245 | 210 | 22 | 315 | 279 | 22 |
| 110 | 119 | 22 | 180 | 166 | 22 | 250 | 214 | 22 | >315 | n/a | — |

Table 2-35: BD: Jeanty (Fetal Age)
 Jeanty: Radiology 143: 513, 1982
 Unit: BD (mm); Age (Days); SD (mm)

| BD | Age | SD | BD | Age | SD | BD | Age | SD | BD | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <15 | n/a | — | 28 | 127 | 0 | 42 | 185 | 0 | 56 | 243 | 0 |
| 15 | 73 | 0 | 29 | 131 | 0 | 43 | 189 | 0 | 57 | 247 | 0 |
| 16 | 77 | 0 | 30 | 135 | 0 | 44 | 193 | 0 | 58 | 251 | 0 |
| 17 | 81 | 0 | 31 | 139 | 0 | 45 | 197 | 0 | 59 | 256 | 0 |
| 18 | 85 | 0 | 32 | 143 | 0 | 46 | 201 | 0 | 60 | 260 | 0 |
| 19 | 89 | 0 | 33 | 147 | 0 | 47 | 206 | 0 | 61 | 264 | 0 |
| 20 | 93 | 0 | 34 | 152 | 0 | 48 | 210 | 0 | 62 | 268 | 0 |
| 21 | 97 | 0 | 35 | 156 | 0 | 49 | 214 | 0 | 63 | 272 | 0 |
| 22 | 102 | 0 | 36 | 160 | 0 | 50 | 218 | 0 | 64 | 276 | 0 |
| 23 | 106 | 0 | 37 | 164 | 0 | 51 | 222 | 0 | 65 | 281 | 0 |
| 24 | 110 | 0 | 38 | 168 | 0 | 52 | 226 | 0 | >65 | n/a | — |
| 25 | 114 | 0 | 39 | 172 | 0 | 53 | 231 | 0 | | | |
| 26 | 118 | 0 | 40 | 177 | 0 | 54 | 235 | 0 | | | |
| 27 | 122 | 0 | 41 | 181 | 0 | 55 | 239 | 0 | | | |

Table 2-36: BPD: Jeanty (Fetal Age)

Jeanty: Radiology 143: 513, 1982

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <10 | n/a | n/a | n/a | 53 | 18w4d | 21w1d | 23w6d |
| 10 | 6w4d | 9w1d | 11w6d | 54 | 18w6d | 21w4d | 24w1d |
| 11 | 6w6d | 9w4d | 12w1d | 55 | 19w1d | 21w6d | 24w4d |
| 12 | 7w0d | 9w5d | 12w3d | 56 | 19w4d | 22w1d | 24w6d |
| 13 | 7w2d | 10w0d | 12w5d | 57 | 19w6d | 22w4d | 25w1d |
| 14 | 7w4d | 10w2d | 12w6d | 58 | 20w1d | 22w6d | 25w4d |
| 15 | 7w6d | 10w4d | 13w1d | 59 | 20w4d | 23w1d | 25w6d |
| 16 | 8w1d | 10w6d | 13w3d | 60 | 20w6d | 23w4d | 26w1d |
| 17 | 8w3d | 11w1d | 13w5d | 61 | 21w1d | 23w6d | 26w4d |
| 18 | 8w4d | 11w2d | 14w0d | 62 | 21w4d | 24w1d | 26w6d |
| 19 | 8w6d | 11w4d | 14w1d | 63 | 21w6d | 24w4d | 27w1d |
| 20 | 9w1d | 11w6d | 14w4d | 64 | 22w1d | 24w6d | 27w4d |
| 21 | 9w3d | 12w1d | 14w6d | 65 | 22w4d | 25w2d | 27w6d |
| 22 | 9w5d | 12w3d | 15w0d | 66 | 22w6d | 25w4d | 28w2d |
| 23 | 9w6d | 12w4d | 15w2d | 67 | 23w2d | 26w0d | 28w4d |
| 24 | 10w1d | 12w6d | 15w4d | 68 | 23w5d | 26w3d | 29w0d |
| 25 | 10w4d | 13w1d | 15w6d | 69 | 24w0d | 26w5d | 29w3d |
| 26 | 10w5d | 13w3d | 16w1d | 70 | 24w3d | 27w1d | 29w6d |
| 27 | 11w0d | 13w5d | 16w3d | 71 | 24w6d | 27w4d | 30w1d |
| 28 | 11w2d | 14w0d | 16w4d | 72 | 25w1d | 27w6d | 30w4d |
| 29 | 11w4d | 14w1d | 16w6d | 73 | 25w4d | 28w2d | 30w6d |
| 30 | 11w6d | 14w4d | 17w1d | 74 | 26w0d | 28w5d | 31w2d |
| 31 | 12w1d | 14w6d | 17w3d | 75 | 26w3d | 29w1d | 31w5d |
| 32 | 12w2d | 15w1d | 17w5d | 76 | 26w6d | 29w4d | 32w1d |
| 33 | 12w4d | 15w2d | 18w0d | 77 | 27w1d | 29w6d | 32w4d |
| 34 | 12w6d | 15w4d | 18w2d | 78 | 27w4d | 30w2d | 33w0d |
| 35 | 13w1d | 15w6d | 18w4d | 79 | 28w0d | 30w5d | 33w3d |
| 36 | 13w4d | 16w1d | 18w6d | 80 | 28w4d | 31w1d | 33w6d |
| 37 | 13w5d | 16w3d | 19w1d | 81 | 28w6d | 31w4d | 34w2d |
| 38 | 14w0d | 16w5d | 19w3d | 82 | 29w2d | 32w0d | 34w5d |
| 39 | 14w2d | 17w0d | 19w5d | 83 | 29w6d | 32w4d | 35w1d |
| 40 | 14w4d | 17w2d | 19w6d | 84 | 30w1d | 32w6d | 35w4d |
| 41 | 14w6d | 17w4d | 20w1d | 85 | 30w5d | 33w3d | 36w0d |
| 42 | 15w1d | 17w6d | 20w4d | 86 | 31w1d | 33w6d | 36w4d |
| 43 | 15w3d | 18w1d | 20w6d | 87 | 31w4d | 34w2d | 37w0d |
| 44 | 15w5d | 18w3d | 21w1d | 88 | 32w1d | 34w6d | 37w3d |
| 45 | 16w0d | 18w5d | 21w3d | 89 | 32w4d | 35w2d | 37w6d |
| 46 | 16w2d | 19w0d | 21w5d | 90 | 33w0d | 35w5d | 38w3d |
| 47 | 16w4d | 19w2d | 22w0d | 91 | 33w4d | 36w1d | 38w6d |
| 48 | 16w6d | 19w4d | 22w2d | 92 | 34w0d | 36w5d | 39w3d |
| 49 | 17w1d | 19w6d | 22w4d | 93 | 34w4d | 37w1d | 39w6d |
| 50 | 17w4d | 20w2d | 22w6d | 94 | 35w0d | 37w5d | 40w3d |
| 51 | 17w6d | 20w4d | 23w1d | 95 | 35w4d | 38w2d | 40w6d |
| 52 | 18w1d | 20w6d | 23w4d | >95 | n/a | n/a | n/a |

Table 2-37: BPD: Jeanty (Fetal Growth)

Jeanty: Radiology 143: 513, 1982

Unit: Age (Weeks/Days); Min/Mean/Max (mm); Table/Graph Range: 5%:95%

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|--------|-----|------|-----|--------|-----|------|-----|
| 10.0+0 | 9 | 14 | 18 | 26.0+0 | 62 | 67 | 71 |
| 11.0+0 | 13 | 17 | 22 | 27.0+0 | 65 | 70 | 74 |
| 12.0+0 | 16 | 21 | 25 | 28.0+0 | 68 | 72 | 77 |
| 13.0+0 | 20 | 24 | 29 | 29.0+0 | 70 | 75 | 79 |
| 14.0+0 | 23 | 28 | 32 | 30.0+0 | 73 | 77 | 82 |
| 15.0+0 | 27 | 31 | 36 | 31.0+0 | 75 | 79 | 84 |
| 16.0+0 | 30 | 35 | 39 | 32.0+0 | 77 | 82 | 86 |
| 17.0+0 | 34 | 38 | 43 | 33.0+0 | 79 | 84 | 88 |
| 18.0+0 | 37 | 42 | 46 | 34.0+0 | 81 | 86 | 90 |
| 19.0+0 | 40 | 45 | 49 | 35.0+0 | 83 | 87 | 92 |
| 20.0+0 | 44 | 48 | 53 | 36.0+0 | 84 | 89 | 93 |
| 21.0+0 | 47 | 51 | 56 | 37.0+0 | 86 | 90 | 95 |
| 22.0+0 | 50 | 55 | 59 | 38.0+0 | 87 | 91 | 96 |
| 23.0+0 | 53 | 58 | 62 | 39.0+0 | 88 | 93 | 97 |
| 24.0+0 | 56 | 61 | 65 | 40.0+0 | 89 | 93 | 98 |
| 25.0+0 | 59 | 64 | 68 | | | | |

Table 2-38: CRL: Jeanty (Fetal Age)

Jeanty: Radiology 143: 513, 1982

Unit: CRL (mm); Age (Days); SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <5 | n/a | — | 17 | 58 | 5 | 30 | 69 | 7 | 43 | 77 | 7 |
| 5 | 44 | 4 | 18 | 59 | 5 | 31 | 70 | 7 | 44 | 78 | 7 |
| 6 | 45 | 4 | 19 | 60 | 5 | 32 | 70 | 7 | 45 | 79 | 7 |
| 7 | 46 | 4 | 20 | 61 | 5 | 33 | 71 | 7 | 46 | 79 | 7 |
| 8 | 48 | 4 | 21 | 62 | 6 | 34 | 72 | 7 | 47 | 80 | 7 |
| 9 | 50 | 4 | 22 | 63 | 6 | 35 | 73 | 7 | 48 | 81 | 7 |
| 10 | 51 | 4 | 23 | 64 | 6 | 36 | 73 | 7 | 49 | 81 | 7 |
| 11 | 52 | 4 | 24 | 65 | 6 | 37 | 74 | 7 | 50 | 82 | 7 |
| 12 | 53 | 4 | 25 | 66 | 6 | 38 | 75 | 7 | 51 | 83 | 7 |
| 13 | 54 | 4 | 26 | 67 | 7 | 39 | 76 | 7 | 52 | 83 | 7 |
| 14 | 55 | 4 | 27 | 67 | 7 | 40 | 76 | 7 | 53 | 84 | 7 |
| 15 | 56 | 5 | 28 | 67 | 7 | 41 | 76 | 7 | 54 | 85 | 7 |
| 16 | 57 | 5 | 29 | 68 | 7 | 42 | 77 | 7 | >54 | n/a | — |

Table 2-39: FIB: Jeanty (Fetal Growth)
 Fetal Limb Bimetry (Letter), Radiology 147:602, 1983
 Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range: 5%:95%

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|-----|------|-----|-----|-----|------|-----|
| 11 | 2 | 2 | 2 | 26 | 32 | 39 | 43 |
| 12 | 5 | 5 | 5 | 27 | 35 | 41 | 47 |
| 13 | 8 | 8 | 8 | 28 | 36 | 43 | 47 |
| 14 | 6 | 11 | 10 | 29 | 40 | 45 | 50 |
| 15 | 10 | 14 | 18 | 30 | 38 | 47 | 52 |
| 16 | 6 | 17 | 22 | 31 | 40 | 48 | 57 |
| 17 | 7 | 19 | 31 | 32 | 40 | 50 | 56 |
| 18 | 10 | 22 | 28 | 33 | 43 | 51 | 59 |
| 19 | 18 | 24 | 30 | 34 | 46 | 52 | 56 |
| 20 | 18 | 27 | 30 | 35 | 51 | 54 | 57 |
| 21 | 24 | 29 | 34 | 36 | 51 | 55 | 56 |
| 22 | 21 | 31 | 37 | 37 | 55 | 56 | 58 |
| 23 | 23 | 33 | 44 | 38 | 54 | 57 | 59 |
| 24 | 26 | 35 | 41 | 39 | 55 | 58 | 62 |
| 25 | 33 | 37 | 42 | 40 | 54 | 59 | 62 |

Table 2-40: FL: Jeanty (Fetal Age)

Jeanty: Radiology 143: 513, 1982

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <14 | n/a | n/a | n/a | 48 | 24w0d | 26w1d | 28w3d |
| 14 | 11w5d | 13w6d | 16w1d | 49 | 24w3d | 26w4d | 28w6d |
| 15 | 12w0d | 14w1d | 16w3d | 50 | 24w6d | 27w0d | 29w1d |
| 16 | 12w3d | 14w4d | 16w6d | 51 | 25w1d | 27w3d | 29w4d |
| 17 | 12w5d | 14w6d | 17w1d | 52 | 25w4d | 27w6d | 30w0d |
| 18 | 13w0d | 15w1d | 17w3d | 53 | 26w0d | 28w1d | 30w3d |
| 19 | 13w3d | 15w4d | 17w6d | 54 | 26w3d | 28w4d | 30w6d |
| 20 | 13w5d | 15w6d | 18w1d | 55 | 26w6d | 29w1d | 31w2d |
| 21 | 14w1d | 16w2d | 18w4d | 56 | 27w2d | 29w4d | 31w5d |
| 22 | 14w3d | 16w4d | 18w6d | 57 | 27w5d | 29w6d | 32w1d |
| 23 | 14w5d | 16w6d | 19w1d | 58 | 28w1d | 30w2d | 32w4d |
| 24 | 15w1d | 17w2d | 19w4d | 59 | 28w4d | 30w5d | 32w6d |
| 25 | 15w3d | 17w4d | 19w6d | 60 | 28w6d | 31w1d | 33w2d |
| 26 | 15w6d | 18w0d | 20w1d | 61 | 29w3d | 31w4d | 33w6d |
| 27 | 16w1d | 18w2d | 20w4d | 62 | 29w6d | 32w0d | 34w1d |
| 28 | 16w4d | 18w5d | 20w6d | 63 | 30w1d | 32w3d | 34w4d |
| 29 | 16w6d | 19w0d | 21w1d | 64 | 30w5d | 32w6d | 35w1d |
| 30 | 17w1d | 19w3d | 21w4d | 65 | 31w1d | 33w2d | 35w4d |
| 31 | 17w4d | 19w6d | 22w0d | 66 | 31w4d | 33w5d | 35w6d |
| 32 | 17w6d | 20w1d | 22w2d | 67 | 32w0d | 34w1d | 36w3d |
| 33 | 18w2d | 20w4d | 22w5d | 68 | 32w3d | 34w4d | 36w6d |
| 34 | 18w5d | 20w6d | 23w1d | 69 | 32w6d | 35w0d | 37w1d |
| 35 | 19w0d | 21w1d | 23w3d | 70 | 33w2d | 35w4d | 37w5d |
| 36 | 19w3d | 21w4d | 23w6d | 71 | 33w5d | 35w6d | 38w1d |
| 37 | 19w6d | 22w0d | 24w1d | 72 | 34w1d | 36w3d | 38w4d |
| 38 | 20w1d | 22w3d | 24w4d | 73 | 34w4d | 36w6d | 39w0d |
| 39 | 20w4d | 22w5d | 24w6d | 74 | 35w1d | 37w2d | 39w4d |
| 40 | 20w6d | 23w1d | 25w2d | 75 | 35w4d | 37w5d | 39w6d |
| 41 | 21w2d | 23w4d | 25w5d | 76 | 36w0d | 38w1d | 40w3d |
| 42 | 21w5d | 23w6d | 26w1d | 77 | 36w3d | 38w4d | 40w6d |
| 43 | 22w1d | 24w2d | 26w4d | 78 | 36w6d | 39w1d | 41w2d |
| 44 | 22w4d | 24w5d | 26w6d | 79 | 37w2d | 39w4d | 41w5d |
| 45 | 22w6d | 25w0d | 27w1d | 80 | 37w6d | 40w0d | 42w1d |
| 46 | 23w1d | 25w3d | 27w4d | >80 | n/a | n/a | n/a |
| 47 | 23w4d | 25w6d | 28w0d | | | | |

Table 2-41: FL: Jeanty (Fetal Growth)

Jeanty: Radiology 143: 513, 1982

Unit: Age (Weeks/Days); Min/Mean/Max (mm); Table/Graph Range: 5%:95%

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|--------|-----|------|-----|--------|-----|------|-----|
| 12.0+0 | 4 | 8 | 13 | 27.0+0 | 45 | 49 | 54 |
| 13.0+0 | 6 | 11 | 16 | 28.0+0 | 47 | 52 | 56 |
| 14.0+0 | 9 | 14 | 18 | 29.0+0 | 50 | 54 | 59 |
| 15.0+0 | 12 | 17 | 21 | 30.0+0 | 52 | 56 | 61 |
| 16.0+0 | 15 | 20 | 24 | 31.0+0 | 54 | 59 | 63 |
| 17.0+0 | 18 | 23 | 27 | 32.0+0 | 56 | 61 | 65 |
| 18.0+0 | 21 | 25 | 30 | 33.0+0 | 58 | 63 | 67 |
| 19.0+0 | 24 | 28 | 33 | 34.0+0 | 60 | 65 | 69 |
| 20.0+0 | 26 | 31 | 36 | 35.0+0 | 62 | 67 | 71 |
| 21.0+0 | 29 | 34 | 38 | 36.0+0 | 64 | 68 | 73 |
| 22.0+0 | 32 | 36 | 41 | 37.0+0 | 65 | 70 | 74 |
| 23.0+0 | 35 | 39 | 44 | 38.0+0 | 67 | 71 | 76 |
| 24.0+0 | 37 | 42 | 46 | 39.0+0 | 68 | 73 | 77 |
| 25.0+0 | 40 | 44 | 49 | 40.0+0 | 70 | 74 | 79 |
| 26.0+0 | 42 | 47 | 51 | | | | |

Table 2-42: HC: Jeanty (Fetal Age)

Jeanty: Radiology 143: 513, 1982

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <80 | n/a | n/a | n/a | 225 | 22w3d | 24w3d | 26w2d |
| 80 | 11w3d | 13w2d | 15w2d | 230 | 22w6d | 24w6d | 26w6d |
| 85 | 11w5d | 13w5d | 15w4d | 235 | 23w3d | 25w3d | 27w2d |
| 90 | 11w7d | 13w7d | 15w6d | 240 | 23w6d | 25w6d | 27w6d |
| 95 | 12w2d | 14w2d | 16w2d | 245 | 24w3d | 26w3d | 28w2d |
| 100 | 12w4d | 14w4d | 16w4d | 250 | 24w7d | 26w6d | 28w6d |
| 105 | 12w7d | 14w6d | 16w6d | 255 | 25w4d | 27w3d | 29w3d |
| 110 | 13w2d | 15w2d | 17w1d | 260 | 26w0d | 28w0d | 29w7d |
| 115 | 13w4d | 15w4d | 17w4d | 265 | 26w4d | 28w4d | 30w4d |
| 120 | 13w6d | 15w6d | 17w6d | 270 | 27w1d | 29w1d | 31w1d |
| 125 | 14w2d | 16w2d | 18w1d | 275 | 27w6d | 29w5d | 31w5d |
| 130 | 14w4d | 16w4d | 18w4d | 280 | 28w3d | 30w2d | 32w2d |
| 135 | 14w7d | 16w6d | 18w6d | 285 | 28w7d | 30w7d | 32w6d |
| 140 | 15w2d | 17w2d | 19w2d | 290 | 29w4d | 31w4d | 33w4d |
| 145 | 15w5d | 17w4d | 19w4d | 295 | 30w2d | 32w1d | 34w1d |
| 150 | 16w0d | 17w7d | 19w7d | 300 | 30w6d | 32w6d | 34w6d |
| 155 | 16w3d | 18w3d | 20w2d | 305 | 31w4d | 33w4d | 35w3d |
| 160 | 16w6d | 18w5d | 20w5d | 310 | 32w2d | 34w1d | 36w1d |
| 165 | 17w1d | 19w1d | 21w1d | 315 | 32w6d | 34w6d | 36w6d |
| 170 | 17w4d | 19w4d | 21w3d | 320 | 33w4d | 35w4d | 37w4d |
| 175 | 17w7d | 19w6d | 21w6d | 325 | 34w2d | 36w2d | 38w2d |
| 180 | 18w3d | 20w2d | 22w2d | 330 | 35w0d | 37w0d | 38w7d |
| 185 | 18w6d | 20w5d | 22w5d | 335 | 35w6d | 37w5d | 39w5d |
| 190 | 19w1d | 21w1d | 23w1d | 340 | 36w4d | 38w4d | 40w3d |
| 195 | 19w4d | 21w4d | 23w4d | 345 | 37w2d | 39w2d | 41w2d |
| 200 | 20w1d | 22w0d | 23w7d | 350 | 38w1d | 40w0d | 42w0d |
| 205 | 20w4d | 22w3d | 24w3d | 355 | 38w6d | 40w6d | 42w6d |
| 210 | 20w7d | 22w7d | 24w6d | 360 | 39w5d | 41w5d | 43w4d |
| 215 | 21w3d | 23w3d | 25w3d | >360 | n/a | n/a | n/a |
| 220 | 21w6d | 23w6d | 25w6d | | | | |

Table 2-43: HC: Jeanty (Fetal Growth)

Jeanty: Radiology 143: 513, 1982

Unit: Age (Weeks/Days); Min/Mean/Max (mm); Table/Graph Range: 5%:95%

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|--------|-----|------|-----|--------|-----|------|-----|
| 12.0+0 | 51 | 75 | 100 | 27.0+0 | 228 | 252 | 277 |
| 13.0+0 | 64 | 88 | 112 | 28.0+0 | 238 | 262 | 286 |
| 14.0+0 | 76 | 101 | 125 | 29.0+0 | 247 | 271 | 296 |
| 15.0+0 | 89 | 113 | 138 | 30.0+0 | 256 | 281 | 305 |
| 16.0+0 | 101 | 126 | 150 | 31.0+0 | 265 | 289 | 313 |
| 17.0+0 | 114 | 138 | 163 | 32.0+0 | 273 | 297 | 322 |
| 18.0+0 | 126 | 151 | 175 | 33.0+0 | 281 | 305 | 329 |
| 19.0+0 | 138 | 163 | 187 | 34.0+0 | 288 | 312 | 336 |
| 20.0+0 | 150 | 175 | 199 | 35.0+0 | 294 | 319 | 343 |
| 21.0+0 | 162 | 187 | 211 | 36.0+0 | 300 | 325 | 349 |
| 22.0+0 | 174 | 198 | 223 | 37.0+0 | 306 | 330 | 355 |
| 23.0+0 | 185 | 210 | 234 | 38.0+0 | 311 | 335 | 359 |
| 24.0+0 | 196 | 221 | 245 | 39.0+0 | 315 | 339 | 364 |
| 25.0+0 | 207 | 232 | 256 | 40.0+0 | 319 | 343 | 367 |
| 26.0+0 | 218 | 242 | 266 | | | | |

Table 2-44: HL: Jeanty (Fetal Age)
Obstetrical Ultrasound, Table 13.9, 1984

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <10 | n/a | n/a | n/a | 40 | 21w4d | 24w2d | 27w1d |
| 10 | 9w6d | 12w4d | 15w2d | 41 | 22w0d | 24w6d | 27w4d |
| 11 | 10w1d | 12w6d | 15w4d | 42 | 22w4d | 25w2d | 28w0d |
| 12 | 10w3d | 13w1d | 15w6d | 43 | 23w0d | 25w5d | 28w4d |
| 13 | 10w6d | 13w4d | 16w1d | 44 | 23w4d | 26w1d | 29w0d |
| 14 | 11w1d | 13w6d | 16w4d | 45 | 24w0d | 26w5d | 29w4d |
| 15 | 11w3d | 14w1d | 16w6d | 46 | 24w4d | 27w1d | 30w0d |
| 16 | 11w6d | 14w4d | 17w2d | 47 | 25w0d | 27w5d | 30w4d |
| 17 | 21w1d | 14w6d | 17w4d | 48 | 25w4d | 28w1d | 31w0d |
| 18 | 12w4d | 15w1d | 18w0d | 49 | 26w0d | 28w6d | 31w4d |
| 19 | 12w6d | 15w4d | 18w2d | 50 | 26w4d | 29w2d | 32w0d |
| 20 | 13w1d | 15w6d | 18w5d | 51 | 27w1d | 29w6d | 32w4d |
| 21 | 13w4d | 16w2d | 19w1d | 52 | 27w4d | 30w2d | 33w1d |
| 22 | 13w6d | 16w5d | 19w3d | 53 | 28w1d | 30w6d | 33w4d |
| 23 | 14w2d | 17w1d | 19w6d | 54 | 28w5d | 31w3d | 34w1d |
| 24 | 14w5d | 17w3d | 20w1d | 55 | 29w1d | 32w0d | 34w5d |
| 25 | 15w1d | 17w6d | 20w4d | 56 | 29w6d | 32w4d | 35w2d |
| 26 | 15w4d | 18w1d | 21w0d | 57 | 30w2d | 33w1d | 35w6d |
| 27 | 15w6d | 18w4d | 21w3d | 58 | 30w6d | 33w4d | 36w3d |
| 28 | 16w2d | 19w0d | 21w6d | 59 | 31w3d | 34w1d | 36w6d |
| 29 | 16w5d | 19w3d | 22w1d | 60 | 32w0d | 34w6d | 37w4d |
| 30 | 17w1d | 19w6d | 22w4d | 61 | 32w4d | 35w2d | 38w1d |
| 31 | 17w4d | 20w2d | 23w0d | 62 | 33w1d | 35w6d | 38w5d |
| 32 | 18w0d | 20w5d | 23w4d | 63 | 33w6d | 36w4d | 39w2d |
| 33 | 18w3d | 21w1d | 23w6d | 64 | 34w3d | 37w1d | 39w6d |
| 34 | 18w6d | 21w4d | 24w2d | 65 | 35w0d | 37w5d | 40w4d |
| 35 | 19w2d | 22w0d | 24w6d | 66 | 35w4d | 38w2d | 41w1d |
| 36 | 19w5d | 22w4d | 25w1d | 67 | 36w1d | 38w6d | 41w5d |
| 37 | 20w1d | 22w6d | 25w5d | 68 | 36w6d | 39w4d | 42w2d |
| 38 | 20w4d | 23w3d | 26w1d | 69 | 37w3d | 40w1d | 42w6d |
| 39 | 21w1d | 23w6d | 26w4d | >69 | n/a | n/a | n/a |

Table 2-45: Radius: Jeanty (Fetal Growth)
 Fetal Limb Bimetry (Letter), Radiology 147:602, 1983
 Unit: Age (weeks); Min/Mean/Max (mm); Table/Graph Range: 5%:95%

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|-----|------|-----|-----|-----|------|-----|
| <11 | n/a | n/a | n/a | 26 | 30 | 37 | 41 |
| 11 | 5 | 5 | 5 | 27 | 33 | 39 | 45 |
| 12 | 7 | 7 | 7 | 28 | 33 | 40 | 45 |
| 13 | 10 | 10 | 10 | 29 | 36 | 42 | 47 |
| 14 | 8 | 13 | 12 | 30 | 34 | 43 | 49 |
| 15 | 12 | 15 | 19 | 31 | 34 | 44 | 53 |
| 16 | 9 | 18 | 21 | 32 | 37 | 45 | 51 |
| 17 | 11 | 20 | 29 | 33 | 41 | 46 | 51 |
| 18 | 14 | 22 | 26 | 34 | 39 | 47 | 53 |
| 19 | 20 | 24 | 29 | 35 | 38 | 48 | 57 |
| 20 | 21 | 27 | 28 | 36 | 41 | 48 | 54 |
| 21 | 25 | 29 | 32 | 37 | 45 | 49 | 53 |
| 22 | 24 | 31 | 34 | 38 | 45 | 49 | 53 |
| 23 | 26 | 32 | 39 | 39 | 46 | 50 | 54 |
| 24 | 27 | 34 | 38 | 40 | 46 | 50 | 54 |
| 25 | 31 | 36 | 40 | >40 | n/a | n/a | n/a |

Table 2-46: TIB: Jeanty (Fetal Age)
Obstetrical Ultrasound, Table 13.9, 1984

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <10 | n/a | n/a | n/a | 40 | 22w3d | 25w2d | 28w1d |
| 10 | 10w4d | 13w3d | 16w2d | 41 | 22w6d | 25w5d | 28w4d |
| 11 | 10w6d | 13w5d | 16w4d | 42 | 23w2d | 26w1d | 29w1d |
| 12 | 11w1d | 14w1d | 17w0d | 43 | 23w5d | 26w4d | 29w4d |
| 13 | 11w4d | 14w3d | 17w2d | 44 | 24w1d | 27w1d | 30w0d |
| 14 | 11w6d | 14w6d | 17w5d | 45 | 24w4d | 27w4d | 30w4d |
| 15 | 12w1d | 15w1d | 18w0d | 46 | 25w1d | 28w0d | 30w6d |
| 16 | 12w4d | 15w4d | 18w3d | 47 | 25w4d | 28w4d | 31w3d |
| 17 | 13w0d | 15w6d | 18w6d | 48 | 26w1d | 29w0d | 31w6d |
| 18 | 13w2d | 16w1d | 19w1d | 49 | 26w4d | 29w3d | 32w2d |
| 19 | 13w5d | 16w4d | 19w4d | 50 | 27w0d | 29w6d | 32w6d |
| 20 | 14w1d | 17w0d | 19w6d | 51 | 27w4d | 30w3d | 33w2d |
| 21 | 14w4d | 17w3d | 20w2d | 52 | 28w0d | 30w6d | 33w6d |
| 22 | 14w6d | 17w6d | 20w5d | 53 | 28w4d | 31w3d | 34w2d |
| 23 | 15w1d | 18w1d | 21w1d | 54 | 29w0d | 31w6d | 34w6d |
| 24 | 15w4d | 18w4d | 21w3d | 55 | 29w4d | 32w3d | 35w2d |
| 25 | 16w0d | 18w6d | 21w6d | 56 | 30w0d | 32w6d | 35w6d |
| 26 | 16w3d | 19w2d | 22w1d | 57 | 30w4d | 33w3d | 36w2d |
| 27 | 16w6d | 19w5d | 22w4d | 58 | 31w0d | 33w6d | 36w6d |
| 28 | 17w1d | 20w1d | 23w0d | 59 | 31w4d | 34w3d | 37w2d |
| 29 | 17w4d | 20w4d | 23w4d | 60 | 32w0d | 34w6d | 37w6d |
| 30 | 18w1d | 21w0d | 23w6d | 61 | 32w4d | 35w3d | 38w2d |
| 31 | 18w4d | 21w3d | 24w2d | 62 | 33w0d | 35w6d | 38w6d |
| 32 | 18w6d | 21w6d | 24w5d | 63 | 33w4d | 36w4d | 39w3d |
| 33 | 19w2d | 22w1d | 25w1d | 64 | 34w1d | 37w0d | 39w6d |
| 34 | 19w5d | 22w4d | 25w4d | 65 | 34w4d | 37w4d | 40w3d |
| 35 | 20w1d | 23w1d | 26w0d | 66 | 35w1d | 38w0d | 41w0d |
| 36 | 20w4d | 23w4d | 26w3d | 67 | 35w5d | 38w4d | 41w4d |
| 37 | 21w0d | 23w6d | 26w6d | 68 | 36w1d | 39w1d | 42w0d |
| 38 | 21w4d | 24w3d | 27w2d | 69 | 36w6d | 39w5d | 42w4d |
| 39 | 21w6d | 24w6d | 27w5d | >69 | n/a | n/a | n/a |

Table 2-47: ULNA: Jeanty (Fetal Age)

Obstetrical Ultrasound, Table 13.9, 1984

Unit: Meas (mm); Min/Mean/Max (Weeks/Days); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <10 | n/a | n/a | n/a | 38 | 22w1d | 25w1d | 28w1d |
| 10 | 10w1d | 13w1d | 16w1d | 39 | 22w4d | 25w4d | 28w5d |
| 11 | 10w4d | 13w4d | 16w4d | 40 | 23w1d | 26w1d | 29w1d |
| 12 | 10w6d | 13w6d | 16w6d | 41 | 23w4d | 26w5d | 29w5d |
| 13 | 11w1d | 14w1d | 17w2d | 42 | 24w1d | 27w1d | 30w2d |
| 14 | 11w4d | 14w4d | 17w5d | 43 | 24w5d | 27w5d | 30w6d |
| 15 | 11w6d | 15w0d | 18w0d | 44 | 25w1d | 28w2d | 31w2d |
| 16 | 12w2d | 15w3d | 18w3d | 45 | 25w6d | 28w6d | 31w6d |
| 17 | 12w5d | 15w5d | 18w6d | 46 | 26w2d | 29w3d | 32w3d |
| 18 | 13w1d | 16w1d | 19w1d | 47 | 26w6d | 29w6d | 33w0d |
| 19 | 13w4d | 16w4d | 19w4d | 48 | 27w3d | 30w4d | 33w4d |
| 20 | 13w6d | 16w6d | 20w0d | 49 | 28w0d | 31w1d | 34w1d |
| 21 | 14w2d | 17w2d | 20w3d | 50 | 28w4d | 31w4d | 34w5d |
| 22 | 14w5d | 17w5d | 20w6d | 51 | 29w1d | 32w1d | 35w2d |
| 23 | 15w1d | 18w1d | 21w1d | 52 | 29w5d | 32w6d | 35w6d |
| 24 | 15w4d | 18w4d | 21w4d | 53 | 30w2d | 33w3d | 36w3d |
| 25 | 16w0d | 19w0d | 22w1d | 54 | 30w6d | 34w0d | 37w0d |
| 26 | 16w3d | 19w3d | 22w4d | 55 | 31w4d | 34w4d | 37w5d |
| 27 | 16w6d | 19w6d | 22w6d | 56 | 32w1d | 35w1d | 38w2d |
| 28 | 17w2d | 20w2d | 23w3d | 57 | 32w6d | 35w6d | 38w6d |
| 29 | 17w5d | 20w6d | 23w6d | 58 | 33w3d | 36w3d | 39w4d |
| 30 | 18w1d | 21w1d | 24w2d | 59 | 34w0d | 37w1d | 40w1d |
| 31 | 18w4d | 21w5d | 24w6d | 60 | 34w4d | 37w5d | 40w6d |
| 32 | 19w1d | 22w1d | 25w1d | 61 | 35w2d | 38w2d | 41w3d |
| 33 | 19w4d | 22w5d | 25w5d | 62 | 35w6d | 39w0d | 42w0d |
| 34 | 20w1d | 23w1d | 26w1d | 63 | 36w4d | 39w4d | 42w5d |
| 35 | 20w4d | 23w4d | 26w5d | 64 | 37w1d | 40w2d | 43w2d |
| 36 | 21w1d | 24w1d | 27w1d | >64 | n/a | n/a | n/a |
| 37 | 21w4d | 24w4d | 27w5d | | | | |

Table 2-48: AC, JSUM, J Med Ultrasound Vol.28 No.5 (2001)

Unit: AC (cm); Age (w+d); SD (cm)

| Age | AC | Age | 1SD | Age | AC | Age | 1SD |
|-----|------|------|------|-----|------|------|------|
| 16 | 10.4 | 16+0 | 0.57 | 30 | 24.2 | 30+0 | 1.24 |
| 17 | 11.4 | 17+0 | 0.62 | 31 | 25.1 | 31+0 | 1.29 |
| 18 | 12.5 | 18+0 | 0.67 | 32 | 25.9 | 32+0 | 1.33 |
| 19 | 13.5 | 19+0 | 0.71 | 33 | 26.8 | 33+0 | 1.38 |
| 20 | 14.5 | 20+0 | 0.76 | 34 | 27.6 | 34+0 | 1.43 |
| 21 | 15.5 | 21+0 | 0.81 | 35 | 28.4 | 35+0 | 1.48 |
| 22 | 16.5 | 22+0 | 0.86 | 36 | 29.2 | 36+0 | 1.52 |
| 23 | 17.5 | 23+0 | 0.90 | 37 | 29.9 | 37+0 | 1.57 |
| 24 | 18.5 | 24+0 | 0.95 | 38 | 30.6 | 38+0 | 1.62 |
| 25 | 19.5 | 25+0 | 1.00 | 39 | 31.3 | 39+0 | 1.67 |
| 26 | 20.5 | 26+0 | 1.05 | 40 | 31.9 | 40+0 | 1.71 |
| 27 | 21.4 | 27+0 | 1.10 | 41 | 32.5 | 41+0 | 1.76 |
| 28 | 22.4 | 28+0 | 1.14 | 42 | 33.1 | 42+0 | 1.81 |
| 29 | 23.3 | 29+0 | 1.19 | | | | |

Table 2-49: BPD, JSUM, J Med Ultrasound Vol.28 No.5 (2001)

Unit: BPD (mm); Age (w+d); SD (mm)

| Age | BPD | Age | 1SD | Age | BPD | Age | 1SD |
|-----|------|------|------|-----|------|------|------|
| 10 | 12.6 | 10+0 | 2.29 | 27 | 67.4 | 27+0 | 3.23 |
| 11 | 15.9 | 11+0 | 2.34 | 28 | 70.1 | 28+0 | 3.29 |
| 12 | 19.3 | 12+0 | 2.40 | 29 | 72.6 | 29+0 | 3.35 |
| 13 | 22.7 | 13+0 | 2.45 | 30 | 75.1 | 30+0 | 3.40 |
| 14 | 26.1 | 14+0 | 2.51 | 31 | 77.4 | 31+0 | 3.46 |
| 15 | 29.5 | 15+0 | 2.57 | 32 | 79.6 | 32+0 | 3.51 |
| 16 | 32.9 | 16+0 | 2.62 | 33 | 81.7 | 33+0 | 3.57 |
| 17 | 36.3 | 17+0 | 2.68 | 34 | 83.6 | 34+0 | 3.62 |
| 18 | 39.6 | 18+0 | 2.73 | 35 | 85.3 | 35+0 | 3.68 |
| 19 | 43.0 | 19+0 | 2.79 | 36 | 86.9 | 36+0 | 3.74 |
| 20 | 46.2 | 20+0 | 2.84 | 37 | 88.3 | 37+0 | 3.79 |
| 21 | 49.5 | 21+0 | 2.90 | 38 | 89.6 | 38+0 | 3.85 |
| 22 | 52.6 | 22+0 | 2.96 | 39 | 90.6 | 39+0 | 3.90 |
| 23 | 55.7 | 23+0 | 3.01 | 40 | 91.5 | 40+0 | 3.96 |
| 24 | 58.8 | 24+0 | 3.07 | 41 | 92.2 | 41+0 | 4.01 |
| 25 | 61.7 | 25+0 | 3.12 | 42 | 92.6 | 42+0 | 4.07 |
| 26 | 64.6 | 26+0 | 3.18 | | | | |

Table 2-50: CRL, JSUM, J Med Ultrasound Vol.28 No.5 (2001)
Unit: GA (week+day); CRL (mm)

| GA | CRL | | | | |
|-------|------|------|------|------|------|
| | 5% | 10% | 50% | 90% | 95% |
| 7W+0 | 5.7 | 6.8 | 10.1 | 16.0 | 17.2 |
| 7W+2 | 6.0 | 7.3 | 10.5 | 15.7 | 16.4 |
| 7W+4 | 6.5 | 8.1 | 11.3 | 16.0 | 16.6 |
| 7W+6 | 7.2 | 9.0 | 12.5 | 17.0 | 17.5 |
| 8W+1 | 8.1 | 10.2 | 14.0 | 18.4 | 19.1 |
| 8W+3 | 9.1 | 11.6 | 15.8 | 20.4 | 21.3 |
| 8W+5 | 10.3 | 13.1 | 17.8 | 22.7 | 24.0 |
| 9W+0 | 11.7 | 14.9 | 20.0 | 25.4 | 27.0 |
| 9W+2 | 13.3 | 16.7 | 22.5 | 28.3 | 30.3 |
| 9W+4 | 15.1 | 18.7 | 25.0 | 31.4 | 33.7 |
| 9W+6 | 17.1 | 20.9 | 27.6 | 34.6 | 37.3 |
| 10W+1 | 19.2 | 23.1 | 30.3 | 37.8 | 40.7 |
| 10W+3 | 21.5 | 25.4 | 33.1 | 41.0 | 44.1 |
| 10W+5 | 24.1 | 27.9 | 35.8 | 44.1 | 47.1 |
| 11W+0 | 26.7 | 30.4 | 38.4 | 47.0 | 49.8 |
| 11W+2 | 29.6 | 32.9 | 40.9 | 49.6 | 52.1 |
| 11W+4 | 32.7 | 35.5 | 43.3 | 51.9 | 53.8 |

Table 2-51: EFW, JSUM, J Med Ultrasound Vol.28 No.5 (2001)
Unit: EFW (g); Age (w+d); 1SD (g)

| Age | EFW | Age | 1SD | Age | EFW | Age | 1SD |
|-----|-------|------|--------|-----|-------|------|--------|
| 18 | 187 | 18+0 | 30.13 | 30 | 1,470 | 30+0 | 185.98 |
| 19 | 247 | 19+0 | 40.47 | 31 | 1,635 | 31+0 | 202.09 |
| 20 | 313 | 20+0 | 51.30 | 32 | 1,805 | 32+0 | 218.68 |
| 21 | 387 | 21+0 | 62.61 | 33 | 1,980 | 33+0 | 235.75 |
| 22 | 469 | 22+0 | 74.39 | 34 | 2,156 | 34+0 | 253.30 |
| 23 | 560 | 23+0 | 86.66 | 35 | 2,333 | 35+0 | 271.33 |
| 24 | 660 | 24+0 | 99.41 | 36 | 2,507 | 36+0 | 289.84 |
| 25 | 771 | 25+0 | 112.64 | 37 | 2,676 | 37+0 | 308.83 |
| 26 | 892 | 26+0 | 126.35 | 38 | 2,838 | 38+0 | 328.30 |
| 27 | 1,023 | 27+0 | 140.53 | 39 | 2,989 | 39+0 | 348.25 |
| 28 | 1,163 | 28+0 | 155.20 | 40 | 3,125 | 40+0 | 368.68 |
| 29 | 1,313 | 29+0 | 170.35 | 41 | 3,244 | 41+0 | 389.59 |

Table 2-52: FL, JSUM, J Med Ultrasound Vol.28 No.5 (2001)
Unit: FL (mm); Age (w+d); SD (mm)

| Age | FL | Age | 1SD | Age | FL | Age | 1SD |
|-----|------|------|------|-----|------|------|------|
| 16 | 20.1 | 16+0 | 2.64 | 30 | 53.8 | 30+0 | 3.11 |
| 17 | 22.7 | 17+0 | 2.67 | 31 | 55.8 | 31+0 | 3.15 |
| 18 | 25.3 | 18+0 | 2.71 | 32 | 57.8 | 32+0 | 3.18 |
| 19 | 27.8 | 19+0 | 2.74 | 33 | 59.6 | 33+0 | 3.21 |
| 20 | 30.4 | 20+0 | 2.77 | 34 | 61.4 | 34+0 | 3.25 |
| 21 | 32.9 | 21+0 | 2.81 | 35 | 63.0 | 35+0 | 3.28 |
| 22 | 35.4 | 22+0 | 2.84 | 36 | 64.6 | 36+0 | 3.31 |
| 23 | 37.9 | 23+0 | 2.88 | 37 | 66.0 | 37+0 | 3.35 |
| 24 | 40.3 | 24+0 | 2.91 | 38 | 67.4 | 38+0 | 3.38 |
| 25 | 42.7 | 25+0 | 2.94 | 39 | 68.6 | 39+0 | 3.42 |
| 26 | 45.0 | 26+0 | 2.98 | 40 | 69.6 | 40+0 | 3.45 |
| 27 | 47.3 | 27+0 | 3.01 | 41 | 70.6 | 41+0 | 3.48 |
| 28 | 49.6 | 28+0 | 3.04 | 42 | 71.4 | 42+0 | 3.52 |
| 29 | 51.7 | 29+0 | 3.08 | | | | |

Table 2-53: MCA PI values with advance in gestation
JSUM, J Med Ultrasound Vol.28 No.5 (2001)
Unit: Age (Weeks)

| Age | 5% | 10% | 50% | 90% | 95% | Age | 5% | 10% | 50% | 90% | 95% |
|-----|-------|-------|-------|-------|-------|-----|-------|-------|-------|-------|-------|
| 20 | 1.271 | 1.270 | 1.440 | 1.880 | 1.990 | 31 | 1.446 | 1.515 | 1.933 | 2.436 | 2.489 |
| 21 | 1.318 | 1.329 | 1.537 | 1.986 | 2.091 | 32 | 1.425 | 1.493 | 1.915 | 2.420 | 2.468 |
| 22 | 1.359 | 1.381 | 1.623 | 2.080 | 2.182 | 33 | 1.397 | 1.464 | 1.887 | 2.394 | 2.435 |
| 23 | 1.393 | 1.426 | 1.699 | 2.164 | 2.261 | 34 | 1.363 | 1.427 | 1.849 | 2.356 | 2.390 |
| 24 | 1.421 | 1.463 | 1.765 | 2.236 | 2.328 | 35 | 1.324 | 1.383 | 1.800 | 2.308 | 2.335 |
| 25 | 1.444 | 1.493 | 1.820 | 2.298 | 2.385 | 36 | 1.277 | 1.331 | 1.741 | 2.248 | 2.268 |
| 26 | 1.459 | 1.515 | 1.865 | 2.348 | 2.430 | 37 | 1.225 | 1.272 | 1.671 | 2.178 | 2.191 |
| 27 | 1.469 | 1.530 | 1.899 | 2.388 | 2.465 | 38 | 1.167 | 1.205 | 1.591 | 2.096 | 2.102 |
| 28 | 1.473 | 1.537 | 1.923 | 2.416 | 2.488 | 39 | 1.102 | 1.131 | 1.501 | 2.004 | 2.001 |
| 29 | 1.470 | 1.537 | 1.937 | 2.434 | 2.499 | 40 | 1.031 | 1.050 | 1.400 | 1.900 | 1.890 |
| 30 | 1.461 | 1.530 | 1.940 | 2.440 | 2.500 | 41 | 0.954 | 0.961 | 1.289 | 1.786 | 1.767 |

Table 2-54: MCA RI values with advance in gestation
 JSUM, J Med Ultrasound Vol.28 No.5 (2001)
 Unit: Age (Weeks)

| Age | 5% | 10% | 50% | 90% | 95% | Age | 5% | 10% | 50% | 90% | 95% |
|-----|-------|-------|-------|-------|-------|-----|-------|-------|-------|-------|-------|
| 20 | 0.717 | 0.718 | 0.775 | 0.842 | 0.871 | 31 | 0.769 | 0.789 | 0.865 | 0.922 | 0.928 |
| 21 | 0.731 | 0.735 | 0.793 | 0.857 | 0.883 | 32 | 0.762 | 0.783 | 0.862 | 0.920 | 0.925 |
| 22 | 0.742 | 0.749 | 0.808 | 0.871 | 0.894 | 33 | 0.755 | 0.775 | 0.857 | 0.916 | 0.920 |
| 23 | 0.753 | 0.761 | 0.821 | 0.883 | 0.903 | 34 | 0.745 | 0.766 | 0.851 | 0.911 | 0.914 |
| 24 | 0.761 | 0.772 | 0.833 | 0.894 | 0.911 | 35 | 0.733 | 0.754 | 0.843 | 0.904 | 0.907 |
| 25 | 0.767 | 0.780 | 0.743 | 0.903 | 0.918 | 36 | 0.720 | 0.740 | 0.833 | 0.895 | 0.898 |
| 26 | 0.772 | 0.787 | 0.851 | 0.910 | 0.923 | 37 | 0.705 | 0.725 | 0.821 | 0.885 | 0.888 |
| 27 | 0.775 | 0.791 | 0.857 | 0.916 | 0.927 | 38 | 0.688 | 0.707 | 0.808 | 0.873 | 0.876 |
| 28 | 0.776 | 0.793 | 0.862 | 0.920 | 0.929 | 39 | 0.669 | 0.688 | 0.793 | 0.859 | 0.863 |
| 29 | 0.775 | 0.794 | 0.865 | 0.922 | 0.930 | 40 | 0.649 | 0.666 | 0.775 | 0.844 | 0.849 |
| 30 | 0.773 | 0.792 | 0.865 | 0.923 | 0.930 | 41 | 0.627 | 0.643 | 0.757 | 0.827 | 0.833 |

Table 2-55: UMA PI values with advance in gestation
 JSUM, J Med Ultrasound Vol.28 No.5 (2001)
 Unit: Age (Weeks)

| Age | 5% | 10% | 50% | 90% | 95% | Age | 5% | 10% | 50% | 90% | 95% |
|-----|-------|-------|-------|-------|-------|-----|-------|-------|-------|-------|-------|
| 20 | 1.118 | 1.144 | 1.390 | 1.620 | 1.688 | 31 | 0.766 | 0.821 | 0.986 | 1.161 | 1.285 |
| 21 | 1.075 | 1.106 | 1.340 | 1.565 | 1.641 | 32 | 0.747 | 0.802 | 0.965 | 1.135 | 1.261 |
| 22 | 1.034 | 1.069 | 1.293 | 1.513 | 1.597 | 33 | 0.731 | 0.785 | 0.947 | 1.112 | 1.238 |
| 23 | 0.996 | 1.034 | 1.249 | 1.464 | 1.554 | 34 | 0.716 | 0.770 | 0.931 | 1.091 | 1.218 |
| 24 | 0.959 | 1.001 | 1.207 | 1.417 | 1.514 | 35 | 0.704 | 0.757 | 0.918 | 1.073 | 1.199 |
| 25 | 0.925 | 0.970 | 1.168 | 1.373 | 1.475 | 36 | 0.694 | 0.746 | 0.907 | 1.057 | 1.182 |
| 26 | 0.893 | 0.941 | 1.131 | 1.331 | 1.438 | 37 | 0.686 | 0.736 | 0.899 | 1.044 | 1.168 |
| 27 | 0.863 | 0.913 | 1.097 | 1.292 | 1.404 | 38 | 0.681 | 0.728 | 0.893 | 1.033 | 1.155 |
| 28 | 0.836 | 0.887 | 1.065 | 1.255 | 1.371 | 39 | 0.677 | 0.722 | 0.890 | 1.025 | 1.145 |
| 29 | 0.810 | 0.863 | 1.036 | 1.221 | 1.341 | 40 | 0.676 | 0.718 | 0.890 | 1.020 | 1.136 |
| 30 | 0.787 | 0.841 | 1.010 | 1.190 | 1.312 | 41 | 0.677 | 0.716 | 0.892 | 1.017 | 1.129 |

Table 2-56: UMA RI values with advance in gestation
 JSUM, J Med Ultrasound Vol.28 No.5 (2001)
 Unit: Age (Weeks)

| Age | 5% | 10% | 50% | 90% | 95% | Age | 5% | 10% | 50% | 90% | 95% |
|-----|-------|-------|-------|-------|-------|-----|-------|-------|-------|-------|-------|
| 20 | 0.698 | 0.722 | 0.778 | 0.820 | 0.846 | 31 | 0.535 | 0.589 | 0.648 | 0.700 | 0.746 |
| 21 | 0.680 | 0.707 | 0.763 | 0.808 | 0.836 | 32 | 0.524 | 0.580 | 0.640 | 0.690 | 0.738 |
| 22 | 0.663 | 0.692 | 0.749 | 0.796 | 0.826 | 33 | 0.513 | 0.573 | 0.632 | 0.681 | 0.730 |
| 23 | 0.646 | 0.679 | 0.735 | 0.785 | 0.816 | 34 | 0.503 | 0.565 | 0.625 | 0.672 | 0.723 |
| 24 | 0.630 | 0.665 | 0.722 | 0.774 | 0.807 | 35 | 0.494 | 0.559 | 0.619 | 0.663 | 0.716 |
| 25 | 0.615 | 0.653 | 0.710 | 0.763 | 0.798 | 36 | 0.485 | 0.552 | 0.613 | 0.654 | 0.708 |
| 26 | 0.600 | 0.640 | 0.698 | 0.752 | 0.788 | 37 | 0.477 | 0.547 | 0.608 | 0.645 | 0.702 |
| 27 | 0.586 | 0.629 | 0.687 | 0.741 | 0.780 | 38 | 0.469 | 0.542 | 0.603 | 0.636 | 0.695 |
| 28 | 0.572 | 0.618 | 0.676 | 0.730 | 0.771 | 39 | 0.462 | 0.538 | 0.599 | 0.628 | 0.688 |
| 29 | 0.559 | 0.608 | 0.666 | 0.720 | 0.762 | 40 | 0.456 | 0.534 | 0.596 | 0.620 | 0.682 |
| 30 | 0.547 | 0.598 | 0.657 | 0.710 | 0.754 | 41 | 0.450 | 0.531 | 0.593 | 0.612 | 0.676 |

Kurtz

Table 2-57: BPD: Kurtz (Fetal Age)
 Journal of Clinical Ultrasound, 8: 319-326, 1980
 Unit: BPD (mm); Age (Days); SD (mm)

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <21 | n/a | — | 40 | 125 | 4 | 60 | 168 | 5 | 80 | 222 | 5 |
| 21 | 84 | 4 | 41 | 127 | 4 | 61 | 170 | 5 | 81 | 225 | 5 |
| 22 | 87 | 4 | 42 | 129 | 4 | 62 | 173 | 5 | 82 | 229 | 5 |
| 23 | 91 | 4 | 43 | 131 | 4 | 63 | 175 | 5 | 83 | 232 | 5 |
| 24 | 93 | 4 | 44 | 133 | 4 | 64 | 178 | 5 | 84 | 235 | 5 |
| 25 | 95 | 4 | 45 | 135 | 4 | 65 | 181 | 5 | 85 | 238 | 5 |
| 26 | 97 | 4 | 46 | 137 | 4 | 66 | 183 | 5 | 86 | 241 | 5 |
| 27 | 99 | 4 | 47 | 139 | 4 | 67 | 186 | 5 | 87 | 244 | 5 |
| 28 | 101 | 4 | 48 | 141 | 4 | 68 | 188 | 5 | 88 | 248 | 5 |
| 29 | 103 | 4 | 49 | 143 | 4 | 69 | 191 | 5 | 89 | 252 | 5 |
| 30 | 105 | 4 | 50 | 145 | 4 | 70 | 194 | 5 | 90 | 257 | 5 |
| 31 | 107 | 4 | 51 | 147 | 4 | 71 | 196 | 5 | 91 | 262 | 5 |
| 32 | 109 | 4 | 52 | 149 | 4 | 72 | 199 | 5 | 92 | 267 | 5 |
| 33 | 111 | 4 | 53 | 151 | 4 | 73 | 201 | 5 | 93 | 272 | 5 |
| 34 | 113 | 4 | 54 | 153 | 4 | 74 | 204 | 5 | 94 | 276 | 5 |
| 35 | 115 | 4 | 55 | 155 | 5 | 75 | 207 | 5 | 95 | 280 | 5 |
| 36 | 117 | 4 | 56 | 157 | 5 | 76 | 210 | 5 | 96 | 284 | 5 |
| 37 | 119 | 4 | 57 | 160 | 5 | 77 | 213 | 5 | 97 | 288 | 5 |
| 38 | 121 | 4 | 58 | 162 | 5 | 78 | 216 | 5 | 98 | 293 | 5 |
| 39 | 123 | 4 | 59 | 165 | 5 | 79 | 219 | 5 | >98 | n/a | — |

Mayden

Table 2-58: IOD: Mayden (Fetal Age)

Am J Obstet Gynecol 144:289, 1982

Unit: Meas (mm); Mean (Weeks)

| Meas | Mean | Meas | Mean | Meas | Mean | Meas | Mean |
|------|------|------|------|------|------|------|------|
| 5 | 11.6 | 11 | 17.9 | 16 | 24.3 | 19 | 32.5 |
| 5 | 11.6 | 12 | 18.4 | 16 | 24.7 | 19 | 33.0 |
| 6 | 12.1 | 12 | 18.9 | 16 | 25.2 | 19 | 33.5 |
| 6 | 12.6 | 12 | 19.4 | 16 | 25.2 | 19 | 34.0 |
| 6 | 12.6 | 13 | 19.4 | 17 | 25.7 | 19 | 34.4 |
| 7 | 13.1 | 13 | 19.9 | 17 | 26.2 | 19 | 35.0 |
| 7 | 13.6 | 13 | 20.4 | 17 | 26.2 | 19 | 35.4 |
| 7 | 13.6 | 13 | 20.4 | 17 | 26.7 | 19 | 35.9 |
| 8 | 14.1 | 14 | 20.9 | 17 | 27.2 | 19 | 36.4 |
| 8 | 14.6 | 14 | 21.3 | 17 | 27.6 | 19 | 36.9 |
| 8 | 14.6 | 14 | 21.3 | 17 | 28.1 | 19 | 37.3 |
| 9 | 15.0 | 14 | 21.8 | 18 | 28.6 | 19 | 37.8 |
| 9 | 15.5 | 14 | 22.3 | 18 | 29.1 | 19 | 38.3 |
| 9 | 15.5 | 15 | 22.3 | 18 | 29.6 | 19 | 38.3 |
| 10 | 16.0 | 15 | 22.8 | 18 | 30.0 | 19 | 39.3 |
| 10 | 16.5 | 15 | 23.3 | 18 | 30.6 | 19 | 39.8 |
| 10 | 16.5 | 15 | 23.3 | 18 | 31.0 | | |
| 10 | 17.0 | 15 | 23.8 | 18 | 31.5 | | |
| 11 | 17.5 | 16 | 24.3 | 18 | 32.0 | | |

Table 2-59: OOD: Mayden (Fetal Age)

Am J Obstet Gynecol 144:289, 1982

Unit: Meas (mm); Mean (Weeks)

| Meas | Mean | Meas | Mean | Meas | Mean | Meas | Mean |
|------|------|------|------|------|------|------|------|
| 13 | 11.6 | 28 | 17.9 | 42 | 24.3 | 52 | 32.5 |
| 14 | 11.6 | 30 | 18.4 | 43 | 24.7 | 53 | 33.0 |
| 15 | 12.1 | 31 | 18.9 | 43 | 25.2 | 54 | 33.5 |
| 16 | 12.6 | 32 | 19.4 | 44 | 25.2 | 54 | 34.0 |
| 17 | 12.6 | 32 | 19.4 | 44 | 25.7 | 54 | 34.4 |
| 17 | 13.1 | 33 | 19.9 | 45 | 26.2 | 55 | 35.0 |
| 18 | 13.6 | 34 | 20.4 | 45 | 26.2 | 55 | 35.4 |
| 19 | 13.6 | 34 | 20.4 | 46 | 26.7 | 56 | 35.9 |
| 20 | 14.1 | 35 | 20.9 | 46 | 27.2 | 56 | 36.4 |
| 21 | 14.6 | 36 | 21.3 | 47 | 27.6 | 57 | 36.9 |
| 21 | 14.6 | 36 | 21.3 | 47 | 28.1 | 57 | 37.3 |
| 22 | 15.0 | 37 | 21.8 | 48 | 28.6 | 58 | 37.8 |
| 23 | 15.5 | 38 | 22.3 | 48 | 29.1 | 58 | 38.3 |
| 24 | 15.5 | 38 | 22.3 | 49 | 29.6 | 58 | 38.3 |
| 25 | 16.0 | 39 | 22.8 | 50 | 30.0 | 59 | 39.3 |
| 25 | 16.5 | 40 | 23.3 | 50 | 30.6 | 59 | 39.8 |
| 26 | 16.5 | 40 | 23.3 | 51 | 31.0 | | |
| 27 | 17.0 | 41 | 23.8 | 51 | 31.5 | | |
| 27 | 17.5 | 41 | 24.3 | 52 | 32.0 | | |

Mercer

Table 2-60: Ft: Mercer (Fetal Age)
Am J Obstet Gynecol, 156: 350-355, 1987
Unit: Meas (mm); Min/Mean/Max (Weeks); Table/Graph Range: 2SD

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|------|------|------|------|------|------|------|
| <10 | n/a | n/a | n/a | 50 | 24.3 | 26.4 | 28.4 |
| 10 | 11.5 | 12.5 | 13.5 | 52 | 24.9 | 27.1 | 29.3 |
| 12 | 12.1 | 13.1 | 14.2 | 54 | 25.7 | 27.9 | 30.1 |
| 14 | 12.7 | 13.8 | 14.9 | 56 | 26.4 | 28.4 | 30.9 |
| 16 | 13.3 | 14.4 | 15.5 | 58 | 27.1 | 29.4 | 31.8 |
| 18 | 13.9 | 15.1 | 16.3 | 60 | 27.8 | 30.2 | 32.6 |
| 20 | 14.5 | 15.7 | 17.0 | 62 | 28.5 | 31.0 | 33.5 |
| 22 | 15.1 | 16.4 | 17.7 | 64 | 29.3 | 31.8 | 34.3 |
| 24 | 15.7 | 17.1 | 18.4 | 66 | 30.0 | 32.6 | 35.2 |
| 26 | 16.3 | 17.7 | 19.1 | 68 | 30.7 | 33.4 | 36.1 |
| 28 | 16.9 | 18.4 | 19.9 | 70 | 31.5 | 34.2 | 36.9 |
| 30 | 17.6 | 19.1 | 20.6 | 72 | 32.2 | 35.0 | 37.8 |
| 32 | 18.2 | 19.8 | 21.4 | 74 | 33.0 | 35.9 | 38.7 |
| 34 | 18.9 | 20.5 | 22.1 | 76 | 33.8 | 36.8 | 39.6 |
| 36 | 19.5 | 21.2 | 22.9 | 78 | 34.5 | 37.5 | 40.5 |
| 38 | 20.2 | 21.9 | 23.7 | 80 | 35.3 | 38.4 | 41.4 |
| 40 | 20.8 | 22.7 | 24.5 | 82 | 36.1 | 39.2 | 42.4 |
| 42 | 21.5 | 23.4 | 25.2 | 84 | 36.9 | 40.1 | 43.3 |
| 44 | 22.2 | 24.1 | 26.0 | 86 | 37.7 | 41.0 | 44.2 |
| 46 | 22.9 | 24.9 | 26.8 | >86 | n/a | n/a | n/a |
| 48 | 23.6 | 25.6 | 27.6 | | | | |

Merz

Table 2-61: AC: Merz (Fetal Age)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Meas (mm); Min/Mean/Max (Weeks); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <56 | n/a | n/a | n/a | 146 | 19w1d | 20w5d | 22w1d |
| 56 | 10w6d | 12w1d | 13w2d | 148 | 19w2d | 20w6d | 22w3d |
| 58 | 11w1d | 12w2d | 13w4d | 150 | 19w4d | 21w1d | 22w4d |
| 60 | 11w2d | 12w4d | 13w5d | 152 | 19w5d | 21w1d | 22w6d |
| 62 | 11w4d | 12w5d | 13w6d | 154 | 19w6d | 21w3d | 23w0d |
| 64 | 11w5d | 12w6d | 14w1d | 156 | 20w1d | 21w4d | 23w1d |
| 66 | 11w6d | 13w1d | 14w2d | 158 | 20w1d | 21w6d | 23w3d |
| 68 | 12w0d | 13w2d | 14w4d | 160 | 20w3d | 22w0d | 23w4d |
| 70 | 12w1d | 13w4d | 14w5d | 162 | 20w4d | 22w1d | 23w6d |
| 72 | 12w3d | 13w4d | 14w6d | 164 | 20w6d | 22w3d | 24w0d |
| 74 | 12w4d | 13w6d | 15w1d | 166 | 21w0d | 22w4d | 24w1d |
| 76 | 12w6d | 14w0d | 15w2d | 168 | 21w1d | 22w6d | 24w3d |
| 78 | 12w6d | 14w1d | 15w4d | 170 | 21w2d | 23w0d | 24w4d |
| 80 | 13w1d | 14w3d | 15w5d | 172 | 21w4d | 23w1d | 24w6d |
| 82 | 13w2d | 14w4d | 15w6d | 174 | 21w5d | 23w2d | 25w0d |
| 84 | 13w4d | 14w6d | 16w1d | 176 | 21w6d | 23w4d | 25w1d |
| 86 | 13w5d | 15w0d | 16w2d | 178 | 22w1d | 23w5d | 25w3d |
| 88 | 13w6d | 15w1d | 16w4d | 180 | 22w1d | 23w6d | 25w4d |
| 90 | 14w0d | 15w3d | 16w5d | 182 | 22w3d | 24w1d | 25w6d |
| 92 | 14w1d | 15w4d | 16w6d | 184 | 22w4d | 24w2d | 26w0d |
| 94 | 14w3d | 15w5d | 17w1d | 186 | 22w6d | 24w4d | 26w1d |
| 96 | 14w4d | 15w6d | 17w2d | 188 | 23w0d | 24w5d | 26w3d |
| 98 | 14w6d | 16w1d | 17w4d | 190 | 23w1d | 24w6d | 26w4d |
| 100 | 14w6d | 16w2d | 17w5d | 192 | 23w2d | 25w0d | 26w6d |
| 102 | 15w1d | 16w4d | 17w6d | 194 | 23w4d | 25w1d | 27w0d |
| 104 | 15w2d | 16w5d | 18w1d | 196 | 23w5d | 25w3d | 27w1d |
| 106 | 15w4d | 16w6d | 18w2d | 198 | 23w6d | 25w4d | 27w3d |
| 108 | 15w5d | 17w1d | 18w3d | 200 | 24w1d | 25w6d | 27w4d |
| 110 | 15w6d | 17w2d | 18w4d | 202 | 24w2d | 26w0d | 27w6d |
| 112 | 16w0d | 17w3d | 18w6d | 204 | 24w3d | 26w1d | 27w6d |
| 114 | 16w1d | 17w4d | 19w0d | 206 | 24w4d | 26w3d | 28w1d |
| 116 | 16w3d | 17w6d | 19w1d | 208 | 24w6d | 26w4d | 28w2d |
| 118 | 16w4d | 18w0d | 19w3d | 210 | 25w0d | 26w6d | 28w4d |
| 120 | 16w6d | 18w1d | 19w4d | 212 | 25w1d | 27w0d | 28w5d |
| 122 | 17w0d | 18w3d | 19w6d | 214 | 25w2d | 27w1d | 28w6d |
| 124 | 17w1d | 18w4d | 20w0d | 216 | 25w4d | 27w2d | 29w1d |
| 126 | 17w2d | 18w6d | 20w1d | 218 | 25w5d | 27w4d | 29w2d |
| 128 | 17w4d | 19w0d | 20w3d | 220 | 25w6d | 27w5d | 29w4d |
| 130 | 17w5d | 19w1d | 20w4d | 222 | 26w1d | 27w6d | 29w5d |
| 132 | 17w6d | 19w2d | 20w6d | 224 | 26w2d | 28w1d | 29w6d |
| 134 | 18w0d | 19w4d | 21w0d | 226 | 26w3d | 28w2d | 30w1d |
| 136 | 18w1d | 19w5d | 21w1d | 228 | 26w4d | 28w4d | 30w2d |
| 138 | 18w3d | 19w6d | 21w3d | 230 | 26w6d | 28w5d | 30w4d |
| 140 | 18w4d | 20w1d | 21w4d | 232 | 27w0d | 28w6d | 30w5d |
| 142 | 18w6d | 20w2d | 21w6d | 234 | 27w1d | 29w0d | 30w6d |
| 144 | 19w0d | 20w4d | 22w0d | 236 | 27w3d | 29w1d | 31w1d |

Table 2-61: AC: Merz (Fetal Age)(Continued)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Meas (mm); Min/Mean/Max (Weeks); Table/Graph Range: 5%:95%

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| 238 | 27w4d | 29w3d | 31w2d | 298 | 33w0d | 35w1d | 37w1d |
| 240 | 27w5d | 29w4d | 31w4d | 300 | 33w1d | 35w2d | 37w3d |
| 242 | 27w6d | 29w6d | 31w5d | 302 | 33w3d | 35w4d | 37w4d |
| 244 | 28w1d | 30w0d | 31w6d | 304 | 33w4d | 35w5d | 37w6d |
| 246 | 28w2d | 30w1d | 32w1d | 306 | 33w5d | 35w6d | 38w0d |
| 248 | 28w3d | 30w3d | 32w2d | 308 | 33w6d | 36w1d | 38w1d |
| 250 | 28w4d | 30w4d | 32w4d | 310 | 34w1d | 36w2d | 38w3d |
| 252 | 28w6d | 30w6d | 32w5d | 312 | 34w2d | 36w4d | 38w4d |
| 254 | 29w0d | 30w6d | 32w6d | 314 | 34w4d | 36w4d | 38w6d |
| 256 | 29w1d | 31w1d | 33w1d | 316 | 34w4d | 36w6d | 39w0d |
| 258 | 29w3d | 31w2d | 33w2d | 318 | 34w6d | 37w0d | 39w1d |
| 260 | 29w4d | 31w4d | 33w4d | 320 | 35w0d | 37w1d | 39w3d |
| 262 | 29w5d | 31w5d | 33w5d | 322 | 35w1d | 37w3d | 39w4d |
| 264 | 29w6d | 31w6d | 33w6d | 324 | 35w3d | 37w4d | 39w6d |
| 266 | 30w1d | 32w1d | 34w1d | 326 | 35w4d | 37w6d | 40w0d |
| 268 | 30w2d | 32w2d | 34w2d | 328 | 35w5d | 38w0d | 40w1d |
| 270 | 30w4d | 32w4d | 34w4d | 330 | 35w6d | 38w1d | 40w3d |
| 272 | 30w4d | 32w5d | 34w5d | 332 | 36w1d | 38w3d | 40w4d |
| 274 | 30w6d | 32w6d | 34w6d | 334 | 36w2d | 38w4d | 40w6d |
| 276 | 31w0d | 33w0d | 35w1d | 336 | 36w4d | 38w5d | 41w0d |
| 278 | 31w1d | 33w1d | 35w2d | 338 | 36w5d | 38w6d | 41w1d |
| 280 | 31w3d | 33w3d | 35w4d | 340 | 36w6d | 39w1d | 41w3d |
| 282 | 31w4d | 33w4d | 35w5d | 342 | 37w0d | 39w2d | 41w4d |
| 284 | 31w5d | 33w6d | 35w6d | 344 | 37w1d | 39w4d | 41w6d |
| 286 | 31w6d | 34w0d | 36w1d | 346 | 37w3d | 39w5d | 42w0d |
| 288 | 32w1d | 34w1d | 36w2d | 348 | 37w4d | 39w6d | 42w1d |
| 290 | 32w2d | 34w3d | 36w4d | >348 | n/a | n/a | n/a |
| 292 | 32w4d | 34w4d | 36w5d | | | | |
| 294 | 32w4d | 34w5d | 36w6d | | | | |
| 296 | 32w6d | 34w6d | 37w1d | | | | |

Table 2-62: AC: Merz (Fetal Growth)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range 5%:95%)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|-----|------|-----|------|-----|------|-----|
| 12.5 | 50 | 62 | 74 | 27.5 | 202 | 222 | 242 |
| 13.0 | 55 | 67 | 80 | 28.0 | 207 | 227 | 247 |
| 13.5 | 60 | 73 | 85 | 28.5 | 212 | 232 | 252 |
| 14.0 | 65 | 78 | 91 | 29.0 | 217 | 237 | 257 |
| 14.5 | 71 | 83 | 96 | 29.5 | 221 | 242 | 263 |
| 15.0 | 76 | 89 | 102 | 30.0 | 226 | 247 | 268 |
| 15.5 | 81 | 94 | 108 | 30.5 | 231 | 252 | 273 |
| 16.0 | 86 | 100 | 114 | 31.0 | 235 | 257 | 278 |
| 16.5 | 91 | 105 | 119 | 31.5 | 240 | 262 | 283 |
| 17.0 | 96 | 111 | 125 | 32.0 | 244 | 266 | 288 |
| 17.5 | 102 | 116 | 131 | 32.5 | 249 | 271 | 293 |
| 18.0 | 107 | 122 | 136 | 33.0 | 253 | 276 | 298 |
| 18.5 | 112 | 127 | 142 | 33.5 | 258 | 280 | 303 |
| 19.0 | 117 | 132 | 148 | 34.0 | 262 | 285 | 308 |
| 19.5 | 122 | 138 | 153 | 34.5 | 266 | 289 | 313 |
| 20.0 | 127 | 143 | 159 | 35.0 | 270 | 294 | 317 |
| 20.5 | 133 | 149 | 165 | 35.5 | 275 | 298 | 322 |
| 21.0 | 138 | 154 | 170 | 36.0 | 279 | 303 | 327 |
| 21.5 | 143 | 159 | 176 | 36.5 | 283 | 307 | 331 |
| 22.0 | 148 | 165 | 181 | 37.0 | 287 | 311 | 336 |
| 22.5 | 153 | 170 | 187 | 37.5 | 290 | 315 | 340 |
| 23.0 | 158 | 175 | 193 | 38.0 | 294 | 319 | 344 |
| 23.5 | 163 | 181 | 198 | 38.5 | 298 | 323 | 348 |
| 24.0 | 168 | 186 | 204 | 39.0 | 301 | 327 | 352 |
| 24.5 | 173 | 191 | 209 | 39.5 | 305 | 331 | 356 |
| 25.0 | 178 | 196 | 215 | 40.0 | 308 | 334 | 360 |
| 25.5 | 183 | 202 | 220 | 40.5 | 311 | 338 | 364 |
| 26.0 | 188 | 207 | 226 | 41.0 | 314 | 341 | 367 |
| 26.5 | 193 | 212 | 231 | 41.5 | 317 | 343 | 370 |
| 27.0 | 198 | 217 | 236 | | | | |

Table 2-63: BPD: Merz (Fetal Age)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: BPD (mm); % Age (Weeks/Days)

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <21 | n/a | n/a | n/a | 62 | 22w1d | 24w1d | 26w1d |
| 21 | 10w5d | 12w1d | 13w5d | 63 | 22w4d | 24w4d | 26w4d |
| 22 | 10w6d | 12w3d | 13w6d | 64 | 22w6d | 24w6d | 26w6d |
| 23 | 11w1d | 12w5d | 14w1d | 65 | 23w1d | 25w1d | 27w1d |
| 24 | 11w4d | 13w0d | 14w4d | 66 | 23w4d | 25w4d | 27w4d |
| 25 | 11w5d | 13w1d | 14w5d | 67 | 23w6d | 25w6d | 27w6d |
| 26 | 12w0d | 13w4d | 15w0d | 68 | 24w1d | 26w1d | 28w2d |
| 27 | 12w1d | 13w6d | 15w3d | 69 | 24w3d | 26w4d | 28w4d |
| 28 | 12w4d | 14w1d | 15w5d | 70 | 24w5d | 26w6d | 28w6d |
| 29 | 12w5d | 14w2d | 15w6d | 71 | 25w1d | 27w1d | 29w2d |
| 30 | 13w0d | 14w4d | 16w1d | 72 | 25w4d | 27w4d | 29w5d |
| 31 | 13w2d | 14w6d | 16w4d | 73 | 25w6d | 27w6d | 30w0d |
| 32 | 13w4d | 15w1d | 16w6d | 74 | 26w1d | 28w2d | 30w3d |
| 33 | 13w6d | 15w3d | 17w0d | 75 | 26w4d | 28w4d | 30w5d |
| 34 | 14w0d | 15w5d | 17w3d | 76 | 26w6d | 29w0d | 31w1d |
| 35 | 14w2d | 16w0d | 17w5d | 77 | 27w1d | 29w3d | 31w4d |
| 36 | 14w4d | 16w2d | 18w0d | 78 | 27w4d | 29w6d | 32w0d |
| 37 | 14w6d | 16w4d | 18w1d | 79 | 27w6d | 30w1d | 32w2d |
| 38 | 15w1d | 16w6d | 18w4d | 80 | 28w2d | 30w4d | 32w5d |
| 39 | 15w3d | 17w1d | 18w6d | 81 | 28w5d | 30w6d | 33w1d |
| 40 | 15w5d | 17w3d | 19w1d | 82 | 29w1d | 31w2d | 33w4d |
| 41 | 15w6d | 17w5d | 19w4d | 83 | 29w4d | 31w5d | 33w6d |
| 42 | 16w1d | 18w0d | 19w6d | 84 | 29w6d | 32w1d | 34w2d |
| 43 | 16w4d | 18w2d | 20w1d | 85 | 30w2d | 32w4d | 34w6d |
| 44 | 16w6d | 18w4d | 20w3d | 86 | 30w5d | 32w6d | 35w1d |
| 45 | 17w1d | 18w6d | 20w5d | 87 | 31w0d | 33w2d | 35w4d |
| 46 | 17w3d | 19w1d | 21w0d | 88 | 31w4d | 33w6d | 36w1d |
| 47 | 17w4d | 19w3d | 21w1d | 89 | 31w6d | 34w1d | 36w4d |
| 48 | 17w6d | 19w5d | 21w4d | 90 | 32w2d | 34w4d | 36w6d |
| 49 | 18w1d | 20w0d | 21w6d | 91 | 32w6d | 35w1d | 37w3d |
| 50 | 18w4d | 20w3d | 22w1d | 92 | 33w1d | 35w4d | 37w6d |
| 51 | 18w6d | 20w5d | 22w4d | 93 | 33w4d | 35w6d | 38w1d |
| 52 | 19w1d | 21w0d | 22w6d | 94 | 34w0d | 36w3d | 38w6d |
| 53 | 19w3d | 21w2d | 23w1d | 95 | 34w4d | 36w6d | 39w2d |
| 54 | 19w5d | 21w4d | 23w4d | 96 | 34w6d | 37w2d | 39w5d |
| 55 | 20w0d | 21w6d | 23w6d | 97 | 35w3d | 37w6d | 40w1d |
| 56 | 20w2d | 22w1d | 24w1d | 98 | 35w6d | 38w2d | 40w5d |
| 57 | 20w4d | 22w4d | 24w3d | 99 | 36w3d | 38w6d | 41w1d |
| 58 | 20w6d | 22w6d | 24w5d | 100 | 36w6d | 39w2d | 41w6d |
| 59 | 21w1d | 23w1d | 25w1d | 101 | 37w2d | 39w6d | 42w2d |
| 60 | 21w4d | 23w4d | 25w4d | 102 | 37w6d | 40w2d | 42w6d |
| 61 | 21w6d | 23w6d | 25w6d | >102 | n/a | n/a | n/a |

Table 2-64: BPD: Merz (Fetal Growth)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range 5%:95%)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|-----|------|-----|------|-----|------|-----|
| 12.5 | 21 | 25 | 29 | 27.5 | 68 | 73 | 78 |
| 13.0 | 23 | 26 | 30 | 28.0 | 69 | 74 | 79 |
| 13.5 | 24 | 28 | 31 | 28.5 | 71 | 76 | 81 |
| 14.0 | 25 | 29 | 33 | 29.0 | 72 | 77 | 82 |
| 14.5 | 27 | 31 | 35 | 29.5 | 73 | 78 | 84 |
| 15.0 | 28 | 32 | 36 | 30.0 | 74 | 80 | 85 |
| 15.5 | 30 | 34 | 38 | 30.5 | 76 | 81 | 86 |
| 16.0 | 31 | 35 | 39 | 31.0 | 77 | 82 | 88 |
| 16.5 | 33 | 37 | 41 | 31.5 | 78 | 83 | 89 |
| 17.0 | 35 | 39 | 43 | 32.0 | 79 | 85 | 90 |
| 17.5 | 36 | 40 | 45 | 32.5 | 80 | 86 | 91 |
| 18.0 | 38 | 42 | 46 | 33.0 | 81 | 87 | 92 |
| 18.5 | 40 | 44 | 48 | 33.5 | 82 | 88 | 93 |
| 19.0 | 41 | 46 | 50 | 34.0 | 83 | 89 | 95 |
| 19.5 | 43 | 47 | 52 | 34.5 | 84 | 90 | 96 |
| 20.0 | 45 | 49 | 53 | 35.0 | 85 | 91 | 97 |
| 20.5 | 46 | 51 | 55 | 35.5 | 86 | 92 | 97 |
| 21.0 | 48 | 52 | 57 | 36.0 | 87 | 92 | 98 |
| 21.5 | 49 | 54 | 59 | 36.5 | 87 | 93 | 99 |
| 22.0 | 51 | 56 | 60 | 37.0 | 88 | 94 | 100 |
| 22.5 | 53 | 57 | 62 | 37.5 | 89 | 95 | 101 |
| 23.0 | 54 | 59 | 64 | 38.0 | 89 | 95 | 101 |
| 23.5 | 56 | 61 | 65 | 38.5 | 90 | 96 | 102 |
| 24.0 | 57 | 62 | 67 | 39.0 | 90 | 96 | 103 |
| 24.5 | 59 | 64 | 69 | 39.5 | 91 | 97 | 103 |
| 25.0 | 61 | 65 | 70 | 40.0 | 91 | 97 | 103 |
| 25.5 | 62 | 67 | 72 | 40.5 | 91 | 97 | 104 |
| 26.0 | 64 | 68 | 73 | 41.0 | 91 | 98 | 104 |
| 26.5 | 65 | 70 | 75 | 41.5 | 92 | 98 | 104 |
| 27.0 | 66 | 71 | 77 | | | | |

Table 2-65: FL: Merz (Fetal Age)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: FL (mm); % Age (Weeks/Days)

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| <10 | n/a | n/a | n/a | 46 | 23w4d | 25w3d | 27w1d |
| 10 | 11w1d | 12w2d | 13w4d | 47 | 24w0d | 25w6d | 27w4d |
| 11 | 11w4d | 12w5d | 13w6d | 48 | 24w3d | 26w1d | 28w0d |
| 12 | 11w6d | 13w0d | 14w1d | 49 | 24w5d | 26w4d | 28w2d |
| 13 | 12w1d | 13w2d | 14w4d | 50 | 25w1d | 26w6d | 28w5d |
| 14 | 12w3d | 13w5d | 15w0d | 51 | 25w4d | 27w2d | 29w1d |
| 15 | 12w5d | 14w0d | 15w2d | 52 | 25w6d | 27w5d | 29w4d |
| 16 | 13w1d | 14w3d | 15w5d | 53 | 26w1d | 28w1d | 30w0d |
| 17 | 13w3d | 14w5d | 16w0d | 54 | 26w4d | 28w4d | 30w4d |
| 18 | 13w6d | 15w1d | 16w3d | 55 | 27w0d | 29w0d | 31w0d |
| 19 | 14w1d | 15w3d | 16w5d | 56 | 27w3d | 29w3d | 31w3d |
| 20 | 14w4d | 15w6d | 17w1d | 57 | 27w6d | 29w6d | 31w6d |
| 21 | 14w6d | 16w1d | 17w3d | 58 | 28w1d | 30w1d | 32w1d |
| 22 | 15w1d | 16w4d | 17w6d | 59 | 28w4d | 30w4d | 32w4d |
| 23 | 15w3d | 16w6d | 18w1d | 60 | 29w0d | 31w0d | 33w0d |
| 24 | 15w6d | 17w1d | 18w4d | 61 | 29w4d | 31w4d | 33w4d |
| 25 | 16w1d | 17w4d | 19w1d | 62 | 29w6d | 31w6d | 33w6d |
| 26 | 16w3d | 17w6d | 19w3d | 63 | 30w2d | 32w2d | 34w2d |
| 27 | 16w6d | 18w2d | 19w6d | 64 | 30w6d | 32w6d | 34w6d |
| 28 | 17w1d | 18w4d | 20w1d | 65 | 31w1d | 33w1d | 35w1d |
| 29 | 17w4d | 19w0d | 20w4d | 66 | 31w4d | 33w4d | 35w4d |
| 30 | 17w6d | 19w3d | 20w6d | 67 | 32w0d | 34w1d | 36w1d |
| 31 | 18w1d | 19w5d | 21w1d | 68 | 32w3d | 34w4d | 36w4d |
| 32 | 18w4d | 20w1d | 21w4d | 69 | 32w6d | 35w0d | 37w1d |
| 33 | 18w6d | 20w4d | 22w1d | 70 | 33w2d | 35w3d | 37w4d |
| 34 | 19w1d | 20w6d | 22w3d | 71 | 33w6d | 35w6d | 38w0d |
| 35 | 19w4d | 21w1d | 22w6d | 72 | 34w1d | 36w2d | 38w3d |
| 36 | 20w0d | 21w4d | 23w1d | 73 | 34w4d | 36w6d | 39w0d |
| 37 | 20w2d | 21w6d | 23w4d | 74 | 35w1d | 37w2d | 39w4d |
| 38 | 20w5d | 22w2d | 23w6d | 75 | 35w4d | 37w5d | 39w6d |
| 39 | 21w0d | 22w5d | 24w3d | 76 | 36w0d | 38w1d | 40w3d |
| 40 | 21w3d | 23w1d | 24w6d | 77 | 36w4d | 38w5d | 40w6d |
| 41 | 21w5d | 23w3d | 25w1d | 78 | 37w0d | 39w1d | 41w3d |
| 42 | 22w1d | 23w6d | 25w4d | 79 | 37w3d | 39w4d | 41w6d |
| 43 | 22w4d | 24w1d | 25w6d | 80 | 37w6d | 40w1d | 42w2d |
| 44 | 22w6d | 24w4d | 26w3d | >80 | n/a | n/a | n/a |
| 45 | 23w1d | 25w0d | 26w6d | | | | |

Table 2-66: FL: Merz (Fetal Growth)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range 5%:95%)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|-----|------|-----|------|-----|------|-----|
| 12.5 | 6 | 9 | 12 | 27.5 | 48 | 52 | 57 |
| 13.0 | 8 | 11 | 14 | 28.0 | 49 | 53 | 58 |
| 13.5 | 10 | 13 | 16 | 28.5 | 50 | 55 | 59 |
| 14.0 | 11 | 15 | 18 | 29.0 | 51 | 56 | 60 |
| 14.5 | 13 | 16 | 20 | 29.5 | 52 | 57 | 61 |
| 15.0 | 15 | 18 | 21 | 30.0 | 53 | 58 | 62 |
| 15.5 | 16 | 20 | 23 | 30.5 | 54 | 59 | 63 |
| 16.0 | 18 | 21 | 25 | 31.0 | 55 | 60 | 64 |
| 16.5 | 19 | 23 | 26 | 31.5 | 56 | 61 | 66 |
| 17.0 | 21 | 24 | 28 | 32.0 | 57 | 62 | 67 |
| 17.5 | 22 | 26 | 29 | 32.5 | 58 | 63 | 68 |
| 18.0 | 24 | 27 | 31 | 33.0 | 59 | 64 | 69 |
| 18.5 | 25 | 29 | 32 | 33.5 | 60 | 65 | 70 |
| 19.0 | 27 | 30 | 34 | 34.0 | 61 | 66 | 71 |
| 19.5 | 28 | 32 | 35 | 34.5 | 62 | 67 | 72 |
| 20.0 | 29 | 33 | 37 | 35.0 | 63 | 68 | 73 |
| 20.5 | 31 | 35 | 38 | 35.5 | 64 | 69 | 74 |
| 21.0 | 32 | 36 | 40 | 36.0 | 65 | 70 | 74 |
| 21.5 | 33 | 37 | 41 | 36.5 | 66 | 70 | 75 |
| 22.0 | 35 | 39 | 42 | 37.0 | 66 | 71 | 76 |
| 22.5 | 36 | 40 | 44 | 37.5 | 67 | 72 | 77 |
| 23.0 | 37 | 41 | 45 | 38.0 | 68 | 73 | 78 |
| 23.5 | 39 | 43 | 46 | 38.5 | 69 | 74 | 79 |
| 24.0 | 40 | 44 | 48 | 39.0 | 69 | 74 | 79 |
| 24.5 | 41 | 45 | 49 | 39.5 | 70 | 75 | 80 |
| 25.0 | 42 | 46 | 50 | 40.0 | 71 | 76 | 81 |
| 25.5 | 43 | 48 | 52 | 40.5 | 71 | 76 | 81 |
| 26.0 | 45 | 49 | 53 | 41.0 | 72 | 77 | 82 |
| 26.5 | 46 | 50 | 54 | 41.5 | 72 | 77 | 83 |
| 27.0 | 47 | 51 | 55 | | | | |

Table 2-67: HC: Merz (Fetal Age)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: HC (mm); % Age (Weeks/Days)

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| >72 | n/a | n/a | n/a | 172 | 17w6d | 19w2d | 20w6d |
| 72 | 11w0d | 12w1d | 13w1d | 174 | 17w6d | 19w3d | 20w6d |
| 74 | 11w1d | 12w2d | 13w4d | 176 | 18w0d | 19w4d | 21w1d |
| 76 | 11w1d | 12w3d | 13w4d | 178 | 18w1d | 19w6d | 21w3d |
| 78 | 11w2d | 12w4d | 13w5d | 180 | 18w2d | 19w6d | 21w4d |
| 80 | 11w4d | 12w5d | 13w6d | 182 | 18w4d | 20w1d | 21w5d |
| 82 | 11w4d | 12w6d | 14w0d | 184 | 18w4d | 20w1d | 21w6d |
| 84 | 11w5d | 12w6d | 14w1d | 186 | 18w6d | 20w3d | 22w0d |
| 86 | 11w6d | 13w1d | 14w2d | 188 | 19w0d | 20w4d | 22w1d |
| 88 | 12w0d | 13w1d | 14w3d | 190 | 19w1d | 20w5d | 22w2d |
| 90 | 12w1d | 13w2d | 14w4d | 192 | 19w2d | 20w6d | 22w4d |
| 92 | 12w2d | 13w4d | 14w5d | 194 | 19w4d | 21w1d | 22w5d |
| 94 | 12w3d | 13w4d | 14w6d | 196 | 19w4d | 21w1d | 22w6d |
| 96 | 12w4d | 13w5d | 14w6d | 198 | 19w5d | 21w3d | 23w0d |
| 98 | 12w5d | 13w6d | 15w1d | 200 | 19w6d | 21w4d | 23w2d |
| 100 | 12w6d | 14w0d | 15w1d | 202 | 20w0d | 21w5d | 23w3d |
| 102 | 12w6d | 14w1d | 15w4d | 204 | 20w1d | 21w6d | 23w4d |
| 104 | 13w0d | 14w2d | 15w4d | 206 | 20w3d | 22w1d | 23w6d |
| 106 | 13w1d | 14w3d | 15w5d | 208 | 20w4d | 22w1d | 23w6d |
| 108 | 13w2d | 14w4d | 15w6d | 210 | 20w5d | 22w3d | 24w1d |
| 110 | 13w3d | 14w5d | 16w0d | 212 | 20w6d | 22w4d | 24w2d |
| 112 | 13w4d | 14w6d | 16w1d | 214 | 21w0d | 22w5d | 24w3d |
| 114 | 13w5d | 15w0d | 16w2d | 216 | 21w1d | 22w6d | 24w4d |
| 116 | 13w6d | 15w1d | 16w3d | 218 | 21w3d | 23w1d | 24w6d |
| 118 | 14w0d | 15w2d | 16w4d | 220 | 21w4d | 23w2d | 25w0d |
| 120 | 14w1d | 15w3d | 16w5d | 222 | 21w6d | 23w4d | 25w1d |
| 122 | 14w1d | 15w4d | 17w0d | 224 | 21w6d | 23w4d | 25w2d |
| 124 | 14w2d | 15w5d | 17w1d | 226 | 22w1d | 23w6d | 25w4d |
| 126 | 14w3d | 15w6d | 17w1d | 228 | 22w1d | 24w0d | 25w6d |
| 128 | 14w4d | 16w0d | 17w3d | 230 | 22w3d | 24w1d | 26w0d |
| 130 | 14w5d | 16w1d | 17w4d | 232 | 22w4d | 24w3d | 26w1d |
| 132 | 14w6d | 16w2d | 17w5d | 234 | 22w5d | 24w4d | 26w2d |
| 134 | 15w0d | 16w3d | 17w6d | 236 | 22w6d | 24w5d | 26w4d |
| 136 | 15w1d | 16w4d | 18w0d | 238 | 23w1d | 24w6d | 26w5d |
| 138 | 15w2d | 16w5d | 18w1d | 240 | 23w2d | 25w1d | 26w6d |
| 140 | 15w4d | 16w6d | 18w2d | 242 | 23w4d | 25w2d | 27w1d |
| 142 | 15w4d | 17w0d | 18w3d | 244 | 23w5d | 25w4d | 27w2d |
| 144 | 15w6d | 17w1d | 18w4d | 246 | 23w6d | 25w5d | 27w4d |
| 146 | 15w6d | 17w2d | 18w5d | 248 | 24w1d | 25w6d | 27w5d |
| 148 | 16w0d | 17w4d | 19w0d | 250 | 24w1d | 26w0d | 27w6d |
| 150 | 16w1d | 17w4d | 19w1d | 252 | 24w3d | 26w1d | 28w0d |
| 152 | 16w2d | 17w6d | 19w2d | 254 | 24w4d | 26w3d | 28w1d |
| 154 | 16w3d | 17w6d | 19w3d | 256 | 24w6d | 26w4d | 28w3d |
| 156 | 16w4d | 18w1d | 19w4d | 258 | 25w0d | 26w6d | 28w4d |
| 158 | 16w5d | 18w1d | 19w5d | 260 | 25w1d | 27w0d | 28w6d |
| 160 | 16w6d | 18w3d | 19w6d | 262 | 25w3d | 27w1d | 29w0d |
| 162 | 17w0d | 18w4d | 20w0d | 264 | 25w4d | 27w3d | 29w1d |
| 164 | 17w1d | 18w5d | 20w1d | 266 | 25w6d | 27w4d | 29w3d |
| 166 | 17w2d | 18w6d | 20w2d | 268 | 26w0d | 27w6d | 29w4d |
| 168 | 17w4d | 19w0d | 20w4d | 270 | 26w1d | 28w1d | 30w0d |
| 170 | 17w4d | 19w1d | 20w4d | 272 | 26w3d | 28w2d | 30w1d |

Table 2-67: HC: Merz (Fetal Age)(Continued)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: HC (mm); % Age (Weeks/Days)

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| 274 | 26w4d | 28w4d | 30w3d | 322 | 32w0d | 34w1d | 36w1d |
| 276 | 26w6d | 28w5d | 30w4d | 324 | 32w2d | 34w3d | 36w4d |
| 278 | 27w0d | 28w6d | 30w6d | 326 | 32w4d | 34w5d | 36w6d |
| 280 | 27w1d | 29w1d | 31w0d | 328 | 32w6d | 34w6d | 37w0d |
| 282 | 27w3d | 29w2d | 31w1d | 330 | 33w1d | 35w1d | 37w2d |
| 284 | 27w5d | 29w4d | 31w4d | 332 | 33w2d | 35w4d | 37w5d |
| 286 | 27w6d | 29w6d | 31w5d | 334 | 33w4d | 35w6d | 38w0d |
| 288 | 28w1d | 30w0d | 31w6d | 336 | 33w6d | 36w1d | 38w2d |
| 290 | 28w2d | 30w1d | 32w1d | 338 | 34w1d | 36w3d | 38w4d |
| 292 | 28w4d | 30w4d | 32w3d | 340 | 34w3d | 36w4d | 38w6d |
| 294 | 28w6d | 30w5d | 32w4d | 342 | 34w5d | 36w6d | 39w1d |
| 296 | 29w0d | 30w6d | 32w6d | 344 | 35w0d | 37w1d | 39w3d |
| 298 | 29w1d | 31w1d | 33w0d | 346 | 35w2d | 37w4d | 39w5d |
| 300 | 29w3d | 31w3d | 33w3d | 348 | 35w4d | 37w6d | 40w1d |
| 302 | 29w4d | 31w4d | 33w4d | 350 | 35w6d | 38w1d | 40w4d |
| 304 | 29w6d | 31w6d | 33w6d | 352 | 36w1d | 38w4d | 40w6d |
| 306 | 30w1d | 32w1d | 34w1d | 354 | 36w4d | 38w6d | 41w1d |
| 308 | 30w2d | 32w2d | 34w2d | 356 | 36w6d | 39w1d | 41w3d |
| 310 | 30w4d | 32w4d | 34w4d | 358 | 37w1d | 39w4d | 41w6d |
| 312 | 30w6d | 32w6d | 34w6d | 360 | 37w4d | 39w6d | 42w1d |
| 314 | 31w1d | 33w1d | 35w1d | 362 | 37w6d | 40w1d | 42w3d |
| 316 | 31w3d | 33w3d | 35w3d | 364 | 38w1d | 40w4d | 42w6d |
| 318 | 31w4d | 33w4d | 35w4d | >364 | n/a | n/a | n/a |
| 320 | 31w6d | 33w6d | 36w0d | | | | |

Table 2-68: HC: Merz (Fetal Growth)
 Habilitationsschrift, Mainz University Women's Hospital, 1988
 Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range 5%:95%)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|-----|------|-----|------|-----|------|-----|
| 12.5 | 80 | 92 | 104 | 27.5 | 253 | 268 | 284 |
| 13.0 | 84 | 96 | 108 | 28.0 | 258 | 273 | 289 |
| 13.5 | 89 | 101 | 113 | 28.5 | 263 | 278 | 294 |
| 14.0 | 94 | 106 | 119 | 29.0 | 268 | 283 | 299 |
| 14.5 | 100 | 112 | 124 | 29.5 | 272 | 288 | 303 |
| 15.0 | 105 | 118 | 130 | 30.0 | 277 | 292 | 308 |
| 15.5 | 111 | 124 | 137 | 30.5 | 281 | 297 | 313 |
| 16.0 | 117 | 130 | 143 | 31.0 | 285 | 301 | 317 |
| 16.5 | 123 | 136 | 149 | 31.5 | 289 | 305 | 321 |
| 17.0 | 130 | 143 | 156 | 32.0 | 293 | 309 | 325 |
| 17.5 | 136 | 149 | 162 | 32.5 | 297 | 313 | 329 |
| 18.0 | 142 | 155 | 168 | 33.0 | 300 | 316 | 333 |
| 18.5 | 148 | 162 | 175 | 33.5 | 303 | 320 | 336 |
| 19.0 | 155 | 168 | 181 | 34.0 | 307 | 323 | 340 |
| 19.5 | 161 | 174 | 188 | 34.5 | 310 | 326 | 343 |
| 20.0 | 167 | 181 | 194 | 35.0 | 313 | 329 | 346 |
| 20.5 | 173 | 187 | 201 | 35.5 | 315 | 332 | 349 |
| 21.0 | 180 | 193 | 207 | 36.0 | 318 | 335 | 352 |
| 21.5 | 186 | 200 | 214 | 36.5 | 320 | 337 | 354 |
| 22.0 | 192 | 206 | 220 | 37.0 | 322 | 339 | 356 |
| 22.5 | 198 | 212 | 226 | 37.5 | 324 | 341 | 359 |
| 23.0 | 204 | 218 | 232 | 38.0 | 326 | 343 | 361 |
| 23.5 | 210 | 224 | 238 | 38.5 | 327 | 345 | 362 |
| 24.0 | 216 | 230 | 244 | 39.0 | 329 | 346 | 364 |
| 24.5 | 221 | 236 | 250 | 39.5 | 330 | 348 | 365 |
| 25.0 | 227 | 241 | 256 | 40.0 | 331 | 349 | 366 |
| 25.5 | 232 | 247 | 262 | 40.5 | 332 | 349 | 367 |
| 26.0 | 238 | 253 | 267 | 41.0 | 332 | 350 | 368 |
| 26.5 | 243 | 258 | 273 | 41.5 | 332 | 350 | 369 |
| 27.0 | 248 | 263 | 278 | | | | |

Moore

Table 2-69: AFI: Moore
Unit: Age (Days); Min/Max (mm); Table/Graph Range (2.5%: 97.5%)

| Age | Min | Max | CRL | Age | SD | CRL | Age | SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| 16 | 73 | 201 | 25 | 89 | 240 | 34 | 72 | 278 |
| 17 | 77 | 211 | 26 | 89 | 242 | 35 | 70 | 279 |
| 18 | 80 | 220 | 27 | 85 | 245 | 36 | 68 | 279 |
| 19 | 83 | 225 | 28 | 86 | 249 | 37 | 66 | 275 |
| 20 | 86 | 230 | 29 | 84 | 254 | 38 | 65 | 269 |
| 21 | 88 | 233 | 30 | 82 | 258 | 39 | 64 | 255 |
| 22 | 89 | 235 | 31 | 79 | 263 | 40 | 63 | 240 |
| 23 | 90 | 237 | 32 | 77 | 269 | 41 | 63 | 216 |
| 24 | 90 | 238 | 33 | 74 | 274 | 42 | 63 | 192 |

Nelson

Table 2-70: CRL: Nelson (Fetal Age)
Journal of Clinical Ultrasound, 9: 67-70, 1981
Unit: CRL (mm); GA (Days)

| CRL | GA | CRL | GA | CRL | GA | CRL | GA | CRL | GA |
|-----|----|-----|----|-----|----|-----|-----|-----|-----|
| 14 | 59 | 34 | 71 | 54 | 83 | 74 | 95 | 94 | 107 |
| 15 | 60 | 35 | 72 | 55 | 84 | 75 | 96 | 95 | 108 |
| 16 | 61 | 36 | 73 | 56 | 85 | 76 | 97 | 96 | 109 |
| 17 | 61 | 37 | 73 | 57 | 85 | 77 | 97 | 97 | 109 |
| 18 | 62 | 38 | 74 | 58 | 86 | 78 | 98 | 98 | 110 |
| 19 | 62 | 39 | 74 | 59 | 86 | 79 | 98 | 99 | 110 |
| 20 | 63 | 40 | 75 | 60 | 87 | 80 | 99 | 100 | 111 |
| 21 | 64 | 41 | 76 | 61 | 88 | 81 | 100 | 101 | 112 |
| 22 | 64 | 42 | 76 | 62 | 88 | 82 | 100 | 102 | 112 |
| 23 | 65 | 43 | 77 | 63 | 89 | 83 | 101 | 103 | 113 |
| 24 | 65 | 44 | 77 | 64 | 89 | 84 | 101 | 104 | 113 |
| 25 | 66 | 45 | 78 | 65 | 90 | 85 | 102 | 105 | 114 |
| 26 | 67 | 46 | 79 | 66 | 91 | 86 | 103 | 106 | 115 |
| 27 | 67 | 47 | 79 | 67 | 91 | 87 | 103 | 107 | 115 |
| 28 | 68 | 48 | 80 | 68 | 92 | 88 | 104 | 108 | 116 |
| 29 | 68 | 49 | 80 | 69 | 92 | 89 | 104 | 109 | 116 |
| 30 | 69 | 50 | 81 | 70 | 93 | 90 | 105 | 110 | 117 |
| 31 | 70 | 51 | 82 | 71 | 94 | 91 | 106 | 111 | 118 |
| 32 | 70 | 52 | 82 | 72 | 94 | 92 | 106 | | |
| 33 | 71 | 53 | 83 | 73 | 95 | 93 | 107 | | |

Osaka

Table 2-71: BPD: Osaka (Fetal Age)
Osaka University Method 1989, 3 by Univ. Osaka

Unit: BPD (mm); Age (Days); SD (mm)

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <13 | n/a | — | 33 | 107 | 2.4 | 54 | 152 | 3.0 | 75 | 203 | 3.5 |
| 13 | 70 | 1.9 | 34 | 109 | 2.5 | 55 | 154 | 3.0 | 76 | 206 | 3.5 |
| 14 | 71 | 1.9 | 35 | 112 | 2.5 | 56 | 157 | 3.0 | 77 | 209 | 3.5 |
| 15 | 73 | 1.9 | 36 | 114 | 2.5 | 57 | 159 | 3.1 | 78 | 212 | 3.5 |
| 16 | 75 | 1.9 | 37 | 116 | 2.5 | 58 | 161 | 3.1 | 79 | 214 | 3.6 |
| 17 | 77 | 2.0 | 38 | 118 | 2.6 | 59 | 164 | 3.1 | 80 | 217 | 3.6 |
| 18 | 78 | 2.0 | 39 | 120 | 2.6 | 60 | 166 | 3.1 | 81 | 220 | 3.6 |
| 19 | 80 | 2.0 | 40 | 122 | 2.6 | 61 | 168 | 3.2 | 82 | 224 | 3.6 |
| 20 | 82 | 2.1 | 41 | 124 | 2.7 | 62 | 171 | 3.2 | 83 | 227 | 3.6 |
| 21 | 84 | 2.1 | 42 | 126 | 2.7 | 63 | 173 | 3.2 | 84 | 230 | 3.7 |
| 22 | 86 | 2.1 | 43 | 128 | 2.7 | 64 | 175 | 3.2 | 85 | 234 | 3.7 |
| 23 | 88 | 2.1 | 44 | 130 | 2.7 | 65 | 178 | 3.3 | 86 | 237 | 3.7 |
| 24 | 90 | 2.2 | 45 | 132 | 2.8 | 66 | 180 | 3.3 | 87 | 238 | 3.7 |
| 25 | 92 | 2.2 | 46 | 135 | 2.8 | 67 | 182 | 3.3 | 88 | 245 | 3.7 |
| 26 | 94 | 2.2 | 47 | 137 | 2.8 | 68 | 185 | 3.3 | 89 | 249 | 3.8 |
| 27 | 96 | 2.3 | 48 | 139 | 2.8 | 69 | 187 | 3.3 | 90 | 254 | 3.8 |
| 28 | 98 | 2.3 | 49 | 141 | 2.9 | 70 | 190 | 3.4 | 91 | 259 | 3.8 |
| 29 | 99 | 2.3 | 50 | 143 | 2.9 | 71 | 193 | 3.4 | 92 | 265 | 3.8 |
| 30 | 101 | 2.3 | 51 | 145 | 2.9 | 72 | 195 | 3.4 | 93 | 273 | 3.9 |
| 31 | 103 | 2.4 | 52 | 148 | 2.9 | 73 | 198 | 3.4 | 94 | 280 | 3.9 |
| 32 | 105 | 2.4 | 53 | 150 | 3.0 | 74 | 200 | 3.5 | >94 | n/a | — |

Table 2-72: CRL: Osaka (Fetal Age)
Osaka University Method 1989, 3 by Univ. Osaka

Unit: CRL (mm); Age (Days); SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <9 | n/a | — | 23 | 65 | 4.0 | 38 | 75 | 5.5 | 53 | 84 | 6.9 |
| 9 | 50 | 1.7 | 24 | 66 | 4.1 | 39 | 76 | 5.7 | 54 | 85 | 7.0 |
| 10 | 52 | 2.0 | 25 | 66 | 4.1 | 40 | 76 | 5.7 | 55 | 85 | 7.0 |
| 11 | 53 | 2.2 | 26 | 67 | 4.3 | 41 | 77 | 5.8 | 56 | 86 | 7.2 |
| 12 | 55 | 2.5 | 27 | 68 | 4.5 | 42 | 77 | 5.8 | 57 | 86 | 7.2 |
| 13 | 56 | 2.6 | 28 | 69 | 4.6 | 43 | 78 | 6.0 | 58 | 87 | 7.3 |
| 14 | 57 | 2.8 | 29 | 69 | 4.6 | 44 | 79 | 6.1 | 59 | 87 | 7.3 |
| 15 | 58 | 2.9 | 30 | 70 | 4.8 | 45 | 79 | 6.1 | 60 | 88 | 7.5 |
| 16 | 59 | 3.1 | 31 | 71 | 4.9 | 46 | 80 | 6.3 | 61 | 89 | 7.6 |
| 17 | 60 | 3.2 | 32 | 71 | 4.9 | 47 | 80 | 6.3 | 62 | 89 | 7.6 |
| 18 | 61 | 3.4 | 33 | 72 | 5.1 | 48 | 81 | 6.4 | 63 | 90 | 7.8 |
| 19 | 62 | 3.5 | 34 | 73 | 5.2 | 49 | 82 | 6.6 | >63 | n/a | — |
| 20 | 63 | 3.7 | 35 | 73 | 5.2 | 50 | 83 | 6.7 | | | |
| 21 | 63 | 3.7 | 36 | 74 | 5.4 | 51 | 83 | 6.7 | | | |
| 22 | 64 | 3.8 | 37 | 74 | 5.4 | 52 | 83 | 6.7 | | | |

Table 2-73: EFW: Osaka (Fetal Age)
 Osaka University Method 1989, 3 by Univ. Osaka
 Unit: EFW (grams); Age (Days); SD (grams)

| EFW | Age | SD | EFW | Age | SD | EFW | Age | SD | EFW | Age | SD |
|------|-----|----|------|-----|-----|------|-----|-----|-------|-----|-----|
| <137 | n/a | — | 590 | 160 | 81 | 1420 | 203 | 171 | 2360 | 242 | 268 |
| 137 | 112 | 29 | 600 | 160 | 81 | 1440 | 204 | 174 | 2380 | 243 | 271 |
| 140 | 113 | 29 | 610 | 161 | 83 | 1460 | 205 | 176 | 2400 | 244 | 274 |
| 150 | 115 | 29 | 620 | 162 | 85 | 1480 | 206 | 178 | 2420 | 245 | 276 |
| 160 | 116 | 30 | 630 | 162 | 85 | 1500 | 207 | 181 | 2440 | 245 | 276 |
| 170 | 118 | 30 | 640 | 163 | 87 | 1520 | 208 | 183 | 2460 | 246 | 279 |
| 180 | 120 | 31 | 650 | 164 | 89 | 1540 | 209 | 185 | 2480 | 247 | 282 |
| 190 | 121 | 32 | 660 | 164 | 89 | 1560 | 210 | 188 | 2500 | 248 | 285 |
| 200 | 123 | 33 | 670 | 165 | 91 | 1580 | 210 | 188 | 2520 | 249 | 288 |
| 210 | 124 | 34 | 680 | 165 | 91 | 1600 | 211 | 190 | 2540 | 249 | 288 |
| 220 | 126 | 35 | 690 | 166 | 92 | 1620 | 212 | 192 | 2560 | 250 | 290 |
| 230 | 127 | 36 | 700 | 167 | 94 | 1640 | 213 | 195 | 2580 | 251 | 293 |
| 240 | 128 | 37 | 720 | 168 | 96 | 1660 | 214 | 197 | 2600 | 252 | 296 |
| 250 | 130 | 39 | 740 | 169 | 98 | 1680 | 215 | 200 | 2620 | 253 | 299 |
| 260 | 131 | 40 | 760 | 170 | 100 | 1700 | 216 | 202 | 2640 | 254 | 302 |
| 270 | 132 | 41 | 780 | 171 | 102 | 1720 | 216 | 202 | 2660 | 254 | 302 |
| 280 | 133 | 42 | 800 | 173 | 106 | 1740 | 217 | 204 | 2680 | 255 | 305 |
| 290 | 134 | 43 | 820 | 174 | 108 | 1760 | 218 | 207 | 2700 | 256 | 308 |
| 300 | 135 | 44 | 840 | 175 | 110 | 1780 | 219 | 209 | 2720 | 257 | 311 |
| 310 | 136 | 45 | 860 | 176 | 112 | 1800 | 220 | 212 | 2740 | 258 | 314 |
| 320 | 137 | 46 | 880 | 177 | 114 | 1820 | 220 | 212 | 2760 | 259 | 317 |
| 330 | 138 | 48 | 900 | 178 | 116 | 1840 | 221 | 214 | 2780 | 259 | 317 |
| 340 | 139 | 49 | 920 | 179 | 118 | 1860 | 222 | 217 | 2800 | 260 | 320 |
| 350 | 140 | 50 | 940 | 180 | 120 | 1880 | 223 | 219 | 2820 | 261 | 323 |
| 360 | 141 | 51 | 960 | 181 | 123 | 1900 | 224 | 222 | 2840 | 262 | 326 |
| 370 | 142 | 53 | 980 | 182 | 125 | 1920 | 224 | 222 | 2860 | 263 | 329 |
| 380 | 143 | 54 | 1000 | 183 | 127 | 1940 | 225 | 224 | 2880 | 264 | 332 |
| 390 | 144 | 56 | 1020 | 185 | 131 | 1960 | 226 | 227 | 2900 | 265 | 335 |
| 400 | 145 | 57 | 1040 | 186 | 133 | 1980 | 227 | 229 | 2920 | 266 | 339 |
| 410 | 146 | 58 | 1060 | 187 | 135 | 2000 | 228 | 232 | 2940 | 266 | 339 |
| 420 | 147 | 60 | 1080 | 188 | 138 | 2020 | 229 | 234 | 2960 | 267 | 342 |
| 430 | 148 | 61 | 1100 | 189 | 140 | 2040 | 229 | 234 | 2980 | 268 | 345 |
| 440 | 149 | 63 | 1120 | 190 | 142 | 2060 | 230 | 237 | 3000 | 269 | 348 |
| 450 | 149 | 63 | 1140 | 191 | 144 | 2080 | 231 | 239 | 3020 | 270 | 352 |
| 460 | 150 | 65 | 1160 | 192 | 146 | 2100 | 232 | 242 | 3040 | 271 | 355 |
| 470 | 151 | 66 | 1180 | 193 | 149 | 2120 | 233 | 244 | 3060 | 272 | 358 |
| 480 | 152 | 68 | 1200 | 194 | 151 | 2140 | 233 | 244 | 3080 | 273 | 362 |
| 490 | 153 | 69 | 1220 | 195 | 153 | 2160 | 234 | 247 | 3100 | 274 | 365 |
| 500 | 153 | 69 | 1240 | 195 | 153 | 2180 | 235 | 250 | 3120 | 275 | 369 |
| 510 | 154 | 71 | 1260 | 196 | 155 | 2200 | 236 | 252 | 3140 | 276 | 372 |
| 520 | 155 | 73 | 1280 | 197 | 158 | 2220 | 236 | 252 | 3160 | 277 | 376 |
| 530 | 155 | 73 | 1300 | 198 | 160 | 2240 | 237 | 255 | 3180 | 278 | 379 |
| 540 | 156 | 74 | 1320 | 199 | 162 | 2260 | 238 | 257 | 3200 | 279 | 383 |
| 550 | 157 | 76 | 1340 | 200 | 164 | 2280 | 239 | 260 | 3220 | 280 | 387 |
| 560 | 157 | 76 | 1360 | 201 | 167 | 2300 | 240 | 263 | >3220 | n/a | — |
| 570 | 158 | 78 | 1380 | 202 | 169 | 2320 | 241 | 265 | | | |
| 580 | 159 | 80 | 1400 | 203 | 171 | 2340 | 241 | 265 | | | |

Table 2-74: FL: Osaka (Fetal Age)
Osaka University Method 1989, 3 by Univ. Osaka
Unit: FL (mm); Age (Days); SD (mm)

| FL | Age | SD | FL | Age | SD | FL | Age | SD | FL | Age | SD |
|----|-----|-----|----|-----|-----|----|-----|-----|-----|-----|-----|
| <9 | n/a | — | 25 | 127 | 2.3 | 42 | 172 | 2.6 | 59 | 227 | 2.9 |
| 9 | 91 | 2.1 | 26 | 130 | 2.3 | 43 | 175 | 2.6 | 60 | 230 | 2.9 |
| 10 | 93 | 2.1 | 27 | 132 | 2.3 | 44 | 178 | 2.6 | 61 | 235 | 2.9 |
| 11 | 95 | 2.1 | 28 | 135 | 2.4 | 45 | 181 | 2.6 | 62 | 239 | 2.9 |
| 12 | 97 | 2.2 | 29 | 137 | 2.4 | 46 | 184 | 2.6 | 63 | 242 | 3.0 |
| 13 | 99 | 2.2 | 30 | 140 | 2.4 | 47 | 186 | 2.6 | 64 | 247 | 3.0 |
| 14 | 102 | 2.2 | 31 | 142 | 2.4 | 48 | 190 | 2.7 | 65 | 250 | 3.0 |
| 15 | 104 | 2.2 | 32 | 145 | 2.4 | 49 | 193 | 2.7 | 66 | 255 | 3.0 |
| 16 | 106 | 2.2 | 33 | 147 | 2.4 | 50 | 196 | 2.7 | 67 | 258 | 3.0 |
| 17 | 108 | 2.2 | 34 | 150 | 2.4 | 51 | 199 | 2.7 | 68 | 260 | 3.1 |
| 18 | 110 | 2.2 | 35 | 152 | 2.5 | 52 | 202 | 2.6 | 69 | 269 | 3.1 |
| 19 | 113 | 2.2 | 36 | 155 | 2.5 | 53 | 205 | 2.8 | 70 | 274 | 3.1 |
| 20 | 115 | 2.3 | 37 | 158 | 2.5 | 54 | 209 | 2.8 | 71 | 279 | 3.2 |
| 21 | 118 | 2.3 | 38 | 162 | 2.5 | 55 | 212 | 2.8 | >71 | n/a | — |
| 22 | 120 | 2.3 | 39 | 163 | 2.5 | 56 | 216 | 2.8 | | | |
| 23 | 122 | 2.3 | 40 | 166 | 2.5 | 57 | 220 | 2.8 | | | |
| 24 | 125 | 2.3 | 41 | 169 | 2.6 | 58 | 223 | 2.9 | | | |

Table 2-75: FTA: Osaka (Fetal Age)
Osaka University Method 1989, 3 by Univ. Osaka
Unit: FTA (mm²); Age (Days); SD (mm²)

| FTA | Age | SD | FTA | Age | SD | FTA | Age | SD | FTA | Age | SD |
|------|-----|-----|------|-----|-----|------|-----|-----|-------|-----|------|
| <560 | n/a | — | 2600 | 159 | 330 | 4800 | 205 | 560 | 7000 | 246 | 800 |
| 560 | 98 | 120 | 2700 | 162 | 340 | 4900 | 207 | 570 | 7100 | 248 | 820 |
| 600 | 100 | 120 | 2800 | 164 | 350 | 5000 | 209 | 580 | 7200 | 250 | 830 |
| 700 | 103 | 130 | 2900 | 166 | 360 | 5100 | 211 | 590 | 7300 | 252 | 840 |
| 800 | 108 | 150 | 3000 | 168 | 370 | 5200 | 213 | 600 | 7400 | 254 | 860 |
| 900 | 113 | 160 | 3100 | 170 | 380 | 5300 | 215 | 610 | 7500 | 256 | 870 |
| 1000 | 115 | 170 | 3200 | 173 | 390 | 5400 | 216 | 620 | 7600 | 258 | 880 |
| 1100 | 117 | 170 | 3300 | 175 | 400 | 5500 | 218 | 630 | 7700 | 260 | 900 |
| 1200 | 122 | 190 | 3400 | 177 | 410 | 5600 | 220 | 640 | 7800 | 262 | 910 |
| 1300 | 125 | 200 | 3500 | 179 | 420 | 5700 | 222 | 650 | 7900 | 264 | 930 |
| 1400 | 128 | 210 | 3600 | 181 | 430 | 5800 | 224 | 670 | 8000 | 265 | 930 |
| 1500 | 130 | 220 | 3700 | 183 | 440 | 5900 | 226 | 680 | 8100 | 268 | 960 |
| 1600 | 134 | 230 | 3800 | 185 | 450 | 6000 | 227 | 680 | 8200 | 270 | 970 |
| 1700 | 137 | 240 | 3900 | 187 | 460 | 6100 | 229 | 700 | 8300 | 273 | 990 |
| 1800 | 139 | 250 | 4000 | 189 | 470 | 6200 | 231 | 710 | 8400 | 274 | 1000 |
| 1900 | 142 | 260 | 4100 | 191 | 480 | 6300 | 233 | 720 | 8500 | 276 | 1010 |
| 2000 | 145 | 270 | 4200 | 193 | 490 | 6400 | 235 | 730 | 8600 | 279 | 1040 |
| 2100 | 147 | 280 | 4300 | 195 | 500 | 6500 | 237 | 750 | 8660 | 280 | 1040 |
| 2200 | 150 | 290 | 4400 | 197 | 510 | 6600 | 238 | 750 | >8660 | n/a | — |
| 2300 | 152 | 300 | 4500 | 199 | 520 | 6700 | 240 | 760 | | | |
| 2400 | 155 | 310 | 4600 | 201 | 530 | 6800 | 242 | 780 | | | |
| 2500 | 157 | 330 | 4700 | 203 | 540 | 6900 | 244 | 790 | | | |

Table 2-76: HL: Osaka (Fetal Age)
 Osaka University Method 1989, 3 by Univ. Osaka
 Unit: HL (mm); Age (Days); SD (mm)

| HL | Age | SD | HL | Age | SD | HL | Age | SD | HL | Age | SD |
|-----|-----|-----|----|-----|-----|----|-----|-----|-----|-----|-----|
| <10 | n/a | — | 23 | 123 | 2.2 | 37 | 164 | 2.4 | 51 | 217 | 2.6 |
| 10 | 91 | 2.0 | 24 | 126 | 2.2 | 38 | 167 | 2.4 | 52 | 222 | 2.6 |
| 11 | 93 | 2.0 | 25 | 129 | 2.2 | 39 | 170 | 2.4 | 53 | 227 | 2.7 |
| 12 | 96 | 2.0 | 26 | 132 | 2.2 | 40 | 174 | 2.4 | 54 | 232 | 2.7 |
| 13 | 98 | 2.1 | 27 | 134 | 2.2 | 41 | 178 | 2.4 | 55 | 237 | 2.7 |
| 14 | 100 | 2.1 | 28 | 137 | 2.2 | 42 | 182 | 2.5 | 56 | 242 | 2.7 |
| 15 | 103 | 2.1 | 29 | 140 | 2.3 | 43 | 185 | 2.5 | 57 | 248 | 2.8 |
| 16 | 105 | 2.1 | 30 | 143 | 2.3 | 44 | 188 | 2.5 | 58 | 254 | 2.8 |
| 17 | 108 | 2.1 | 31 | 145 | 2.3 | 45 | 192 | 2.5 | 59 | 260 | 2.8 |
| 18 | 110 | 2.1 | 32 | 149 | 2.3 | 46 | 196 | 2.5 | 60 | 267 | 2.9 |
| 19 | 113 | 2.1 | 33 | 151 | 2.3 | 47 | 200 | 2.5 | 61 | 275 | 2.9 |
| 20 | 115 | 2.1 | 34 | 155 | 2.3 | 48 | 204 | 2.6 | 62 | 280 | 2.9 |
| 21 | 117 | 2.1 | 35 | 158 | 2.3 | 49 | 208 | 2.6 | >62 | n/a | — |
| 22 | 121 | 2.2 | 36 | 161 | 2.4 | 50 | 213 | 2.6 | | | |

Paris

Table 2-77: BPD: Paris (Fetal Age)
Unit: BPD (mm); Age (Days); SD (mm)

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <13 | n/a | — | 33 | 110 | 3 | 54 | 158 | 4 | 75 | 210 | 5 |
| 13 | 77 | 3 | 34 | 113 | 3 | 55 | 161 | 4 | 76 | 213 | 5 |
| 14 | 78 | 3 | 35 | 115 | 3 | 56 | 163 | 4 | 77 | 217 | 5 |
| 15 | 79 | 3 | 36 | 117 | 3 | 57 | 165 | 4 | 78 | 220 | 5 |
| 16 | 80 | 3 | 37 | 119 | 3 | 58 | 168 | 4 | 79 | 224 | 5 |
| 17 | 81 | 3 | 38 | 121 | 3 | 59 | 170 | 4 | 80 | 227 | 5 |
| 18 | 82 | 3 | 39 | 123 | 3 | 60 | 172 | 4 | 81 | 231 | 5 |
| 19 | 83 | 3 | 40 | 126 | 4 | 61 | 175 | 4 | 82 | 234 | 5 |
| 20 | 84 | 3 | 41 | 128 | 4 | 62 | 177 | 4 | 83 | 238 | 5 |
| 21 | 85 | 3 | 42 | 130 | 4 | 63 | 179 | 4 | 84 | 242 | 5 |
| 22 | 87 | 3 | 43 | 133 | 4 | 64 | 182 | 4 | 85 | 247 | 5 |
| 23 | 89 | 3 | 44 | 135 | 4 | 65 | 184 | 4 | 86 | 252 | 5 |
| 24 | 91 | 3 | 45 | 137 | 4 | 66 | 187 | 4 | 87 | 256 | 5 |
| 25 | 93 | 3 | 46 | 140 | 4 | 67 | 189 | 4 | 88 | 261 | 5 |
| 26 | 95 | 3 | 47 | 142 | 4 | 68 | 192 | 4 | 89 | 266 | 5 |
| 27 | 97 | 3 | 48 | 144 | 4 | 69 | 194 | 4 | 90 | 287 | 5 |
| 28 | 100 | 3 | 49 | 147 | 4 | 70 | 197 | 4 | >90 | n/a | — |
| 29 | 102 | 3 | 50 | 149 | 4 | 71 | 199 | 4 | | | |
| 30 | 104 | 3 | 51 | 151 | 4 | 72 | 202 | 4 | | | |
| 31 | 106 | 3 | 52 | 154 | 4 | 73 | 204 | 4 | | | |
| 32 | 108 | 3 | 53 | 156 | 4 | 74 | 207 | 4 | | | |

Table 2-78: CRL: Paris (Fetal Age)
Unit: CRL (mm); Age (Days); SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <5 | n/a | — | 25 | 64 | 7 | 46 | 78 | 7 | 67 | 90 | 7 |
| 5 | 42 | 4 | 26 | 65 | 7 | 47 | 79 | 7 | 68 | 90 | 7 |
| 6 | 43 | 4 | 27 | 66 | 7 | 48 | 79 | 7 | 69 | 91 | 7 |
| 7 | 44 | 4 | 28 | 66 | 7 | 49 | 80 | 7 | 70 | 91 | 7 |
| 8 | 46 | 4 | 29 | 67 | 7 | 50 | 80 | 7 | 71 | 91 | 7 |
| 9 | 47 | 4 | 30 | 68 | 7 | 51 | 81 | 7 | 72 | 92 | 7 |
| 10 | 49 | 4 | 31 | 69 | 7 | 52 | 82 | 7 | 73 | 92 | 7 |
| 11 | 50 | 4 | 32 | 70 | 7 | 53 | 82 | 7 | 74 | 93 | 7 |
| 12 | 51 | 4 | 33 | 70 | 7 | 54 | 83 | 7 | 75 | 93 | 7 |
| 13 | 52 | 4 | 34 | 71 | 7 | 55 | 84 | 7 | 76 | 94 | 7 |
| 14 | 53 | 4 | 35 | 71 | 7 | 56 | 84 | 7 | 77 | 94 | 7 |
| 15 | 54 | 4 | 36 | 72 | 7 | 57 | 85 | 7 | 78 | 94 | 7 |
| 16 | 55 | 5 | 37 | 73 | 7 | 58 | 85 | 7 | 79 | 95 | 7 |
| 17 | 56 | 5 | 38 | 73 | 7 | 59 | 86 | 7 | 80 | 95 | 7 |
| 18 | 57 | 5 | 39 | 74 | 7 | 60 | 86 | 7 | 81 | 96 | 7 |
| 19 | 58 | 6 | 40 | 74 | 7 | 61 | 87 | 7 | 82 | 96 | 7 |
| 20 | 59 | 6 | 41 | 75 | 7 | 62 | 87 | 7 | 83 | 97 | 7 |
| 21 | 60 | 6 | 42 | 76 | 7 | 63 | 88 | 7 | 84 | 97 | 7 |
| 22 | 61 | 6 | 43 | 76 | 7 | 64 | 88 | 7 | 85 | 98 | 7 |
| 23 | 63 | 6 | 44 | 77 | 7 | 65 | 89 | 7 | >85 | n/a | — |
| 24 | 63 | 7 | 45 | 77 | 7 | 66 | 89 | 7 | | | |

Table 2-79: FL: Paris (Fetal Age)
Unit: FL (mm); Age (Days); SD (mm)

| FL | Age | SD | FL | Age | SD | FL | Age | SD | FL | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <15 | n/a | — | 31 | 137 | 5 | 48 | 183 | 5 | 65 | 238 | 5 |
| 15 | 98 | 4 | 32 | 139 | 5 | 49 | 186 | 5 | 66 | 241 | 5 |
| 16 | 100 | 4 | 33 | 142 | 5 | 50 | 189 | 5 | 67 | 245 | 5 |
| 17 | 102 | 4 | 34 | 145 | 5 | 51 | 192 | 5 | 68 | 248 | 5 |
| 18 | 105 | 4 | 35 | 148 | 5 | 52 | 194 | 5 | 69 | 252 | 5 |
| 19 | 107 | 4 | 36 | 150 | 5 | 53 | 197 | 5 | 70 | 255 | 5 |
| 20 | 109 | 4 | 37 | 153 | 5 | 54 | 200 | 5 | 71 | 259 | 5 |
| 21 | 112 | 4 | 38 | 156 | 5 | 55 | 203 | 5 | 72 | 262 | 5 |
| 22 | 114 | 4 | 39 | 159 | 5 | 56 | 206 | 5 | 73 | 266 | 5 |
| 23 | 116 | 4 | 40 | 161 | 5 | 57 | 210 | 5 | 74 | 271 | 5 |
| 24 | 119 | 4 | 41 | 164 | 5 | 58 | 213 | 5 | 75 | 276 | 5 |
| 25 | 121 | 4 | 42 | 167 | 5 | 59 | 217 | 5 | 76 | 281 | 5 |
| 26 | 123 | 4 | 43 | 170 | 5 | 60 | 219 | 5 | 77 | 287 | 5 |
| 27 | 126 | 5 | 44 | 172 | 5 | 61 | 221 | 5 | >77 | n/a | — |
| 28 | 128 | 5 | 45 | 175 | 5 | 62 | 224 | 5 | | | |
| 29 | 131 | 5 | 46 | 178 | 5 | 63 | 231 | 5 | | | |
| 30 | 134 | 5 | 47 | 181 | 5 | 64 | 234 | 5 | | | |

Table 2-80: Ft: Paris (Fetal Age)
Unit: Ft (mm); Age (Days); SD (mm)

| Ft | Age | SD | Ft | Age | SD | Ft | Age | SD | Ft | Age | SD |
|-----|-----|----|----|-----|----|----|-----|----|-----|-----|----|
| <13 | n/a | — | 29 | 133 | 4 | 46 | 173 | 4 | 63 | 221 | 4 |
| 13 | 91 | 2 | 30 | 135 | 4 | 47 | 175 | 4 | 64 | 224 | 4 |
| 14 | 94 | 2 | 31 | 137 | 4 | 48 | 178 | 4 | 65 | 227 | 4 |
| 15 | 97 | 2 | 32 | 140 | 4 | 49 | 180 | 4 | 66 | 231 | 5 |
| 16 | 100 | 2 | 33 | 142 | 4 | 50 | 183 | 4 | 67 | 234 | 5 |
| 17 | 103 | 3 | 34 | 144 | 4 | 51 | 185 | 4 | 68 | 238 | 5 |
| 18 | 106 | 3 | 35 | 147 | 4 | 52 | 188 | 4 | 69 | 242 | 5 |
| 19 | 109 | 3 | 36 | 149 | 4 | 53 | 190 | 4 | 70 | 246 | 5 |
| 20 | 112 | 4 | 37 | 151 | 4 | 54 | 193 | 4 | 71 | 250 | 5 |
| 21 | 114 | 4 | 38 | 154 | 4 | 55 | 196 | 4 | 72 | 254 | 5 |
| 22 | 116 | 4 | 39 | 156 | 4 | 56 | 199 | 4 | 73 | 258 | 5 |
| 23 | 119 | 4 | 40 | 158 | 4 | 57 | 202 | 4 | 74 | 262 | 5 |
| 24 | 121 | 4 | 41 | 161 | 4 | 58 | 205 | 4 | 75 | 266 | 6 |
| 25 | 123 | 4 | 42 | 163 | 4 | 59 | 208 | 4 | >75 | n/a | — |
| 26 | 126 | 4 | 43 | 165 | 4 | 60 | 211 | 4 | | | |
| 27 | 128 | 4 | 44 | 168 | 4 | 61 | 215 | 4 | | | |
| 28 | 130 | 4 | 45 | 170 | 4 | 62 | 218 | 4 | | | |

Table 2-81: TAD: Paris (Fetal Age)
Unit: TAD (mm); Age (Days); SD (mm)

| TAD | Age | SD | TAD | Age | SD | TAD | Age | SD | TAD | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <10 | n/a | — | 32 | 116 | 0 | 55 | 171 | 0 | 78 | 229 | 0 |
| 10 | 84 | 0 | 33 | 118 | 0 | 56 | 174 | 0 | 79 | 232 | 0 |
| 11 | 84 | 0 | 34 | 120 | 0 | 57 | 176 | 0 | 80 | 234 | 0 |
| 12 | 85 | 0 | 35 | 122 | 0 | 58 | 179 | 0 | 81 | 237 | 0 |
| 13 | 86 | 0 | 36 | 124 | 0 | 59 | 181 | 0 | 82 | 239 | 0 |
| 14 | 87 | 0 | 37 | 126 | 0 | 60 | 184 | 0 | 83 | 242 | 0 |
| 15 | 87 | 0 | 38 | 128 | 0 | 61 | 186 | 0 | 84 | 245 | 0 |
| 16 | 88 | 0 | 39 | 131 | 0 | 62 | 189 | 0 | 85 | 248 | 0 |
| 17 | 89 | 0 | 40 | 133 | 0 | 63 | 191 | 0 | 86 | 252 | 0 |
| 18 | 90 | 0 | 41 | 136 | 0 | 64 | 194 | 0 | 87 | 255 | 0 |
| 19 | 91 | 0 | 42 | 138 | 0 | 65 | 196 | 0 | 88 | 259 | 0 |
| 20 | 92 | 0 | 43 | 141 | 0 | 66 | 199 | 0 | 89 | 262 | 0 |
| 21 | 94 | 0 | 44 | 143 | 0 | 67 | 201 | 0 | 90 | 266 | 0 |
| 22 | 96 | 0 | 45 | 146 | 0 | 68 | 204 | 0 | 91 | 269 | 0 |
| 23 | 98 | 0 | 46 | 148 | 0 | 69 | 207 | 0 | 92 | 273 | 0 |
| 24 | 100 | 0 | 47 | 151 | 0 | 70 | 209 | 0 | 93 | 276 | 0 |
| 25 | 102 | 0 | 48 | 153 | 0 | 71 | 212 | 0 | 94 | 280 | 0 |
| 26 | 104 | 0 | 49 | 156 | 0 | 72 | 214 | 0 | 95 | 283 | 0 |
| 27 | 106 | 0 | 50 | 158 | 0 | 73 | 217 | 0 | 96 | 287 | 0 |
| 28 | 108 | 0 | 51 | 161 | 0 | 74 | 219 | 0 | >96 | n/a | — |
| 29 | 110 | 0 | 52 | 163 | 0 | 75 | 222 | 0 | | | |
| 30 | 112 | 0 | 53 | 166 | 0 | 76 | 224 | 0 | | | |
| 31 | 114 | 0 | 54 | 169 | 0 | 77 | 227 | 0 | | | |

Rempen

Table 2-82: BPD: Rempen (Fetal Age)

Der Frauenarzt 32, 4 (1991) 425-30

Known LMP (left)—Unknown LMP (right)

Unit: BPD (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD | BPD | Age | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|-----|-------|-----|
| <2 | n/a | — | 15 | 10w2d | 4 | <3 | n/a | — | 16 | 10w4d | 8 |
| 2 | 6w2d | 4 | 16 | 10w5d | 4 | 3 | 6w6d | 8 | 17 | 10w6d | 8 |
| 3 | 6w4d | 4 | 17 | 11w0d | 4 | 4 | 7w1d | 8 | 18 | 11w1d | 8 |
| 4 | 6w6d | 4 | 18 | 11w2d | 4 | 5 | 7w3d | 8 | 19 | 11w3d | 8 |
| 5 | 7w1d | 4 | 19 | 11w5d | 4 | 6 | 7w5d | 8 | 20 | 11w5d | 8 |
| 6 | 7w4d | 4 | 20 | 12w0d | 4 | 7 | 8w0d | 8 | 21 | 12w0d | 8 |
| 7 | 7w6d | 4 | 21 | 12w2d | 4 | 8 | 8w2d | 8 | 22 | 12w2d | 8 |
| 8 | 8w1d | 4 | 22 | 12w4d | 4 | 9 | 8w4d | 8 | 23 | 12w4d | 8 |
| 9 | 8w3d | 4 | 23 | 13w0d | 4 | 10 | 8w6d | 8 | 24 | 12w6d | 8 |
| 10 | 8w5d | 4 | 24 | 13w2d | 4 | 11 | 9w1d | 8 | 25 | 13w1d | 8 |
| 11 | 9w1d | 4 | >24 | n/a | — | 12 | 9w3d | 8 | 26 | 13w3d | 8 |
| 12 | 9w3d | 4 | | | | 13 | 9w5d | 8 | 27 | 13w5d | 8 |
| 13 | 9w5d | 4 | | | | 14 | 10w0d | 8 | >27 | n/a | — |
| 14 | 10w0d | 4 | | | | 15 | 10w2d | 8 | | | |

Table 2-83: BPD: Rempen (Fetal Growth)

Der Frauenarzt 32, 4 (1991) 425-30

Unit: Age (Weeks/Days); Mean (mm); 2SD (mm); Table/Graph Range (5%:95%)

| Age | Mean | SD | Age | Mean | SD | Age | Mean | SD |
|------|------|-----|-------|------|-----|-------|------|-----|
| 6w2d | 2.0 | 3.7 | 8w5d | 9.8 | 3.7 | 11w1d | 17.4 | 3.7 |
| 6w3d | 2.5 | 3.7 | 8w6d | 10.3 | 3.7 | 11w2d | 17.9 | 3.7 |
| 6w4d | 3.0 | 3.7 | 9w0d | 10.7 | 3.7 | 11w3d | 18.3 | 3.7 |
| 6w5d | 3.4 | 3.7 | 9w1d | 11.2 | 3.7 | 11w4d | 18.7 | 3.7 |
| 6w6d | 3.9 | 3.7 | 9w2d | 11.6 | 3.7 | 11w5d | 19.2 | 3.7 |
| 7w0d | 4.3 | 3.7 | 9w3d | 12.1 | 3.7 | 11w6d | 19.6 | 3.7 |
| 7w1d | 4.8 | 3.7 | 9w4d | 12.5 | 3.7 | 12w0d | 20.0 | 3.7 |
| 7w2d | 5.3 | 3.7 | 9w5d | 13.0 | 3.7 | 12w1d | 20.5 | 3.7 |
| 7w3d | 5.7 | 3.7 | 9w6d | 13.4 | 3.7 | 12w2d | 20.9 | 3.7 |
| 7w4d | 6.2 | 3.7 | 10w0d | 13.9 | 3.7 | 12w3d | 21.3 | 3.7 |
| 7w5d | 6.7 | 3.7 | 10w1d | 14.3 | 3.7 | 12w4d | 21.8 | 3.7 |
| 7w6d | 7.1 | 3.7 | 10w2d | 14.8 | 3.7 | 12w5d | 22.2 | 3.7 |
| 8w0d | 7.6 | 3.7 | 10w3d | 15.2 | 3.7 | 12w6d | 22.6 | 3.7 |
| 8w1d | 8.0 | 3.7 | 10w4d | 15.7 | 3.7 | 13w0d | 23.1 | 3.7 |
| 8w2d | 8.5 | 3.7 | 10w5d | 16.1 | 3.7 | 13w1d | 23.5 | 3.7 |
| 8w3d | 8.9 | 3.7 | 10w6d | 16.5 | 3.7 | 13w2d | 23.9 | 3.7 |
| 8w4d | 9.4 | 3.7 | 11w0d | 17.0 | 3.7 | | | |

Table 2-84: CRL: Rempen (Fetal Age)

Der Frauenarzt 32, 4 (1991) 425-30

Known LMP (left)—Unknown LMP (right)

Unit: CRL (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD | CRL | Age | 2SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|-----|-------|-----|
| <1 | n/a | — | 40 | 10w5d | 8 | <2 | n/a | — | 41 | 10w5d | 7 |
| 1 | 5w5d | 8 | 41 | 10w6d | 8 | 2 | 6w0d | 6 | 42 | 10w6d | 6 |
| 2 | 5w6d | 8 | 42 | 11w0d | 8 | 3 | 6w1d | 6 | 43 | 11w0d | 7 |
| 3 | 6w0d | 8 | 43 | 11w0d | 8 | 4 | 6w2d | 6 | 44 | 11w0d | 7 |
| 4 | 6w1d | 8 | 44 | 11w1d | 8 | 5 | 6w3d | 6 | 45 | 11w1d | 6 |
| 5 | 6w2d | 8 | 45 | 11w2d | 8 | 6 | 6w4d | 6 | 46 | 11w2d | 7 |
| 6 | 6w3d | 8 | 46 | 11w2d | 8 | 7 | 6w5d | 6 | 47 | 11w2d | 7 |
| 7 | 6w4d | 8 | 47 | 11w3d | 8 | 8 | 6w6d | 6 | 48 | 11w3d | 6 |
| 8 | 6w6d | 8 | 48 | 11w4d | 8 | 9 | 7w0d | 6 | 49 | 11w4d | 7 |
| 9 | 6w6d | 8 | 49 | 11w4d | 8 | 10 | 7w1d | 6 | 50 | 11w4d | 6 |
| 10 | 7w0d | 8 | 50 | 11w5d | 8 | 11 | 7w2d | 6 | 51 | 11w5d | 6 |
| 11 | 7w2d | 8 | 51 | 11w6d | 8 | 12 | 7w3d | 6 | 52 | 11w5d | 7 |
| 12 | 7w2d | 8 | 52 | 12w0d | 8 | 13 | 7w4d | 7 | 53 | 11w6d | 6 |
| 13 | 7w4d | 8 | 53 | 12w0d | 8 | 14 | 7w5d | 7 | 54 | 12w0d | 7 |
| 14 | 7w4d | 8 | 54 | 12w1d | 8 | 15 | 7w6d | 7 | 55 | 12w0d | 7 |
| 15 | 7w5d | 8 | 55 | 12w2d | 8 | 16 | 7w6d | 7 | 56 | 12w1d | 6 |
| 16 | 7w6d | 8 | 56 | 12w2d | 8 | 17 | 8w0d | 7 | 57 | 12w1d | 7 |
| 17 | 8w0d | 8 | 57 | 12w3d | 8 | 18 | 8w1d | 6 | 58 | 12w2d | 6 |
| 18 | 8w1d | 8 | 58 | 12w3d | 8 | 19 | 8w2d | 6 | 59 | 12w3d | 7 |
| 19 | 8w2d | 8 | 59 | 12w4d | 8 | 20 | 8w3d | 6 | 60 | 12w3d | 6 |
| 20 | 8w3d | 8 | 60 | 12w5d | 8 | 21 | 8w4d | 7 | 61 | 12w4d | 7 |
| 21 | 8w4d | 8 | 61 | 12w5d | 8 | 22 | 8w5d | 7 | 62 | 12w4d | 6 |
| 22 | 8w5d | 8 | 62 | 12w6d | 8 | 23 | 8w5d | 7 | 63 | 12w5d | 7 |
| 23 | 8w6d | 8 | 63 | 13w0d | 8 | 24 | 8w6d | 7 | 64 | 12w5d | 7 |
| 24 | 8w6d | 8 | 64 | 13w0d | 8 | 25 | 9w0d | 6 | 65 | 12w6d | 6 |
| 25 | 9w0d | 8 | 65 | 13w1d | 8 | 26 | 9w1d | 6 | 66 | 12w6d | 7 |
| 26 | 9w1d | 8 | 66 | 13w2d | 8 | 27 | 9w2d | 7 | 67 | 13w0d | 6 |
| 27 | 9w2d | 8 | >66 | n/a | — | 28 | 9w3d | 7 | 68 | 13w0d | 7 |
| 28 | 9w3d | 8 | | | | 29 | 9w3d | 7 | 69 | 13w1d | 6 |
| 29 | 9w4d | 8 | | | | 30 | 9w4d | 7 | 70 | 13w1d | 7 |
| 30 | 9w4d | 8 | | | | 31 | 9w5d | 7 | 71 | 13w2d | 7 |
| 31 | 9w5d | 8 | | | | 32 | 9w6d | 7 | 72 | 13w2d | 6 |
| 32 | 9w6d | 8 | | | | 33 | 9w6d | 7 | 73 | 13w3d | 7 |
| 33 | 10w0d | 8 | | | | 34 | 10w0d | 6 | 74 | 13w3d | 6 |
| 34 | 10w1d | 8 | | | | 35 | 10w1d | 6 | 75 | 13w4d | 7 |
| 35 | 10w1d | 8 | | | | 36 | 10w2d | 7 | 76 | 13w4d | 6 |
| 36 | 10w2d | 8 | | | | 37 | 10w2d | 7 | 77 | 13w4d | 7 |
| 37 | 10w3d | 8 | | | | 38 | 10w3d | 6 | 78 | 13w5d | 6 |
| 38 | 10w4d | 8 | | | | 39 | 10w4d | 6 | >78 | n/a | — |
| 39 | 10w4d | 8 | | | | 40 | 10w5d | 7 | | | |

Table 2-85: CRL: Rempen (Fetal Growth)
 Der Frauenarzt 32, 4 (1991) 425-30
 Unit: Age (Weeks/Days); Mean (mm); 2SD (mm); Table/Graph Range (5%:95%)

| Age | Mean | SD | Age | Mean | SD | Age | Mean | SD |
|------|------|-----|-------|------|-----|-------|------|-----|
| 5w5d | 1.2 | 7.8 | 8w2d | 18.9 | 7.8 | 10w6d | 41.3 | 7.8 |
| 5w6d | 2.1 | 7.8 | 8w3d | 20.0 | 7.8 | 11w0d | 42.6 | 7.8 |
| 6w0d | 3.0 | 7.8 | 8w4d | 21.1 | 7.8 | 11w1d | 44.0 | 7.8 |
| 6w1d | 3.8 | 7.8 | 8w5d | 22.3 | 7.8 | 11w2d | 45.4 | 7.8 |
| 6w2d | 4.7 | 7.8 | 8w6d | 23.5 | 7.8 | 11w3d | 46.9 | 7.8 |
| 6w3d | 5.7 | 7.8 | 9w0d | 24.6 | 7.8 | 11w4d | 48.3 | 7.8 |
| 6w4d | 6.6 | 7.8 | 9w1d | 25.8 | 7.8 | 11w5d | 49.8 | 7.8 |
| 6w5d | 7.5 | 7.8 | 9w2d | 27.0 | 7.8 | 11w6d | 51.2 | 7.8 |
| 6w6d | 8.5 | 7.8 | 9w3d | 28.3 | 7.8 | 12w0d | 52.7 | 7.8 |
| 7w0d | 9.5 | 7.8 | 9w4d | 29.5 | 7.8 | 12w1d | 54.2 | 7.8 |
| 7w1d | 10.5 | 7.8 | 9w5d | 30.7 | 7.8 | 12w2d | 55.7 | 7.8 |
| 7w2d | 11.5 | 7.8 | 9w6d | 32.0 | 7.8 | 12w3d | 57.3 | 7.8 |
| 7w3d | 12.5 | 7.8 | 10w0d | 33.3 | 7.8 | 12w4d | 58.8 | 7.8 |
| 7w4d | 13.5 | 7.8 | 10w1d | 34.6 | 7.8 | 12w5d | 60.3 | 7.8 |
| 7w5d | 14.6 | 7.8 | 10w2d | 35.9 | 7.8 | 12w6d | 61.9 | 7.8 |
| 7w6d | 15.6 | 7.8 | 10w3d | 37.2 | 7.8 | 13w0d | 63.5 | 7.8 |
| 8w0d | 16.7 | 7.8 | 10w4d | 38.5 | 7.8 | 13w1d | 65.1 | 7.8 |
| 8w1d | 17.8 | 7.8 | 10w5d | 39.9 | 7.8 | 13w2d | 66.7 | 7.8 |

Table 2-86: GS: Rempen (Fetal Age)

Der Frauenarzt 32, 4 (1991) 425-30

Known LMP (left)—Unknown LMP (right)

Unit: GS (mm); Age (Weeks/Days); 2SD (mm [Known LMP] or day [Unknown LMP])

| GS | Age | 2SD | GS | Age | 2SD | GS | Age | 2SD | GS | Age | 2SD |
|----|------|-----|-----|-------|-----|----|------|-----|-----|-------|-----|
| <1 | n/a | — | 38 | 9w1d | 11 | <1 | n/a | — | 38 | 9w1d | 10 |
| 1 | 4w4d | 11 | 39 | 9w2d | 11 | 1 | 4w5d | 10 | 39 | 9w2d | 10 |
| 2 | 4w5d | 11 | 40 | 9w4d | 11 | 2 | 4w6d | 10 | 40 | 9w3d | 10 |
| 3 | 4w6d | 11 | 41 | 9w5d | 11 | 3 | 5w0d | 10 | 41 | 9w4d | 10 |
| 4 | 5w0d | 11 | 42 | 9w6d | 11 | 4 | 5w1d | 10 | 42 | 9w5d | 10 |
| 5 | 5w0d | 11 | 43 | 10w0d | 11 | 5 | 5w2d | 10 | 43 | 9w6d | 10 |
| 6 | 5w1d | 11 | 44 | 10w1d | 11 | 6 | 5w2d | 10 | 44 | 9w6d | 10 |
| 7 | 5w2d | 11 | 45 | 10w2d | 11 | 7 | 5w3d | 10 | 45 | 10w0d | 10 |
| 8 | 5w3d | 11 | 46 | 10w3d | 11 | 8 | 5w4d | 10 | 46 | 10w1d | 10 |
| 9 | 5w3d | 11 | 47 | 10w4d | 11 | 9 | 5w5d | 10 | 47 | 10w2d | 10 |
| 10 | 5w4d | 11 | 48 | 10w6d | 11 | 10 | 5w5d | 10 | 48 | 10w3d | 10 |
| 11 | 5w5d | 11 | 49 | 11w0d | 11 | 11 | 5w6d | 10 | 49 | 10w4d | 10 |
| 12 | 5w6d | 11 | 50 | 11w1d | 11 | 12 | 6w0d | 10 | 50 | 10w5d | 10 |
| 13 | 6w0d | 11 | 51 | 11w2d | 11 | 13 | 6w1d | 10 | 51 | 10w6d | 10 |
| 14 | 6w0d | 11 | 52 | 11w4d | 11 | 14 | 6w2d | 10 | 52 | 11w0d | 10 |
| 15 | 6w1d | 11 | 53 | 11w5d | 11 | 15 | 6w2d | 10 | 53 | 11w1d | 10 |
| 16 | 6w2d | 11 | 54 | 12w0d | 11 | 16 | 6w3d | 10 | 54 | 11w2d | 10 |
| 17 | 6w3d | 11 | 55 | 12w1d | 11 | 17 | 6w4d | 10 | 55 | 11w3d | 10 |
| 18 | 6w4d | 11 | 56 | 12w2d | 11 | 18 | 6w5d | 10 | 56 | 11w4d | 10 |
| 19 | 6w5d | 11 | 57 | 12w4d | 11 | 19 | 6w6d | 10 | 57 | 11w5d | 10 |
| 20 | 6w6d | 11 | 58 | 12w5d | 11 | 20 | 6w6d | 10 | 58 | 11w6d | 10 |
| 21 | 6w6d | 11 | 59 | 13w0d | 11 | 21 | 7w0d | 10 | 59 | 12w0d | 10 |
| 22 | 7w0d | 11 | 60 | 13w1d | 11 | 22 | 7w1d | 10 | 60 | 12w1d | 10 |
| 23 | 7w1d | 11 | >60 | n/a | — | 23 | 7w2d | 10 | 61 | 12w2d | 10 |
| 24 | 7w2d | 11 | | | | 24 | 7w3d | 10 | 62 | 12w3d | 10 |
| 25 | 7w3d | 11 | | | | 25 | 7w4d | 10 | 63 | 12w4d | 10 |
| 26 | 7w4d | 11 | | | | 26 | 7w4d | 10 | 64 | 12w5d | 10 |
| 27 | 7w5d | 11 | | | | 27 | 7w5d | 10 | 65 | 12w6d | 10 |
| 28 | 7w6d | 11 | | | | 28 | 7w6d | 10 | 66 | 13w0d | 10 |
| 29 | 8w0d | 11 | | | | 29 | 8w0d | 10 | 67 | 13w1d | 10 |
| 30 | 8w0d | 11 | | | | 30 | 8w1d | 10 | 68 | 13w2d | 10 |
| 31 | 8w1d | 11 | | | | 31 | 8w2d | 10 | 69 | 13w3d | 10 |
| 32 | 8w2d | 11 | | | | 32 | 8w3d | 10 | 70 | 13w4d | 10 |
| 33 | 8w3d | 11 | | | | 33 | 8w3d | 10 | 71 | 13w5d | 10 |
| 34 | 8w4d | 11 | | | | 34 | 8w4d | 10 | 72 | 14w0d | 10 |
| 35 | 8w5d | 11 | | | | 35 | 8w5d | 10 | 73 | 14w1d | 10 |
| 36 | 8w6d | 11 | | | | 36 | 8w6d | 10 | >73 | n/a | — |
| 37 | 9w0d | 11 | | | | 37 | 9w0d | 10 | | | |

Table 2-87: GS: Rempen (Fetal Growth)

Der Frauenarzt 32, 4 (1991) 425-30

Unit: Age (Weeks/Days); Mean (mm); 2SD (mm); Table/Graph Range (5%:95%)

| Age | Mean | SD | Age | Mean | SD | Age | Mean | SD |
|------|------|------|-------|------|------|-------|------|------|
| 4w4d | 0.5 | 10.5 | 7w4d | 26.2 | 10.5 | 10w4d | 46.6 | 10.5 |
| 4w5d | 1.8 | 10.5 | 7w5d | 27.3 | 10.5 | 10w5d | 47.4 | 10.5 |
| 4w6d | 3.2 | 10.5 | 7w6d | 28.4 | 10.5 | 10w6d | 48.2 | 10.5 |
| 5w0d | 4.5 | 10.5 | 8w0d | 29.5 | 10.5 | 11w0d | 49.0 | 10.5 |
| 5w1d | 5.8 | 10.5 | 8w1d | 30.5 | 10.5 | 11w1d | 49.8 | 10.5 |
| 5w2d | 7.1 | 10.5 | 8w2d | 31.6 | 10.5 | 11w2d | 50.6 | 10.5 |
| 5w3d | 8.4 | 10.5 | 8w3d | 32.6 | 10.5 | 11w3d | 51.4 | 10.5 |
| 5w4d | 9.7 | 10.5 | 8w4d | 33.6 | 10.5 | 11w4d | 52.1 | 10.5 |
| 5w5d | 10.9 | 10.5 | 8w5d | 34.6 | 10.5 | 11w5d | 52.9 | 10.5 |
| 5w6d | 12.2 | 10.5 | 8w6d | 35.6 | 10.5 | 11w6d | 53.6 | 10.5 |
| 6w0d | 13.4 | 10.5 | 9w0d | 36.6 | 10.5 | 12w0d | 54.3 | 10.5 |
| 6w1d | 14.6 | 10.5 | 9w1d | 37.6 | 10.5 | 12w1d | 55.1 | 10.5 |
| 6w2d | 15.9 | 10.5 | 9w2d | 38.5 | 10.5 | 12w2d | 55.8 | 10.5 |
| 6w3d | 17.1 | 10.5 | 9w3d | 39.5 | 10.5 | 12w3d | 56.4 | 10.5 |
| 6w4d | 18.3 | 10.5 | 9w4d | 40.4 | 10.5 | 12w4d | 57.1 | 10.5 |
| 6w5d | 19.4 | 10.5 | 9w5d | 41.3 | 10.5 | 12w5d | 57.8 | 10.5 |
| 6w6d | 20.6 | 10.5 | 9w6d | 42.2 | 10.5 | 12w6d | 58.4 | 10.5 |
| 7w0d | 21.7 | 10.5 | 10w0d | 43.1 | 10.5 | 13w0d | 59.1 | 10.5 |
| 7w1d | 22.9 | 10.5 | 10w1d | 44.0 | 10.5 | 13w1d | 59.7 | 10.5 |
| 7w2d | 24.0 | 10.5 | 10w2d | 44.9 | 10.5 | 13w2d | 60.3 | 10.5 |
| 7w3d | 25.1 | 10.5 | 10w3d | 45.7 | 10.5 | | | |

Robinson

Table 2-88: CRL: Robinson (Fetal Age)

Br J Gynecol, 82: 702, 1975

Unit: CRL (mm); Age (Days); SD (mm)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|----|-----|-----|----|-----|-----|----|-----|-----|----|
| <7 | n/a | — | 26 | 64 | 5 | 46 | 78 | 7 | 66 | 90 | 7 |
| 7 | 45 | 4 | 27 | 65 | 5 | 47 | 79 | 7 | 67 | 90 | 7 |
| 8 | 46 | 4 | 28 | 66 | 6 | 48 | 79 | 7 | 68 | 91 | 7 |
| 9 | 47 | 4 | 29 | 67 | 6 | 49 | 80 | 7 | 69 | 91 | 7 |
| 10 | 48 | 4 | 30 | 68 | 6 | 50 | 81 | 7 | 70 | 91 | 7 |
| 11 | 50 | 4 | 31 | 69 | 7 | 51 | 82 | 7 | 71 | 92 | 7 |
| 12 | 52 | 4 | 32 | 69 | 7 | 52 | 83 | 7 | 72 | 92 | 7 |
| 13 | 53 | 4 | 33 | 70 | 7 | 53 | 83 | 7 | 73 | 93 | 7 |
| 14 | 54 | 4 | 34 | 70 | 7 | 54 | 83 | 7 | 74 | 93 | 7 |
| 15 | 55 | 4 | 35 | 71 | 7 | 55 | 84 | 7 | 75 | 93 | 7 |
| 16 | 56 | 4 | 36 | 72 | 7 | 56 | 84 | 7 | 76 | 94 | 7 |
| 17 | 57 | 4 | 37 | 72 | 7 | 57 | 84 | 7 | 77 | 94 | 7 |
| 18 | 58 | 4 | 38 | 73 | 7 | 58 | 85 | 7 | 78 | 95 | 7 |
| 19 | 59 | 4 | 39 | 74 | 7 | 59 | 85 | 7 | 79 | 95 | 7 |
| 20 | 60 | 4 | 40 | 74 | 7 | 60 | 86 | 7 | 80 | 96 | 7 |
| 21 | 60 | 4 | 41 | 75 | 7 | 61 | 86 | 7 | 81 | 97 | 7 |
| 22 | 61 | 4 | 42 | 75 | 7 | 62 | 87 | 7 | 82 | 98 | 7 |
| 23 | 62 | 4 | 43 | 76 | 7 | 63 | 88 | 7 | >82 | n/a | — |
| 24 | 63 | 5 | 44 | 77 | 7 | 64 | 89 | 7 | | | |
| 25 | 64 | 5 | 45 | 77 | 7 | 65 | 90 | 7 | | | |

Tokyo

Table 2-89: APTDxTTD: Tokyo (Fetal Age)

Tokyo University Method 1986, 6 by University Tokyo

Unit: Meas (cm²); Age (Weeks/Days); SD (Days)

| Meas | Age | SD | Meas | Age | SD | Meas | Age | SD |
|------|-------|-------|------|-------|-------|------|-------|-------|
| <10 | n/a | — | 38 | 25w6d | ± 10d | 68 | 33w3d | ± 15d |
| 10 | 16w1d | ± 8d | 40 | 26w3d | ± 11d | 70 | 33w6d | ± 16d |
| 12 | 17w0d | ± 8d | 42 | 27w0d | ± 11d | 72 | 34w2d | ± 16d |
| 14 | 17w6d | ± 8d | 44 | 27w3d | ± 11d | 74 | 34w6d | ± 17d |
| 16 | 18w4d | ± 8d | 46 | 28w0d | ± 12d | 76 | 35w3d | ± 17d |
| 18 | 19w3d | ± 8d | 48 | 28w4d | ± 12d | 78 | 35w6d | ± 17d |
| 20 | 20w1d | ± 8d | 50 | 29w0d | ± 12d | 80 | 36w3d | ± 18d |
| 22 | 20w6d | ± 9d | 52 | 29w3d | ± 13d | 82 | 37w0d | ± 18d |
| 24 | 21w4d | ± 9d | 54 | 30w0d | ± 13d | 84 | 37w4d | ± 18d |
| 26 | 22w2d | ± 9d | 56 | 30w3d | ± 13d | 86 | 38w1d | ± 18d |
| 28 | 22w6d | ± 9d | 58 | 31w0d | ± 14d | 88 | 38w5d | ± 19d |
| 30 | 23w4d | ± 9d | 60 | 31w3d | ± 14d | 90 | 39w2d | ± 19d |
| 32 | 24w1d | ± 10d | 62 | 31w6d | ± 14d | >90 | n/a | — |
| 34 | 24w5d | ± 10d | 64 | 32w3d | ± 15d | | | |
| 36 | 25w2d | ± 10d | 66 | 32w6d | ± 15d | | | |

Table 2-90: APTDxTTD by Gestational Age: Tokyo
Tokyo University Method 1986, 6 by University Tokyo

| Weeks | -1.64 SD | -1.5 SD | -1.28 SD | Mean | +1.28 SD | +1.5 SD | +1.64 SD |
|-------|----------|---------|----------|-------|----------|---------|----------|
| 16.0 | 7.0 | 7.4 | 7.9 | 11.2 | 14.6 | 15.1 | 15.5 |
| 17.0 | 8.7 | 9.0 | 9.7 | 13.3 | 17.0 | 17.6 | 18.0 |
| 18.0 | 10.5 | 10.9 | 11.6 | 15.6 | 19.6 | 20.3 | 20.7 |
| 19.0 | 12.5 | 13.0 | 13.7 | 18.1 | 22.4 | 23.2 | 23.6 |
| 20.0 | 14.7 | 15.2 | 16.1 | 20.8 | 25.5 | 26.3 | 26.8 |
| 21.0 | 17.1 | 17.6 | 18.5 | 23.6 | 28.8 | 29.6 | 30.2 |
| 22.0 | 19.6 | 20.2 | 21.2 | 26.7 | 32.2 | 33.2 | 33.8 |
| 23.0 | 22.2 | 22.9 | 23.9 | 29.9 | 35.9 | 36.9 | 37.5 |
| 24.0 | 25.0 | 25.7 | 26.8 | 33.2 | 39.7 | 40.8 | 41.5 |
| 25.0 | 27.9 | 28.6 | 29.8 | 36.7 | 43.6 | 44.8 | 45.6 |
| 26.0 | 30.9 | 31.7 | 33.0 | 40.3 | 47.7 | 49.0 | 49.8 |
| 27.0 | 33.9 | 34.8 | 36.2 | 44.1 | 52.0 | 53.3 | 54.2 |
| 28.0 | 37.1 | 38.0 | 39.4 | 47.9 | 56.3 | 57.8 | 58.7 |
| 29.0 | 40.3 | 41.3 | 42.8 | 51.8 | 60.8 | 62.3 | 63.3 |
| 30.0 | 43.5 | 44.5 | 46.2 | 55.7 | 65.3 | 66.9 | 68.0 |
| 31.0 | 46.8 | 47.9 | 49.6 | 59.7 | 69.9 | 71.6 | 72.7 |
| 32.0 | 50.0 | 51.2 | 53.0 | 63.8 | 74.5 | 76.4 | 77.6 |
| 33.0 | 53.3 | 54.5 | 56.5 | 67.8 | 79.2 | 81.2 | 82.4 |
| 34.0 | 56.5 | 57.8 | 59.9 | 71.9 | 83.9 | 86.0 | 87.3 |
| 35.0 | 59.7 | 61.1 | 63.3 | 75.9 | 88.6 | 90.8 | 92.2 |
| 36.0 | 62.8 | 64.3 | 66.6 | 79.9 | 93.3 | 95.6 | 97.0 |
| 37.0 | 65.9 | 67.4 | 69.8 | 83.9 | 97.9 | 100.3 | 101.9 |
| 38.0 | 68.8 | 70.4 | 72.9 | 87.7 | 102.5 | 105.0 | 106.7 |
| 39.0 | 71.6 | 73.3 | 76.0 | 91.5 | 107.0 | 109.7 | 111.4 |
| 40.0 | 74.3 | 76.1 | 78.9 | 95.1 | 111.4 | 114.2 | 116.0 |
| 41.0 | 76.8 | 78.6 | 81.6 | 98.6 | 115.7 | 118.6 | 120.5 |
| 42.0 | 79.1 | 81.1 | 84.1 | 102.0 | 119.8 | 122.9 | 124.8 |

Table 2-91: BPD: Tokyo (Fetal Age)
Tokyo University Method 1986, 6 by University Tokyo
Unit: BPD (mm); Age (Days); SD (Days)

| BPD | Age | SD | BPD | Age | SD | BPD | Age | SD | BPD | Age | SD |
|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <20 | n/a | — | 38 | 123 | ± 5 | 57 | 164 | ± 6 | 76 | 213 | ± 8 |
| 20 | 85 | ± 6 | 39 | 125 | ± 5 | 58 | 167 | ± 6 | 77 | 216 | ± 8 |
| 21 | 87 | ± 6 | 40 | 127 | ± 5 | 59 | 169 | ± 6 | 78 | 218 | ± 8 |
| 22 | 89 | ± 6 | 41 | 129 | ± 5 | 60 | 171 | ± 6 | 79 | 221 | ± 8 |
| 23 | 92 | ± 6 | 42 | 131 | ± 5 | 61 | 174 | ± 7 | 80 | 225 | ± 8 |
| 24 | 94 | ± 6 | 43 | 133 | ± 5 | 62 | 176 | ± 7 | 81 | 228 | ± 8 |
| 25 | 96 | ± 6 | 44 | 135 | ± 5 | 63 | 179 | ± 7 | 82 | 231 | ± 8 |
| 26 | 98 | ± 6 | 45 | 138 | ± 6 | 64 | 181 | ± 7 | 83 | 234 | ± 9 |
| 27 | 100 | ± 6 | 46 | 140 | ± 6 | 65 | 183 | ± 7 | 84 | 238 | ± 9 |
| 28 | 102 | ± 6 | 47 | 142 | ± 6 | 66 | 186 | ± 7 | 85 | 241 | ± 9 |
| 29 | 102 | ± 6 | 48 | 144 | ± 6 | 67 | 188 | ± 7 | 86 | 245 | ± 9 |
| 30 | 106 | ± 5 | 49 | 146 | ± 6 | 68 | 191 | ± 7 | 87 | 249 | ± 9 |
| 31 | 108 | ± 5 | 50 | 148 | ± 6 | 69 | 194 | ± 7 | 88 | 253 | ± 9 |
| 32 | 110 | ± 5 | 51 | 151 | ± 6 | 70 | 196 | ± 7 | 89 | 258 | ± 9 |
| 33 | 112 | ± 5 | 52 | 153 | ± 6 | 71 | 199 | ± 8 | 90 | 262 | ± 9 |
| 34 | 114 | ± 5 | 53 | 154 | ± 6 | 72 | 201 | ± 8 | >90 | n/a | — |
| 35 | 116 | ± 5 | 54 | 157 | ± 6 | 73 | 204 | ± 8 | | | |
| 36 | 118 | ± 5 | 55 | 160 | ± 6 | 74 | 207 | ± 8 | | | |
| 37 | 120 | ± 5 | 56 | 162 | ± 6 | 75 | 210 | ± 8 | | | |

Table 2-92: CRL: Tokyo (Fetal Age)
Tokyo University Method 1986, 6 by University Tokyo
Unit: CRL (mm); Age (Days); SD (Days)

| CRL | Age | SD | CRL | Age | SD | CRL | Age | SD | CRL | Age | SD |
|-----|-----|------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| <13 | n/a | — | 22 | 64 | ± 7 | 32 | 73 | ± 7 | 42 | 81 | ± 7 |
| 13 | 55 | ± 8 | 23 | 65 | ± 7 | 33 | 74 | ± 7 | 43 | 81 | ± 7 |
| 14 | 56 | ± 9 | 24 | 66 | ± 7 | 34 | 74 | ± 7 | 44 | 82 | ± 7 |
| 15 | 57 | ± 10 | 25 | 67 | ± 7 | 35 | 75 | ± 7 | 45 | 83 | ± 7 |
| 16 | 58 | ± 8 | 26 | 68 | ± 7 | 36 | 76 | ± 7 | 46 | 84 | ± 7 |
| 17 | 59 | ± 9 | 27 | 68 | ± 7 | 37 | 77 | ± 7 | 47 | 84 | ± 7 |
| 18 | 60 | ± 10 | 28 | 69 | ± 7 | 38 | 78 | ± 7 | 48 | 85 | ± 7 |
| 19 | 61 | ± 8 | 29 | 70 | ± 7 | 39 | 78 | ± 7 | 49 | 86 | ± 7 |
| 20 | 62 | ± 9 | 30 | 71 | ± 7 | 40 | 79 | ± 7 | 50 | 86 | ± 7 |
| 21 | 63 | ± 7 | 31 | 72 | ± 7 | 41 | 80 | ± 7 | >50 | n/a | — |

Table 2-93: FL: Tokyo (Fetal Age)
Tokyo University Method 1986, 6 by University Tokyo
Unit: FL (mm); Age (Days); SD (mm)

| FL | Age | SD | FL | Age | SD | FL | Age | SD | FL | Age | SD |
|-----|-----|-----|----|-----|-----|----|-----|-----|-----|-----|-----|
| <33 | n/a | — | 43 | 175 | ± 6 | 54 | 210 | ± 7 | 65 | 251 | ± 8 |
| 33 | 143 | ± 6 | 44 | 178 | ± 6 | 55 | 214 | ± 7 | 66 | 256 | ± 8 |
| 34 | 146 | ± 6 | 45 | 181 | ± 6 | 56 | 217 | ± 7 | 67 | 260 | ± 8 |
| 35 | 149 | ± 6 | 46 | 185 | ± 7 | 57 | 220 | ± 7 | 68 | 266 | ± 7 |
| 36 | 153 | ± 6 | 47 | 188 | ± 7 | 58 | 224 | ± 7 | 69 | 271 | ± 7 |
| 37 | 156 | ± 6 | 48 | 191 | ± 7 | 59 | 228 | ± 8 | 70 | 278 | ± 7 |
| 38 | 159 | ± 6 | 49 | 194 | ± 7 | 60 | 231 | ± 8 | 71 | 286 | ± 6 |
| 39 | 162 | ± 6 | 50 | 197 | ± 7 | 61 | 235 | ± 8 | >71 | n/a | — |
| 40 | 166 | ± 6 | 51 | 200 | ± 7 | 62 | 239 | ± 8 | | | |
| 41 | 169 | ± 6 | 52 | 204 | ± 7 | 63 | 243 | ± 8 | | | |
| 42 | 172 | ± 6 | 53 | 207 | ± 7 | 64 | 247 | ± 8 | | | |

Table 2-94: GS: Tokyo (Fetal Age)
Tokyo University Method 1986, 6 by University Tokyo
Unit: GS (mm); Age (Days); SD (Days)

| GS | Age | SD | GS | Age | SD | GS | Age | SD | GS | Age | SD |
|-----|-----|-----|----|-----|------|----|-----|-----|-----|-----|-----|
| <12 | n/a | — | 22 | 43 | ± 7 | 33 | 56 | ± 0 | 44 | 66 | ± 0 |
| 12 | 31 | ± 7 | 23 | 44 | ± 7 | 34 | 57 | ± 0 | 45 | 67 | ± 0 |
| 13 | 32 | ± 7 | 24 | 46 | ± 7 | 35 | 58 | ± 0 | 46 | 68 | ± 0 |
| 14 | 33 | ± 7 | 25 | 47 | ± 7 | 36 | 59 | ± 0 | 47 | 69 | ± 0 |
| 15 | 34 | ± 7 | 26 | 48 | ± 8 | 37 | 60 | ± 0 | 48 | 70 | ± 0 |
| 16 | 36 | ± 7 | 27 | 49 | ± 9 | 38 | 61 | ± 0 | 49 | 71 | ± 0 |
| 17 | 37 | ± 7 | 28 | 50 | ± 10 | 39 | 62 | ± 0 | 50 | 72 | ± 0 |
| 18 | 38 | ± 7 | 29 | 51 | ± 0 | 40 | 63 | ± 0 | >50 | n/a | — |
| 19 | 40 | ± 7 | 30 | 52 | ± 0 | 41 | 64 | ± 0 | | | |
| 20 | 41 | ± 7 | 31 | 53 | ± 0 | 42 | 65 | ± 0 | | | |
| 21 | 42 | ± 7 | 32 | 55 | ± 0 | 43 | 65 | ± 0 | | | |

Table 2-95: LV: Tokyo (Fetal Age)
Tokyo University Method 1986, 6 by University Tokyo
Unit: LV (mm); Age (Days); SD (Days)

| LV | Age | SD | LV | Age | SD | LV | Age | SD | LV | Age | SD |
|-----|-----|-----|----|-----|------|----|-----|------|-----|-----|------|
| <44 | n/a | — | 55 | 181 | ± 7 | 67 | 217 | ± 10 | 79 | 260 | ± 10 |
| 44 | 154 | ± 5 | 56 | 183 | ± 8 | 68 | 220 | ± 10 | 80 | 264 | ± 10 |
| 45 | 157 | ± 5 | 57 | 186 | ± 8 | 69 | 224 | ± 10 | 81 | 267 | ± 10 |
| 46 | 159 | ± 5 | 58 | 189 | ± 8 | 70 | 227 | ± 11 | 82 | 271 | ± 10 |
| 47 | 161 | ± 5 | 59 | 192 | ± 8 | 71 | 231 | ± 11 | 83 | 275 | ± 10 |
| 48 | 163 | ± 5 | 60 | 195 | ± 9 | 72 | 234 | ± 11 | 84 | 278 | ± 10 |
| 49 | 166 | ± 6 | 61 | 198 | ± 9 | 73 | 238 | ± 11 | 85 | 282 | ± 10 |
| 50 | 168 | ± 6 | 62 | 201 | ± 9 | 74 | 241 | ± 11 | 86 | 285 | ± 10 |
| 51 | 171 | ± 6 | 63 | 204 | ± 9 | 75 | 245 | ± 11 | >86 | n/a | — |
| 52 | 173 | ± 6 | 64 | 207 | ± 10 | 76 | 249 | ± 11 | | | |
| 53 | 176 | ± 7 | 65 | 210 | ± 10 | 77 | 252 | ± 11 | | | |
| 54 | 178 | ± 7 | 66 | 213 | ± 10 | 78 | 256 | ± 11 | | | |

Tokyo Shinozuka

Table 2-96: AC: Tokyo Shinozuka (Fetal Age)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: AC (cm); Age (Weeks/Days); SD (cm)

| AC | Age | 1SD | AC | Age | 1SD | AC | Age | 1SD |
|-----|-------|-----|----|-------|-----|-----|-------|-----|
| <10 | n/a | — | 18 | 23w3d | 0.9 | 27 | 33w1d | 1.4 |
| 10 | 15w3d | 0.5 | 19 | 24w3d | 1.0 | 28 | 34w2d | 1.4 |
| 11 | 16w4d | 0.6 | 20 | 25w3d | 1.0 | 29 | 35w4d | 1.5 |
| 12 | 17w4d | 0.6 | 21 | 26w3d | 1.1 | 30 | 37w0d | 1.6 |
| 13 | 18w4d | 0.7 | 22 | 27w3d | 1.1 | 31 | 38w2d | 1.6 |
| 14 | 19w4d | 0.7 | 23 | 28w4d | 1.2 | 32 | 39w6d | 1.7 |
| 15 | 20w3d | 0.8 | 24 | 29w4d | 1.2 | 33 | 41w2d | 1.8 |
| 16 | 21w3d | 0.8 | 25 | 30w5d | 1.3 | >33 | n/a | — |
| 17 | 22w3d | 0.9 | 26 | 31w6d | 1.3 | | | |

Table 2-97: AC: Tokyo Shinozuka (Fetal Growth)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: Age (Weeks/Days); Min/Mean/Max (cm); Table/Graph Range: 1.64SD

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|------|------|------|-----|------|------|------|
| 16 | 9.3 | 10.9 | 12.5 | 30 | 22.0 | 24.7 | 27.3 |
| 17 | 10.3 | 12.0 | 13.6 | 31 | 22.8 | 25.6 | 28.3 |
| 18 | 11.2 | 13.0 | 14.7 | 32 | 23.5 | 26.5 | 29.2 |
| 19 | 12.2 | 14.0 | 15.8 | 33 | 24.3 | 27.3 | 30.1 |
| 20 | 13.1 | 15.1 | 16.9 | 34 | 25.0 | 28.1 | 31.0 |
| 21 | 14.0 | 16.1 | 18.0 | 35 | 25.7 | 28.9 | 31.9 |
| 22 | 15.0 | 17.1 | 19.1 | 36 | 26.4 | 29.7 | 32.7 |
| 23 | 15.9 | 18.1 | 20.2 | 37 | 27.0 | 30.4 | 33.5 |
| 24 | 16.8 | 19.1 | 21.2 | 38 | 27.6 | 31.1 | 34.3 |
| 25 | 17.7 | 20.1 | 22.3 | 39 | 28.2 | 31.8 | 35.0 |
| 26 | 18.6 | 21.0 | 23.3 | 40 | 28.8 | 32.4 | 35.7 |
| 27 | 19.5 | 22.0 | 24.4 | 41 | 29.3 | 33.0 | 36.4 |
| 28 | 20.3 | 22.9 | 25.4 | 42 | 29.7 | 33.6 | 37.0 |
| 29 | 21.1 | 23.8 | 26.4 | | | | |

Table 2-98: AxT (APTDxTTD): Tokyo Shinozuka (Fetal Age)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: AxT (mm); Age (Weeks/Days); SD (cm²)

| AxT | Age | 1SD | AxT | Age | 1SD | AxT | Age | 1SD |
|-----|-------|-----|-----|-------|-----|-----|-------|------|
| <10 | n/a | — | 38 | 25w6d | 5.5 | 68 | 33w3d | 8.8 |
| 10 | 16w1d | 2.5 | 40 | 26w3d | 5.7 | 70 | 33w6d | 9.1 |
| 12 | 17w0d | 2.7 | 42 | 27w0d | 6.0 | 72 | 34w2d | 9.3 |
| 14 | 17w6d | 2.9 | 44 | 27w3d | 6.1 | 74 | 34w6d | 9.6 |
| 16 | 18w4d | 3.1 | 46 | 28w0d | 6.4 | 76 | 35w3d | 9.9 |
| 18 | 19w3d | 3.4 | 48 | 28w4d | 6.6 | 78 | 35w6d | 10.1 |
| 20 | 20w1d | 3.6 | 50 | 29w0d | 6.8 | 80 | 36w3d | 10.2 |
| 22 | 20w6d | 3.8 | 52 | 29w3d | 7.0 | 82 | 37w0d | 10.7 |
| 24 | 21w4d | 4.0 | 54 | 30w0d | 7.2 | 84 | 37w4d | 11.0 |
| 26 | 22w2d | 4.3 | 56 | 30w3d | 7.4 | 86 | 38w1d | 11.3 |
| 28 | 22w6d | 4.4 | 58 | 31w0d | 7.7 | 88 | 38w5d | 11.7 |
| 30 | 23w4d | 4.7 | 60 | 31w3d | 7.9 | 90 | 39w2d | 12.0 |
| 32 | 24w1d | 4.9 | 62 | 31w6d | 8.1 | >90 | n/a | — |
| 34 | 24w5d | 5.1 | 64 | 32w3d | 8.4 | | | |
| 36 | 25w2d | 5.3 | 66 | 32w6d | 8.6 | | | |

Table 2-99: AxT (APTDxTTD): Tokyo Shinozuka (Fetal Growth)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: Age (Weeks); Min/Mean/Max (cm²); Table/Graph Range: 1.64SD

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|------|------|------|-----|------|-------|-------|
| 16 | 7.0 | 11.2 | 15.5 | 30 | 43.5 | 55.7 | 68.0 |
| 17 | 8.7 | 13.3 | 18.0 | 31 | 46.8 | 59.7 | 72.7 |
| 18 | 10.5 | 15.6 | 20.7 | 32 | 50.0 | 63.8 | 77.6 |
| 19 | 12.5 | 18.1 | 23.6 | 33 | 53.3 | 67.8 | 82.4 |
| 20 | 14.7 | 20.8 | 26.8 | 34 | 56.5 | 71.9 | 87.3 |
| 21 | 17.1 | 23.6 | 30.2 | 35 | 59.7 | 75.9 | 92.2 |
| 22 | 19.6 | 26.7 | 33.8 | 36 | 62.8 | 79.9 | 97.0 |
| 23 | 22.2 | 29.9 | 37.5 | 37 | 65.9 | 83.9 | 101.9 |
| 24 | 25.0 | 33.2 | 41.5 | 38 | 68.8 | 87.7 | 106.7 |
| 25 | 27.9 | 36.7 | 45.6 | 39 | 71.6 | 91.5 | 111.4 |
| 26 | 30.9 | 40.3 | 49.8 | 40 | 74.3 | 95.1 | 116.0 |
| 27 | 33.9 | 44.1 | 54.2 | 41 | 76.8 | 98.6 | 120.5 |
| 28 | 37.1 | 47.9 | 58.7 | 42 | 79.1 | 102.0 | 124.8 |
| 29 | 40.3 | 51.8 | 63.3 | | | | |

Table 2-100: BPD: Tokyo Shinozuka (Fetal Age)
 Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996
 Unit: BPD (mm); Age (Weeks/Days); SD (mm)

| BPD | Age | 1SD | BPD | Age | 1SD | BPD | Age | 1SD |
|-----|-------|-----|-----|-------|-----|-----|-------|-----|
| <13 | n/a | — | 39 | 17w6d | 2.7 | 66 | 26w3d | 3.2 |
| 13 | 10w1d | 2.3 | 40 | 18w1d | 2.7 | 67 | 26w6d | 3.2 |
| 14 | 10w3d | 2.3 | 41 | 18w3d | 2.8 | 68 | 27w2d | 3.3 |
| 15 | 10w5d | 2.3 | 42 | 18w5d | 2.8 | 69 | 27w4d | 3.3 |
| 16 | 11w0d | 2.3 | 43 | 19w0d | 2.8 | 70 | 28w0d | 3.3 |
| 17 | 11w2d | 2.4 | 44 | 19w2d | 2.8 | 71 | 28w3d | 3.3 |
| 18 | 11w4d | 2.4 | 45 | 19w4d | 2.8 | 72 | 28w5d | 3.3 |
| 19 | 11w6d | 2.4 | 46 | 20w0d | 2.8 | 73 | 29w1d | 3.4 |
| 20 | 12w1d | 2.4 | 47 | 20w2d | 2.9 | 74 | 29w4d | 3.4 |
| 21 | 12w3d | 2.4 | 48 | 20w4d | 2.9 | 75 | 30w0d | 3.4 |
| 22 | 12w6d | 2.4 | 49 | 20w6d | 2.9 | 76 | 30w3d | 3.4 |
| 23 | 13w1d | 2.5 | 50 | 21w1d | 2.9 | 77 | 30w6d | 3.4 |
| 24 | 13w3d | 2.5 | 51 | 21w3d | 2.9 | 78 | 31w2d | 3.5 |
| 25 | 13w5d | 2.5 | 52 | 21w6d | 2.9 | 79 | 31w5d | 3.5 |
| 26 | 14w0d | 2.5 | 53 | 22w1d | 3.0 | 80 | 32w1d | 3.5 |
| 27 | 14w2d | 2.5 | 54 | 22w3d | 3.0 | 81 | 32w5d | 3.6 |
| 28 | 14w4d | 2.5 | 55 | 22w5d | 3.0 | 82 | 33w1d | 3.6 |
| 29 | 14w6d | 2.6 | 56 | 23w1d | 3.0 | 83 | 33w5d | 3.6 |
| 30 | 15w1d | 2.6 | 57 | 23w3d | 3.0 | 84 | 34w2d | 3.6 |
| 31 | 15w3d | 2.6 | 58 | 23w5d | 3.1 | 85 | 34w6d | 3.7 |
| 32 | 15w5d | 2.6 | 59 | 24w1d | 3.1 | 86 | 35w3d | 3.7 |
| 33 | 16w0d | 2.6 | 60 | 24w3d | 3.1 | 87 | 36w0d | 3.7 |
| 34 | 16w2d | 2.6 | 61 | 24w5d | 3.1 | 88 | 36w5d | 3.8 |
| 35 | 16w4d | 2.7 | 62 | 25w1d | 3.1 | 89 | 37w4d | 3.8 |
| 36 | 16w6d | 2.7 | 63 | 25w3d | 3.1 | 90 | 38w3d | 3.9 |
| 37 | 17w1d | 2.7 | 64 | 25w5d | 3.2 | >90 | n/a | — |
| 38 | 17w4d | 2.7 | 65 | 26w1d | 3.2 | | | |

Table 2-101: BPD: Tokyo Shinozuka (Fetal Growth)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range: 1.64SD

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|------|------|------|-----|------|------|------|
| 10 | 10.5 | 14.3 | 18.1 | 27 | 63.4 | 68.7 | 74.1 |
| 11 | 13.7 | 17.6 | 21.5 | 28 | 65.9 | 71.4 | 76.8 |
| 12 | 17.0 | 21.0 | 25.0 | 29 | 68.3 | 73.9 | 79.4 |
| 13 | 20.4 | 24.4 | 28.5 | 30 | 70.6 | 76.3 | 81.9 |
| 14 | 23.7 | 27.8 | 32.0 | 31 | 72.8 | 78.5 | 84.2 |
| 15 | 27.0 | 31.2 | 35.5 | 32 | 74.8 | 80.6 | 86.5 |
| 16 | 30.3 | 34.6 | 39.0 | 33 | 76.7 | 82.6 | 88.5 |
| 17 | 33.5 | 38.0 | 42.4 | 34 | 78.5 | 84.5 | 90.4 |
| 18 | 36.8 | 41.3 | 45.8 | 35 | 80.1 | 86.1 | 92.2 |
| 19 | 40.0 | 44.6 | 49.2 | 36 | 81.5 | 87.6 | 93.8 |
| 20 | 43.2 | 47.9 | 52.6 | 37 | 82.7 | 89.0 | 95.2 |
| 21 | 46.3 | 51.1 | 55.9 | 38 | 83.8 | 90.1 | 96.5 |
| 22 | 49.3 | 54.2 | 59.1 | 39 | 84.6 | 91.1 | 97.5 |
| 23 | 52.3 | 57.3 | 62.3 | 40 | 85.3 | 91.8 | 98.4 |
| 24 | 55.2 | 60.3 | 65.3 | 41 | 85.8 | 92.4 | 99.0 |
| 25 | 58.0 | 63.2 | 68.4 | 42 | 86.0 | 92.8 | 99.5 |
| 26 | 60.8 | 66.0 | 71.3 | | | | |

Table 2-102: CRL: Tokyo Shinozuka (Fetal Age)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: CRL (mm); Age (Weeks/Days); SD (mm)

| CRL | Age | 1SD | CRL | Age | 1SD | CRL | Age | 1SD |
|-----|------|-----|-----|-------|-----|-----|-------|-----|
| <5 | n/a | — | 20 | 8w6d | 3.7 | 36 | 10w6d | 5.9 |
| 5 | 6w3d | 1.1 | 21 | 9w0d | 3.9 | 37 | 11w0d | 6.0 |
| 6 | 6w4d | 1.3 | 22 | 9w1d | 4.0 | 38 | 11w0d | 6.0 |
| 7 | 6w6d | 1.6 | 23 | 9w2d | 4.2 | 39 | 11w1d | 6.2 |
| 8 | 7w0d | 1.7 | 24 | 9w3d | 4.3 | 40 | 11w2d | 6.3 |
| 9 | 7w1d | 1.9 | 25 | 9w4d | 4.5 | 41 | 11w3d | 6.5 |
| 10 | 7w2d | 2.0 | 26 | 9w4d | 4.5 | 42 | 11w3d | 6.5 |
| 11 | 7w3d | 2.2 | 27 | 9w5d | 4.6 | 43 | 11w4d | 6.6 |
| 12 | 7w4d | 2.3 | 28 | 9w6d | 4.8 | 44 | 11w5d | 6.8 |
| 13 | 7w5d | 2.5 | 29 | 10w0d | 4.9 | 45 | 11w6d | 6.9 |
| 14 | 7w6d | 2.6 | 30 | 10w1d | 5.1 | 46 | 11w6d | 6.9 |
| 15 | 8w1d | 2.9 | 31 | 10w2d | 5.2 | 47 | 12w0d | 7.1 |
| 16 | 8w2d | 3.1 | 32 | 10w3d | 5.4 | 48 | 12w1d | 7.2 |
| 17 | 8w3d | 3.3 | 33 | 10w4d | 5.5 | 49 | 12w1d | 7.2 |
| 18 | 8w4d | 3.4 | 34 | 10w5d | 5.7 | 50 | 12w2d | 7.4 |
| 19 | 8w5d | 3.6 | 35 | 10w6d | 5.9 | >50 | n/a | — |

Tokyo Shinozuka

Table 2-103: CRL: Tokyo Shinozuka (Fetal Growth)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: Age (Weeks/Days); Mean (mm); SD (mm); Table/Graph Range: 1.64SD

| Age | Mean | SD | Age | Mean | SD | Age | Mean | SD |
|------|------|-----|-------|------|-----|-------|------|-----|
| 7w0d | 7.9 | 1.7 | 9w1d | 21.2 | 4.0 | 11w2d | 40.0 | 6.3 |
| 7w1d | 8.6 | 1.9 | 9w2d | 22.3 | 4.2 | 11w3d | 41.4 | 6.5 |
| 7w2d | 9.3 | 2.0 | 9w3d | 23.4 | 4.3 | 11w4d | 42.9 | 6.6 |
| 7w3d | 10.1 | 2.2 | 9w4d | 24.5 | 4.5 | 11w5d | 44.4 | 6.8 |
| 7w4d | 10.9 | 2.3 | 9w5d | 25.7 | 4.6 | 11w6d | 45.9 | 6.9 |
| 7w5d | 11.7 | 2.5 | 9w6d | 26.8 | 4.8 | 12w0d | 47.4 | 7.1 |
| 7w6d | 12.5 | 2.6 | 10w0d | 28.0 | 4.9 | 12w1d | 49.0 | 7.2 |
| 8w0d | 13.4 | 2.8 | 10w1d | 29.3 | 5.1 | 12w2d | 50.6 | 7.4 |
| 8w1d | 14.3 | 2.9 | 10w2d | 30.5 | 5.2 | 12w3d | 52.2 | 7.5 |
| 8w2d | 15.2 | 3.1 | 10w3d | 31.8 | 5.4 | 12w4d | 53.9 | 7.7 |
| 8w3d | 16.1 | 3.3 | 10w4d | 33.1 | 5.5 | 12w5d | 55.5 | 7.8 |
| 8w4d | 17.1 | 3.4 | 10w5d | 34.4 | 5.7 | 12w6d | 57.2 | 8.0 |
| 8w5d | 18.1 | 3.6 | 10w6d | 35.8 | 5.9 | 13w0d | 58.9 | 8.2 |
| 8w6d | 19.1 | 3.7 | 11w0d | 37.1 | 6.0 | | | |
| 9w0d | 20.1 | 3.9 | 11w1d | 38.5 | 6.2 | | | |

Table 2-104: EFW: Tokyo Shinozuka (Fetal Age)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: EFW (grams); Age (Weeks/Days); SD (grams)

| EFW | Age | 1SD | EFW | Age | 1SD | EFW | Age | 1SD |
|------|-------|-----|------|-------|-----|-------|-------|-----|
| <250 | n/a | — | 1200 | 28w3d | 162 | 2200 | 34w2d | 258 |
| 250 | 19w3d | 45 | 1250 | 28w5d | 166 | 2250 | 34w4d | 264 |
| 300 | 20w0d | 51 | 1300 | 29w1d | 173 | 2300 | 34w6d | 269 |
| 350 | 20w4d | 58 | 1350 | 29w3d | 177 | 2350 | 35w1d | 274 |
| 400 | 21w2d | 66 | 1400 | 29w5d | 181 | 2400 | 35w3d | 279 |
| 450 | 21w5d | 71 | 1450 | 30w0d | 186 | 2450 | 35w5d | 284 |
| 500 | 22w2d | 78 | 1500 | 30w2d | 191 | 2500 | 35w7d | 290 |
| 550 | 22w6d | 85 | 1550 | 30w5d | 197 | 2550 | 36w2d | 295 |
| 600 | 23w2d | 90 | 1600 | 30w7d | 202 | 2600 | 36w4d | 301 |
| 650 | 23w6d | 98 | 1650 | 31w2d | 207 | 2650 | 36w6d | 306 |
| 700 | 24w2d | 103 | 1700 | 31w4d | 211 | 2700 | 37w2d | 314 |
| 750 | 24w5d | 109 | 1750 | 31w6d | 216 | 2750 | 37w4d | 320 |
| 800 | 25w2d | 116 | 1800 | 32w1d | 221 | 2800 | 37w6d | 325 |
| 850 | 25w5d | 122 | 1850 | 32w3d | 226 | 2850 | 38w1d | 331 |
| 900 | 26w1d | 128 | 1900 | 32w5d | 231 | 2900 | 38w4d | 340 |
| 950 | 26w4d | 134 | 1950 | 32w7d | 236 | 2950 | 38w6d | 345 |
| 1000 | 26w6d | 138 | 2000 | 33w1d | 238 | 3000 | 39w2d | 354 |
| 1050 | 27w2d | 145 | 2050 | 33w3d | 243 | >3000 | n/a | — |
| 1100 | 27w5d | 151 | 2100 | 33w5d | 248 | | | |
| 1150 | 28w0d | 155 | 2150 | 34w0d | 253 | | | |

Table 2-105: EFW: Tokyo Shinozuka (Fetal Growth)
 Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996
 Unit: Age (Weeks); Min/Mean/Max (grams); Table/Graph Range: 1.64SD

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|------|------|------|-----|------|------|------|
| 18 | 158 | 216 | 274 | 30 | 1234 | 1552 | 1870 |
| 19 | 204 | 279 | 355 | 31 | 1375 | 1720 | 2064 |
| 20 | 256 | 349 | 442 | 32 | 1520 | 1892 | 2265 |
| 21 | 314 | 427 | 539 | 33 | 1667 | 2068 | 2469 |
| 22 | 381 | 513 | 645 | 34 | 1814 | 2244 | 2675 |
| 23 | 456 | 609 | 761 | 35 | 1960 | 2420 | 2880 |
| 24 | 541 | 714 | 888 | 36 | 2102 | 2592 | 3083 |
| 25 | 634 | 830 | 1026 | 37 | 2236 | 2758 | 3280 |
| 26 | 737 | 956 | 1175 | 38 | 2360 | 2915 | 3469 |
| 27 | 849 | 1092 | 1334 | 39 | 2471 | 3059 | 3647 |
| 28 | 970 | 1237 | 1504 | 40 | 2565 | 3187 | 3809 |
| 29 | 1099 | 1391 | 1683 | 41 | 2639 | 3296 | 3952 |

Table 2-106: FL: Tokyo Shinozuka (Fetal Age)
 Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996
 Unit: FL (mm); Age (Weeks/Days); SD (mm)

| FL | Age | 1SD | FL | Age | 1SD | FL | Age | 1SD |
|-----|-------|-----|----|-------|-----|-----|-------|-----|
| <20 | n/a | — | 37 | 22w2d | 2.9 | 55 | 30w5d | 3.1 |
| 20 | 16w1d | 2.6 | 38 | 22w5d | 2.9 | 56 | 31w2d | 3.2 |
| 21 | 16w3d | 2.7 | 39 | 23w1d | 2.9 | 57 | 31w6d | 3.2 |
| 22 | 16w6d | 2.7 | 40 | 23w4d | 2.9 | 58 | 32w3d | 3.2 |
| 23 | 17w1d | 2.7 | 41 | 24w0d | 2.9 | 59 | 33w0d | 3.2 |
| 24 | 17w3d | 2.7 | 42 | 24w3d | 2.9 | 60 | 33w3d | 3.2 |
| 25 | 17w6d | 2.7 | 43 | 24w6d | 2.9 | 61 | 34w0d | 3.2 |
| 26 | 18w1d | 2.7 | 44 | 25w3d | 3.0 | 62 | 34w4d | 3.3 |
| 27 | 18w3d | 2.7 | 45 | 25w6d | 3.0 | 63 | 35w1d | 3.3 |
| 28 | 18w6d | 2.7 | 46 | 26w2d | 3.0 | 64 | 35w5d | 3.3 |
| 29 | 19w1d | 2.7 | 47 | 26w5d | 3.0 | 65 | 36w2d | 3.3 |
| 30 | 19w4d | 2.8 | 48 | 27w2d | 3.0 | 66 | 37w0d | 3.3 |
| 31 | 20w0d | 2.8 | 49 | 27w5d | 3.0 | 67 | 37w4d | 3.4 |
| 32 | 20w2d | 2.8 | 50 | 28w2d | 3.1 | 68 | 38w1d | 3.4 |
| 33 | 20w5d | 2.8 | 51 | 28w5d | 3.1 | 69 | 38w5d | 3.4 |
| 34 | 21w1d | 2.8 | 52 | 29w2d | 3.1 | 70 | 39w3d | 3.4 |
| 35 | 21w3d | 2.8 | 53 | 29w5d | 3.1 | >70 | n/a | — |
| 36 | 21w6d | 2.8 | 54 | 30w2d | 3.1 | | | |

Tokyo Shinozuka

Table 2-107: FL: Tokyo Shinozuka (Fetal Growth)

Shinozuka Jpn J Med Ultrasonics vol 23: 12 1996

Unit: Age (Weeks); Min/Mean/Max (mm); Table/Graph Range: 1.64SD

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|-----|------|------|------|-----|------|------|------|
| 16 | 17.1 | 21.4 | 25.8 | 30 | 49.7 | 54.8 | 60.0 |
| 17 | 19.6 | 24.0 | 28.4 | 31 | 51.6 | 56.8 | 62.0 |
| 18 | 22.1 | 26.5 | 31.0 | 32 | 53.5 | 58.7 | 64.0 |
| 19 | 24.6 | 29.1 | 33.6 | 33 | 55.2 | 60.5 | 65.8 |
| 20 | 27.1 | 31.6 | 36.2 | 34 | 56.9 | 62.2 | 67.6 |
| 21 | 29.5 | 34.1 | 38.8 | 35 | 58.4 | 63.8 | 69.2 |
| 22 | 31.9 | 36.6 | 41.3 | 36 | 59.9 | 65.3 | 70.8 |
| 23 | 34.3 | 39.1 | 43.8 | 37 | 61.2 | 66.7 | 72.2 |
| 24 | 36.7 | 41.5 | 46.3 | 38 | 62.4 | 68.0 | 73.6 |
| 25 | 39.0 | 43.9 | 48.7 | 39 | 63.5 | 69.1 | 74.7 |
| 26 | 41.3 | 46.2 | 51.1 | 40 | 64.4 | 70.1 | 75.8 |
| 27 | 43.5 | 48.4 | 53.4 | 41 | 65.3 | 71.0 | 76.7 |
| 28 | 45.6 | 50.6 | 55.7 | 42 | 65.9 | 71.7 | 77.5 |
| 29 | 47.7 | 52.8 | 57.9 | | | | |

Williams

Table 2-108: EFW: Williams (Fetal Growth)
Unit: Age (Weeks); Min/Mean/Max (grams)

| Age | Min | Mean | Max | Age | Min | Mean | Max |
|------|------|------|------|------|------|------|------|
| 22.0 | 320 | 513 | 746 | 34.0 | 1728 | 2394 | 3132 |
| 23.0 | 365 | 589 | 861 | 35.0 | 1974 | 2628 | 3333 |
| 24.0 | 417 | 675 | 989 | 36.0 | 2224 | 2849 | 3521 |
| 25.0 | 477 | 773 | 1132 | 37.0 | 2455 | 3052 | 3706 |
| 26.0 | 546 | 882 | 1289 | 38.0 | 2642 | 3227 | 3867 |
| 27.0 | 627 | 1005 | 1463 | 39.0 | 2790 | 3364 | 3994 |
| 28.0 | 720 | 1143 | 1653 | 40.0 | 2881 | 3462 | 4080 |
| 29.0 | 829 | 1298 | 1859 | 41.0 | 2946 | 3524 | 4127 |
| 30.0 | 955 | 1484 | 2136 | 42.0 | 3011 | 3589 | 4185 |
| 31.0 | 1100 | 1695 | 2402 | 43.0 | 3044 | 3626 | 4221 |
| 32.0 | 1284 | 1920 | 2673 | 44.0 | 3043 | 3633 | 4233 |
| 33.0 | 1499 | 2155 | 2910 | | | | |

Yarkoni

Table 2-109: CLA:Yarkoni S, Journal of Ultrasound in Medicine, 4:467-470, 1985
(Fetal Age)

Unit: Meas (mm); Min/Mean/Max (Weeks/Days)

| Meas | Min | Mean | Max | Meas | Min | Mean | Max |
|------|-------|-------|-------|------|-------|-------|-------|
| 11 | 8w3d | 13w6d | 17w2d | 29 | 23w2d | 28w5d | 32w1d |
| 12 | 9w1d | 14w4d | 18w1d | 30 | 24w0d | 29w4d | 34w0d |
| 13 | 10w0d | 14w3d | 19w6d | 31 | 25w6d | 29w2d | 34w6d |
| 14 | 11w6d | 15w2d | 20w5d | 32 | 26w5d | 30w1d | 35w4d |
| 15 | 12w5d | 16w1d | 21w4d | 33 | 27w4d | 31w0d | 35w3d |
| 16 | 12w3d | 18w0d | 21w3d | 34 | 27w3d | 32w6d | 36w2d |
| 17 | 13w2d | 18w5d | 22w2d | 35 | 28w1d | 33w5d | 37w1d |
| 18 | 14w1d | 19w4d | 23w0d | 36 | 29w0d | 33w3d | 39w0d |
| 19 | 16w0d | 19w3d | 24w6d | 37 | 30w6d | 34w2d | 39w5d |
| 20 | 16w6d | 20w2d | 25w5d | 38 | 31w5d | 35w1d | 40w4d |
| 21 | 17w4d | 21w1d | 26w4d | 39 | 32w4d | 37w0d | 40w3d |
| 22 | 17w3d | 22w6d | 26w2d | 40 | 32w2d | 37w6d | 41w2d |
| 23 | 18w2d | 23w5d | 27w1d | 41 | 33w1d | 38w4d | 42w0d |
| 24 | 19w1d | 24w4d | 28w0d | 42 | 35w0d | 38w3d | 43w6d |
| 25 | 21w0d | 24w3d | 29w6d | 43 | 35w6d | 39w2d | 44w5d |
| 26 | 21w5d | 25w1d | 30w5d | 44 | 36w5d | 40w1d | 45w4d |
| 27 | 22w4d | 26w0d | 30w3d | 45 | 36w3d | 41w6d | 45w3d |
| 28 | 22w3d | 27w6d | 31w2d | | | | |

Chapter 3

Acoustic information

This chapter includes the following information:

| | |
|--|------------|
| • The real-time display of acoustic output indices | 153 |
| • Thermal Index | 153 |
| • Mechanical Index: | 154 |
| • Concerns Surrounding the Use of Diagnostic Ultrasound | 155 |
| • Default Settings and Output Levels | 155 |
| • Controls Affecting Acoustic Output | 156 |
| • Track 3 Summary Table | 158 |
| • Acoustic Parameters as Measured in Water | 161 |
| • Definitions, symbols and abbreviations | 161 |
| • Translations of definitions, symbols, and abbreviations | 162 |
| • Acoustic Output Reporting Tables for Track 3/IEC 60601-2-37 | 186 |
| • Explanation of Footnotes | 186 |
| • Multiple focal-zones | 186 |
| • Operating Conditions | 187 |
| • Transducer Model: 3Sc-RS | 188 |
| • Transducer Model: M4S-RS | 194 |
| • Transducer Model: 5S-RS | 200 |
| • Transducer Model: 6S-RS | 206 |
| • Transducer Model: 7S-RS | 212 |
| • Transducer Model: 10S-RS | 218 |
| • Transducer Model: 12S-RS | 224 |
| • Transducer Model: e8C-RS | 230 |
| • Transducer Model: 3C-RS | 235 |
| • Transducer Model: 4C-RS | 240 |
| • Transducer Model: 8C-RS | 245 |

| | |
|---|-----|
| • Transducer Model: 8L-RS..... | 250 |
| • Transducer Model: 9L-RS..... | 255 |
| • Transducer Model: 12L-RS..... | 260 |
| • Transducer Model: i12L-RS..... | 265 |
| • Transducer Model: 6T-RS | 270 |
| • Transducer Model: 6Tc-RS..... | 276 |
| • Transducer Model: 9T-RS | 282 |
| • Transducer Model: P2D-RS..... | 288 |
| • Transducer Model: P6D-RS..... | 289 |
| • Transducer Model: AcuNav8 | 290 |
| • Transducer Model: AcuNav™ 10/SoundStar™ 3D 10FG / eco 10FG | 296 |

The real-time display of acoustic output indices

The Vivid q N has real-time display features according to Track 3 in the FDA 510(k) Guidance. It displays both a thermal (TI) and a mechanical (MI) index in all operating modes. These two indices are intended to estimate the potential for thermal and mechanical bioeffects induced by ultrasound. Both TI and MI are displayed with increments of 0.1. Neither are displayed if the value is below 0.4. The displayed (estimated) TI and MI are nominal values.

Thermal Index

TI is defined as:

$$TI = \frac{W_0}{W_{deg}}$$

where: W_0 is the time-averaged acoustic power and W_{deg} is the estimated power necessary to raise the target tissue one degree C.

The displayed TI is an estimate of temperature increase of soft tissue or bone, presented to make it easier for the operator to implement the ALARA (As Low As Reasonably Achievable) principle. There are three thermal index categories:

- **TIS:** Soft tissue thermal index. The main TI category. Used for applications that do not image bone.
- **TIB:** Bone thermal index (bone located in a focal region). Used for fetal application.
- **TIC:** Cranial bone thermal index (bone located close to the surface). Used for transcranial application.

The Vivid q N chooses the correct category based on mode of operation and chosen application, and presents only one TI to the operator. It is therefore important that the operator chooses the right application.

Mechanical Index:

MI is the estimated likelihood of tissue damage due to cavitation. MI is defined as:

$$MI = \frac{p_{r.3}(z_{sp})}{\sqrt{f_c}}$$

where $p_{r.3}$ is the derated (attenuated) peak rarefactional (negative) pressure (MPa) and f_c is the center frequency (MHz).

The MI will not exceed a value of 1.9 according to Track 3 in the FDA 510(k) guidance of 2008.

Display Accuracy and Acoustic Measurement Uncertainties

The display accuracy and measurement precision of the output display are summarized in the table below.

Accuracy of the output display (TI, MI) parameters depends on the measurement system precision, the acoustic model used to calculate the parameters and variation in the acoustic output of probes and systems. The measurement precision and overall accuracy of the measurements have been assessed by determining both the random and the systematic uncertainties and given in percent at 95% confidence level.

| Parameter | Estimated accuracy ^a | Measurement precision | |
|--------------|---------------------------------|-----------------------|-------------|
| | | M/Color/PW/CW | 2D/CFM mode |
| Pressure, MI | ± 25% | ± 15% | |
| Power, TI | ± 50% | ± 30% | ± 40% |
| Frequency | ± 1% ^b | ± 1% | |

a. Accuracy = (Measured value - displayed value)/displayed value * 100%

b. The displayed Frequency value shown on the screen may differ from the actual value of f_c which is used to calculate the estimated MI parameter

Concerns Surrounding the Use of Diagnostic Ultrasound

During a diagnostic ultrasound examination, high frequency sound penetrates and interacts with tissue in and around the area of anatomy to be imaged. Only a small portion of this sound energy is reflected back to the probe for use in constructing the image while the remainder is dissipated within the tissue. The interaction of sound energy with tissue at sufficiently high levels can produce biological effects (aka bioeffects) of either a mechanical or thermal nature. Although the generation of bioeffect is intentional with therapeutic ultrasound, it is generally undesired in diagnostic applications and may be harmful in some conditions.

Note: The American Institute of Ultrasound in Medicine (AIUM) has published a document entitled "Medical Ultrasound Safety". This three part document covers Bioeffects and Biophysics, Prudent Use and Implementing ALARA.

Ultrasound users should read the AIUM documents to become more familiar with Ultrasound safety. A copy of this document is included as part of the DOC-CD (GE P/N 5390423).

To contact the AIUM concerning their publications:

- In the USA, by telephone at 1-800-638-5352.
- To write them, use the following address:

AIUM
14750 Sweitzer Lane
Suite 100
Laurel, MD, USA 20707-5906

Default Settings and Output Levels

The default acoustic output level will not exceed a thermal index (TI) of 3.0 or a mechanical index (MI) of 1.5.

The maximum default TI is 50% of the maximum possible TI (6.0) and the maximum default MI is 80% of the maximum possible MI (1.9).

The output level will not exceed the default level until the user intentionally increases the power level by adjusting the power control on the system.

The output level will return to default each time

- a new probe is chosen
- a new application is chosen
- a new patient is chosen.

Controls Affecting Acoustic Output

The initial means by which the user can affect acoustic output are by 1) selecting a probe, 2) selecting an application (exam category) and then 3) selecting the imaging mode or particular imaging characteristics. After these selections are made, the only user control that can affect the output is the acoustic output control. This is achieved through an acoustic output control scheme in which all parameters that directly or indirectly affect acoustic output are fed to the control algorithm. The algorithm estimates all relevant parameters and compares them to the FDA limits.

Output levels remain below the limits with a 90% confidence margin. The absolute maximum allowable output for all applications is:

- $ISPTA \leq 720 \text{ mW/cm}^2$
- $MI \leq 1.9$
- $TI \leq 6$




Figure 3-1 The Vivid q N Acoustic Output Control Scheme

Track 3 Summary Table

The following table summarizes the mode/probe combinations for which the global maximum displayed MI or TI may be greater than 1.0. For each probe/mode combination checked, a Track 3 acoustic output table exists.

Not all probes listed may be supported worldwide. Please refer to your local language User Manual for an overview of the probes that are supported in your country.

| Operating Mode | Transducer Model | | | |
|----------------|------------------|--------|--------|-------|
| | | 3Sc-RS | M4S-RS | 5S-RS |
| 2D | | Yes | Yes | Yes |
| M-Mode | | Yes | Yes | Yes |
| PW Doppler | | Yes | Yes | Yes |
| CW Doppler | | Yes | Yes | Yes |
| CFM | | Yes | Yes | Yes |
| CMM | | Yes | Yes | Yes |

| Operating Mode | Transducer Model | | | |
|----------------|------------------|-------|--------|--------|
| | 6S-RS | 7S-RS | 10S-RS | 12S-RS |
| 2D | Yes | Yes | Yes | Yes |
| M-Mode | Yes | Yes | Yes | Yes |
| PW Doppler | Yes | Yes | Yes | Yes |
| CW Doppler | Yes | Yes | Yes | Yes |
| CFM | Yes | Yes | Yes | Yes |
| CMM | Yes | Yes | Yes | Yes |

| Operating Mode | Transducer Model | | | |
|-------------------|------------------|-------|--------|---------|
| | 8L-RS | 9L-RS | 12L-RS | i12L-RS |
| 2D | Yes | Yes | Yes | Yes |
| M-Mode | Yes | Yes | Yes | Yes |
| PW Doppler | Yes | Yes | Yes | Yes |
| CW Doppler | - | - | - | - |
| CFM | Yes | Yes | Yes | Yes |
| CMM | Yes | Yes | Yes | Yes |

| Operating Mode | Transducer Model | | |
|-------------------|------------------|-------|--------|
| | 4C-RS | 8C-RS | e8C-RS |
| 2D | Yes | Yes | Yes |
| M-Mode | Yes | Yes | Yes |
| PW Doppler | Yes | Yes | Yes |
| CW Doppler | - | - | - |
| CFM | Yes | Yes | Yes |
| CMM | Yes | Yes | Yes |

| Operating Mode | Transducer Model | | |
|-------------------|------------------|--------|-------|
| | 6T-RS | 6Tc-RS | 9T-RS |
| 2D | Yes | Yes | Yes |
| M-Mode | Yes | Yes | Yes |
| PW Doppler | Yes | Yes | Yes |
| CW Doppler | Yes | Yes | Yes |
| CFM | Yes | Yes | Yes |
| CMM | Yes | Yes | Yes |

| Operating Mode | Transducer Model | |
|----------------|------------------|--------|
| | P2D-RS | P6D-RS |
| 2D | - | - |
| M-Mode | - | - |
| PW Doppler | - | - |
| CW Doppler | Yes | Yes |
| CFM | - | - |
| CMM | - | - |

| Operating Mode | Transducer Model | |
|----------------|------------------|---|
| | AcuNav8™ | AcuNav™ 10/ SoundStar™ 3D 10FG / eco 10FG |
| 2D | Yes | Yes |
| M-Mode | Yes | Yes |
| PW Doppler | Yes | Yes |
| CW Doppler | Yes | Yes |
| CFM | Yes | Yes |
| CMM | Yes | Yes |

Acoustic Parameters as Measured in Water

Definitions, symbols and abbreviations

The following definitions, symbols and abbreviations are used in the acoustic output reporting tables in this chapter:

| FDA | IEC | Meaning—IEC 60601-2-37 / FDA & NEMA UD2, UD3 |
|------------------------|-----------------------------|---|
| a | α | <i>Acoustic Attenuation Coefficient / Derating factor (usually 0.3 dB/cm-MHz)</i> |
| A_{aprt} | A_{aprt} | <i>-12db Output Beam Area / Active aperture area</i> |
| | C_{MI} | <i>Normalizing Coefficient</i> |
| D_{eq} | D_{eq} | <i>Equivalent Aperture Diameter / (same)</i> |
| d_{-6} | d_{-6} | <i>Pulse Beam Width / Beam diameter at -6 dB</i> |
| d_{eq} | d_{eq} | <i>Equivalent Beam Diameter</i> |
| f_c | f_{awf} | <i>Acoustic Working Frequency / Center frequency</i> |
| I_{pa} | I_{pa} | <i>Pulse-Average Intensity</i> |
| $I_{\text{pa},3}$ | $I_{\text{pa},\alpha}$ | <i>Attenuated Pulse-Average Intensity</i> |
| PII | I_{pi} | <i>Pulse-Intensity Integral</i> |
| $PII_{.3}$ | $I_{\text{pi},\alpha}$ | <i>Attenuated Pulse-Intensity Integral</i> |
| I_{TA} | $I_{\text{ta}}(z)$ | <i>Temporal-Average Intensity</i> |
| $I_{\text{TA},3}(Z)$ | $I_{\text{ta},\alpha}(z)$ | <i>Attenuated Temporal-Average Intensity / (at depth z)</i> |
| $I_{\text{SPTA}}(Z)$ | $I_{\text{zpta}}(z)$ | <i>Spatial-Peak Temporal-Average Intensity</i> |
| $I_{\text{SPTA},3}(Z)$ | $I_{\text{zpta},\alpha}(z)$ | <i>Attenuated Spatial-Peak Temporal-Average Intensity</i> |
| MI | MI | <i>Mechanical Index</i> |
| W_o | P | <i>Output Power / Time average acoustic power at the source</i> |
| $W_{.3}(Z)$ | P_{α} | <i>Attenuated Output Power / Time average acoustic power derated to depth z</i> |
| W_{o1} | P_1 | <i>Bounded Output Power / Power emitted from the central 1cm of aperture</i> |
| PII | p_i | <i>Pulse Pressure Squared Integral / Pulse intensity integral</i> |
| p_r | p_r | <i>Peak-Rarefactional Acoustic Pressure / (same)</i> |
| $p_{r,3}$ | $p_{r,\alpha}$ | <i>Attenuated Peak-Rarefactional Acoustic Pressure / (same)</i> |
| PRF | p_{rr} | <i>Pulse Repetition Rate / Pulse repetition frequency</i> |

Acoustic information

| FDA | IEC | Meaning—IEC 60601-2-37 / FDA & NEMA UD2, UD3 |
|------------------------------------|-----------------------|---|
| TI | <i>TI</i> | <i>Thermal Index / (same)</i> |
| TIB | <i>TIB</i> | <i>Bone Thermal Index / (same)</i> |
| TIC | <i>TIC</i> | <i>Cranial-Bone Thermal Index / (same)</i> |
| TIS | <i>TIS</i> | <i>Soft-Tissue Thermal Index / (same)</i> |
| PD | t_d | <i>Pulse Duration / (same)</i> |
| X ₋₁₂ :Y ₋₁₂ | X, Y | <i>-12 dB Output Beam Dimensions / (same)</i> |
| Z | <i>z</i> | <i>Distance from the Source to a Specified Point / (same)</i> |
| Z _{sp} | <i>z_b</i> | <i>Depth for TIB / Depth at which the relevant index is maximum</i> |
| Z _{bp} | <i>z_{bp}</i> | <i>Break-Point Depth / (same)</i> |
| Z _{sp} | <i>z_s</i> | <i>Depth for TIS / Depth at which the relevant index is maximum</i> |

Translations of definitions, symbols, and abbreviations

The following definitions, symbols and abbreviations, used in the acoustic output reporting tables, have been translated into several languages. See the following tables:

| Language | Page |
|-----------------------|------|
| English | 163 |
| Chinese | 165 |
| Danish | 166 |
| French | 168 |
| German | 170 |
| Italian | 172 |
| Japanese | 174 |
| Latin America Spanish | 176 |
| Brazilian Portuguese | 178 |
| Swedish | 180 |
| Norwegian | 182 |
| Finnish | 184 |

English

| Symbol | Unit | Definition |
|-------------------------------|--------------------|---|
| MI | n/a | Mechanical Index |
| TIS_{scan} | n/a | Soft Tissue Thermal Index in auto-scanning mode |
| TIS_{non-scan} | n/a | Soft Tissue Thermal Index in non-auto-scanning mode |
| TIB | n/a | Bone Thermal Index |
| TIC | n/a | Cranial Thermal Index |
| A_{aprt} | cm ² | Area of the active aperture |
| P_{r.3} | MPa | Derated peak rarefactional pressure (MPa) associated with the transmit pattern giving rise to the value reported under MI |
| W₀ | mW | Ultrasonic power, except for TIS_{scan} , in which case it is the ultrasonic power passing through a one centimeter window |
| W_{.3(z1)} | mW | Derated ultrasonic power at axial distance z₁ |
| I_{TA.3(z1)} | mW/cm ² | Derated spatial-peak, temporal-average intensity at axial distance z₁ |
| Z₁ | cm | Axial distance corresponding to the location of max [min(W_{.3(z)} , I_{TA.3(z)} × 1 cm ²)], where $z \geq z_{bp}$ |
| Z_{bp} | cm | $1.69(A_{aprt})^{1/2}$ |
| Z_{sp} | cm | For MI , the axial distance at which pr.3 is measured For TIB , the axial distance at which TIB is a maximum (i.e. $Z_{sp} = Z_{B.3}$) |
| d_{eq(z)} | cm | Equivalent beam diameter as a function of axial distance z , and is equal to $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, where I_{TA(z)} is the temporal-average intensity as a function of z |

Acoustic information

| Symbol | Unit | Definition |
|--|-------------------|--|
| f_c | Mhz | Center frequency For MI , f_c is the center frequency associated with the transmit pattern giving rise to the maximum reported value of MI . For TI , for combined modes involving transmit patterns of unequal center frequency, f_c is defined as the overall range of center frequencies of the respective transmit patterns |
| Dim. of A_{aprt} | cm | Active aperture dimensions for the azimuthal and elevational planes |
| PD | μs | Pulse duration associated with the transmit pattern giving rise to the reported value of MI |
| PRF | Hz | Pulse repetition frequency associated with the transmit pattern giving rise to the reported value of MI |
| P_r @ PII_{max} | MPa | Peak rarefactional pressure at the point where the freefield, spatial-peak pulse intensity integral is a maximum |
| D_{eq} @ PII_{max} | cm | Equivalent beam diameter at the point where the freefield, spatial-peak pulse intensity integral is a maximum |
| FL | cm | Focal length, or azimuthal and elevational lengths, if different |
| I_{PA.3} @ MI_{max} | W/cm ² | Derated pulse average intensity at the point of maximum reported MI |
| ROI | n/a | Region Of Interest |
| TB | n/a | Trackball |
| CF | n/a | Color Flow Mode |
| CM | n/a | Color M Mode |
| PW/CW | n/a | Pulsed Wave/Continuous Wave Doppler |

Chinese

| 符号 | 单位 | 定义 |
|--------------------------------------|--------------------|--|
| MI | 不适用 | 机械指数 |
| TIS _{scan} | 不适用 | 自动扫描模式下的软组织热敏指数 |
| TIS _{non-scan} | 不适用 | 非自动扫描模式下的软组织热敏指数 |
| TIB | 不适用 | 骨组织热敏指数 |
| TIC | 不适用 | 头盖骨热敏指数 |
| A _{aprt} | cm ² | 有效孔径区 |
| P _{r,3} | MPa | 与可以产生 MI 报告值的传送模式关联的减额最大稀薄压 (MPa) |
| W ₀ | mW | 超声功率, 除了在使用 TIS _{scan} 的情况下, 此时, 它是指一厘米窗口内通过的超声功率 |
| W _{.3(z₁)} | mW | 减额超声功率 (轴长 z ₁) |
| I _{TA,3(z₁)} | mW/cm ² | 减额空间峰值、时间平均强度 (轴长 z ₁) |
| z ₁ | cm | 与 max [min(W _{.3(z)} , I _{TA,3(z)} x 1 cm ²)] 位置相对应的轴长, 其中 z ≥ z _{bp} |
| z _{bp} | | 1.69(A _{aprt}) ^{1/2} |
| d _{eq(z)} | cm | 等效波束直径, 轴长 z 的函数, 等于 [(4/π)(W ₀ /I _{TA(z)})] ^{1/2} , 其中 I _{TA(z)} 为时间平均强度, z 的函数 |
| f _c | MHz | 中心频率 对于 MI, f _c 为与可以产生 MI 最大报告值的传送模式关联的中心频率。 对于 TI, 用于包括不同中心频率传送模式的组合模式时, f _c 定义为各个传送模式的所有中心频率 |
| Dim. of A _{aprt} | cm | 用于水平面和垂直平面的有效孔径 |
| PD | μs | 与可以产生 MI 报告值的传送模式关联的脉冲持续时间 |
| PRF | Hz | 与可以产生 MI 报告值的传送模式关联的脉冲重复频率 |
| P _{r @ PII_{max}} | MPa | 自由场、空间峰值脉冲强度积分最大处的峰值稀薄压 |
| d _{eq @ PII_{max}} | cm | 自由场、空间峰值脉冲强度积分最大处的等效波束直径 |
| FL | cm | 焦距, 或者水平长度和垂直高度 (如果不同) |
| I _{PA,3 @ MI_{max}} | W/cm ² | MI 最大报告值处的减额脉冲平均强度 |
| ROI | 不适用 | 兴趣区 |
| TB | 不适用 | 轨迹球 |
| CF | 不适用 | 彩色模式 |
| CM | 不适用 | 彩色 M 模式 |
| PW/CW | 不适用 | 脉冲波 / 连续波多普勒 |

Danish

| Symboler | Enhed | Definition |
|--|--------------------|--|
| MI | n/a | Mekanisk indeks |
| TIS_{scan} | n/a | Termisk indeks for blødt væv i automatisk scanningsmode |
| TIS_{non-scan} | n/a | Termisk indeks over blødt væv i ikke-automatisk scanningsmode |
| TIB | n/a | Termisk indeks for knogler |
| TIC | n/a | Termisk indeks for kranieknogle |
| A_{aptr} | cm ² | Område af den aktive blænde |
| P_{r.3} | MPa | Belastningsreduceret, fortyndet maksimumtryk (MPa), der er knyttet til det sendemønster, der giver værdien, som er angivet for MI |
| W_o | mW | Ultralydeffekt, undtagen for TIS_{scan} hvor ultralydeffekten passerer gennem et vindue på 1 cm |
| W_{.3(z₁)} | mW | Belastningsreduceret ultralydeffekt ved aksialafstand z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Belastningsreduceret, tidsmæssigt gennemsnitsintensitet med rumligt maksimum ved aksialafstand z₁ |
| z₁ | cm | Aksialafstanden svarer til placeringen af maks. [min. (W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], hvor $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{blænde})^{1/2}$ |
| d_{eq(z)} | cm | Tilsvarende strålediameter som funktion af aksialafstanden z , og lig med $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, hvor I_{TA(z)} er den tidsmæssige gennemsnitsintensitet som funktion af z |
| f_c | MHz | Centerfrekvens Vedr. MI er f_c centerfrekvensen, der er knyttet til det sendemønster, der giver den maksimalt rapporterede værdi af MI . Vedr. TI i kombinerede mode, der involverer sendemønstre med ulige centerfrekvens, defineres f_c som det overordnede område af centerfrekvenser for de pågældende sendemønstre |
| Dim. of A_{aptr} | cm | Mål for aktiv blænde for azimuth- og elevationsplan |
| PD | µs | Pulsvarighed, der er knyttet til det sendemønster, der giver værdien, som er angivet for MI |
| PRF | Hz | Pulsvarighedsfrekvensen, der er knyttet til det sendemønster, der giver værdien, som er angivet for MI |
| P_r @ PII_{max} | MPa | Maksimalt fortyndet tryk ved det punkt, hvor det frie, rumlige maksimum for pulsintensitetsintegralet er størst |

| Symboler | Enhed | Definition |
|-----------------------|-------------------|---|
| $d_{eq} @ PII_{max}$ | cm | Tilsvarende strålediameter ved det punkt, hvor det frie, rumlige maksimum for pulsintensitetsintegralet er størst |
| FL | cm | Fokuslængde eller azimuth- og elevationslængde, hvis de er forskellige |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | Belastningsreduceret gennemsnitspulsintensitet ved det maksimalt rapporterede punkt MI |
| ROI | n/a | Interesseområde |
| TB | n/a | Trackball |
| CF | n/a | Farve-Flow-Mode |
| CM | n/a | Farve-M-Mode |
| PW/CW | n/a | Pulsed Wave/Continuous Wave Doppler |

French

| Symboles | Unité | Définition |
|--|--------------------|---|
| MI | n/d | Indice mécanique |
| TIS_{scan} | n/d | Indice thermique pour les tissus mous en mode d'auto-examen |
| TIS_{non-scan} | n/d | Indice thermique pour les tissus mous en mode de non-auto-examen |
| TIB | n/d | Indice thermique pour les os |
| TIC | n/d | Indice thermique crânien |
| A_{aprt} | cm ² | Zone d'ouverture active |
| P_{r.3} | MPa | Tensions rares de pic non notées (MPa) associées au schéma de transmission et donnant lieu à la valeur indiquée sous MI |
| W_o | mW | Puissance échographique, sauf pour l'examen ITS_{acq} , auquel cas il s'agit de la puissance échographique passant par une fenêtre d'un cm. |
| W_{.3(z₁)} | mW | Puissance échographique non cotée à distance axiale z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Pic spatial non coté, intensité temporelle moyenne à distance axiale z₁ |
| z₁ | cm | Distance axiale correspondant à l'emplacement du max [min(W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], où $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Diamètre de faisceau équivalent comme fonction de la distance axiale z , et égal à $[(4/\pi)(W_o/I_{TA}(z))]^{1/2}$, où I_{TA(z)} est l'intensité moyenne temporelle fonction de z |
| f_c | MHz | Fréquence centrale Pour MI , f_c est la fréquence centrale associée au schéma de transmission qui donne lieu à la valeur rapportée maximale de MI . Pour TI , pour des modes combinés impliquant des schémas de transmission de fréquence centrale inégale, f_c est défini comme la gamme totale des fréquences centrales des schémas de transmission respectifs |
| Dim. of A_{aprt} | cm | Dimensions d'ouverture active pour les plans azimutaux et d'élévation |
| PD | μs | Durée de pulsation associée au schéma de transmission donnant lieu à la valeur rapportée de MI |
| PRF | Hz | Fréquence de répétition associée au schéma de transmission donnant lieu à la valeur rapportée de MI |
| P_r @ PII_{max} | MPa | Tension de pic rare au point maximal de champ libre, d'intégrale d'intensité de pic spatial |

| Symboles | Unité | Définition |
|-----------------------|-------------------|---|
| $d_{eq} @ PII_{max}$ | cm | Diamètre de faisceau équivalent au point maximal de champ libre, d'intégrale d'intensité de pic spatial |
| FL | cm | Longueur focale ou longueurs azimutales et d'élévation, si elles sont différentes |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | Intensité moyenne de pulsation non cotée au point maximum reporté MI |
| ROI | n/d | Région d'intérêt |
| TB | n/d | Trackball |
| CF | n/d | Mode de flux de couleurs |
| CM | n/d | Mode M Couleur |
| PW/CW | n/d | Doppler à ondes pulsées/continues |

German

| Symbole | Einheit | Bedeutung |
|--|---------------------|---|
| MI | nicht zutreffend | Mechanischer Index |
| TIS_{scan} | nicht zutreffend | Soft Tissue Thermal Index im Auto-Scanning-Modus |
| TIS_{non-scan} | nicht zutreffend | Soft Tissue Thermal Index im Nicht-Auto-Scanning-Modus |
| TIB | nicht zutreffend | Bone Thermal Index |
| TIC | nicht zutreffend | Cranial Thermal Index |
| A_{aprt} | cm ² | Fläche der aktiven Apertur |
| P_{r.3} | MPa | Freigesetzter maximaler Verdünnungs-Druck (MPa) bei dem verwendeten Sendemuster, das zu dem unter MI angegebenen Wert führt |
| W₀ | mW | Ultraschallleistung, außer beim TIS_{scan} , bei dem es sich um die Ultraschallleistung durch ein Ein-Zentimeter-Fenster handelt |
| W_{.3(z₁)} | mW | Freigesetzte Ultraschallleistung bei Axialabstand z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Freigesetzter räumlicher Spitzen- und zeitlicher Mittelwert der Intensität im Axialabstand z₁ |
| z₁ | cm | Axialabstand entsprechend der Position von max. [min.(W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], wobei $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Äquivalenter Strahldurchmesser als Funktion des Axialabstands z und gleich $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, wobei I_{TA}(z) die zeitlich gemittelte Intensität als Funktion von z ist. |
| f_c | MHz | Mittenfrequenz Für MI ist f_c die Mittenfrequenz bei dem Sendemuster, das zum Maximalwert von MI führt. Für TI bei kombinierten Betriebsarten mit Sendemustern von ungleicher Mittenfrequenz ist f_c definiert als der Gesamtbereich der Mittenfrequenzen der jeweiligen Sendemuster |
| Dim. of A_{aprt} | cm | Maße der aktiven Apertur für die Azimutal- und die Elevationsebene |
| PD | µs | Pulsdauer des Sendemusters, das zum angegebenen Wert von MI führt |
| PRF | Hz | Pulswiederholungsfrequenz des Sendemusters, das zum angegebenen Wert von MI führt |

| Symbole | Einheit | Bedeutung |
|------------------------|-------------------|---|
| $P_r @ PII_{max.}$ | MPa | Maximaler Verdünnungs-Druck an dem Punkt, an dem der räumliche Spitzenwert des Pulsintensitätsintegrals im freien Feld ein Maximum ist |
| $d_{eq} @ PII_{max.}$ | cm | Äquivalenter Strahldurchmesser an dem Punkt, an dem der räumliche Spitzenwert des Pulsintensitätsintegrals im freien Feld ein Maximum ist |
| FL | cm | Fokuslänge bzw. Azimutal- und Elevationslänge, falls unterschiedlich |
| $I_{PA,3} @ MI_{max.}$ | W/cm ² | Reduzierter Pulsmittelwert der Intensität am Punkt des maximalen angegebenen MI |
| ROI | nicht zutreffend | Einstellbare Ausschnittsgröße |
| TB | nicht zutreffend | Trackball |
| CF | nicht zutreffend | Farbfluss-Modus |
| CM | nicht zutreffend | Farb-M-Modus |
| PW/CW | nicht zutreffend | Pulsed-Wave-/Continuous-Wave-Doppler |

Italian

| Simboli | Unità | Definizione |
|--|--------------------|--|
| MI | n/a | Indice Meccanico |
| TIS_{scan} | n/a | Indice termico tessuti molli in modalità di scansione automatica |
| TIS_{non-scan} | n/a | Indice termico tessuti molli in modalità di scansione non automatica |
| TIB | n/a | Indice termico delle ossa |
| TIC | n/a | Indice termico cranico |
| A_{aprt} | cm ² | Area dell'apertura attiva |
| P_{r.3} | MPa | Pressione di rarefazione di picco a prestazioni ridotte (MPa) associata allo schema di trasmissione che genera il valore riportato alla voce MI |
| W_o | mW | Potenza ultrasuoni, tranne per TIS_{scansione} nel qual caso corrisponde alla potenza degli ultrasuoni che passa attraverso una finestra di un centimetro |
| W_{.3(z₁)} | mW | Potenza ultrasuoni a prestazioni ridotte in corrispondenza della distanza assiale z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Intensità media temporale, picco spaziale a prestazioni ridotte in corrispondenza della distanza assiale z₁ |
| z₁ | cm | Distanza assiale corrispondente alla posizione di max [min(W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], dove $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Diametro raggio equivalente in funzione della distanza assiale z e pari a $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, dove I_{TA}(z) è l'intensità media temporale in funzione di z |
| f_c | MHz | Frequenza della parte centrale Per MI , f_c è la frequenza della parte centrale associata allo schema di trasmissione che genera il valore massimo riportato di MI Per TI , per le modalità che comportano schemi di trasmissione con frequenza della parte centrale ineguale, f_c è definito come la gamma totale delle frequenze della parte centrale dei rispettivi schemi di trasmissione |
| Dim. of A_{aprt} | cm | Le dimensioni dell'apertura attiva per i piani azimutali e verticali |
| PD | μs | Durata degli impulsi associata allo schema di trasmissione che genera il valore riportato di MI |
| PRF | Hz | Frequenza di ripetizione degli impulsi associata allo schema di trasmissione che genera il valore riportato di MI |

| Simboli | Unità | Definizione |
|-----------------------|-------------------|---|
| $P_r @ PII_{max}$ | MPa | Pressione di rarefazione di picco in corrispondenza del punto in cui l'integrale dell'intensità degli impulsi di picco spaziale a campo libero è al massimo |
| $d_{eq} @ PII_{max}$ | cm | Diametro del raggio equivalente in corrispondenza del punto in cui l'integrale dell'intensità degli impulsi di picco spaziale a campo libero è al massimo |
| FL | cm | Lunghezza focale o lunghezze azimutali e verticali, se diverse |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | Intensità della media degli impulsi a prestazioni ridotte in corrispondenza del punto del valore massimo riportato MI |
| ROI | n/a | Regione di interesse |
| TB | n/a | Trackball |
| CF | n/a | Color Flow |
| CM | n/a | Color M-Mode |
| PW/CW | n/a | Doppler PW/CW |

Japanese

| 記 号 | 単 位 | 定 義 |
|--|--------------------|--|
| MI | n/a | 機械的指数。 |
| TIS_{scan} | n/a | 軟組織熱的指数。自動スキャンモードで使用。 |
| TIS_{non-scan} | n/a | 軟組織 熱的指数。自動スキャンモード以外で使用。 |
| TIB | n/a | 軟骨熱的指数。 |
| TIC | n/a | 頭蓋熱的指数。 |
| A_{aprt} | cm ² | アクティブアパーチャーの面積。 |
| P_{r,3} | MPa | 緩和ピーク疎密圧力 (MPa)。計測値 MI を発生する透過パターンで使用。 |
| W₀ | mW | 超音波出力。 TIS_{scan} では使用せず。TIS スキャンでは、1 センチメートルの窓を通過する超音波出力です。 |
| W_{.3(z₁)} | mW | 軸距離が z₁ のときの緩和超音波出力。 |
| I_{TA,3(z₁)} | mW/cm ² | 緩和空間ピーク。軸距離が z₁ のときの時間平均密度。 |
| z₁ | cm | 最大位置に対する軸距離 [min(W_{.3(z)} , I_{TA,3(z)}) × 1 cm ²] ここで、 $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | 軸距離 z に換算した等価ビーム直径。[(4/π)(W ₀ /I _{TA(z)})] ^{1/2} と同等。 ここで、I _{TA(z)} は z に換算した時間平均密度。 |
| f_c | (MHz) | 中心周波数。 MI では、 f_c は、最大計測値 MI を発生する透過パターンで使用する中心周波数。 TI の場合、不等中心周波数が発生する組み合わせモードでは、 f_c は、各透過パターンの中心周波数の全範囲であると定義できます。 |
| Dim. of A_{aprt} | cm | 方位平面と高度平面におけるアクティブアパーチャーの規模。 |
| PD | μ s | パルス持続時間。計測値 MI を発生する透過パターンで使用。 |
| PRF | Hz | パルス繰り返し周波数。計測値 MI を発生する透過パターンで使用。 |
| P_r @ PII_{max} | MPa | フリーフィールド、空間 ピークパルス密度積分が最大のときのピーク希薄圧力。 |
| d_{eq} @ PII_{max} | cm | フリーフィールド、空間 ピークパルス密度積分が最大になる地点での等価ビーム直径。 |
| FL | cm | 焦点距離、方位距離、または高度距離 (異なる場合)。 |
| I_{PA,3} @ MI_{max} | W/cm ² | 最大計測値 MI を発生する地点における緩和パルス平均密度。 |

| 記 号 | 単 位 | 定 義 |
|-------|-----|----------------|
| 関心領域 | n/a | 関心領域 |
| TB | n/a | トラックボール |
| CF | n/a | カラーフローモード |
| CM | n/a | カラー M モード |
| PW/CW | n/a | パルス波 / 連続波ドプラー |

Latin America Spanish

| Símbolos | Unidad | Definición |
|--|--------------------|---|
| MI | n/a | Índice Mecánico |
| TIS_{scan} | n/a | Índice Termal del Tejido Suave en el modo de auto-examinación |
| TIS_{non-scan} | n/a | Índice Termal del Tejido Suave en el modo de no-auto-examinación |
| TIB | n/a | Índice Termal del hueso |
| TIC | n/a | Índice Termal Craneal |
| A_{aprt} | cm ² | Área de la abertura activa |
| P_{r.3} | MPa | Presión rarefaccional máxima desratiza (MPa) asociada con el patrón transmitido aumentando el valor reportado bajo MI |
| W₀ | mW | Potencia ultrasónica, con excepción para la Examinación TIS en la cual el caso es que la potencia ultrasónica pasando a través de una ventana de centímetro |
| W_{.3(z₁)} | mW | Potencia ultrasónica desratizada a una distancia axial z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Pico espacial desratizado, intensidad del promedio-temporal en la intensidad axial z₁ |
| z₁ | cm | Distancia axial correspondiente a la ubicación de máx [mín(W _{.3(z)} , I _{TA.3(z)} x 1 cm ²)], donde $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Diámetro del haz equivalente como una función de distancia axial z , y es igual a $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, donde I _{TA} (z) es la intensidad del promedio temporal como una función de z |
| fc | MHz | Centro de Frecuencia Para MI , fc es el centro de frecuencia asociado con el patrón de transmisión aumentando al máximo el valor reportado de MI . para TI , los modos combinados incluyendo los patrones transmitidos del centro de la frecuencia desigual, fc es definido como el rango general del centro de frecuencias de los patrones respectivos transmitidos |
| Dim. of A_{aprt} | cm | Dimensiones de abertura activa para los planos "azimuthal" y elevacionales |
| PD | μs | Duración del Pulso asociado con el patrón transmitido aumentando el valor reportado de MI |
| PRF | Hz | Frecuencia de reproducción del pulso asociado con el patrón transmitido aumentado en el valor reportado de MI |
| P_r @ PII_{max} | MPa | Presión rarefaccional máxima al punto del campo libre, intensidad integral pico del pulso espacial es un máximo |

| Símbolos | Unidad | Definición |
|-----------------------|-------------------|---|
| $d_{eq} @ PII_{max}$ | cm | El diámetro del haz equivalente al punto donde el campo libre, el pico espacial, intensidad integral del pulso es un máximo |
| FL | cm | Longitud focal o longitudes "azimutal" y elevacional, si es diferente |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | La intensidad del promedio del pulso desratizado al punto máximo reportado de MI |
| ROI | n/a | Región de Interés |
| TB | n/a | "Trackball" |
| CF | n/a | Modo del Flujo de Color |
| CM | n/a | Modo de Color |
| PW/CW | n/a | Onda Pulsada/Onda Doppler Continua |

Brazilian Portuguese

| Símbolos | Unidade | Definição |
|--|--------------------|---|
| MI | n/d | Índice mecânico |
| TIS_{scan} | n/d | Índice térmico do tecido mole no modo de varredura automática |
| TIS_{non-scan} | n/d | Índice térmico do tecido mole no modo de varredura não automática |
| TIB | n/d | Índice térmico do osso |
| TIC | n/d | Índice térmico craniano |
| A_{aprt} | cm ² | Área da abertura ativa |
| P_{r.3} | MPa | Pressão de rarefação de pico reduzido (MPa) associada com a elevação fornecida do padrão de transmissão para o valor relatado sob MI |
| W_o | mW | A energia ultra-sônica, exceto por TIS_{varr.} no qual a energia ultra-sônica passa por uma janela de um centímetro |
| W_{.3(z₁)} | mW | Energia ultra-sônica reduzida na distância axial z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Pico espacial reduzido, intensidade média temporal na distância axial z₁ |
| z₁ | cm | Distância axial correspondente ao local de máx [W_{.3(z)} , I_{TA.3(z)} x 1 cm ²], onde $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Diâmetro equivalente do feixe como uma função da distância axial z é igual a $[(4/\pi)(W_o/I_{TA}(z))]^{1/2}$, onde I_{TA(z)} é a intensidade média temporal como uma função de z |
| f_c | Mhz | Frequência central Para MI , f_c é a frequência central associada com a elevação fornecida do padrão de transmissão para o valor máximo relatado de MI . Para TI , para os modos combinados envolvendo padrões de transmissão de frequência central desigual, f_c é definido como o intervalo total de frequências centrais dos respectivos padrões de transmissão |
| Dim. of A_{aprt} | cm | Dimensões da abertura ativa para os planos azimutais e de elevação |
| PD | μs | Duração do pulso associado à elevação fornecida do padrão de transmissão para o valor relatado de MI |
| PRF | Hz | Frequência de repetição do pulso associado à elevação fornecida do padrão de transmissão para o valor relatado de MI |
| P_r @ PII_{max} | MPa | Pressão de rarefação do pico no ponto onde o campo livre, o integral de pulso do pico espacial é um máximo |

| Símbolos | Unidade | Definição |
|-----------------------|-------------------|---|
| $d_{eq} @ PII_{max}$ | cm | Diâmetro de feixe equivalente no ponto onde o campo livre, o integral de pulso do pico espacial é um máximo |
| CF | cm | Comprimento focal ou comprimentos de azimuth e elevação, se forem diferentes |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | Intensidade média do pulso reduzida no ponto do MI máximo relatado |
| ROI | n/d | Região de interesse |
| TB | n/d | Trackball |
| CF | n/d | Modo de fluxo colorido |
| CM | n/d | Modo M colorido |
| PW/CW | n/d | Doppler de onda pulsada/onda contínua |

Swedish

| Symboler | Enhet | Definition |
|---|--------------------|---|
| MI | n/a | Mekaniskt index |
| TIS_{scan} | n/a | Termiskt index för mjuk vävnad i automatiskt skanningsmode |
| TIS_{non-scan} | n/a | Termiskt index för mjuk vävnad i icke-automatiskt skanningsmode |
| TIB | n/a | Termiskt index för benvävnad |
| TIC | n/a | Termiskt index för kranialt |
| A_{aprt} | cm ² | Område för aktiv bländare |
| P_{r.3} | MPa | Undervärderat topstryck (MPa) associerat med rörelsemönstret som resulterar i värdet som rapporteras under MI |
| W_o | mW | Ultraljudskraft med undantag för TIS-skanning då i vilket fall ultraljudskraften passerar genom ett en centimeter tjock fönster |
| W_{.3(z₁)} | mW | Undervärderad ultraljudskraft vid axiell distans z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Undervärderad spatial topp, temporal genomsnittsinensitet via axiell diskans z₁ |
| z₁ | cm | Axiell distans korresponderande mot lokaliseringen av max [min(W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], där $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Ekvivalent stråldiameter som en funktion av axiell distans z och är lika med $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, där I_{TA(z)} är den temporal genomsnittsinensiteten som en funktion av z |
| f_c | MHz | Centrumfrekvens För MI , f_c är centrumfrekvensen associerad med överföringsmönstret som ger upphov till det maximala rapportvärdet av MI . För TI , för kombinerade inställningar (mode) som involverar överföringsmönster av olika centrumfrekvens f_c sår definierad som genomsnittsintervallet av centrumfrekvenser av respektive överförelsemönster |
| Dim. of A_{aprt} | cm | Aktiva bländardimensioner för azimuthal- och lutande plan |
| PD | μs | Pulstryck associerat med överförelsemönstret som ger upphov till det rapporterade värdet i MI |
| PRF | Hz | Pulsrepetitionsfrekvens associerat med överförelsemönstret som ger upphov till det rapporterade värdet i MI |
| P_r @ PII_{max} | MPa | Ovanligt topstryck när frifältet, spatiala toppvärdet för pulsintensitetens integral är på max |
| d_{eq} @ PII_{max} | cm | Ekvivalent stråldiameter när frifältet, spatiala toppvärdet för pulsintensitetens integral är på max |

| Symboler | Enhet | Definition |
|---|-------------------|--|
| FL | cm | Fokal längd eller azimuthal- och lutande längder är olika |
| $I_{PA.3} @ MI_{max}$ | W/cm ² | Undervärderad pulsgenomsnittsintensitet vid maximalt rapporterad MI |
| ROI | n/a | Studerat område |
| TB | n/a | Styrkula |
| CF | n/a | Färgflödesläge |
| CM | n/a | Färg-M-mode |
| PW/CW | n/a | Pulsed Wave (PW-)/Continuous Wave (CW)-doppler |

Norwegian

| Symboler | Enhet | Definisjon |
|--|--------------------|--|
| MI | n/a | Mekanisk Indeks |
| TIS_{scan} | n/a | Bløtdel Thermal Index i auto-skanning modus |
| TIS_{non-scan} | n/a | Bløtdel Thermal Index i non-auto-skanning modus |
| TIB | n/a | Bone Thermal Index |
| TIC | n/a | Kraniell Thermal Index |
| A_{aprt} | cm ² | Område for den aktive åpningen |
| P_{r.3} | MPa | Redusert maksimalt trykk (MPa) assosiert med sendemønsteret som gir grunnlag for verdien som angis under MI |
| W_o | mW | Ultralydeffekt, bortsett fra TIS_{scan} hvor det er ultralydeffekten som passerer gjennom et 1 centimeter vindu. |
| W_{.3(z₁)} | mW | Redusert ultralydeffekt i aksial avstand z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Redusert romlig-peak, temporal-gjennomsnitt intensitet ved aksial avstand z₁ |
| z₁ | cm | Aksial distanse svarende til plasseringen av maks [min(W_{.3(z)} , I_{TA.3(z)} × 1 cm ²)], hvor z ≥ z_{bp} |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Ekvivalent strålediameter som en funksjon av aksial distanse z , er lik $[(4/\pi)(W_o/I_{TA}(z))]^{1/2}$, hvor I_{TA(z)} er den temporale-gjennomsnitt intensiteten som en funksjon av z |
| f_c | MHz | Senterfrekvensen for MI , f_c er senterfrekvensen som er forbundet med sendemønsteret som er bakgrunnen for den maksimale rapporterte verdien av MI . For TI , for kombinerte moduser som involverer sendemønstre av ulik senterfrekvens, f_c er definert som det samlede området av senterfrekvenser for de respektive sendemønstrene |
| Dim. of A_{aprt} | cm | Aktive apertur dimensjoner for de azimuthale og opphevede planene |
| PD | μs | Pulsvarighet assosiert med sendemønsteret gir grunnlag for den rapporterte verdien av MI |
| PRF | Hz | Puls repetisjonsfrekvens assosiert med sendemønsteret som gir grunnlag for den rapporterte verdien av MI |
| P_r @ PII_{max} | MPa | Peak trykket ved det punkt hvor, romlig-peak pulshintensitet integralet er ved maksimum |

| | | |
|--|-------------------|--|
| d_{eq} @ PII_{max} | cm | Ekvivalent strålediameter ved det punktet hvor romlig-peak pulsintensitet integralet er ved maksimum |
| FL | cm | Fokal lengde, eller azimutale og høydelengder, er forskjellige |
| I_{PA.3} @ MI_{max} | W/cm ² | Redusert puls gjennomsnitt intensitet ved punktet for maksimum MI |
| ROI | n/a | Interesseområde |
| TB | n/a | Trackball |
| CF | n/a | Fargedoppler modus |
| CM | n/a | Farge M Mode |
| PW/CW | n/a | Pulset/Kontinuerlig Doppler |

Finnish

| Symbolit | Laite | Kuvaus |
|---|--------------------|--|
| MI | e/k | Mekaaninen indeksi |
| TIS_{scan} | e/k | Kudoksen lämpöindeksi automaattisessa skannaustilassa |
| TIS_{non-scan} | e/k | Kudoksen lämpöindeksi manuaalisessa skannaustilassa |
| TIB | e/k | Luun lämpöindeksi |
| TIC | e/k | Kalloluun lämpöindeksi |
| A_{aprt} | cm ² | Aktiivisen apertuurin alue |
| P_{r.3} | MPa | Alennettu huippuvaimenemisen paine (MPa), joka liittyy siirtotapaan ja nostaa kohdassa MI ilmoitettua arvoa. |
| W_o | mW | Ultraääniteho, lukuun ottamatta TIS_{scan} , jolloin se on yhden senttimetrin levyisen ikkunan kautta kulkeva ultraääniteho. |
| W_{.3(z₁)} | mW | Alennettu ultraääniteho aksiaalisella etäisyydellä z₁ |
| I_{TA.3(z₁)} | mW/cm ² | Alennettu spatiaalihiippu, väliaikainen tiheyskeskiarvo aksiaalisella etäisyydellä z₁ |
| z₁ | cm | Maksimin sijaintia vastaava aksiaalinen etäisyys [minimi (W_{.3(z)} , I_{TA.3(z)} x 1 cm ²)], jossa $z \geq z_{bp}$ |
| z_{bp} | | $1.69(A_{aprt})^{1/2}$ |
| d_{eq(z)} | cm | Vastaava säteen halkaisija aksiaalisen etäisyyden z toimintona, joka vastaa $[(4/\pi)(W_0/I_{TA}(z))]^{1/2}$, jossa I_{TA(z)} on z :n toiminnon lämpökeskiarvon tiheys. |
| f_c | MHz | Keskustaajuus Kohdan MI , f_c keskustaajuus liittyy siirtotapaan ja nostaa kohdassa MI ilmoitettua maksimiarvoa. TI yhdistelmätiloille, jotka liittyvät erilaisten keskustaajuuksien siirtokuvioihin, f_c määritetään vastaavien keskuskuvioiden kokonaisalueena. |
| Dim. of A_{aprt} | cm | Aktiivisen apertuurin mitat atsimutaalisille ja kohotetuille tasoille. |
| PD | μs | Siirtokuvioon liittyvä pulssin kesto, joka nostaa kohdassa MI ilmoitettua arvoa. |
| PRF | Hz | Siirtokuvioon liittyvä pulssin toistotaajuus, joka nostaa kohdassa MI ilmoitettua arvoa. |
| P_r @ PII_{max} | MPa | Huippuohentumisen paine pisteessä, jossa vapaa-kenttä, spatiaalisen huippupulssin tiheyden integraali, on maksimiarvossa. |
| d_{eq} @ PII_{max} | cm | Vastaava säteen halkaisija pisteessä, jossa vapaa-kenttä, spatiaalisen huippupulssin tiheyden integraali, on maksimiarvossa. |

| | | |
|---|-------------------|---|
| FL | cm | Tarkennuspituus tai atsimutaalinen ja kohotettu pituus (jos arvot eroavat). |
| I_{PA.3} @ M_I_{max} | W/cm ² | Alennetun pulssikeskiarvon tiheys maksimipisteessä, joka ilmoitetaan kohdassa MI |
| ROI | e/k | Kiinnostusalueet |
| TB | e/k | Ohjauspallo |
| CF | e/k | Värvirtaustila |
| CM | e/k | Väriellinen M-tila |
| PW/CW | e/k | Pulssi-/jatkuva doppler |

Acoustic Output Reporting Tables for Track 3/IEC 60601-2-37

Not all probes listed may be supported worldwide. Please refer to your local language User Manual for an overview of the probes that are supported in your country.

Explanation of Footnotes

The mechanical and thermal indices may be replaced by one of the following footnotes because of the reasons listed:

- p. Display of this index is not required for this operating mode.
- q. This probe is not intended for transcranial or neonatal cephalic uses.
- r. This formulation for TIS is less than that for an alternate formulation in this mode.

If so, the table entries are replaced by a “#”, meaning: no data are provided for this operating condition since the maximum reported value is not reported for the reason listed.

If neither an index or a footnote is given, this means that the index is irrelevant for this transducer/mode combination.

Multiple focal-zones

When using multiple focal-zones on Vivid q N, the time in one frame is divided between the different focal-zones. When measuring this, the MI is found as the maximum MI of all zones:

$$MI = \max_{\text{all zones}} (MI)$$

while the TI and W_0 is found as the time-weighted sum of all zones:

$$TI = \sum_{\text{all zones}} TI_{\text{zone}} \cdot t_{\text{zone}}$$
$$W_0 = \sum_{\text{all zones}} W_{0\text{zone}} \cdot t_{\text{zone}}$$

t_{zone} is the time fraction used per zone in a frame.

Some of the parameters in the acoustic output report tables will have one value per zone. In this case, the range of the parameter values is reported. The number of zones and which zone has the greater MI is also given in the tables.

Operating Conditions

All table entries are with the operating conditions specified at the end of the table.

Transducer Model: 3Sc-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|--|---|---|---|-------|--------|----------|---|--------------------------|--------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.51 | 0.65 | - | - | - | 1.93 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.13 | | | | | |
| | P | | W_0 (mW) | | 137.84 | - | | - | 137.84 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.60 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 2.15 | 1.90 | - | - | - | 1.90 |
| | Dim of A_{aprt} | X (cm) | | | 1.92 | - | - | - | 1.92 |
| | | Y (cm) | | | 1.30 | - | - | - | 1.30 |
| Other Info | t_d | | PD (μsec) | 0.68 | | | | | |
| | p_{rr} | | PRF (Hz) | 24 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 2.40 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.74 | - | - | | 0.74 |
| | | FL_y (cm) | | | 0.48 | - | - | | 0.48 |
| $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 162 | | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 2 | 3 | - | - | - | 3 |
| | Frequency | | (MHz) | 2.42 | 2.00 | - | - | - | 2.00 |
| | Width | | (deg or ratio to max width) | 75.00 | 60.00 | - | - | - | 60.00 |
| | Depth | | (mm) | 300 | 300 | - | - | - | 300 |
| | Focus | | (mm) | 18 | 118 | - | - | - | 118 |

Transducer Model: 3Sc-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|------|-----------------------|-----------------------|-----------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | A _{aprt} ≤ 1 | A _{aprt} > 1 | | |
| Maximum Index Value | | | | 1.49 | - | - | 0.37 | 1.18 | 0.82 |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | | |
| | P _{r,α} | P _{r,3} | (Mpa) | 1.95 | | | | | |
| | P | W ₀ | (mW) | | - | - | | 51.14 | 58.41 |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | Min of [W _{.3} (Z ₁), I _{TA,3} (Z ₁)] | (mW) | | | | 24.26 | | |
| | Z _s | Z ₁ | (cm) | | | | 6.70 | | |
| | Z _{bp} | Z _{bp} | (cm) | | | | 2.58 | | |
| | Z _b | Z _{sp} | (cm) | | | | | 6.20 | |
| | z at max I _{pi α} | Z _{sp} | (cm) | 5.80 | | | | | |
| | d _{eq} (Z _b) | d _{eq} (Z _{sp}) | (cm) | | | | | 0.43 | |
| | f _{awf} | f _c | (MHz) | 1.90 | - | - | 1.90 | 1.90 | 1.90 |
| | Dim of A _{aprt} | X | (cm) | | - | - | 1.92 | 1.92 | 1.92 |
| | | Y | (cm) | | - | - | 1.30 | 1.30 | 1.30 |
| Other Info | t _d | PD | (μsec) | 0.90 | | | | | |
| | p _{rr} | PRF | (Hz) | 1000 | | | | | |
| | p _r at max I _{pi} | P _r @ P _{II} _{max} | (Mpa) | 2.86 | | | | | |
| | d _{eq} at max I _{pi} | d _{eq} @ P _{II} _{max} | (cm) | | | | | 0.43 | |
| | Focal Length | FL _x | (cm) | | - | - | 0.72 | | 0.72 |
| | | FL _y | (cm) | | - | - | 0.48 | | 0.48 |
| | I _{pi α} at max MI | I _{PA,3} @ MI _{max} | (W/cm ²) | 164 | | | | | |
| Operator Control | Power | (dB) | 0 | - | - | 0 | 0 | 0 | |
| | Beam angle | (deg) | 0 | - | - | 0 | 0 | 0 | |
| | Frequency | (MHz) | 2.00 | - | - | 2.00 | 2.00 | 2.00 | |
| | Depth | (mm) | 300 | - | - | 300 | 300 | 300 | |
| | Focus | (mm) | 93 | - | - | 118 | 93 | 118 | |

Transducer Model: 3Sc-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|----------------------|--|------|------|----------------------|-----------------------|-----------------|--------|
| | | | | | scan | non-scan | | | |
| | | | | | | A _{aprt} ≤1 | A _{aprt} > 1 | | |
| Maximum Index Value | | | | 1.48 | - | - | 1.43 | 2.57 | 2.43 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | P _{r,α} | | P _{r,3} (Mpa) | 1.72 | | | | | |
| | P | | W ₀ (mW) | | - | - | | 103.09 | 135.51 |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA.3} (Z ₁)] (mW) | | | | 25.02 | | |
| | Z _s | | Z ₁ (cm) | | | | 6.80 | | |
| | Z _{bp} | | Z _{bp} (cm) | | | | 2.34 | | |
| | Z _b | | Z _{sp} (cm) | | | | | 6.20 | |
| | z at max I _{pi α} | | Z _{sp} (cm) | 5.60 | | | | | |
| | d _{eq} (Z _b) | | d _{eq} (Z _{sp}) (cm) | | | | | 0.41 | |
| | f _{awf} | | f _c (MHz) | 2.45 | - | - | 3.60 | 2.25 | 3.60 |
| | Dim of A _{aprt} | X (cm) | | | - | - | 1.59 | 1.59 | 1.59 |
| | | Y (cm) | | | - | - | 1.30 | 1.30 | 1.30 |
| Other Info | t _d | | PD (μsec) | 1.03 | | | | | |
| | prf | | PRF (Hz) | 247 | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} (Mpa) | 2.76 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} (cm) | | | | | 0.41 | |
| | Focal Length | FL _x (cm) | | | - | - | 0.56 | | 0.56 |
| | | FL _y (cm) | | | - | - | 0.28 | | 0.28 |
| | I _{pi α} at max MI | | I _{PA.3} @ MI _{max} (W/cm ²) | 128 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | - | 0 | 0 | 0 |
| | PRF | | (Hz) | 247 | - | - | 5618 | 5618 | 5618 |
| | ROI span | | (mm) | 155 | - | - | 10 | 10 | 10 |
| | ROI center | | (mm) | 30 | - | - | 110 | 110 | 110 |
| | Sample Volume | | (mm) | 0.86 | - | - | 1.45 | 1.52 | 1.45 |
| | Frequency | | (MHz) | 2.50 | - | - | 3.64 | 2.22 | 3.64 |

Transducer Model: 3Sc-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|--|--------|---|----------------------|-------|----------|---|-------------------|-------|----------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | | $A_{aprt} > 1$ |
| Maximum Index Value | | | | 1.50 | 1.89 | - | - | - | 3.96 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 2.10 | | | | | |
| | P | | W_0 | (mW) | | 187.83 | - | | - | 187.83 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{TA.3}(Z_1)]$ | (mW) | | | | - | | |
| | z_s | | z_1 | (cm) | | | | - | | |
| | z_{bp} | | z_{bp} | (cm) | | | | - | | |
| | z_b | | z_{sp} | (cm) | | | | | - | |
| | z at max $I_{pi\ \alpha}$ | | z_{sp} | (cm) | 1.60 | | | | | |
| | $d_{eq}(z_b)$ | | $d_{eq}(Z_{sp})$ | (cm) | | | | | - | |
| | f_{awf} | | f_c | (MHz) | 2.10 | 2.48 | - | - | - | 2.48 |
| | Dim of A_{aprt} | X | | (cm) | | 1.38 | - | - | - | 1.38 |
| | | Y | | (cm) | | 1.30 | - | - | - | 1.30 |
| Other Info | t_d | | PD | (μsec) | 0.67 | | | | | |
| | prr | | PRF | (Hz) | 34 | | | | | |
| | p_r at max I_{pi} | | $P_r@ P_{II_{max}}$ | (Mpa) | 2.36 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq}@$ $P_{II_{max}}$ | (cm) | | | | | - | |
| | Focal Length | FL_x | | (cm) | | 0.20 | - | - | | 0.20 |
| | | FL_y | | (cm) | | 0.64 | - | - | | 0.64 |
| | $I_{pi\ \alpha}$ at max MI | | $I_{PA.3}@$ MI_{max} | (W/cm ²) | 163 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 | |
| | Framerate | | (index) | 0 | 0 | - | - | - | 0 | |
| | PRF | | (Hz) | 750 | 750 | - | - | - | 750 | |
| | ROI Span | | (mm) | 20 | 10 | - | - | - | 10 | |
| | ROI Center | | (mm) | 70 | 30 | - | - | - | 30 | |
| | Sample Volume | | (mm) | 1.00 | 1.50 | - | - | - | 1.50 | |
| | ROI Width | | (deg or ratio to max width) | 15.00 | 90.00 | - | - | - | 90.00 | |
| | Frequency | | (MHz) | 3.64 | 2.50 | - | - | - | 2.50 | |

Transducer Model: 3Sc-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|--|---|---|---|----------------------|------|----------|--------|--------------------------|--------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.37 | - | 1.84 | - | 2.52 | 2.49 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 1.79 | | | | | |
| | P | | W_0 | (mW) | | - | 106.43 | | 106.43 | |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | - | | | |
| | Z_s | | Z_1 | (cm) | | | - | | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | - | | | |
| | Z_b | | Z_{sp} | (cm) | | | | 2.60 | | |
| | z at max $I_{\text{pi } \alpha}$ | | Z_{sp} | (cm) | 5.50 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | 0.50 | | |
| | f_{awf} | | f_c | (MHz) | 1.90 | - | 3.63 | - | 3.64 | 3.63 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.69 | - | 0.69 | 0.69 |
| | | Y | | (cm) | | - | 1.30 | - | 1.30 | 1.30 |
| Other Info | t_d | | PD | (μsec) | 0.96 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1046 | | | | | |
| | p_r at max I_{pi} | | $P_r@$ PII_{max} | (Mpa) | 2.57 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}}@$ PII_{max} | (cm) | | | | 0.50 | | |
| | Focal Length | FL_x | | (cm) | | - | 0.28 | - | | 0.28 |
| | | FL_y | | (cm) | | - | 0.92 | - | | 0.92 |
| $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3}@$ MI_{max} | | (W/cm ²) | 155 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Sample vol. position | | (mm) | 51 | - | 20 | - | 20 | 20 | |
| | Sample Volume | | (mm) | 1.04 | - | 6.02 | - | 6.02 | 6.02 | |
| | Scale | | (m/s) | 0.40 | - | 2.91 | - | 2.91 | 2.91 | |
| | Frequency | | (MHz) | 1.82 | - | 3.64 | - | 3.64 | 3.64 | |

Transducer Model: 3Sc-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|-----|------|----------------------------------|----------------------------|------------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | (a) | - | $A_{\text{aprt}} \leq 1$ 1.47 | $A_{\text{aprt}} > 1$ - | 4.70 3.89 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 139.34 | | 166.13 166.13 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},.3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 3.50 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.51 |
| | f_{awf} | f_c (MHz) | # | - | 2.22 | - | 1.82 1.82 |
| | Dim of A_{aprt} | X (cm) | | - | 0.69 | - | 0.69 0.69 |
| | | Y (cm) | | - | 1.30 | - | 1.30 1.30 |
| Other Info | t_d | PD (μsec) | # | | | | |
| | p_{rr} | PRF (Hz) | # | | | | |
| | p_r at max $I_{\text{pi},\alpha}$ | $P_r @ P_{\text{II,max}}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | | 0.51 |
| | Focal Length | FL_x (cm) | | - | 0.46 | - | 0.52 |
| | | FL_y (cm) | | - | 0.64 | - | 0.68 |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},.3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. Position | (mm) | - | - | 93 | - | 72 72 |
| | Frequency | (MHz) | - | - | 2.22 | - | 1.82 1.82 |

Note:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this **operating condition** since the **global maximum** index value is not reported for the reason listed.

Transducer Model: M4S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|--|---|---|---|-------|--------|----------|---|--------------------------|--------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.43 | 0.85 | - | - | - | 2.07 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.16 | | | | | |
| | P | | W_0 (mW) | | 146.28 | - | | - | 146.28 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 5.68 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 2.65 | 2.25 | - | - | - | 2.25 |
| | Dim of A_{aprt} | X (cm) | | | 1.89 | - | - | - | 1.89 |
| | | Y (cm) | | | 1.30 | - | - | - | 1.30 |
| Other Info | t_d | | PD (ms) | 0.57 | | | | | |
| | p_{rr} | | PRF (Hz) | 31 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 3.63 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 1.60 | - | - | | 1.60 |
| | | FL_y (cm) | | | 0.54 | - | - | | 0.54 |
| $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} (W/cm ²) | 240 | | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 1 | 1 | - | - | - | 3 |
| | Frequency | | (MHz) | 3.33 | 1.60 | - | - | - | 2.50 |
| | Width | | (deg or ratio to max width) | 65.00 | 10.00 | - | - | - | 75.00 |
| | Depth | | (mm) | 300 | 300 | - | - | - | 300 |
| | Focus | | (mm) | 70 | 150 | - | - | - | 300 |

Transducer Model: M4S-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|---|------|------|--------------------------|-----------------------|-----------------|-------|------|
| | | | | | scan | non-scan | | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | | |
| Maximum Index Value | | | | 1.45 | - | - | 0.40 | 1.47 | 0.86 | |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.36 | | | | | | |
| | P | | W_0 (mW) | | - | - | | 61.87 | 65.81 | |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | 33.39 | | | |
| | Z_s | | Z_1 (cm) | | | | 4.80 | | | |
| | z_{bp} | | z_{bp} (cm) | | | | 2.79 | | | |
| | z_{b} | | z_{sp} (cm) | | | | | 6.80 | | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 6.00 | | | | | | |
| | $d_{\text{eq}}(z_{\text{b}})$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.44 | | |
| | f_{awf} | | f_{c} (MHz) | 2.80 | - | - | 2.05 | 1.65 | 2.05 | |
| | Dim of A_{aprt} | X | | (cm) | | - | - | 2.20 | 2.20 | 2.20 |
| | | Y | | (cm) | | - | - | 1.30 | 1.30 | 1.30 |
| Other Info | t_{d} | | PD (ms) | 0.56 | | | | | | |
| | p_{rr} | | PRF (Hz) | 1000 | | | | | | |
| | p_{r} at max I_{pi} | | $P_{\text{r}}@ P_{\text{II}_{\text{max}}}$ (MPa) | 4.21 | | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}}@ P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.44 | | |
| | Focal Length | FL _x | | (cm) | | - | - | 1.60 | | 1.60 |
| | | FL _y | | (cm) | | - | - | 0.40 | | 0.40 |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3}@ MI_{\text{max}}$ (W/cm ²) | 225 | | | | | | |
| | Power | | (dB) | 0 | - | - | 0 | 0 | 0 | |
| | Beam Angle | | (deg) | 0 | - | - | 0 | 0 | 0 | |
| | Frequency | | (MHz) | 3.64 | - | - | 2.50 | 1.60 | 2.50 | |
| | Depth | | (mm) | 300 | - | - | 300 | 300 | 300 | |
| | Focus | | (mm) | 70 | - | - | 250 | 95 | 250 | |

Transducer Model: M4S-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.49 | - | - 0.81 | 2.32 | 1.86 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 1.80 | | | | |
| | P | W_0 (mW) | | - | - | 96.53 | 99.59 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 30.94 | | |
| | Z_s | Z_1 (cm) | | | 6.60 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 2.39 | | |
| | Z_b | Z_{sp} (cm) | | | | 6.60 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 5.69 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.36 | |
| | f_{awf} | f_c (MHz) | 2.50 | - | - 2.50 | 2.50 | 2.50 |
| | Dim of A_{aprt} | X (cm) | | - | - 1.61 | 1.61 | 1.61 |
| | | Y (cm) | | - | - 1.30 | 1.30 | 1.30 |
| Other Info | t_d | PD (ms) | 0.62 | | | | |
| | prr | PRF (Hz) | 250 | | | | |
| | p_r at max $I_{pi\alpha}$ | $P_r@ PII_{max}$ (MPa) | 2.94 | | | | |
| | d_{eq} at max $I_{pi\alpha}$ | $d_{eq}@$ PII_{max} (cm) | | | | 0.36 | |
| | Focal Length | FL_x (cm) | | - | - 0.36 | | 0.36 |
| | | FL_y (cm) | | - | - 0.32 | | 0.34 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3}@$ MI_{max} (W/cm ²) | 184 | | | | |
| | Power | (dB) | 0 | - | - 0 | 0 | 0 |
| | PRF | (Hz) | 250 | - | - 5618 | 5618 | 5618 |
| | ROI Span | (mm) | 10 | - | - 10 | 10 | 10 |
| | ROI Center | (mm) | 80 | - | - 110 | 110 | 110 |
| | Sample Volume | (mm) | 0.76 | - | - 1.51 | 1.51 | 1.51 |
| Operator Control | Frequency | (MHz) | 2.50 | - | - 2.50 | 2.50 | 2.50 |

Transducer Model: M4S-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | | |
|------------------------------|---|--------|--|----------------------|-------|----------|---|--------------------------|------|-----------------------|------|
| | | | | | scan | non-scan | | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | | 1.38 | 1.64 | - | - | - | 3.77 | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (MPa) | 0.78 | | | | | | |
| | P | | W_0 | (mW) | | 110.79 | - | | - | 112.56 | |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | | |
| | Z_s | | Z_1 | (cm) | | | | - | | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | | |
| | Z_b | | Z_{sp} | (cm) | | | | | - | | |
| | z at max $I_{\text{pi } \alpha}$ | | Z_{sp} | (cm) | 7.89 | | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - | | |
| | f_{awf} | | f_c | (MHz) | 2.68 | 2.70 | - | - | - | - | 2.65 |
| | Dim of A_{aprt} | X | | (cm) | | 0.96 | - | - | - | - | 0.96 |
| | | Y | | (cm) | | 0.65 | - | - | - | - | 0.65 |
| Other Info | t_d | | PD | (ms) | 1.85 | | | | | | |
| | p_{rr} | | PRF | (Hz) | 250 | | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 1.62 | | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - | | |
| | Focal Length | FL_x | | (cm) | | 0.14 | - | - | | 0.14 | |
| | | FL_y | | (cm) | | 0.92 | - | - | | 0.94 | |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 27 | | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - | 0 | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - | 0 | |
| | Framerate | | (index) | 0 | 0 | - | - | - | - | 0 | |
| | PRF | | (Hz) | 250 | 250 | - | - | - | - | 250 | |
| | ROI Span | | (mm) | 10 | 10 | - | - | - | - | 10 | |
| | ROI Center | | (mm) | 120 | 20 | - | - | - | - | 20 | |
| | Sample Voume | | (mm) | 1.43 | 1.43 | - | - | - | - | 1.43 | |
| | ROI Width | | (deg or ratio to max width) | 15.00 | 15.00 | - | - | - | - | 15.00 | |
| | Frequency | | (MHz) | 2.67 | 2.67 | - | - | - | - | 2.67 | |

Transducer Model: M4S-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|---|----------------------|------|------|----------|---|--------------------------|--------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.46 | - | 1.69 | - | 3.56 | 3.73 |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 2.16 | | | | | |
| | P | W_0 | (mW) | | - | 97.97 | | 126.14 | 133.15 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | | |
| | Z_b | Z_{sp} | (cm) | | | | | 4.49 | |
| | z at max $I_{\text{pi } \alpha}$ | Z_{sp} | (cm) | 4.86 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.37 | |
| | f_{awf} | f_c | (MHz) | 2.40 | - | 3.63 | - | 2.50 | 2.50 |
| | Dim of A_{aprt} | X | (cm) | | - | 0.82 | - | 0.96 | 0.96 |
| | | Y | (cm) | | - | 0.65 | - | 0.65 | 0.65 |
| Other Info | t_d | PD | (ms) | 0.95 | | | | | |
| | p_{rr} | PRF | (Hz) | 1266 | | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 3.24 | | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.37 | |
| | Focal Length | FL_x | (cm) | | - | 0.30 | - | | 0.36 |
| | | FL_y | (cm) | | - | 0.44 | - | | 0.56 |
| | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 170 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 |
| | Beam Angle | | (cm) | 0 | - | 0 | - | 0 | 0 |
| | Sample Volume Position | | (mm) | 51 | - | 31 | - | 31 | 31 |
| | Sample Volume | | (mm) | 1.00 | - | 8.55 | - | 8.58 | 8.58 |
| | Scale | | (m/s) | 0.40 | - | 2.40 | - | 2.40 | 2.40 |
| | Frequency | | (MHz) | 2.50 | - | 3.64 | - | 2.50 | 2.50 |

Transducer Model: M4S-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|-----|------|--|-----------------|-----------------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 1.42 - | 4.05 | 3.79 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | # | | | | |
| | P | W_0 (mW) | | - | 118.94 | | 129.09 135.03 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{r,3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 3.78 |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.45 |
| | f_{awf} | f_c (MHz) | # | - | 2.50 | - | 1.82 1.82 |
| | Dim of A_{aprt} | X (cm) | | - | 0.96 | - | 0.96 0.96 |
| | | Y (cm) | | - | 0.65 | - | 0.65 0.65 |
| Other Info | t_d | PD (ms) | # | | | | |
| | p_{rr} | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (MPa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | | 0.45 |
| | Focal Length | FL_x (cm) | | - | 0.38 | - | 0.50 |
| | | FL_y (cm) | | - | 0.32 | - | 0.36 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | # | | | | |
| | Sample Volume Position (mm) | | - | - | 31 | - | 31 31 |
| | Frequency (MHz) | | - | - | 2.50 | - | 1.82 1.82 |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 5S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.51 | 0.61 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.15 | | | | | |
| | P | | W_0 (mW) | | 70.71 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 4.58 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 2.20 | 2.20 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 1.24 | - | - | - | # |
| | | Y (cm) | | | 1.20 | - | - | - | # |
| Other Info | t_d | | PD (μsec) | 0.73 | | | | | |
| | p_{rr} | | PRF (Hz) | 400 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 3.05 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.40 | - | - | | # |
| | | FL_y (cm) | | | 0.28 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 224 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 3 | 3 | - | - | - | - |
| | Frequency | | (MHz) | 2.00 | 2.00 | - | - | - | - |
| | ROI Width | | (deg or ratio to max width) | 10.00 | 30.00 | - | - | - | - |
| | Depth | | (mm) | 300 | 300 | - | - | - | - |
| | Focus | | (mm) | 90 | 90 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 5S-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.51 | - | - 0.23 | 0.97 | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.62 | | | | |
| | P | W_0 (mW) | | - | - | 18.03 | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 8.50 | | |
| | Z_s | Z_1 (cm) | | | 4.58 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 1.97 | | |
| | Z_b | Z_{sp} (cm) | | | | 4.58 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 4.58 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.34 | |
| | f_{awf} | f_c (MHz) | 2.38 | - | - | 2.38 | # |
| | Dim of A_{aprt} | X (cm) | | - | - | 1.24 | # |
| | | Y (cm) | | - | - | 1.20 | # |
| Other Info | t_d | PD (μ sec) | 0.66 | | | | |
| | prr | PRF (Hz) | 2431 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.35 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | 0.34 | |
| | Focal Length | FL_x (cm) | | - | - | 0.46 | # |
| | | FL_y (cm) | | - | - | 0.28 | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 143 | | | | |
| | Power | (dB) | 0 | - | - | 0 | - |
| | Beam angle | (deg) | 0 | - | - | 0 | - |
| | Frequency | (MHz) | 2.22 | - | - | 2.22 | - |
| | Depth | (mm) | 300 | - | - | 300 | - |
| | Focus | (mm) | 90 | - | - | 90 | - |

Acoustic information

Transducer Model: 5S-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.51 | - | - 0.43 | 1.41 | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.63 | | | | |
| | P | W_0 (mW) | | - | - | 4.28 | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 3.24 | | |
| | Z_s | Z_1 (cm) | | | 4.68 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 1.97 | | |
| | Z_b | Z_{sp} (cm) | | | | 4.68 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 4.48 | | | | |
| | $d_{eq}(z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.39 | |
| | f_{awf} | f_c (MHz) | 2.38 | - | - 2.20 | 2.38 | # |
| | Dim of A_{aprt} | X (cm) | | - | - 1.24 | 1.24 | # |
| | | Y (cm) | | - | - 1.20 | 1.20 | # |
| Other Info | t_d | PD (μ sec) | 0.66 | | | | |
| | prr | PRF (Hz) | 2431 | | | | |
| | p_r at max $I_{pi\alpha}$ | $P_r @ PII_{max}$ (Mpa) | 2.35 | | | | |
| | d_{eq} at max $I_{pi\alpha}$ | $d_{eq} @$ PII_{max} (cm) | | | | 0.39 | |
| | Focal Length | FL_x (cm) | | - | - 0.38 | | # |
| | | FL_y (cm) | | - | - 0.32 | | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 139 | | | | |
| | Power | (dB) | 0 | - | - 0 | 0 | - |
| | PDF | (Hz) | 4505 | - | - 249 | 2825 | - |
| | ROI Span | (mm) | 115 | - | - 10 | 245 | - |
| | ROI Center | (mm) | 60 | - | - 20 | 60 | - |
| | Sample Volume | (mm) | 1.44 | - | - 0.96 | 1.44 | - |
| Operator Control | Frequency | (MHz) | 2.67 | - | - 2.67 | 2.67 | - |

Transducer Model: 5S-RS

Operating Mode: CFM

| Index Label | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--------------------|-------|--|-----------------|-----|
| | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | 1.51 | 0.56 | - | - | (b) |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (Mpa) | 2.24 | | | | |
| | P W_0 (mW) | | 43.44 | - | | # |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s Z_1 (cm) | | | | - | |
| | Z_{bp} Z_{bp} (cm) | | | | - | |
| | Z_b Z_{sp} (cm) | | | | | - |
| | z at max $I_{\text{pi } \alpha}$ Z_{sp} (cm) | 4.58 | | | | |
| | $d_{\text{eq}}(Z_b)$ $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - |
| | f_{awf} f_c (MHz) | 2.20 | 2.50 | - | - | - |
| | Dim of A_{aprt} | X (cm) | 1.24 | - | - | - |
| | | Y (cm) | 1.20 | - | - | - |
| Other Info | t_d PD (μsec) | 0.73 | | | | |
| | p_{rr} PRF (Hz) | 400 | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 3.17 | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | - |
| | Focal Length | FL_x (cm) | 0.64 | - | - | - |
| | | FL_y (cm) | 0.28 | - | - | - |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 214 | | | | |
| | Power (dB) | 0 | 0 | - | - | - |
| | Tilt (deg) | 0 | 0 | - | - | - |
| | Framerate (index) | 1 | 1 | - | - | - |
| | PDF (Hz) | 1500 | 1500 | - | - | - |
| | ROI Span (mm) | 70 | 70 | - | - | - |
| | ROI Center (mm) | 100 | 100 | - | - | - |
| | Sample Volume (mm) | 1.26 | 1.26 | - | - | - |
| | ROI Width (deg or ratio to max width) | 15.00 | 15.00 | - | - | - |
| | Frequency (MHz) | 2.50 | 2.50 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 5S-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.10 | - | 1.16 1.12 | 2.78 | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.36 | | | | |
| | P | W_0 (mW) | | - | 70.96 | 131.70 | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 51.62 | | |
| | Z_s | Z_1 (cm) | | | 4.99 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 1.97 | | |
| | Z_b | Z_{sp} (cm) | | | | 5.09 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 5.09 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.49 | |
| | f_{awf} | f_c (MHz) | 2.65 | - | 2.70 2.63 | 2.65 | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.77 1.24 | 1.24 | # |
| | | Y (cm) | | - | 1.20 1.20 | 1.20 | # |
| Other Info | t_d | PD (μ sec) | 2.49 | | | | |
| | prr | PRF (Hz) | 1374 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.17 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | 0.49 | |
| | Focal Length | FL_x (cm) | | - | 0.32 1.04 | | # |
| | | FL_y (cm) | | - | 0.28 0.26 | | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 81 | | | | |
| | Power | (dB) | 0 | - | 0 0 | 0 | - |
| | Beam angle | (deg) | 0 | - | 0 0 | 0 | - |
| | Sample vol. position | (mm) | 279 | - | 20 279 | 279 | - |
| | Sample volume | (mm) | 2.44 | - | 2.44 2.44 | 2.44 | - |
| | Scale | (m/sec) | 0.40 | - | 0.40 0.40 | 0.40 | - |
| | Frequency | (MHz) | 2.67 | - | 2.67 2.67 | 2.67 | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 5S-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|-----|------|--|-----------------|----------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 1.46 - | 3.69 | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 110.41 | 110.41 | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 3.58 |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.41 |
| | f_{awf} | f_c (MHz) | # | - | 2.50 | - | 2.50 # |
| | Dim of A_{aprt} | X (cm) | | - | 0.54 | - | 0.54 # |
| | | Y (cm) | | - | 1.20 | - | 1.20 # |
| Other Info | t_d | PD (μ sec) | # | | | | |
| | prr | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | | 0.41 |
| | Focal Length | FL_x (cm) | | - | 0.46 | - | # |
| | | FL_y (cm) | | - | 0.44 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. position (mm) | | - | - | 279 | - | 279 - |
| | Frequency (MHz) | | - | - | 2.50 | - | 2.50 - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.47 | 1.17 | - | - | - | 1.85 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.33 | | | | | |
| | P | | W_0 (mW) | | 63.01 | - | | - | 80.37 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 2.89 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 3.00 | 3.95 | - | - | - | 3.00 |
| | Dim of A_{aprt} | X (cm) | | | 1.02 | - | - | - | 1.02 |
| | | Y (cm) | | | 0.90 | - | - | - | 0.90 |
| Other Info | t_d | | PD (ms) | 0.43 | | | | | |
| | p_{rr} | | PRF (Hz) | 240 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 3.17 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.20 | - | - | | 0.50 |
| | | FL_y (cm) | | | 0.28 | - | - | | 0.22 |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} (W/cm ²) | 261 | | | | | |
| | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 3 | 2 | - | - | - | 2 |
| | Frequency | | (MHz) | 3.33 | 3.08 | - | - | - | 3.08 |
| | Width | | (deg or ratio to max width) | 20.00 | 10.00 | - | - | - | 10.00 |
| Depth | | (mm) | 160 | 160 | - | - | - | 160 | |
| Focus | | (mm) | 50 | 80 | - | - | - | 80 | |

Transducer Model: 6S-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|----------------------|------|----------|-------|--------------------------|-------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.28 | - | 0.76 | - | 0.57 | 0.73 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (MPa) | 2.48 | | | | | |
| | P | | W_0 | (mW) | | - | 31.81 | | 31.81 | |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{,3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | z_s | | z_1 | (cm) | | | | - | | |
| | z_{bp} | | z_{bp} | (cm) | | | | - | | |
| | z_b | | z_{sp} | (cm) | | | | | 1.56 | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | 2.56 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(z_{\text{sp}})$ | (cm) | | | | | 0.74 | |
| | f_{awf} | | f_c | (MHz) | 3.95 | - | 5.00 | - | 5.00 | 5.00 |
| | Dim of A_{aprt} | X | | (cm) | | - | 1.02 | - | 1.02 | 1.02 |
| | | Y | | (cm) | | - | 0.90 | - | 0.90 | 0.90 |
| Other Info | t_d | | PD | (ms) | 0.50 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 3.51 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.74 | |
| | Focal Length | FL _x | | (cm) | | - | 0.86 | - | | 0.86 |
| | | FL _y | | (cm) | | - | 0.64 | - | | 0.64 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI _{max} | (W/cm ²) | 272 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam Angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Frequency | | (MHz) | 4.00 | - | 5.00 | - | 5.00 | 5.00 | |
| | Depth | | (mm) | 160 | - | 160 | - | 160 | 160 | |
| | Focus | | (mm) | 25 | - | 100 | - | 100 | 100 | |

Transducer Model: 6S-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.31 | - | 1.05 | - | 1.82 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 1.79 | | | | |
| | P | W_0 (mW) | | - | 51.55 | | 51.55 |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 3.88 |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 3.08 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.32 |
| | f_{awf} | f_c (MHz) | 3.38 | - | 3.10 | - | 3.20 |
| | Dim of A_{aprt} | X (cm) | | - | 1.02 | - | 1.02 |
| | | Y (cm) | | - | 0.90 | - | 0.90 |
| Other Info | t_d | PD (μs) | 1.48 | | | | |
| | prr | PRF (Hz) | 2500 | | | | |
| | p_r at max $I_{pi\alpha}$ | $P_r @ PII_{max}$ (MPa) | 2.56 | | | | |
| | d_{eq} at max $I_{pi\alpha}$ | $d_{eq} @ PII_{max}$ (cm) | | | | | 0.32 |
| | Focal Length | FL_x (cm) | | - | 0.58 | - | 0.58 |
| | | FL_y (cm) | | - | 0.24 | - | 0.24 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @ MI_{max}$ (W/cm ²) | 137 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 |
| | PRF | (Hz) | 2500 | - | 4505 | - | 4505 |
| | ROI Span | (mm) | 60 | - | 10 | - | 10 |
| | ROI Center | (mm) | 20 | - | 110 | - | 110 |
| | Sample Volume | (mm) | 1.10 | - | 0.64 | - | 0.64 |
| Operator Control | Frequency | (MHz) | 3.33 | - | 3.08 | - | 3.08 |

Transducer Model: 6S-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|-------------|--|-------|-------|----------|---|-------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | |
| Maximum Index Value | | | | 1.47 | 1.25 | - | - | - | 1.86 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 1.64 | | | | | |
| | P | | W_0 (mW) | | 73.82 | - | | - | 73.82 |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | | Min of $[W_{,3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} (cm) | | | | - | | |
| | Z_b | | Z_{sp} (cm) | | | | | - | |
| | z at max $I_{pi\alpha}$ | | Z_{sp} (cm) | 3.99 | | | | | |
| | $d_{eq}(Z_b)$ | | $d_{eq}(Z_{sp})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 3.25 | 3.35 | - | - | - | 3.35 |
| | Dim of A_{aprt} | X (cm) | | | 1.02 | - | - | - | 1.02 |
| | | Y (cm) | | | 0.90 | - | - | - | 0.90 |
| Other Info | t_d | | PD (μs) | 0.71 | | | | | |
| | p_{rr} | | PRF (Hz) | 502 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{II_{\max}}$ (MPa) | 2.56 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq} @$ $P_{II_{\max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.54 | - | - | | 0.54 |
| | | FL_y (cm) | | | 0.22 | - | - | | 0.22 |
| | $I_{pi\alpha}$ at max MI | | $I_{PA,3} @$ MI_{\max} (W/cm ²) | 103 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 0 | 2 | - | - | - | 0 |
| | PRF | | (Hz) | 2000 | 3500 | - | - | - | 2000 |
| | ROI Span | | (mm) | 10 | 30 | - | - | - | 10 |
| | ROI Center | | (mm) | 100 | 40 | - | - | - | 100 |
| | Sample Voume | | (mm) | 0.59 | 1.07 | - | - | - | 0.59 |
| | ROI Width | | (deg or ratio to max width) | 20.00 | 15.00 | - | - | - | 20.00 |
| | Frequency | | (MHz) | 3.33 | 3.08 | - | - | - | 3.33 |

Acoustic information

Transducer Model: 6S-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|---|----------------------|------|------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.20 | - | 1.53 | - | 2.23 | 1.93 |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 1.98 | | | | | |
| | P | W_0 | (mW) | | - | 80.89 | | 50.80 | 55.39 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | | |
| | Z_b | Z_{sp} | (cm) | | | | | 2.89 | |
| | z at max $I_{\text{pi } \alpha}$ | Z_{sp} | (cm) | 2.88 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.27 | |
| | f_{awf} | f_c | (MHz) | 3.20 | - | 3.98 | - | 3.20 | 3.93 |
| | Dim of A_{aprt} | X | (cm) | | - | 1.02 | - | 0.51 | 0.45 |
| | | Y | (cm) | | - | 0.90 | - | 0.90 | 0.90 |
| Other Info | t_d | PD | (μs) | 1.08 | | | | | |
| | p_{rr} | PRF | (Hz) | 1684 | | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 2.72 | | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.27 | |
| | Focal Length | FL_x | (cm) | | - | 0.88 | - | | 0.28 |
| | | FL_y | (cm) | | - | 0.38 | - | | 0.34 |
| | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 185 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 |
| | Beam Angle | | (cm) | 0 | - | 0 | - | 0 | 0 |
| | Sample Volume Position | | (mm) | 20 | - | 114 | - | 10 | 10 |
| | Sample Volume | | (mm) | 0.97 | - | 0.99 | - | 0.97 | 0.99 |
| | Scale | | (m/s) | 0.40 | - | 0.40 | - | 0.80 | 0.80 |
| | Frequency | | (MHz) | 3.08 | - | 4.00 | - | 3.08 | 4.00 |

Transducer Model: 6S-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|-----|------|--|-----------------|-------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 0.75 - | 1.99 | 1.65 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | # | | | | |
| | P | W_0 (mW) | | - | 39.45 | 47.36 | 47.36 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | 2.64 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.31 | |
| | f_{awf} | f_c (MHz) | # | - | 4.00 | 3.08 | 3.08 |
| | Dim of A_{aprt} | X (cm) | | - | 0.45 | 0.45 | 0.45 |
| | | Y (cm) | | - | 0.90 | 0.90 | 0.90 |
| Other Info | t_d | PD (μs) | # | | | | |
| | prr | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (MPa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | 0.31 | |
| | Focal Length | FL_x (cm) | | - | 0.30 | 0.34 | 0.34 |
| | | FL_y (cm) | | - | 0.34 | 0.40 | 0.40 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample Volume Position | (mm) | - | - | 103 | 155 | 155 |
| | Frequency | (MHz) | - | - | 4.00 | 3.08 | 3.08 |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 7S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.29 | 1.31 | - | - | - | 1.59 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.36 | | | | | |
| | P | | W_0 (mW) | | 57.79 | - | | - | 57.79 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 2.80 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 3.40 | 4.75 | - | - | - | 4.75 |
| | Dim of A_{aprt} | X (cm) | | | 0.93 | - | - | - | 0.93 |
| | | Y (cm) | | | 0.70 | - | - | - | 0.70 |
| Other Info | t_d | | PD (μsec) | 0.42 | | | | | |
| | p_{rr} | | PRF (Hz) | 35 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II} \text{ max}}$ (Mpa) | 3.28 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II} \text{ max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.66 | - | - | | 0.66 |
| | | FL_y (cm) | | | 0.22 | - | - | | 0.22 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 303 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 1 | 1 | - | - | - | 1 |
| | Frequency | | (MHz) | 3.33 | 4.00 | - | - | - | 5.00 |
| | Width | | (deg or ratio to max width) | 75.00 | 75.00 | - | - | - | 90.00 |
| | Depth | | (mm) | 160 | 160 | - | - | - | 160 |
| | Focus | | (mm) | 50 | 45 | - | - | - | 80 |

Transducer Model: 7S-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|------|------|----------------------------------|----------------------------|----------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | 1.24 | - | $A_{\text{aprt}} \leq 1$ 0.36 | $A_{\text{aprt}} > 1$ - | 0.73 0.44 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 2.19 | | | | |
| | P | W_0 (mW) | | - | 15.83 | | 12.84 15.83 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 2.80 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 2.90 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.19 |
| | f_{awf} | f_c (MHz) | 3.40 | - | 4.80 | - | 3.90 4.80 |
| | Dim of A_{aprt} | X (cm) | | - | 0.93 | - | 0.93 0.93 |
| | | Y (cm) | | - | 0.70 | - | 0.70 0.70 |
| Other Info | t_d | PD (μsec) | 0.42 | | | | |
| | p_{rr} | PRF (Hz) | 1000 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 3.07 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.19 |
| | Focal Length | FL_x (cm) | | - | 0.76 | - | 0.76 |
| | | FL_y (cm) | | - | 0.22 | - | 0.22 |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm^2) | 281 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 0 |
| | Beam angle | (deg) | 0 | - | 0 | - | 0 0 |
| | Frequency | (MHz) | 3.33 | - | 5.00 | - | 4.00 5.00 |
| | Depth | (mm) | 160 | - | 160 | - | 160 160 |
| | Focus | (mm) | 50 | - | 140 | - | 40 140 |

Transducer Model: 7S-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|--------|--|----------------------|------|----------|-------|-----------------|-------|--------------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | | | $A_{\text{aprt}} \leq 1$ |
| Maximum Index Value | | | | 1.24 | - | 1.02 | - | 1.61 | 1.32 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 1.99 | | | | | |
| | P | | W_0 | (mW) | | - | 43.59 | | 30.68 | 43.59 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | |
| | Z_b | | Z_{sp} | (cm) | | | | | 3.40 | |
| | z at max $I_{\text{pi}} \propto$ | | Z_{sp} | (cm) | 3.10 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.21 | |
| | f_{awf} | | f_c | (MHz) | 3.55 | - | 4.43 | - | 3.68 | 4.43 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.93 | - | 0.93 | 0.93 |
| | | Y | | (cm) | | - | 0.70 | - | 0.70 | 0.70 |
| Other Info | t_d | | PD | (μsec) | 0.72 | | | | | |
| | p_{rr} | | PRF | (Hz) | 250 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II,max}}$ | (Mpa) | 2.91 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | | 0.21 | |
| | Focal Length | FL_x | | (cm) | | - | 0.40 | - | | 0.40 |
| | | FL_y | | (cm) | | - | 0.16 | - | | 0.16 |
| | $I_{\text{pi}} \propto$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} | (W/cm ²) | 177 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | PRF | | (Hz) | 250 | - | 4000 | - | 3497 | 4000 | |
| | ROI Span | | (mm) | 10 | - | 300 | - | 60 | 300 | |
| | ROI Center | | (mm) | 80 | - | 60 | - | 60 | 60 | |
| | Sample Volume | | (mm) | 0.56 | - | 1.08 | - | 1.02 | 1.08 | |
| | Frequency | | (MHz) | 3.64 | - | 4.44 | - | 3.64 | 4.44 | |

Transducer Model: 7S-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | | TIB non-scan | TIC |
|------------------------------|--|-----------------|---|----------------------|-------|-----------------------|-----------------------|---|-----------------|-------|
| | | | | | scan | non-scan | | | | |
| | | | | | | A _{aprt} ≤ 1 | A _{aprt} > 1 | | | |
| Maximum Index Value | | | | 1.25 | 1.50 | - | - | - | - | 1.77 |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | P _{r,α} | | P _{r,3} | (Mpa) | 1.36 | | | | | |
| | P | | W ₀ | (mW) | | 59.84 | - | | - | 49.22 |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA,3} (Z ₁)] | (mW) | | | | - | | |
| | Z _s | | Z ₁ | (cm) | | | | - | | |
| | Z _{bp} | | Z _{bp} | (cm) | | | | - | | |
| | Z _b | | Z _{sp} | (cm) | | | | | - | |
| | z at max I _{pi α} | | Z _{sp} | (cm) | 3.20 | | | | | |
| | d _{eq} (Z _b) | | d _{eq} (Z _{sp}) | (cm) | | | | | - | |
| | f _{awf} | | f _c | (MHz) | 3.80 | 4.98 | - | - | - | 4.45 |
| | Dim of A _{aprt} | X | | (cm) | | 0.93 | - | - | - | 0.93 |
| | | Y | | (cm) | | 0.70 | - | - | - | 0.70 |
| Other Info | t _d | | PD | (μsec) | 0.68 | | | | | |
| | p _{rr} | | PRF | (Hz) | 248 | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} | (Mpa) | 2.07 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} | (cm) | | | | | - | |
| | Focal Length | FL _x | | (cm) | | 0.76 | - | - | | 0.78 |
| | | FL _y | | (cm) | | 0.20 | - | - | | 0.28 |
| | I _{pi α} at max MI | | I _{PA,3} @ MI _{max} | (W/cm ²) | 99 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 | |
| | Framerate | | (index) | 0 | 0 | - | - | - | 0 | |
| | PRF | | (Hz) | 250 | 250 | - | - | - | 4000 | |
| | ROI Span | | (mm) | 10 | 10 | - | - | - | 300 | |
| | ROI Center | | (mm) | 40 | 40 | - | - | - | 60 | |
| | Sample Volume | | (mm) | 0.56 | 0.56 | - | - | - | 1.08 | |
| | ROI Width | | (deg or ratio to max width) | 25.00 | 25.00 | - | - | - | 15.00 | |
| | Frequency | | (MHz) | 3.64 | 3.64 | - | - | - | 4.44 | |

Transducer Model: 7S-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|--|--|----------------------------|------|----------|-------|--------------------------|-------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.13 | - | 1.80 | - | 2.85 | 2.34 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 2.07 | | | | | |
| | P | | W_0 | (mW) | | - | 85.07 | | 79.87 | 85.07 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | |
| | Z_b | | Z_{sp} | (cm) | | | | | 2.10 | |
| | z at max $I_{\text{pi } \alpha}$ | | Z_{sp} | (cm) | 2.80 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.33 | |
| | f_{awf} | | f_c | (MHz) | 3.63 | - | 4.45 | - | 4.43 | 4.45 |
| | Dim of A_{aprt} | | X | (cm) | | - | 0.93 | - | 0.93 | 0.93 |
| | | | Y | (cm) | | - | 0.70 | - | 0.70 | 0.70 |
| Other Info | t_d | | PD | (μsec) | 1.07 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1873 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (Mpa) | 2.93 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.33 | |
| | Focal Length | | FL_x | (cm) | | - | 0.78 | - | | 0.78 |
| | | | FL_y | (cm) | | - | 0.28 | - | | 0.28 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm^2) | 217 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Sample vol. position | | (mm) | 31 | - | 83 | - | 83 | 83 | |
| | Sample Volume | | (mm) | 0.98 | - | 0.99 | - | 0.99 | 0.99 | |
| | Scale | | (m/sec) | 0.40 | - | 3.85 | - | 3.85 | 3.85 | |
| | Frequency | | (MHz) | 3.64 | - | 4.44 | - | 4.44 | 4.44 | |

Transducer Model: 7S-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|-----|------|----------------------------------|----------------------------|----------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | (a) | - | $A_{\text{aprt}} \leq 1$ 0.43 | $A_{\text{aprt}} > 1$ - | 1.22 0.93 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 22.40 | | 21.59 22.40 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 2.00 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.23 |
| | f_{awf} | f_c (MHz) | # | - | 4.00 | - | 4.00 4.00 |
| | Dim of A_{aprt} | X (cm) | | - | 0.41 | - | 0.41 0.41 |
| | | Y (cm) | | - | 0.70 | - | 0.70 0.70 |
| Other Info | t_d | PD (μsec) | # | | | | |
| | p_{rr} | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.23 |
| | Focal Length | FL_x (cm) | | - | 0.24 | - | 0.24 |
| | | FL_y (cm) | | - | 0.28 | - | 0.28 |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. Position | (mm) | - | - | 51 | - | 51 51 |
| Control | Frequency | (MHz) | - | - | 4.00 | - | 4.00 4.00 |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this **operating condition** since the **global maximum** index value is not reported for the reason listed.

Transducer Model: 10S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|--|--|--|---|-------|-------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 0.82 | 0.31 | - | - | - | 0.45 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 1.78 | | | | | |
| | P | | W_0 (mW) | | 10.94 | - | | - | 11.25 |
| | Min of [$P_\alpha(Z_s)$, $I_{\text{ta},\alpha}(Z_s)$] | | Min of [$W_{.3}(Z_1)$, $I_{\text{TA},3}(Z_1)$] (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.90 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 5.50 | 5.88 | - | - | - | 5.38 |
| | Dim of A_{aprt} | X (cm) | | | 0.62 | - | - | - | 0.62 |
| | | Y (cm) | | | 0.50 | - | - | - | 0.50 |
| Other Info | t_d | | PD (μsec) | 0.40 | | | | | |
| | p_{rr} | | PRF (Hz) | 44 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 2.56 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.48 | - | - | | 0.44 |
| | | FL_y (cm) | | | 0.14 | - | - | | 0.08 |
| $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 141 | | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 1 | 1 | - | - | - | 3 |
| | Frequency | | (MHz) | 6.15 | 5.00 | - | - | - | 6.15 |
| | Width | | (deg or ratio to max width) | 90.00 | 90.00 | - | - | - | 90.00 |
| | Depth | | (mm) | 120 | 120 | - | - | - | 120 |
| | Focus | | (mm) | 40 | 7 | - | - | - | 100 |

Transducer Model: 10S-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|----------------------|------|----------|------|--------------------------|------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.22 | | 0.05 | - | 0.26 | 0.08 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 2.73 | | | | | |
| | P | | W_0 | (mW) | | - | 1.87 | | 1.97 | 2.01 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | z_{bp} | | z_{bp} | (cm) | | | | - | | |
| | z_b | | z_{sp} | (cm) | | | | | 1.50 | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | 1.40 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.10 | |
| | f_{awf} | | f_c | (MHz) | 5.00 | - | 5.50 | - | 5.00 | 5.00 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.62 | - | 0.62 | 0.62 |
| | | Y | | (cm) | | - | 0.50 | - | 0.50 | 0.50 |
| Other Info | t_d | | PD | (μsec) | 0.32 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (Mpa) | 3.48 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.10 | |
| | Focal Length | FL _x | | (cm) | | - | 0.40 | - | | 0.08 |
| | | FL _y | | (cm) | | - | 0.10 | - | | 0.08 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 430 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Frequency | | (MHz) | 5.00 | - | 6.15 | - | 5.00 | 5.00 | |
| | Depth | | (mm) | 120 | - | 120 | - | 120 | 120 | |
| | Focus | | (mm) | 15 | - | 100 | - | 20 | 20 | |

Transducer Model: 10S-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|--------------------|--|------|------|----------|------|--------------------------|------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.16 | - | 0.20 | - | 0.76 | 0.35 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 1.52 | | | | | |
| | P | | W_0 (mW) | | - | 7.86 | | 7.86 | 7.86 |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} (cm) | | | | - | | |
| | Z_b | | Z_{sp} (cm) | | | | | 2.10 | |
| | z at max $I_{\text{pi} \alpha}$ | | Z_{sp} (cm) | 2.00 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.14 | |
| | f_{awf} | | f_c (MHz) | 4.45 | - | 4.80 | - | 4.80 | 4.80 |
| | Dim of A_{aprt} | X (cm) | | | - | 0.62 | - | 0.62 | 0.62 |
| | | Y (cm) | | | - | 0.50 | - | 0.50 | 0.50 |
| Other Info | t_d | | PD (μsec) | 3.03 | | | | | |
| | p_{rr} | | PRF (Hz) | 5000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II} \text{max}}$ (Mpa) | 2.06 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II} \text{max}}$ (cm) | | | | | 0.14 | |
| | Focal Length | FL_x (cm) | | | - | 0.16 | - | | 0.16 |
| | | FL_y (cm) | | | - | 0.12 | - | | 0.12 |
| | $I_{\text{pi} \alpha}$ at max MI | | $I_{\text{PA}.3} @ MI_{\text{max}}$ (W/cm ²) | 23 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | 0 | 0 | |
| | PRF (Hz) | | 5000 | - | 5000 | - | 5000 | 5000 | |
| | ROI Span (mm) | | 45 | - | 45 | | 45 | 45 | |
| | ROI Center (mm) | | 30 | - | 40 | - | 40 | 40 | |
| | Sample Vol. (mm) | | 0.50 | | 0.35 | - | 0.35 | 0.35 | |
| | Frequency (MHz) | | 4.00 | - | 4.00 | - | 4.00 | 4.00 | |

Transducer Model: 10S-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | | TIB non-scan | TIC |
|------------------------------|---|--|--|----------------------------|-------|--------------------------|-----------------------|---|-----------------|-------|
| | | | | | scan | non-scan | | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | | |
| Maximum Index Value | | | | 0.82 | 0.46 | - | - | - | 0.80 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 0.86 | | | | | |
| | P | | W_0 | (mW) | | 15.93 | - | | - | 15.93 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | z_s | | z_1 | (cm) | | | | - | | |
| | z_{bp} | | z_{bp} | (cm) | | | | - | | |
| | z_b | | z_{sp} | (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | 2.20 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(z_{\text{sp}})$ | (cm) | | | | | - | |
| | f_{awf} | | f_c | (MHz) | 5.68 | 5.05 | - | - | - | 5.05 |
| | Dim of A_{aprt} | | X | (cm) | | 0.62 | - | - | - | 0.62 |
| | | | Y | (cm) | | 0.50 | - | - | - | 0.50 |
| Other Info | t_d | | PD | (μsec) | 0.94 | | | | | |
| | p_{rr} | | PRF | (Hz) | 248 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (Mpa) | 1.33 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - | |
| | Focal Length | | FL_x | (cm) | | 0.28 | - | - | | 0.28 |
| | | | FL_y | (cm) | | 0.12 | - | - | | 0.12 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm^2) | 38 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 | |
| | Framerate | | (index) | 0 | 0 | - | - | - | 0 | |
| | PRF | | (Hz) | 250 | 248 | - | - | - | 248 | |
| | ROI Span | | (mm) | 70 | 25 | - | - | - | 25 | |
| | ROI Center | | (mm) | 35 | 50 | - | - | - | 50 | |
| | Sample vol. | | (mm) | 0.78 | 0.74 | - | - | - | 0.74 | |
| | ROI Width | | (deg or ratio to max width) | 15.00 | 20.00 | - | - | - | 20.00 | |
| | Frequency | | (MHz) | 5.71 | 5.00 | - | - | - | 5.00 | |

Transducer Model: 10S-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|--------|--|----------------------------|------|----------|--------------------------|-----------------------|-------|-------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | | 0.67 | - | 0.47 | - | 0.84 | 0.59 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 1.34 | | | | | |
| | P | | W_0 | (mW) | | - | 14.89 | | 13.70 | 14.89 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | |
| | Z_b | | Z_{sp} | (cm) | | | | | 1.40 | |
| | z at max $I_{\text{pi} \alpha}$ | | Z_{sp} | (cm) | 1.80 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | 0.20 | |
| | f_{awf} | | f_c | (MHz) | 4.08 | - | 6.65 | - | 6.63 | 6.65 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.62 | - | 0.62 | 0.62 |
| | | Y | | (cm) | | - | 0.50 | - | 0.50 | 0.50 |
| Other Info | t_d | | PD | (μsec) | 1.74 | | | | | |
| | p_{rr} | | PRF | (Hz) | 2075 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (Mpa) | 1.72 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | 0.20 | |
| | Focal Length | FL_x | | (cm) | | - | 0.52 | - | | 0.52 |
| | | FL_y | | (cm) | | - | 0.08 | - | | 0.08 |
| Operator Control | $I_{\text{pi} \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm^2) | 74 | | | | | |
| | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Sample vol. position | | (mm) | 20 | - | 155 | - | 155 | 155 | |
| | Sample Volume | | (mm) | 1.49 | - | 1.53 | - | 1.53 | 1.53 | |
| | Scale | | (m/sec) | 0.40 | - | 0.60 | - | 0.60 | 0.60 | |
| Frequency | | (MHz) | 4.00 | - | 6.67 | - | 6.67 | 6.67 | | |

Transducer Model: 10S-RS

Operating Mode: CW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|----------------------|------|----------|------|--------------------------|------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | (a) | - | 0.09 | - | 0.45 | 0.29 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | # | | | | | |
| | P | | W_0 | (mW) | | - | 4.88 | | 4.88 | |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | - | | | |
| | Z_s | | Z_1 | (cm) | | | - | | | |
| | z_{bp} | | z_{bp} | (cm) | | | - | | | |
| | z_b | | z_{sp} | (cm) | | | | 1.20 | | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | # | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | 0.18 | | |
| | f_{awf} | | f_c | (MHz) | # | - | 4.00 | - | 4.00 | 4.00 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.27 | - | 0.27 | 0.27 |
| | | Y | | (cm) | | - | 0.50 | - | 0.50 | 0.50 |
| Other Info | t_d | | PD | (μsec) | # | | | | | |
| | p_{rr} | | PRF | (Hz) | # | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II,max}}$ | (Mpa) | # | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | 0.18 | | |
| | Focal Length | FL _x | | (cm) | | - | 0.20 | - | | 0.20 |
| | | FL _y | | (cm) | | - | 0.14 | - | | 0.14 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI _{max} | (W/cm ²) | # | | | | | |
| Operator Control | Sample vol. position (mm) | | | - | - | 103 | - | 103 | 103 | |
| | Frequency (MHz) | | | - | - | 4.00 | - | 4.00 | 4.00 | |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this **operating condition** since the **global maximum** index value is not reported for the reason listed.

Transducer Model: 12S-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.33 | 0.80 | - | - | - | 1.10 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.71 | | | | | |
| | P | | W_0 (mW) | | 22.48 | - | | - | 28.98 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.84 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 4.88 | 7.25 | - | - | - | 4.88 |
| | Dim of A_{aprt} | X (cm) | | | 0.62 | - | - | - | 0.62 |
| | | Y (cm) | | | 0.55 | - | - | - | 0.55 |
| Other Info | t_d | | PD (μsec) | 0.31 | | | | | |
| | p_{rr} | | PRF (Hz) | 128 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 3.70 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.49 | - | - | | 0.42 |
| | | FL_y (cm) | | | 0.23 | - | - | | 0.13 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 284 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 3 | 2 | - | - | - | 2 |
| | Frequency | | (MHz) | 5.00 | 8.00 | - | - | - | 5.00 |
| | Width | | (deg or ratio to max width) | 20.00 | 10 | - | - | - | 20.00 |
| | Depth | | (mm) | 120 | 120 | - | - | - | 120 |
| | Focus | | (mm) | 44 | 89 | - | - | - | 99 |

Transducer Model: 12S-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|----------------------|------|----------|------|--------------------------|------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.31 | | 0.12 | - | 0.31 | 0.20 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 2.92 | | | | | |
| | P | | W_0 | (mW) | | - | 5.22 | | 5.00 | 5.22 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ | (mW) | | | - | | | |
| | Z_s | | Z_1 | (cm) | | | - | | | |
| | z_{bp} | | z_{bp} | (cm) | | | - | | | |
| | z_b | | z_{sp} | (cm) | | | | 2.05 | | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | 2.05 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | 0.18 | | |
| | f_{awf} | | f_c | (MHz) | 5.00 | - | 4.88 | - | 5.00 | 4.88 |
| | Dim of A_{aprt} | X | | (cm) | | - | 0.62 | - | 0.62 | 0.62 |
| | | Y | | (cm) | | - | 0.55 | - | 0.55 | 0.55 |
| Other Info | t_d | | PD | (μsec) | 0.30 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II,max}}$ | (Mpa) | 4.16 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | 0.18 | | |
| | Focal Length | FL _x | | (cm) | | - | 0.40 | - | | 0.40 |
| | | FL _y | | (cm) | | - | 0.14 | - | | 0.14 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI _{max} | (W/cm ²) | 354 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | 0 | - | 0 | 0 | |
| | Beam angle | | (deg) | 0 | - | 0 | - | 0 | 0 | |
| | Frequency | | (MHz) | 5.00 | - | 5.00 | - | 5.00 | 5.00 | |
| | Depth | | (mm) | 120 | - | 120 | - | 120 | 120 | |
| | Focus | | (mm) | 44 | - | 89 | - | 44 | 89 | |

Transducer Model: 12S-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|--|------|------|----------|------|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.19 | - | 0.58 | - | 1.50 | 0.91 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 1.91 | | | | | |
| | P | | W_0 (mW) | | - | 22.72 | | 18.61 | 22.72 |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{,3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | z_s | | z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | 2.05 | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 2.15 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.15 | |
| | f_{awf} | | f_c (MHz) | 5.05 | - | 5.10 | - | 5.10 | 5.10 |
| | Dim of A_{aprt} | X (cm) | | | - | 0.62 | - | 0.62 | 0.62 |
| | | Y (cm) | | | - | 0.55 | - | 0.55 | 0.55 |
| Other Info | t_d | | PD (μsec) | 0.56 | | | | | |
| | p_{rr} | | PRF (Hz) | 5000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 2.78 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.15 | |
| | Focal Length | FL_x (cm) | | | - | 0.28 | - | | 0.28 |
| | | FL_y (cm) | | | - | 0.12 | - | | 0.12 |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 207 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | 0 | 0 | |
| | PRF (Hz) | | 5000 | - | 3497 | - | 5000 | 3497 | |
| | ROI Span (mm) | | 75 | - | 45 | - | 75 | 45 | |
| | ROI Center (mm) | | 30 | - | 55 | - | 30 | 55 | |
| | Sample Vol. (mm) | | 0.47 | - | 0.47 | - | 0.47 | 0.47 | |
| | Frequency (MHz) | | 5.00 | - | 5.00 | - | 5.00 | 5.00 | |

Transducer Model: 12S-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | | |
|------------------------------|--|--|--|--------|----------------------|----------------------|-----------------------|-----------------|------|------|-------|
| | | | | | scan | non-scan | | | | | |
| | | | | | | A _{aprt} ≤1 | A _{aprt} > 1 | | | | |
| Maximum Index Value | | | | 1.23 | 0.24 | - | - | - | 0.42 | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | P _{r,α} | | P _{r,3} | (Mpa) | 0.63 | | | | | | |
| | P | | W ₀ | (mW) | | 7.59 | - | | 7.59 | | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W ₃ (Z ₁), I _{TA,3} (Z ₁)] | | | | - | | | | |
| | Z _s | | Z ₁ | | | | - | | | | |
| | Z _{bp} | | Z _{bp} | | | | - | | | | |
| | Z _b | | Z _{sp} | | | | | - | | | |
| | z at max I _{pi α} | | Z _{sp} | | 2.15 | | | | | | |
| | d _{eq} (z _b) | | d _{eq} (Z _{sp}) | | | | | - | | | |
| | f _{awf} | | f _c | | 5.10 | 4.30 | - | - | - | 4.30 | |
| | Dim of A _{aprt} | | X | | | 0.62 | - | - | - | 0.62 | |
| | | | Y | | | 0.55 | - | - | - | 0.55 | |
| Other Info | t _d | | PD | (μsec) | 0.84 | | | | | | |
| | prr | | PRF | (Hz) | 248 | | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II_{max}} | (Mpa) | 0.92 | | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II_{max}} | (cm) | | | | - | | | |
| | Focal Length | | FL _x | | | 0.30 | - | - | | 0.30 | |
| | | | FL _y | | | 0.13 | - | - | | 0.13 | |
| | I _{pi α} at max MI | | I _{PA,3} @ MI _{max} | | (W/cm ²) | 18 | | | | | |
| Operator Control | Power | | | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | | | (index) | 0 | 0 | - | - | - | 0 |
| | PRF | | | | (Hz) | 250 | 250 | - | - | - | 280 |
| | ROI Span | | | | (mm) | 40 | 40 | - | - | - | 40 |
| | ROI Center | | | | (mm) | 35 | 25 | - | - | - | 25 |
| | Sample vol. | | | | (mm) | 0.60 | 0.57 | - | - | - | 0.57 |
| | ROI Width | | (deg or ratio to max width) | | | 15.00 | 15.00 | - | - | - | 15.00 |
| | Frequency | | | | (MHz) | 5.00 | 4.00 | - | - | - | 4.00 |

Transducer Model: 12S-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|--|-------------|------|------|----------|--------------------------|-----------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | | 0.91 | - | 1.81 | - | 2.32 | 2.16 |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (Mpa) | | | 1.95 | | | | | |
| | P W_0 (mW) | | | | - | 56.88 | | 53.84 | 56.88 |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | | | - | | |
| | Z_s Z_1 (cm) | | | | | | - | | |
| | Z_{bp} Z_{bp} (cm) | | | | | | - | | |
| | Z_b Z_{sp} (cm) | | | | | | | 1.53 | |
| | z at max $I_{\text{pi} \alpha}$ Z_{sp} (cm) | | | 2.05 | | | | | |
| | $d_{\text{eq}}(z_b)$ $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | | | 0.26 | |
| | f_{awf} f_c (MHz) | | | 4.98 | - | 6.68 | - | 6.68 | 6.68 |
| | Dim of A_{aprt} | | X (cm) | | - | 0.62 | - | 0.62 | 0.62 |
| | | | Y (cm) | | - | 0.55 | - | 0.55 | 0.55 |
| Other Info | t_d PD (μsec) | | | 1.21 | | | | | |
| | p_{rr} PRF (Hz) | | | 2577 | | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | | | 2.77 | | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | | | 0.26 | |
| | Focal Length | | FL_x (cm) | | - | 0.49 | - | | 0.49 |
| | | | FL_y (cm) | | - | 0.29 | - | | 0.29 |
| | $I_{\text{pi} \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | | | 147 | | | | | |
| Operator Control | Power (dB) | | | 0 | - | 0 | - | 0 | 0 |
| | Beam angle (deg) | | | 0 | - | 0 | - | 0 | 0 |
| | Sample vol. position (mm) | | | 238 | - | 124 | - | 124 | 124 |
| | Sample Volume (mm) | | | 1.03 | - | 1.03 | - | 1.03 | 1.03 |
| | Scale (m/sec) | | | 0.40 | - | 0.40 | - | 0.40 | 0.40 |
| | Frequency (MHz) | | | 5.00 | - | 6.67 | - | 6.67 | 6.67 |

Transducer Model: 12S-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|-----|------|----------------------------------|----------------------------|----------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | (a) | - | $A_{\text{aprt}} \leq 1$ 0.61 | $A_{\text{aprt}} > 1$ - | 1.45 1.39 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 22.52 | | 21.85 22.52 |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 1.66 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.18 |
| | f_{awf} | f_c (MHz) | # | - | 5.71 | - | 5.71 5.71 |
| | Dim of A_{aprt} | X (cm) | | - | 0.23 | - | 0.23 0.23 |
| | | Y (cm) | | - | 0.55 | - | 0.55 0.55 |
| Other Info | t_d | PD (μsec) | # | | | | |
| | p_{rr} | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.18 |
| | Focal Length | FL_x (cm) | | - | 0.26 | - | 0.26 |
| | | FL_y (cm) | | - | 0.11 | - | 0.11 |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. position | (mm) | - | - | 114 | - | 31 114 |
| | Frequency | (MHz) | - | - | 5.71 | - | 5.71 5.71 |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this **operating condition** since the **global maximum** index value is not reported for the reason listed.

Transducer Model: e8C-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|-------------|---|--------|--------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.14 | 0.31 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.38 | | | | | |
| | P | | W_0 (mW) | | 12.21 | - | | - | # |
| | Min of [$P_\alpha(Z_s)$, $I_{\text{ta},\alpha}(Z_s)$] | | Min of [$W_{.3}(Z_1)$, $I_{\text{TA}.3}(Z_1)$] (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.79 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 4.75 | 5.25 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 0.58 | - | - | - | # |
| | | Y (cm) | | | 0.50 | - | - | - | # |
| Other Info | t_d | | PD (μsec) | 0.32 | | | | | |
| | p_{rr} | | PRF (Hz) | 13 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II} \text{ max}}$ (Mpa) | 3.20 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II} \text{ max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.32 | - | - | | # |
| | | FL_y (cm) | | | 0.20 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 254 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 1 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 5.00 | 5.00 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 120.00 | 120.00 | - | - | - | - |
| | Depth | | (mm) | 140 | 140 | - | - | - | - |
| | Focus | | (mm) | 28 | 28 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: e8C-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|------|------|----------------------------------|----------------------------|-------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | 0.96 | | $A_{\text{aprt}} \leq 1$ 0.10 | $A_{\text{aprt}} > 1$ - | 0.38 (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 2.17 | | | | |
| | P | W_0 (mW) | | - | 2.57 | | 2.57 # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 1.27 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 1.47 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.20 |
| | f_{awf} | f_c (MHz) | 4.75 | - | 4.13 | - | 4.13 # |
| | Dim of A_{aprt} | X (cm) | | - | 0.41 | - | 0.41 # |
| | | Y (cm) | | - | 0.50 | - | 0.50 # |
| Other Info | t_d | PD (μsec) | 0.33 | | | | |
| | p_{rr} | PRF (Hz) | 1000 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ (Mpa) | 2.76 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | | 0.20 |
| | Focal Length | FL_x (cm) | | - | 0.14 | - | # |
| | | FL_y (cm) | | - | 0.30 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 199 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 - |
| | Beam angle | (deg) | 0 | - | 0 | - | 0 - |
| | Frequency | (MHz) | 5.00 | - | 4.00 | - | 4.00 - |
| | Depth | (mm) | 140 | - | 140 | - | 140 - |
| | Focus | (mm) | 21 | - | 21 | - | 21 - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: e8C-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|--|------|------|----------|------|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 0.93 | - | 0.13 | - | 0.47 | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 0.86 | | | | | |
| | P | | W_0 (mW) | | - | 3.85 | | 3.85 | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | z_s | | z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | 1.76 | |
| | z at max $I_{\text{pi} \alpha}$ | | z_{sp} (cm) | 1.97 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(z_{\text{sp}})$ (cm) | | | | | 0.14 | |
| | f_{awf} | | f_c (MHz) | 5.50 | - | 5.05 | - | 5.05 | # |
| | Dim of A_{aprt} | X (cm) | | | - | 0.75 | - | 0.75 | # |
| | | Y (cm) | | | - | 0.50 | - | 0.50 | # |
| Other Info | t_d | | PD (μsec) | 0.29 | | | | | |
| | p_{rr} | | PRF (Hz) | 250 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 1.25 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | 0.14 | |
| | Focal Length | FL_x (cm) | | | - | 0.08 | - | | # |
| | | FL_y (cm) | | | - | 0.24 | - | | # |
| | $I_{\text{pi} \alpha}$ at max MI | | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 23 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | 0 | - | |
| | PDF (Hz) | | 250 | - | 5000 | - | 5000 | - | |
| | ROI Span (mm) | | 15 | - | 20 | - | 20 | - | |
| | ROI Center (mm) | | 20 | - | 20 | - | 20 | - | |
| | Sample vol. (mm) | | 0.37 | - | 0.88 | - | 0.88 | - | |
| | Frequency (MHz) | | 5.00 | - | 5.00 | - | 5.00 | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: e8C-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | | |
|------------------------------|--|-----------------|---|--------|----------------------|----------|------|-----------------------|-----|-----------------------|---|
| | | | | | scan | non-scan | | | | | |
| | | | | | | | | A _{aprt} ≤ 1 | | A _{aprt} > 1 | |
| Maximum Index Value | | | | 1.14 | 0.35 | - | - | - | (b) | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | P _{r,α} | | P _{r,3} | (Mpa) | 0.75 | | | | | | |
| | P | | W ₀ | (mW) | | 16.11 | - | | - | # | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA.3} (Z ₁)] | | (mW) | | | - | | | |
| | Z _s | | Z ₁ | | (cm) | | | - | | | |
| | Z _{bp} | | Z _{bp} | | (cm) | | | - | | | |
| | Z _b | | Z _{sp} | | (cm) | | | | - | | |
| | z at max I _{pi α} | | Z _{sp} | | (cm) | 1.96 | | | | | |
| | d _{eq} (Z _b) | | d _{eq} (Z _{sp}) | | (cm) | | | | - | | |
| | f _{awf} | | f _c | | (MHz) | 4.15 | 4.95 | - | - | - | # |
| | Dim of A _{aprt} | X | | (cm) | | 1.24 | - | - | - | - | # |
| | | Y | | (cm) | | 0.50 | - | - | - | - | # |
| Other Info | t _d | | PD | (μsec) | 1.02 | | | | | | |
| | prr | | PRF | (Hz) | 385 | | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} | (Mpa) | 0.99 | | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} | | (cm) | | | | - | | |
| | Focal Length | FL _x | | (cm) | | 0.92 | - | - | | # | |
| | | FL _y | | (cm) | | 0.28 | - | - | | # | |
| | I _{pi α} at max MI | | I _{PA.3} @ MI _{max} | | (W/cm ²) | 17 | | | | | |
| Operator Control | Power | | (dB) | | 0 | 0 | - | - | - | - | |
| | Tilt | | (deg) | | 0 | 0 | - | - | - | - | |
| | Framerate | | (index) | | 0 | 0 | - | - | - | - | |
| | PRF | | (Hz) | | 5000 | 5000 | - | - | - | - | |
| | ROI Span | | (mm) | | 20 | 20 | - | - | - | - | |
| | ROI Center | | (mm) | | 20 | 20 | - | - | - | - | |
| | Sample vol. | | (mm) | | 0.89 | 0.89 | - | - | - | - | |
| | ROI Width | | (deg or ratio to max width) | | 30.00 | 30.00 | - | - | - | - | |
| | Frequency | | (MHz) | | 4.00 | 4.00 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: e8C-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|------|------|----------------------------------|----------------------------|-------------|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | 1.19 | - | $A_{\text{aprt}} \leq 1$ 0.38 | $A_{\text{aprt}} > 1$ - | 0.93 (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 2.65 | | | | |
| | P | W_0 (mW) | | - | 13.85 | | 9.99 # |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 1.85 |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 1.85 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | 0.13 |
| | f_{awf} | f_c (MHz) | 5.00 | - | 5.70 | - | 5.00 # |
| | Dim of A_{aprt} | X (cm) | | - | 1.24 | - | 0.87 # |
| | | Y (cm) | | - | 0.50 | - | 0.50 # |
| Other Info | t_d | PD (μsec) | 2.31 | | | | |
| | p_{rr} | PRF (Hz) | 449 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ (Mpa) | 3.65 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | | 0.13 |
| | Focal Length | FL_x (cm) | | - | 0.96 | - | # |
| | | FL_y (cm) | | - | 0.28 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 376 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 - |
| | Beam angle | (deg) | 0 | - | 0 | - | 0 - |
| | Sample vol. position | (mm) | 20 | - | 269 | - | 20 - |
| | Sample volume | (mm) | 2.03 | - | 2.00 | - | 2.03 - |
| | Scale | (m/sec) | 0.07 | - | 0.05 | - | 0.07 - |
| | Frequency | (MHz) | 5.00 | - | 5.71 | - | 5.00 - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 3C-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|------|------|----------|---|--------------------------|-----|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.0 | (a) | - | - | - | (b) | |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 1.8 | | | | | | |
| | P | | W_0 (mW) | | # | - | | - | # | |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | | |
| | Z_s | | Z_1 (cm) | | | | - | | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | | |
| | z_b | | z_{sp} (cm) | 4.2 | | | | - | | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | | | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | | |
| | f_{awf} | | f_c (MHz) | 3.5 | # | - | - | - | # | |
| | Dim of A_{aprt} | X | | (cm) | | # | - | - | - | # |
| | | Y | | (cm) | | # | - | - | - | # |
| Other Info | t_d | | PD (ms) | 0.6 | | | | | | |
| | p_{rr} | | PRF (Hz) | 7252 | | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II,max}}$ (MPa) | 3.0 | | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | | - | | |
| | Focal Length | FL _x | | (cm) | | # | - | - | | # |
| | | FL _y | | (cm) | | # | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 113 | | | | | | |
| Operator Control | | | Power (dB) | 0 | # | - | - | - | # | |
| | | | Focus (cm) | 4.6 | # | - | - | - | # | |
| | | | Range (cm) | 8 | # | - | - | - | # | |
| | | | 2D Angle (deg) | 30.0 | # | - | - | - | # | |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed

Transducer Model: 3C-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|-------------|--|------|------|----------|------|-------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | |
| Maximum Index Value | | | | 1.1 | | - | 0.2 | 0.6 | (a) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.0 | | | | | |
| | P | | W_0 (mW) | | - | - | | 18.8 | # |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ | | Min of $[W_3(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | | 12.3 | | |
| | z_s | | z_1 (cm) | | | | 1.8 | | |
| | z_{bp} | | z_{bp} (cm) | | | | 1.8 | | |
| | z_b | | z_{sp} (cm) | 4.2 | | | | 4.0 | |
| | z at max $I_{pi\alpha}$ | | z_{sp} (cm) | | | | | | |
| | $d_{eq}(z_b)$ | | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.27 | |
| | f_{awf} | | f_c (MHz) | 3.5 | - | - | 3.5 | 3.5 | # |
| | Dim of A_{aprt} | X (cm) | | | - | - | 1.1 | 1.1 | # |
| | | Y (cm) | | | - | - | 1.0 | 1.0 | # |
| Other Info | t_d | | PD (ms) | 0.7 | | | | | |
| | p_{rr} | | PRF (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{II_{max}}$ (MPa) | 2.8 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq} @ P_{II_{max}}$ (cm) | | | | | 2.8 | |
| | Focal Length | FL_x (cm) | | | - | - | 4.7 | | # |
| | | FL_y (cm) | | | - | - | 5.0 | | # |
| | $I_{pi\alpha}$ at max MI | | $I_{PA,3} @ MI_{max}$ (W/cm ²) | 166 | | | | | |
| Operator Control | Power | | (dB) | 0 | - | - | 0 | 0 | # |
| | Focus | | (cm) | 4.6 | | | 4.6 | 4.6 | # |
| | Range | | (cm) | 20.0 | - | - | 20.0 | 20.0 | # |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed

Transducer Model: 3C-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|--|-------------|------|------|--------------------------|-----------------------|-----------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | | 0.9 | - | - | 0.2 | 0.5 | (a) |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (MPa) | | | 1.7 | | | | | |
| | P W_0 (mW) | | | | - | - | | 14.5 | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | | | 90.5 | | |
| | Z_s Z_1 (cm) | | | | | | 1.8 | | |
| | z_{bp} z_{bp} (cm) | | | | | | 1.8 | | |
| | z_b z_{sp} (cm) | | | 4.2 | | | | 3.7 | |
| | z at max $I_{\text{pi } \alpha}$ z_{sp} (cm) | | | | | | | | |
| | $d_{\text{eq}}(z_b)$ $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | | | 0.24 | |
| | f_{awf} f_c (MHz) | | | 3.5 | - | - | 3.1 | 3.1 | # |
| | Dim of A_{aprt} | | X (cm) | | - | - | 1.6 | 1.6 | # |
| | | | Y (cm) | | - | - | 1.0 | 1.0 | # |
| Other Info | t_d PD (ms) | | | 0.7 | | | | | |
| | p_{rr} PRF (Hz) | | | 1000 | | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II,max}}$ (MPa) | | | 2.8 | | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | | | | 0.20 | |
| | Focal Length | | FL_x (cm) | | - | - | 3.3 | | # |
| | | | FL_y (cm) | | - | - | 5.0 | | # |
| | $I_{\text{pi } \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | | | 166 | | | | | |
| Operator Control | Power (dB) | | | 0 | - | - | 0 | 0 | # |
| | Range (cm) | | | 8.0 | - | - | 8.0 | 8.0 | # |
| | ROI depth (cm) | | | 3.3 | - | - | 3.3 | 3.3 | # |
| | Velocity (m/sec) | | | Max. | - | - | Max. | Max. | # |
| | ROI length (mm) | | | 45 | - | - | 45 | 45 | # |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed

Transducer Model: 3C-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | | |
|------------------------------|---|--|--|----------------------|------|----------|-----------------|-----|--------------------------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 0.9 | 1.5 | - | - | - | (a) | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (MPa) | 1.7 | | | | | |
| | P | | W_0 | (mW) | | 98.5 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | |
| | Z_b | | Z_{sp} | (cm) | 3.8 | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | Z_{sp} | (cm) | | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - | |
| | f_{awf} | | f_c | (MHz) | 3.1 | 3.1 | - | - | - | # |
| | Dim of A_{aprt} | | X | (cm) | | 1.6 | - | - | - | # |
| | | | Y | (cm) | | 1.0 | - | - | - | # |
| Other Info | t_d | | PD | (ms) | 1.8 | | | | | |
| | p_{rr} | | PRF | (Hz) | 9615 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 2.6 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - | |
| | Focal Length | | FL_x | (cm) | | 3.2 | 3.3 | - | | # |
| | | | FL_y | (cm) | | 5.0 | - | - | | # |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 126 | | | | | |
| | | | Power | (dB) | 0 | 0 | - | - | - | # |
| | | | Range | (cm) | 8.0 | 8.0 | - | - | - | # |
| | | | ROI start depth | (cm) | 1.0 | 1.0 | - | - | - | # |
| | | | Dopp. PRF | (KHz) | 9.2 | 9.2 | - | - | - | # |
| | | | ROI length | (mm) | 45 | 45 | - | - | - | # |
| | | | ROI width | (deg) | 14 | 14 | - | - | - | # |
| | | | 2D Angle | (deg) | 30 | 30 | - | - | - | # |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed

Transducer Model: 3C-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|---|-----------------|--|----------------------|------|--------------------------|-----------------------|-----------------|-----|---|
| | | | | | scan | non-scan | | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | | |
| Maximum Index Value | | | | 0.8 | - | - | 0.4 | 1.8 | (a) | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (MPa) | 1.4 | | | | | |
| | P | | W_0 | (mW) | | - | - | 42.0 | # | |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | 27.8 | | | |
| | Z_s | | Z_1 | (cm) | | | 1.9 | | | |
| | z_{bp} | | z_{bp} | (cm) | | | 1.9 | | | |
| | z_b | | z_{sp} | (cm) | 3.7 | | | 3.7 | | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | 0.24 | | |
| | f_{awf} | | f_c | (MHz) | 3.1 | - | - | 3.1 | 3.1 | # |
| | Dim of A_{aprt} | X | | (cm) | | - | - | 1.3 | 1.3 | # |
| | | Y | | (cm) | | - | - | 1.0 | 1.0 | # |
| Other Info | t_d | | PD | (ms) | 0.9 | | | | | |
| | p_{rr} | | PRF | (Hz) | 2900 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II,max}}$ | (MPa) | 2.0 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | 0.24 | | |
| | Focal Length | FL _x | | (cm) | | - | - | 3.2 | | # |
| | | FL _y | | (cm) | | - | - | 5.0 | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 97 | | | | | |
| Operator Control | | | Power | (dB) | 0 | - | - | 0 | 0 | # |
| | | | Range | (cm) | 8.0 | - | - | 8.0 | 8.0 | # |
| | | | SV depth | (cm) | 3.0 | - | - | 3.0 | 3.0 | # |
| | | | Velocity | (m/sec) | 0.6 | - | - | 3.4 | 3.4 | # |
| | | | SV | (mm) | 1 | - | - | 2 | 2 | # |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed

Transducer Model: 4C-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | |
|------------------------------|--|--|--|-------|-------|----------|-----------------|-----|----------------------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | | A _{aprt} ≤1 |
| Maximum Index Value | | | | 1.53 | 0.47 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | P _{r,α} | | P _{r,3} (Mpa) | 1.15 | | | | | |
| | P | | W ₀ (mW) | | 29.16 | - | | - | # |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA.3} (Z ₁)] (mW) | | | | - | | |
| | Z _s | | Z ₁ (cm) | | | | - | | |
| | z _{bp} | | z _{bp} (cm) | | | | - | | |
| | z _b | | z _{sp} (cm) | | | | | - | |
| | z at max I _{pi α} | | z _{sp} (cm) | 4.63 | | | | | |
| | d _{eq} (z _b) | | d _{eq} (Z _{sp}) (cm) | | | | | - | |
| | f _{awf} | | f _c (MHz) | 1.80 | 2.95 | - | - | - | # |
| | Dim of A _{aprt} | | X (cm) | | 1.69 | - | - | - | # |
| | | | Y (cm) | | 1.30 | - | - | - | # |
| Other Info | t _d | | PD (μsec) | 1.03 | | | | | |
| | p _{rr} | | PRF (Hz) | 25 | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} (Mpa) | 1.53 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} (cm) | | | | | - | |
| | Focal Length | | FL _x (cm) | | 0.10 | - | - | | # |
| | | | FL _y (cm) | | 0.60 | - | - | | # |
| | I _{pi α} at max MI | | I _{PA.3} @ MI _{max} (W/cm ²) | 79 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 0 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 1.51 | 3.08 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 45.00 | 58.00 | - | - | - | - |
| | Depth | | (mm) | 300 | 300 | - | - | - | - |
| | Focus | | (mm) | 45 | 41 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 4C-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|--|-----------------|-------------|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 1.39 | - | - | 0.30 | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.24 | | | |
| | P | W_0 | (mW) | | - | - | 36.36 # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | 13.48 | |
| | Z_s | Z_1 | (cm) | | | 5.42 | |
| | Z_{bp} | Z_{bp} | (cm) | | | 2.37 | |
| | Z_b | Z_{sp} | (cm) | | | 5.42 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 4.47 | | | |
| | $d_{eq}(z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | 0.42 | |
| | f_{awf} | f_c | (MHz) | 3.10 | - | - | 2.65 2.65 # |
| | Dim of A_{aprt} | X | (cm) | | - | - | 1.64 1.64 # |
| | | Y | (cm) | | - | - | 1.30 1.30 # |
| Other Info | t_d | PD | (μsec) | 0.47 | | | |
| | prr | PRF | (Hz) | 1000 | | | |
| | p_r at max I_{pi} | $P_r @ PII_{\text{max}}$ | (Mpa) | 3.62 | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | 0.42 | |
| | Focal Length | FL_x | (cm) | | - | - | 0.55 # |
| | | FL_y | (cm) | | - | - | 0.30 # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 196 | | | |
| | Power | (dB) | | 0 | - | - | 0 0 - |
| | Beam angle | (deg) | | 0 | - | - | 0 0 - |
| | Frequency | (MHz) | | 3.64 | - | - | 2.76 2.76 - |
| | Depth | (mm) | | 300 | - | - | 300 300 - |
| | Focus | (mm) | | 45 | - | - | 90 90 - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 4C-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.44 | - | - | 1.04 | 2.06 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.11 | | | | |
| | P | W_0 (mW) | | - | - | 141.65 | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 53.78 | | |
| | Z_s | Z_1 (cm) | | | 5.88 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 3.26 | | |
| | Z_b | Z_{sp} (cm) | | | | 5.88 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 5.18 | | | | |
| | $d_{eq}(z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.61 | |
| | f_{awf} | f_c (MHz) | 3.60 | - | - | 2.53 | 2.53 |
| | Dim of A_{aprt} | X (cm) | | - | - | 3.04 | 3.04 |
| | | Y (cm) | | - | - | 1.30 | 1.30 |
| Other Info | t_d | PD (μ sec) | 1.26 | | | | |
| | prr | PRF (Hz) | 5814 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.12 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | 0.61 | |
| | Focal Length | FL_x (cm) | | - | - | 1.35 | # |
| | | FL_y (cm) | | - | - | 0.30 | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 69 | | | | |
| | Power | (dB) | 0 | - | - | 0 | 0 |
| | PRF | (Hz) | 5814 | - | - | 4000 | 4000 |
| | ROI Span | (mm) | 100 | - | - | 105 | 105 |
| | ROI Center | (mm) | 40 | - | - | 110 | 110 |
| | Sample vol. | (mm) | 1.13 | - | - | 1.11 | 1.11 |
| | Frequency | (MHz) | 3.64 | - | - | 2.50 | 2.50 |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 4C-RS

Operating Mode: CFM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|-------|--------|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | |
| | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | 1.56 | 1.16 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.78 | | | | |
| | P | W_0 (mW) | | 197.14 | - | | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | Z_{sp} (cm) | 4.98 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 2.25 | 2.65 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | 2.46 | - | - | # |
| | | Y (cm) | | 1.30 | - | - | # |
| Other Info | t_d | PD (μsec) | 2.12 | | | | |
| | p_{rr} | PRF (Hz) | 251 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 2.62 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | 0.20 | - | - | # |
| | | FL_y (cm) | | 0.30 | - | - | # |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 144 | | | | |
| | Power | (dB) | 0 | 0 | - | - | - |
| | Tilt | (deg) | 0 | 0 | - | - | - |
| | Framerate | (index) | 0 | 0 | - | - | - |
| | PRF | (Hz) | 250 | 500 | - | - | - |
| | ROI Span | (mm) | 100 | 135 | - | - | - |
| | ROI Center | (mm) | 20 | 50 | - | - | - |
| | Sample vol. | (mm) | 1.76 | 1.64 | - | - | - |
| | ROI Width | (deg or ratio to max width) | 20.00 | 15.00 | - | - | - |
| | Frequency | (MHz) | 2.22 | 2.67 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 4C-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.15 | - | - 0.70 | 2.79 | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.88 | | | | |
| | P | W_0 (mW) | | - | - | 68.95 | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | 44.20 | | |
| | Z_s | Z_1 (cm) | | | 2.10 | | |
| | Z_{bp} | Z_{bp} (cm) | | | 1.69 | | |
| | Z_b | Z_{sp} (cm) | | | | 2.10 | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 5.39 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | 0.36 | |
| | f_{awf} | f_c (MHz) | 3.03 | - | - 3.08 | 3.08 | # |
| | Dim of A_{aprt} | X (cm) | | - | - 0.87 | 0.87 | # |
| | | Y (cm) | | - | - 1.30 | 1.30 | # |
| Other Info | t_d | PD (μ sec) | 1.90 | | | | |
| | prr | PRF (Hz) | 779 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 3.29 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | 0.36 | |
| | Focal Length | FL_x (cm) | | - | - 0.10 | | # |
| | | FL_y (cm) | | - | - 0.90 | | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 162 | | | | |
| | Power | (dB) | 0 | - | - 0 | 0 | |
| | Beam angle | (deg) | 0 | - | - 0 | 0 | |
| | Sample vol. position | (mm) | 62 | - | - 20 | 20 | |
| | Sample vol. | (mm) | 1.57 | - | - 4.52 | 4.52 | |
| | Scale | (m/sec) | 0.20 | - | - 1.49 | 1.49 | |
| | Frequency | (MHz) | 3.08 | - | - 3.08 | 3.08 | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 8C-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|--|-----------------|---|----------------------|-------|----------|---|-----------------------|-------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | A _{aprt} ≤ 1 | | A _{aprt} > 1 |
| Maximum Index Value | | | | 1.39 | 0.74 | - | - | - | 1.09 | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | P _{r,α} | | P _{r,3} | (Mpa) | 3.65 | | | | | |
| | P | | W ₀ | (mW) | | 37.62 | - | | 37.62 | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA,3} (Z ₁)] | (mW) | | | - | | | |
| | Z _s | | Z ₁ | (cm) | | | - | | | |
| | z _{bp} | | z _{bp} | (cm) | | | - | | | |
| | z _b | | z _{sp} | (cm) | | | | - | | |
| | z at max I _{pi α} | | z _{sp} | (cm) | 0.73 | | | | | |
| | d _{eq} (z _b) | | d _{eq} (Z _{sp}) | (cm) | | | | - | | |
| | f _{awf} | | f _c | (MHz) | 6.88 | 4.88 | - | - | - | 4.88 |
| | Dim of A _{aprt} | X | | (cm) | | 1.18 | - | - | - | 1.18 |
| | | Y | | (cm) | | 0.50 | - | - | - | 0.50 |
| Other Info | t _d | | PD | (μsec) | 0.32 | | | | | |
| | p _{rr} | | PRF | (Hz) | 51 | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} | (Mpa) | 4.34 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} | (cm) | | | | - | | |
| | Focal Length | FL _x | | (cm) | | 0.94 | - | - | | 0.94 |
| | | FL _y | | (cm) | | 0.30 | - | - | | 0.30 |
| | I _{pi α} at max MI | | I _{PA,3} @ MI _{max} | (W/cm ²) | 334 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 | |
| | Framerate | | (index) | 0 | 0 | - | - | - | 0 | |
| | Frequency | | (MHz) | 8.00 | 8.00 | - | - | - | 5.71 | |
| | Width | | (deg or ratio to max width) | 40.00 | 40.00 | - | - | - | 30.00 | |
| | Depth | | (mm) | 140 | 140 | - | - | - | 140 | |
| | Focus | | (mm) | 7 | 7 | - | - | - | 133 | |

Acoustic information

Transducer Model: 8C-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-------------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 0.98 | - | 0.08 | - | 0.15 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 2.55 | | | | |
| | P | W_0 (mW) | | - | 2.41 | | 2.85 3.23 |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 1.47 |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 0.50 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.21 |
| | f_{awf} | f_c (MHz) | 7.00 | - | 6.63 | - | 4.75 4.75 |
| | Dim of A_{aprt} | X (cm) | | - | 0.46 | - | 0.41 0.46 |
| | | Y (cm) | | - | 0.50 | - | 0.50 0.50 |
| Other Info | t_d | PD (μ sec) | 0.23 | | | | |
| | prr | PRF (Hz) | 1000 | | | | |
| | p_r at max $I_{pi\alpha}$ | $P_r@ PII_{max}$ (Mpa) | 2.88 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq}@ PII_{max}$ (cm) | | | | | 0.21 |
| | Focal Length | FL_x (cm) | | - | 0.32 | - | 0.32 |
| | | FL_y (cm) | | - | 0.28 | - | 0.28 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3}@ MI_{max}$ (W/cm ²) | 201 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 0 |
| | Beam angle | (deg) | 0 | - | 0 | - | 0 0 |
| | Frequency | (MHz) | 8.00 | - | 8.00 | - | 5.00 5.00 |
| | Depth | (mm) | 140 | - | 140 | - | 140 140 |
| | Focus | (mm) | 7 | - | 98 | - | 21 133 |

Transducer Model: 8C-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|--|-----------------|------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 0.80 | - | 0.10 | - | 0.39 |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 0.66 | | | |
| | P | W_0 | (mW) | | - | 4.86 | 4.39 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | - | |
| | Z_s | Z_1 | (cm) | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | 1.84 |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 1.95 | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | | 0.16 |
| | f_{awf} | f_c | (MHz) | 5.70 | - | 4.20 | 4.15 |
| | Dim of A_{aprt} | X | (cm) | | - | 0.71 | 0.71 |
| | | Y | (cm) | | - | 0.50 | 0.50 |
| Other Info | t_d | PD | (μ sec) | 0.89 | | | |
| | prr | PRF | (Hz) | 3497 | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 0.97 | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | | 0.16 |
| | Focal Length | FL_x | (cm) | | - | 0.12 | 0.12 |
| | | FL_y | (cm) | | - | 0.22 | 0.22 |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 19 | | | |
| | Power | (dB) | | 0 | - | 0 | 0 |
| | PRF | (Hz) | | 3497 | - | 5000 | 5000 |
| | ROI Span | (mm) | | 20 | - | 20 | 20 |
| | ROI Center | (mm) | | 20 | - | 20 | 20 |
| | Sample vol. | (mm) | | 0.85 | - | 0.89 | 0.89 |
| | Frequency | (MHz) | | 5.71 | - | 4.00 | 4.00 |

Acoustic information

Transducer Model: 8C-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|--|-------|-------|----------|---|--------------------------|-------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.06 | 0.45 | - | - | - | 0.68 |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 1.61 | | | | | |
| | P | | W_0 (mW) | | 16.71 | - | | - | 16.71 |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{\text{TA.3}}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.95 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 5.10 | 5.10 | - | - | - | 5.10 |
| | Dim of A_{aprt} | | X (cm) | | 0.71 | - | - | - | 0.71 |
| | | | Y (cm) | | 0.50 | - | - | - | 0.50 |
| Other Info | t_d | | PD (μsec) | 0.49 | | | | | |
| | p_{rr} | | PRF (Hz) | 313 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 2.26 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | 0.08 | - | - | | 0.08 | |
| | | FL_y (cm) | | 0.28 | - | - | | 0.28 | |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA.3}} @ MI_{\text{max}}$ (W/cm ²) | 116 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | 0 |
| | Tilt | | (deg) | 0 | 0 | - | - | - | 0 |
| | Framerate | | (index) | 1 | 1 | - | - | - | 1 |
| | PDF | | (Hz) | 5000 | 5000 | - | - | - | 5000 |
| | ROI Span | | (mm) | 35 | 35 | - | - | - | 35 |
| | ROI Center | | (mm) | 20 | 20 | - | - | - | 20 |
| | Sample Vol. | | (mm) | 0.50 | 0.50 | - | - | - | 0.50 |
| | ROI Width | | (deg or ratio to max width) | 30.00 | 30.00 | - | - | - | 30.00 |
| | Frequency | | (MHz) | 5.00 | 5.00 | - | - | - | 5.00 |

Transducer Model: 8C-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-------|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.03 | - | 0.34 | - | 0.60 |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 2.05 | | | | |
| | P | W_0 (mW) | | - | 12.66 | | 12.77 |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | | 1.20 |
| | z at max $I_{pi \alpha}$ | Z_{sp} (cm) | 1.95 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | | 0.19 |
| | f_{awf} | f_c (MHz) | 4.00 | - | 5.70 | - | 4.03 |
| | Dim of A_{aprt} | X (cm) | | - | 1.24 | - | 0.44 |
| | | Y (cm) | | - | 0.50 | - | 0.50 |
| Other Info | t_d | PD (μ sec) | 7.27 | | | | |
| | prr | PRF (Hz) | 257 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.69 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | | 0.19 |
| | Focal Length | FL_x (cm) | | - | 1.12 | - | 0.12 |
| | | FL_y (cm) | | - | 0.28 | - | 0.30 |
| Operator Control | $I_{pi \alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 190 | | | | |
| | Power | (dB) | 0 | - | 0 | - | 0 |
| | Beam angle | (deg) | 0 | - | 0 | - | 0 |
| | Sample vol. position | (mm) | 20 | - | 228 | - | 10 |
| | Sample vol. | (mm) | 5.98 | - | 2.00 | - | 2.01 |
| | Scale | (m/sec) | 0.05 | - | 0.05 | - | 0.20 |
| | Frequency | (MHz) | 4.00 | - | 5.71 | - | 4.00 |

Transducer Model: 8L-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|-------------|--|------|-------|----------|---|-------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | |
| Maximum Index Value | | | | 1.59 | 0.99 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 3.48 | | | | | |
| | P | | W_0 (mW) | | 43.87 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{TA.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{pi\ \alpha}$ | | z_{sp} (cm) | 1.74 | | | | | |
| | $d_{eq}(z_b)$ | | $d_{eq}(Z_{sp})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 5.38 | 4.20 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 1.83 | - | - | - | # |
| | | Y (cm) | | | 0.40 | - | - | - | # |
| Other Info | t_d | | PD (μ sec) | 0.31 | | | | | |
| | p_{rr} | | PRF (Hz) | 40 | | | | | |
| | p_r at max I_{pi} | | $P_r@ P_{II_{max}}$ (Mpa) | 4.80 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq}@$ $P_{II_{max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 1.20 | - | - | | # |
| | | FL_y (cm) | | | 0.14 | - | - | | # |
| | $I_{pi\ \alpha}$ at max MI | | $I_{PA.3}@$ MI_{max} (W/cm ²) | 599 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 0 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 5.33 | 4.00 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 1.00 | 0.70 | - | - | - | - |
| | Depth | | (mm) | 100 | 100 | - | - | - | - |
| | Focus | | (mm) | 20 | 45 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 8L-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|------|------|----------------------------------|----------------------------|-----|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | 1.61 | - | $A_{\text{aprt}} \leq 1$ 0.59 | $A_{\text{aprt}} > 1$ - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 3.52 | | | | |
| | P | W_0 (mW) | | - | 14.41 | - | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 1.67 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 5.38 | - | 4.15 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 1.83 | - | # |
| | | Y (cm) | | - | 0.40 | - | # |
| Other Info | t_d | PD (μsec) | 0.31 | | | | |
| | p_{rr} | PRF (Hz) | 1000 | | | | |
| | p_r at max $I_{\text{pi},\alpha}$ | $P_r @ P_{\text{II,max}}$ (Mpa) | 4.80 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 1.20 | - | # |
| | | FL_y (cm) | | - | 0.14 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 5.73 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam angle | (deg) | 0 | - | 0 | - | - |
| | Frequency | (MHz) | 5.33 | - | 4.00 | - | - |
| | Depth | (mm) | 100 | - | 100 | - | - |
| | Focus | (mm) | 20 | - | 45 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 8L-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|-------------|------|------|----------|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | 1.58 | - | 0.15 | - | - | (b) |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (Mpa) | | 1.66 | | | | | |
| | P W_0 (mW) | | | - | 4.12 | | - | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | | - | | |
| | Z_s Z_1 (cm) | | | | | - | | |
| | Z_{bp} Z_{bp} (cm) | | | | | - | | |
| | Z_b Z_{sp} (cm) | | | | | | - | |
| | z at max $I_{\text{pi} \alpha}$ Z_{sp} (cm) | | 1.36 | | | | | |
| | d_{eq} (Z_b) d_{eq} (Z_{sp}) (cm) | | | | | | - | |
| | f_{awf} f_c (MHz) | | 5.00 | - | 6.68 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.51 | - | - | # |
| | | Y (cm) | | - | 0.40 | - | - | # |
| Other Info | t_d PD (μsec) | | 0.95 | | | | | |
| | p_{rr} PRF (Hz) | | 3497 | | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | | 2.10 | | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.04 | - | | # |
| | | FL_y (cm) | | - | 0.20 | - | | # |
| | $I_{\text{pi} \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | | 150 | | | | | |
| Operator Control | Power (dB0) | | 0 | - | 0 | - | - | - |
| | PRF (Hz) | | 3497 | - | 3497 | - | - | - |
| | ROI Span (mm) | | 30 | - | 10 | - | - | - |
| | ROI Center (mm) | | 15 | - | 10 | - | - | - |
| | Sample vol. (mm) | | 0.88 | - | 1.10 | - | - | - |
| | Frequency (MHz) | | 5.00 | - | 6.67 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 8L-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | |
|------------------------------|---|--|-------------|------|-------|----------|-----------------|-----|--------------------------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | | $A_{\text{aprt}} \leq 1$ |
| Maximum Index Value | | | | 1.59 | 0.91 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (Mpa) | | | 2.27 | | | | | |
| | P W_0 (mW) | | | | 33.80 | - | | - | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | | | - | | |
| | Z_s Z_1 (cm) | | | | | | - | | |
| | Z_{bp} Z_{bp} (cm) | | | | | | - | | |
| | Z_b Z_{sp} (cm) | | | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ Z_{sp} (cm) | | | 1.46 | | | | | |
| | $d_{\text{eq}}(z_b)$ $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | | | - | |
| | f_{awf} f_c (MHz) | | | 6.70 | 5.05 | - | - | - | # |
| | Dim of A_{aprt} | | X (cm) | | 1.92 | - | - | - | # |
| | | | Y (cm) | | 0.40 | - | - | - | # |
| Other Info | t_d PD (μsec) | | | 1.28 | | | | | |
| | p_{rr} PRF (Hz) | | | 245 | | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | | | 3.19 | | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | | | - | |
| | Focal Length | | FL_x (cm) | | 1.20 | - | - | | # |
| | | | FL_y (cm) | | 0.12 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | | | 234 | | | | | |
| Operator Control | Power (dB) | | | 0 | 0 | - | - | - | - |
| | Tilt (deg) | | | 0 | 0 | - | - | - | - |
| | Framerate (index) | | | 0 | 2 | - | - | - | - |
| | PRF (Hz) | | | 247 | 750 | - | - | - | - |
| | ROI Span (mm) | | | 10 | 160 | - | - | - | - |
| | ROI Center (mm) | | | 10 | 10 | - | - | - | - |
| | Sample vol. (mm) | | | 1.10 | 1.12 | - | - | - | - |
| | ROI width (deg or ratio to max width) | | | 0.51 | 0.20 | - | - | - | - |
| | Frequency (MHz) | | | 6.67 | 5.00 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 8L-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | | 1.51 | - | 0.17 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.98 | | | | |
| | P | W_0 | (mW) | | - | 7.86 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 1.02 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | | - | |
| | f_{awf} | f_c | (MHz) | 4.95 | - | 4.45 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.42 | - | # |
| | | Y | (cm) | | - | 0.40 | - | # |
| Other Info | t_d | PD | (μ sec) | 1.06 | | | | |
| | prr | PRF | (Hz) | 450 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 3.55 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | | - | |
| | Focal Length | FL_x | (cm) | | - | 1.10 | - | # |
| | | FL_y | (cm) | | - | 0.12 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 528 | | | | |
| | Power | (dB) | | 0 | - | 0 | - | - |
| | Beam angle | (deg) | | 0 | - | 0 | - | - |
| | Sample vol. position | (mm) | | 16 | - | 16 | - | - |
| | Sample vol. | (mm) | | 0.96 | - | 3.97 | - | - |
| | Scale | (m/sec) | | 0.07 | - | 1.99 | - | - |
| | Frequency | (MHz) | | 5.00 | - | 4.44 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9L-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | | | |
|------------------------------|--|-----------------|---|--------|----------------------|----------|-----------------|-----|-----------------------|-----------------------|---|
| | | | | | scan | non-scan | | | | | |
| | | | | | | | | | A _{aprt} ≤ 1 | A _{aprt} > 1 | |
| Maximum Index Value | | | | 1.39 | 0.76 | - | - | - | (b) | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | P _{r,α} | | P _{r,3} | (Mpa) | 1.53 | | | | | | |
| | P | | W ₀ | (mW) | | 47.55 | - | | - | # | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA.3} (Z ₁)] | | (mW) | | | - | | | |
| | Z _s | | Z ₁ | | (cm) | | | - | | | |
| | Z _{bp} | | Z _{bp} | | (cm) | | | - | | | |
| | Z _b | | Z _{sp} | | (cm) | | | | - | | |
| | z at max I _{pi α} | | Z _{sp} | | (cm) | 3.00 | | | | | |
| | d _{eq} (z _b) | | d _{eq} (Z _{sp}) | | (cm) | | | | - | | |
| | f _{awf} | | f _c | | (MHz) | 3.35 | 4.88 | - | - | - | # |
| | Dim of A _{aprt} | X | | (cm) | | 1.45 | - | - | - | - | # |
| | | Y | | (cm) | | 0.60 | - | - | - | - | # |
| Other Info | t _d | | PD | (μsec) | 0.57 | | | | | | |
| | prr | | PRF | (Hz) | 61 | | | | | | |
| | p _r at max I _{pi} | | P _r @ PII _{max} | (Mpa) | 2.17 | | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ PII _{max} | | (cm) | | | | - | | |
| | Focal Length | FL _x | | (cm) | | 0.92 | - | - | | # | |
| | | FL _y | | (cm) | | 0.14 | - | - | | # | |
| | I _{pi α} at max MI | | I _{PA.3} @ MI _{max} | | (W/cm ²) | 158 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - | - | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - | - | |
| | Framerate | | (index) | 0 | 1 | - | - | - | - | - | |
| | Frequency | | (MHz) | 3.08 | 5.71 | - | - | - | - | - | |
| | Width | | (deg or ratio to max width) | 0.70 | 1.00 | - | - | - | - | - | |
| | Depth | | (mm) | 120 | 120 | - | - | - | - | - | |
| | Focus | | (mm) | 30 | 54 | - | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 9L-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|--|----------------------|---|------|------|----------|---|-------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | |
| Maximum Index Value | | | | 1.38 | - | 0.31 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.97 | | | | | |
| | P | | W_0 (mW) | | - | 13.07 | | - | # |
| | Min of $[P_{\alpha}(Z_s), I_{ta,\alpha}(Z_s)]$ | | Min of $[W_3(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} (cm) | | | | - | | |
| | Z_b | | Z_{sp} (cm) | | | | | - | |
| | z at max $I_{pi\ \alpha}$ | | Z_{sp} (cm) | 2.65 | | | | | |
| | $d_{eq}(Z_b)$ | | $d_{eq}(Z_{sp})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 4.88 | - | 5.00 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | - | 1.47 | - | - | # |
| | | Y (cm) | | | - | 0.60 | - | - | # |
| Other Info | t_d | | PD (μ sec) | 0.31 | | | | | |
| | p_{rr} | | PRF (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r@ P_{II_{max}}$ (Mpa) | 4.64 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq}@ P_{II_{max}}$ (cm) | | | | | - | |
| | Focal Length | FL _x (cm) | | | - | 1.20 | - | | # |
| | | FL _y (cm) | | | - | 0.14 | - | | # |
| | $I_{pi\ \alpha}$ at max MI | | $I_{PA,3}@ MI_{max}$ (W/cm ²) | 433 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | - | - | |
| | Beam angle (deg) | | 0 | - | 0 | - | - | - | |
| | Frequency (MHz) | | 5.71 | - | 5.71 | - | - | - | |
| | Depth (mm) | | 120 | - | 120 | - | - | - | |
| | Focus (mm) | | 30 | - | 114 | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9L-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 1.38 | - | 0.12 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.47 | | | | |
| | P | W_0 (mW) | | - | 4.98 | - | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1), I_{\text{TA},.3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 3.38 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 6.55 | - | 4.55 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 1.01 | - | # |
| | | Y (cm) | | - | 0.60 | - | # |
| Other Info | t_d | PD (μsec) | 0.78 | | | | |
| | p_{rr} | PRF (Hz) | 250 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ (Mpa) | 3.14 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.08 | - | # |
| | | FL_y (cm) | | - | 0.14 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},.3} @ MI_{\text{max}}$ (W/cm^2) | 116 | | | | |
| | Power | (dB0) | 0 | - | 0 | - | - |
| | PRF | (Hz) | 250 | - | 3497 | - | - |
| | ROI Span | (mm) | 70 | - | 10 | - | - |
| | ROI Center | (mm) | 30 | - | 20 | - | - |
| | Sample vol. | (mm) | 0.70 | - | 0.74 | - | - |
| | Frequency | (MHz) | 6.67 | - | 4.44 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 9L-RS

Operating Mode: CFM

| Index Label | | | MI | TIS | | | | TIB non-scan | TIC |
|------------------------------|--|-------------|------|--------|----------|---|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | | | |
| | | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | 1.32 | 1.76 | - | - | - | (b) | |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (Mpa) | | 2.05 | | | | | | |
| | P W_0 (mW) | | | 120.90 | - | | - | # | |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | | - | | | |
| | Z_s Z_1 (cm) | | | | | - | | | |
| | Z_{bp} Z_{bp} (cm) | | | | | - | | | |
| | Z_b Z_{sp} (cm) | | | | | | - | | |
| | z at max $I_{\text{pi } \alpha}$ Z_{sp} (cm) | | 2.95 | | | | | | |
| | $d_{\text{eq}}(Z_b)$ $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | | - | | |
| | f_{awf} f_c (MHz) | | 4.50 | 4.45 | - | - | - | # | |
| | Dim of A_{aprt} | X (cm) | | 1.52 | - | - | - | # | |
| | | Y (cm) | | 0.60 | - | - | - | # | |
| Other Info | t_d PD (μsec) | | 0.79 | | | | | | |
| | p_{rr} PRF (Hz) | | 567 | | | | | | |
| | p_r at max I_{pi} $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | | 3.25 | | | | | | |
| | d_{eq} at max I_{pi} $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | | - | | |
| | Focal Length | FL_x (cm) | | 0.08 | - | - | | # | |
| | | FL_y (cm) | | 0.16 | - | - | | # | |
| | $I_{\text{pi } \alpha}$ at max MI $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | | 255 | | | | | | |
| Operator Control | Power (dB) | | 0 | 0 | - | - | - | - | |
| | Tilt (deg) | | 0 | 0 | - | - | - | - | |
| | Framerate (index) | | 0 | 0 | - | - | - | - | |
| | PRF (Hz) | | 8000 | 2500 | - | - | - | - | |
| | ROI Span (mm) | | 10 | 10 | - | - | - | - | |
| | ROI Center (mm) | | 25 | 25 | - | - | - | - | |
| | Sample vol. (mm) | | 0.74 | 1.18 | - | - | - | - | |
| | ROI width (deg or ratio to max width) | | 0.20 | 1.00 | - | - | - | - | |
| | Frequency (MHz) | | 4.44 | 4.44 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9L-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | | 1.09 | - | 0.33 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.22 | | | | |
| | P | W_0 | (mW) | | - | 15.79 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 2.29 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | | - | |
| | f_{awf} | f_c | (MHz) | 4.43 | - | 4.45 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.92 | - | # |
| | | Y | (cm) | | - | 0.60 | - | # |
| Other Info | t_d | PD | (μ sec) | 1.73 | | | | |
| | prr | PRF | (Hz) | 285 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 3.14 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | | - | |
| | Focal Length | FL_x | (cm) | | - | 0.14 | - | # |
| | | FL_y | (cm) | | - | 0.16 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 279 | | | | |
| | Power | (dB) | | 0 | - | 0 | - | - |
| | Beam angle | (deg) | | 0 | - | 0 | - | - |
| | Sample vol. position | (mm) | | 32 | - | 36 | - | - |
| | Sample vol. | (mm) | | 1.54 | - | 5.48 | - | - |
| | Scale | (m/sec) | | 0.05 | - | 0.05 | - | - |
| | Frequency | (MHz) | | 4.44 | - | 4.44 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 12L-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.39 | 1.87 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 3.71 | | | | | |
| | P | | W_0 (mW) | | 79.01 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.27 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 7.50 | 6.38 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 1.26 | - | - | - | # |
| | | Y (cm) | | | 0.40 | - | - | - | # |
| Other Info | t_d | | PD (μsec) | 0.30 | | | | | |
| | p_{rr} | | PRF (Hz) | 77 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | 5.16 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.92 | - | - | | # |
| | | FL_y (cm) | | | 0.08 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI_{max} (W/cm ²) | 475 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 0 | 0 | - | - | - | - |
| | Frequency | | (MHz) | 8.00 | 6.15 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 0.90 | 0.70 | - | - | - | - |
| | Depth | | (mm) | 60 | 60 | - | - | - | - |
| | Focus | | (mm) | 18 | 57 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 12L-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.41 | - | 0.27 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 3.56 | | | | |
| | P | W_0 (mW) | | - | 9.41 | - | # |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ | Min of $[W_3(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 1.54 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 7.50 | - | 6.13 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 1.18 | - | # |
| | | Y (cm) | | - | 0.40 | - | # |
| Other Info | t_d | PD (μ sec) | 0.31 | | | | |
| | prr | PRF (Hz) | 1000 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 5.31 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @ PII_{max}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.66 | - | # |
| | | FL_y (cm) | | - | 0.10 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @ MI_{max}$ (W/cm ²) | 490 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam angle | (deg) | 0 | - | 0 | - | - |
| | Frequency | (MHz) | 8.00 | - | 6.15 | - | - |
| | Depth | (mm) | 60 | - | 60 | - | - |
| | Focus | (mm) | 18 | - | 30 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 12L-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.35 | - | 0.08 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 3.29 | | | | |
| | P | W_0 (mW) | | - | 1.96 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 1.36 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 6.63 | - | 6.68 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.90 | - | # |
| | | Y (cm) | | - | 0.40 | - | # |
| Other Info | t_d | PD (μ sec) | 0.35 | | | | |
| | prr | PRF (Hz) | 247 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 4.50 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.04 | - | # |
| | | FL_y (cm) | | - | 0.10 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 549 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | PRF | (Hz) | 247 | - | 3497 | - | - |
| | ROI Span | (mm) | 10 | - | 35 | - | - |
| | ROI Center | (mm) | 15 | - | 15 | - | - |
| | Sample vol. | (mm) | 0.32 | - | 0.92 | - | - |
| | Frequency | (MHz) | 6.67 | - | 6.67 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 12L-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | | |
|------------------------------|--|--|---------------------------------------|-------|----------------------|----------|------|-------------------|-----|----------------|---|
| | | | | | scan | non-scan | | | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | | $A_{aprt} > 1$ | |
| Maximum Index Value | | | | 1.46 | 1.06 | - | - | - | (b) | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 3.42 | | | | | | |
| | P | | W_0 | (mW) | | 31.57 | - | | - | # | |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{TA.3}(Z_1)]$ | | (mW) | | | - | | | |
| | Z_s | | Z_1 | | (cm) | | | - | | | |
| | Z_{bp} | | Z_{bp} | | (cm) | | | - | | | |
| | Z_b | | Z_{sp} | | (cm) | | | | - | | |
| | z at max $I_{pi\alpha}$ | | Z_{sp} | | (cm) | 1.36 | | | | | |
| | $d_{eq}(Z_b)$ | | $d_{eq}(Z_{sp})$ | | (cm) | | | | - | | |
| | f_{awf} | | f_c | | (MHz) | 5.88 | 6.70 | - | - | - | # |
| | Dim of A_{aprt} | | X | | (cm) | | 0.78 | - | - | - | # |
| | | | Y | | (cm) | | 0.40 | - | - | - | # |
| Other Info | t_d | | PD | | (μ sec) | 0.43 | | | | | |
| | prr | | PRF | | (Hz) | 245 | | | | | |
| | p_r at max I_{pi} | | $P_r@PII_{max}$ | | (Mpa) | 4.51 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq}@PII_{max}$ | | (cm) | | | | - | | |
| | Focal Length | | FL _x | | (cm) | | 0.04 | - | - | | # |
| | | | FL _y | | (cm) | | 0.08 | - | - | | # |
| | $I_{pi\alpha}$ at max MI | | $I_{PA.3}@MI_{max}$ | | (W/cm ²) | 607 | | | | | |
| Operator Control | Power | | (dB) | | 0 | 0 | - | - | - | - | |
| | Tilt | | (deg) | | 0 | 0 | - | - | - | - | |
| | Framerate | | (index) | | 0 | 0 | - | - | - | - | |
| | PRF | | (Hz) | | 247 | 12000 | - | - | - | - | |
| | ROI Span | | (mm) | | 10 | 10 | - | - | - | - | |
| | ROI Center | | (mm) | | 10 | 10 | - | - | - | - | |
| | Sample vol. | | (mm) | | 0.36 | 0.92 | - | - | - | - | |
| | ROI Width | | (deg or ratio to max width) | | 0.20 | 0.20 | - | - | - | - | |
| | Frequency | | (MHz) | | 5.71 | 6.67 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 12L-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.41 | - | 0.63 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 3.53 | | | | |
| | P | W_0 (mW) | | - | 19.85 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 1.28 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 6.63 | - | 6.66 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 1.02 | - | # |
| | | Y (cm) | | - | 0.40 | - | # |
| Other Info | t_d | PD (μ sec) | 1.17 | | | | |
| | prr | PRF (Hz) | 427 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 4.73 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.60 | - | # |
| | | FL_y (cm) | | - | 0.08 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 557 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam angle | (deg) | 0 | - | 0 | - | - |
| | Sample vol. position | (mm) | 20 | - | 40 | - | - |
| | Sample vol. | (mm) | 1.03 | - | 6.48 | - | - |
| | Scale | (m/sec) | 0.05 | - | 1.39 | - | - |
| | Frequency | (MHz) | 6.67 | - | 6.67 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: i12L-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | | | |
|------------------------------|--|-----------------|---|-------|----------------------|----------|-----------------|-----|-----------------------|-----------------------|---|
| | | | | | scan | non-scan | | | | | |
| | | | | | | | | | A _{aprt} ≤ 1 | A _{aprt} > 1 | |
| Maximum Index Value | | | | 1.00 | 0.20 | - | - | - | (b) | | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | P _{r,α} | | P _{r,3} | (MPa) | 2.84 | | | | | | |
| | P | | W ₀ | (mW) | | 6.40 | - | | - | # | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA.3} (Z ₁)] | | (mW) | | | - | | | |
| | Z _s | | Z ₁ | | (cm) | | | - | | | |
| | Z _{bp} | | Z _{bp} | | (cm) | | | - | | | |
| | Z _b | | Z _{sp} | | (cm) | | | | - | | |
| | z at max I _{pi α} | | Z _{sp} | | (cm) | 1.48 | | | | | |
| | d _{eq} (z _b) | | d _{eq} (Z _{sp}) | | (cm) | | | | - | | |
| | f _{awf} | | f _c | | (MHz) | 8.00 | 7.75 | - | - | - | # |
| | Dim of A _{aprt} | X | | (cm) | | 1.12 | - | - | - | - | # |
| | | Y | | (cm) | | 0.70 | - | - | - | - | # |
| Other Info | t _d | | PD | (ms) | 0.18 | | | | | | |
| | prr | | PRF | (Hz) | 42 | | | | | | |
| | p _r at max I _{pi α} | | P _r @ P _{II} _{max} | | (MPa) | 4.27 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} | | (cm) | | | | - | | |
| | Focal Length | FL _x | | (cm) | | 0.04 | - | - | | # | |
| | | FL _y | | (cm) | | 0.20 | - | - | | # | |
| | I _{pi α} at max MI | | I _{PA.3} @ MI _{max} | | (W/cm ²) | 457 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - | - | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - | - | |
| | Framerate | | (index) | 0 | 0 | - | - | - | - | - | |
| | Frequency | | (MHz) | 10.00 | 10.00 | - | - | - | - | - | |
| | Width | | (deg or ratio to max width) | 0.40 | 0.40 | - | - | - | - | - | |
| | Depth | | (mm) | 60 | 60 | - | - | - | - | - | |
| | Focus | | (mm) | 15 | 18 | - | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: i12L-RS

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|--|------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 0.97 | - | 0.02 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.74 | | | | | |
| | P | | W_0 (mW) | | - | 0.49 | | - | # |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} (cm) | | | | - | | |
| | Z_b | | Z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi} \alpha}$ | | Z_{sp} (cm) | 1.21 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 8.00 | - | 8.00 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | - | 0.60 | - | - | # |
| | | Y (cm) | | | - | 0.70 | - | - | # |
| Other Info | t_d | | PD (ms) | 0.17 | | | | | |
| | p_{rr} | | PRF (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 3.83 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @ P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | - | 0.04 | - | | # |
| | | FL_y (cm) | | | - | 0.04 | - | | # |
| | $I_{\text{pi} \alpha}$ at max MI | | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 318 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | - | - | |
| | Beam Angle (deg) | | 0 | - | 0 | - | - | - | |
| | Frequency (MHz) | | 10.00 | - | 10.00 | - | - | - | |
| | Depth (mm) | | 60 | - | 60 | - | - | - | |
| | Focus (mm) | | 15 | - | 15 | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: i12L-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|----------------------|------|------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | | 1.33 | - | 0.08 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 3.42 | | | | |
| | P | W_0 | (mW) | | - | 2.97 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} | (cm) | 0.86 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | - | |
| | f_{awf} | f_c | (MHz) | 6.55 | - | 5.15 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.73 | - | # |
| | | Y | (cm) | | - | 0.70 | - | # |
| Other Info | t_d | PD | (ms) | 0.67 | | | | |
| | p_{rr} | PRF | (Hz) | 500 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ | (MPa) | 4.15 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | - | |
| | Focal Length | FL_x | (cm) | | - | 0.08 | - | # |
| | | FL_y | (cm) | | - | 0.12 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 655 | | | | |
| | Power | (dB) | | 0 | - | 0 | - | - |
| | PRF | (Hz) | | 500 | - | 11905 | - | - |
| | ROI Span | (mm) | | 19 | - | 38 | - | - |
| | ROI Center | (mm) | | 9 | - | 18 | - | - |
| | Sample Volume | (mm) | | 0.61 | - | 0.64 | - | - |
| Operator Control | Frequency | (MHz) | | 6.67 | - | 5.00 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: i12L-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|----------------------|-------|------|--------------------------|-----------------------|-----|
| | | | | | scan | non-scan | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | | 0.86 | 0.21 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 1.75 | | | | |
| | P | W_0 | (mW) | | 7.86 | - | | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | | - |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} | (cm) | 0.82 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - |
| | f_{awf} | f_c | (MHz) | 5.15 | 5.15 | - | - | - |
| | Dim of A_{aprt} | X | (cm) | | 0.18 | - | - | - |
| | | Y | (cm) | | 0.70 | - | - | - |
| Other Info | t_d | PD | 0.52 | 0.82 | | | | |
| | p_{rr} | PRF | (Hz) | 1667 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 2.03 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - |
| | Focal Length | FL_x | (cm) | | 0.14 | - | - | # |
| | | FL_y | (cm) | | 0.06 | - | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 104 | | | | |
| | Power | (dB) | | 0 | 0 | - | - | - |
| | Tilt | (deg) | | 0 | 0 | - | - | - |
| | Framerate | (index) | | 0 | 0 | - | - | - |
| | PRF | (Hz) | | 10000 | 248 | - | - | - |
| | ROI Span | (mm) | | 11 | 3 | - | - | - |
| | ROI Center | (mm) | | 3 | 3 | - | - | - |
| | Sample Voume | (mm) | | 0.64 | 0.64 | - | - | - |
| | ROI Width | (deg or ratio to max width) | | 0.10 | 0.10 | - | - | - |
| | Frequency | (MHz) | | 5.00 | 5.00 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: i12L-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 1.30 | - | 0.03 - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 2.64 | | | | |
| | P | W_0 (mW) | | - | 1.77 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{3,3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 0.70 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 4.13 | - | 4.13 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.31 | - | # |
| | | Y (cm) | | - | 0.70 | - | # |
| Other Info | t_d | PD (ms) | 1.03 | | | | |
| | p_{rr} | PRF (Hz) | 545 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{II_{\text{max}}}$ (MPa) | 2.92 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ $P_{II_{\text{max}}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.12 | - | # |
| | | FL_y (cm) | | - | 0.06 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 294 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | (cm) | 0 | - | 0 | - | - |
| | Sample Volume Position | (mm) | 12 | - | 12 | - | - |
| | Sample Volume | (mm) | 1.04 | - | 1.04 | - | - |
| | Scale | (m/s) | 0.10 | - | 0.10 | - | - |
| | Frequency | (MHz) | 4.00 | - | 4.00 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6T-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | |
|------------------------------|---|--|----------------------|-------|-------|--------------------------|-----------------------|-----|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | | 1.16 | 0.45 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.37 | | | | | |
| | P | W_0 | (mW) | | 21.15 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | | |
| | Z_b | Z_{sp} | (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | Z_{sp} | (cm) | 2.80 | | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - | |
| | f_{awf} | f_c | (MHz) | 4.50 | 4.50 | - | - | - | # |
| | Dim of A_{aprt} | X | (cm) | | 0.93 | - | - | - | # |
| | | Y | (cm) | | 0.80 | - | - | - | # |
| Other Info | t_d | PD | (μs) | 0.47 | | | | | |
| | p_{rr} | PRF | (Hz) | 41 | | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ | (Mpa) | 3.65 | | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - | |
| | Focal Length | FL_x | (cm) | | 0.14 | - | - | | # |
| | | FL_y | (cm) | | 0.32 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA}.3} @$ MI_{max} | (W/cm ²) | 140 | | | | | |
| Operator Control | Power | (dB) | 0 | 0 | - | - | - | - | |
| | Tilt | (deg) | 0 | 0 | - | - | - | - | |
| | Framerate | (index) | 1 | 1 | - | - | - | - | |
| | Frequency | (MHz) | 6.15 | 6.15 | - | - | - | - | |
| | Width | (deg or ratio to max width) | 75.00 | 75.00 | - | - | - | - | |
| | Depth | (mm) | 200 | 200 | - | - | - | - | |
| | Focus | (mm) | 50 | 50 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6T-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.01 | | 0.13 - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.81 | | | | |
| | P | W_0 (mW) | | - | 8.26 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 2.98 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.20 | - | 3.25 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.77 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (μs) | 0.66 | | | | |
| | prr | PRF (Hz) | 1000 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.52 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.22 | - | # |
| | | FL_y (cm) | | - | 0.42 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 121 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam angle | (deg) | 0 | - | 0 | - | - |
| | Frequency | (MHz) | 2.86 | - | 2.86 | - | - |
| | Depth | (mm) | 200 | - | 200 | - | - |
| | Focus | (mm) | 70 | - | 70 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 6T-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.08 | - | 0.65 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.20 | | | | |
| | P | W_0 (mW) | | - | 26.19 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 3.30 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.70 | - | 4.93 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.93 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (μs) | 1.07 | | | | |
| | prr | PRF (Hz) | 4505 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 1.84 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.30 | - | # |
| | | FL_y (cm) | | - | 0.20 | - | # |
| | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 91 | | | | |
| Operator Control | Power | (dB) | 0 | - | 0 | - | - |
| | PRF | (Hz) | 4505 | - | 3425 | - | - |
| | ROI Span | (mm) | 45 | - | 300 | - | - |
| | ROI Center | (mm) | 90 | - | 70 | - | - |
| | Sample Vol. | (mm) | 0.96 | - | 1.00 | - | - |
| | Frequency | (MHz) | 3.63 | - | 5.00 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6T-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|-------|-------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | | 1.11 | 0.41 | - - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 0.77 | | | | |
| | P | W_0 | (mW) | | 14.46 | - | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 3.20 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | | - | |
| | f_{awf} | f_c | (MHz) | 3.75 | 3.45 | - | - | # |
| | Dim of A_{aprt} | X | (cm) | | 0.93 | - | - | # |
| | | Y | (cm) | | 0.80 | - | - | # |
| Other Info | t_d | PD | (μs) | 1.09 | | | | |
| | prr | PRF | (Hz) | 751 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 1.17 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | | - | |
| | Focal Length | FL_x | (cm) | | 0.18 | - | - | # |
| | | FL_y | (cm) | | 0.32 | - | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 26 | | | | |
| | Power | (dB) | | 0 | 0 | - | - | - |
| | Tilt | (deg) | | 0 | 0 | - | - | - |
| | Framerate | (index) | | 3 | 0 | - | - | - |
| | PRF | (Hz) | | 750 | 250 | - | - | - |
| | ROI Span | (mm) | | 95 | 20 | - | - | - |
| | ROI Center | (mm) | | 90 | 40 | - | - | - |
| | Sample vol. | (mm) | | 0.96 | 0.60 | - | - | - |
| | ROI Width | (deg or ratio to max width) | | 15.00 | 30.00 | - | - | - |
| | Frequency | (MHz) | | 3.64 | 3.08 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 6T-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.15 | - | 0.94 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 1.98 | | | | |
| | P | W_0 (mW) | | - | 53.21 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 3.10 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.25 | - | 3.70 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.93 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (μs) | 1.11 | | | | |
| | prr | PRF (Hz) | 1712 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | 2.80 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.28 | - | # |
| | | FL_y (cm) | | - | 0.22 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 158 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam angle | (deg) | 0 | - | 0 | - | - |
| | Sample vol. position | (mm) | 93 | - | 93 | - | - |
| | Sample vol. | (mm) | 1.02 | - | 1.02 | - | - |
| | Scale | (m/sec) | 0.40 | - | 5.08 | - | - |
| | Frequency | (MHz) | 3.08 | - | 3.63 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6T-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|-----|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 0.85 - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 30.40 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | # | - | 4.00 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.41 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (μs) | # | | | | |
| | prr | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.34 | - | # |
| | | FL_y (cm) | | - | 0.60 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. position (mm) | | - | - | 51 | - | - |
| | Frequency (MHz) | | - | - | 4.00 | - | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6Tc-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.23 | 0.40 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.17 | | | | | |
| | P | | W_0 (mW) | | 16.15 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.60 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 3.10 | 5.25 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 0.93 | - | - | - | # |
| | | Y (cm) | | | 0.80 | - | - | - | # |
| Other Info | t_d | | PD (ms) | 0.59 | | | | | |
| | p_{rr} | | PRF (Hz) | 62 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II} \text{ max}}$ (MPa) | 2.57 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II} \text{ max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.04 | - | - | | # |
| | | FL_y (cm) | | | 0.58 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 132 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 3 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 2.86 | 6.15 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 75.00 | 30.00 | - | - | - | - |
| | Depth | | (mm) | 200 | 200 | - | - | - | - |
| | Focus | | (mm) | 20 | 20 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6Tc-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 1.18 | - | 0.11 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 2.08 | | | | |
| | P | W_0 (mW) | | - | 7.59 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 1.67 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.15 | - | 3.15 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.77 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (ms) | 0.61 | | | | |
| | p_{rr} | PRF (Hz) | 1000 | | | | |
| | p_r at max $I_{pi\alpha}$ | $P_r @ P_{II_{\text{max}}}$ (MPa) | 2.49 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ $P_{II_{\text{max}}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.12 | - | # |
| | | FL_y (cm) | | - | 0.62 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 120 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | (deg) | 0 | - | 0 | - | - |
| | Frequency | (MHz) | 2.86 | - | 2.86 | - | - |
| | Depth | (mm) | 200 | - | 200 | - | - |
| | Focus | (mm) | 20 | - | 20 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 6Tc-RS

Operating Mode: CMM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.21 | - | 0.37 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 1.53 | | | | |
| | P | W_0 (mW) | | - | 18.43 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | 3.20 | | | | |
| | $d_{eq}(z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.68 | - | 3.88 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.93 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (ms) | 1.10 | | | | |
| | prr | PRF (Hz) | 1000 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (MPa) | 2.29 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.08 | - | # |
| | | FL_y (cm) | | - | 0.58 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | 123 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | PRF | (Hz) | 1000 | - | 3497 | - | - |
| | ROI Span | (mm) | 100 | - | 20 | - | - |
| | ROI Center | (mm) | 90 | - | 20 | - | - |
| | Sample Volume | (mm) | 0.95 | - | 0.60 | - | - |
| | Frequency | (MHz) | 3.63 | - | 3.63 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6Tc-RS

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|--|-----------------|---|----------------------|-------|----------|---|-----------------------|-----|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | A _{aprt} ≤ 1 | | A _{aprt} > 1 |
| Maximum Index Value | | | | 1.25 | 0.61 | - | - | - | (b) | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | P _{r,α} | | P _{r,3} | (MPa) | 0.96 | | | | | |
| | P | | W ₀ | (mW) | | 28.26 | - | | # | |
| | Min of [P _α (Z _s), I _{ta,α} (Z _s)] | | Min of [W _{.3} (Z ₁), I _{TA,3} (Z ₁)] | (mW) | | | - | | | |
| | Z _s | | Z ₁ | (cm) | | | - | | | |
| | Z _{bp} | | Z _{bp} | (cm) | | | - | | | |
| | Z _b | | Z _{sp} | (cm) | | | | - | | |
| | z at max I _{pi α} | | Z _{sp} | (cm) | 3.50 | | | | | |
| | d _{eq} (Z _b) | | d _{eq} (Z _{sp}) | (cm) | | | | - | | |
| | f _{awf} | | f _c | (MHz) | 3.70 | 3.65 | - | - | - | # |
| | Dim of A _{aprt} | X | | (cm) | | 0.93 | - | - | - | # |
| | | Y | | (cm) | | 0.80 | - | - | - | # |
| Other Info | t _d | | PD | (ms) | 1.09 | | | | | |
| | p _{rr} | | PRF | (Hz) | 251 | | | | | |
| | p _r at max I _{pi} | | P _r @ P _{II} _{max} | (MPa) | 1.50 | | | | | |
| | d _{eq} at max I _{pi} | | d _{eq} @ P _{II} _{max} | (cm) | | | | - | | |
| | Focal Length | FL _x | | (cm) | | 0.12 | - | - | | # |
| | | FL _y | | (cm) | | 0.48 | - | - | | # |
| | I _{pi α} at max MI | | I _{PA,3} @ MI _{max} | (W/cm ²) | 51 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - | |
| | Framerate | | (index) | 0 | 0 | - | - | - | - | |
| | PRF | | (Hz) | 250 | 248 | - | - | - | - | |
| | ROI Span | | (mm) | 100 | 20 | - | - | - | - | |
| | ROI Center | | (mm) | 90 | 20 | - | - | - | - | |
| | Sample Voume | | (mm) | 0.95 | 0.95 | - | - | - | - | |
| | ROI Width | | (deg or ratio to max width) | 30.00 | 15.00 | - | - | - | - | |
| | Frequency | | (MHz) | 3.64 | 3.64 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6Tc-RS

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|------|------|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | |
| | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | 1.33 | - | 0.62 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 2.49 | | | | |
| | P | W_0 (mW) | | - | 37.25 | | # |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 1.21 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 3.50 | - | 3.50 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.49 | - | # |
| | | Y (cm) | | - | 0.80 | - | # |
| Other Info | t_d | PD (ms) | 0.87 | | | | |
| | p_{rr} | PRF (Hz) | 1873 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ (MPa) | 2.88 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.12 | - | # |
| | | FL_y (cm) | | - | 0.72 | - | # |
| | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm ²) | 259 | | | | |
| Operator Control | Power | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | (cm) | 0 | - | 0 | - | - |
| | Sample Volume Position | (mm) | 10 | - | 10 | - | - |
| | Sample Volume | (mm) | 0.94 | - | 0.94 | - | - |
| | Scale | (m/s) | 0.40 | - | 0.40 | - | - |
| | Frequency | (MHz) | 3.30 | - | 3.30 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 6Tc-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 0.56 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | # | | | |
| | P | W_0 | (mW) | | - | 35.86 | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | - | |
| | Z_s | Z_1 | (cm) | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | - | |
| | Z_b | Z_{sp} | (cm) | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | # | | | |
| | $d_{eq}(z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | - | |
| | f_{awf} | f_c | (MHz) | # | - | 3.30 | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.41 | # |
| | | Y | (cm) | | - | 0.80 | # |
| Other Info | t_d | PD | (ms) | # | | | |
| | prr | PRF | (Hz) | # | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (MPa) | # | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | - | |
| | Focal Length | FL_x | (cm) | | - | 0.26 | # |
| | | FL_y | (cm) | | - | 0.66 | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | # | | | |
| | Sample Volume Position | | (mm) | - | - | 31 | - |
| | Frequency | | (MHz) | - | - | 3.30 | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9T-RS

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | | |
|------------------------------|---|---------------|--|----------------------|-------|----------|-----------------|-----|--------------------------|-----------------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |
| Maximum Index Value | | | | 1.19 | 0.44 | - | - | - | (b) | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | 2.80 | | | | | |
| | P | | W_0 | (mW) | | 13.57 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | z_{bp} | | z_{bp} | (cm) | | | | - | | |
| | z_b | | z_{sp} | (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} | (cm) | 0.90 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - | |
| | f_{awf} | | f_c | (MHz) | 5.50 | 6.88 | - | - | - | # |
| | Dim of A_{aprt} | X | | (cm) | | 0.49 | - | - | - | # |
| | | Y | | (cm) | | 0.65 | - | - | - | # |
| Other Info | t_d | | PD | (μsec) | 0.31 | | | | | |
| | p_{rr} | | PRF | (Hz) | 82 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}}_{\text{max}}$ | (Mpa) | 3.31 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}}_{\text{max}}$ | (cm) | | | | | - | |
| | Focal Length | FL_x | | (cm) | | 0.04 | - | - | | # |
| | | FL_y | | (cm) | | 0.56 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} | (W/cm ²) | 233 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - | |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - | |
| | Framerate | | (index) | 1 | 3 | - | - | - | - | |
| | Frequency | | (MHz) | 7.27 | 7.27 | - | - | - | - | |
| | Width | | (deg or ratio to max width) | 60.00 | 30.00 | - | - | - | - | |
| | Depth | | (mm) | 140 | 140 | - | - | - | - | |
| | Focus | | (mm) | 9 | 9 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9T-RS

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | 1.23 | - | 0.10 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.89 | | | |
| | P | W_0 | (mW) | | - | 3.71 | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | - | |
| | Z_s | Z_1 | (cm) | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | - | |
| | Z_b | Z_{sp} | (cm) | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 0.90 | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | - | |
| | f_{awf} | f_c | (MHz) | 5.50 | - | 5.75 | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.49 | # |
| | | Y | (cm) | | - | 0.65 | # |
| Other Info | t_d | PD | (μ sec) | 0.33 | | | |
| | prr | PRF | (Hz) | 1000 | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 3.42 | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | - | |
| | Focal Length | FL_x | (cm) | | - | 0.04 | # |
| | | FL_y | (cm) | | - | 0.56 | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 229 | | | |
| | Power | (dB) | | 0 | - | 0 | - |
| | Beam angle | (deg) | | 0 | - | 0 | - |
| | Frequency | (MHz) | | 7.27 | - | 7.27 | - |
| | Depth | (mm) | | 140 | - | 140 | - |
| | Focus | (mm) | | 9 | - | 9 | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9T-RS

Operating Mode: CMM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|--|------|------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.32 | - | 0.54 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (Mpa) | 2.83 | | | | | |
| | P | | W_0 (mW) | | - | 21.73 | | - | # |
| | Min of $[P_{\alpha}(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | | |
| | z_s | | z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 1.69 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 4.95 | - | 5.05 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | - | 0.59 | - | - | # |
| | | Y (cm) | | | - | 0.65 | - | - | # |
| Other Info | t_d | | PD (μsec) | 0.83 | | | | | |
| | p_{rr} | | PRF (Hz) | 347 | | | | | |
| | p_r at max I_{pi} | | $P_{r@PII_{\text{max}}}$ (Mpa) | 3.77 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}@PII_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | - | 0.08 | - | | # |
| | | FL_y (cm) | | | - | 0.50 | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA},3@MI_{\text{max}}}$ (W/cm^2) | 252 | | | | | |
| Operator Control | Power (dB) | | 0 | - | 0 | - | - | - | |
| | PRF (Hz) | | 347 | - | 2000 | - | - | - | |
| | ROI Span (mm) | | 10 | - | 10 | - | - | - | |
| | ROI Center (mm) | | 20 | - | 20 | - | - | - | |
| | Sample vol. (mm) | | 0.68 | - | 0.68 | - | - | - | |
| | Frequency (MHz) | | 5.00 | - | 5.00 | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9T-RS

Operating Mode: CFM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|-------|-------|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | |
| | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | 1.19 | 0.59 | - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | 0.60 | | | | |
| | P | W_0 (mW) | | 21.44 | - | | # |
| | Min of $[P_\alpha(Z_s), I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1), I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 2.20 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 5.00 | 5.10 | - | - | # |
| | Dim of A_{aprt} | X (cm) | | 0.59 | - | - | # |
| | | Y (cm) | | 0.65 | - | - | # |
| Other Info | t_d | PD (μsec) | 0.81 | | | | |
| | p_{rr} | PRF (Hz) | 349 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ (Mpa) | 0.88 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @ P_{\text{II,max}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | 0.14 | - | - | # |
| | | FL_y (cm) | | 0.44 | - | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @ MI_{\text{max}}$ (W/cm^2) | 13 | | | | |
| | Power | (dB) | 0 | 0 | - | - | - |
| | Tilt | (deg) | 0 | 0 | - | - | - |
| | Framerate | (index) | 0 | 0 | - | - | - |
| | PRF | (Hz) | 350 | 350 | - | - | - |
| | ROI Span | (mm) | 10 | 10 | - | - | - |
| | ROI Center | (mm) | 20 | 35 | - | - | - |
| | Sample vol. | (mm) | 0.68 | 0.68 | - | - | - |
| | ROI Width | (deg or ratio to max width) | 45.00 | 15.00 | - | - | - |
| | Frequency | (MHz) | 5.00 | 5.00 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: 9T-RS

Operating Mode: PW

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|---|----------------------|------|------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | | 1.08 | - | 0.78 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (Mpa) | 2.33 | | | | |
| | P | W_0 | (mW) | | - | 32.71 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} | (cm) | 1.00 | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ | (cm) | | | | - | |
| | f_{awf} | f_c | (MHz) | 4.95 | - | 5.03 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.49 | - | # |
| | | Y | (cm) | | - | 0.65 | - | # |
| Other Info | t_d | PD | (μ sec) | 1.00 | | | | |
| | prr | PRF | (Hz) | 2604 | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ | (Mpa) | 2.76 | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} | (cm) | | | | - | |
| | Focal Length | FL_x | (cm) | | - | 0.08 | - | # |
| | | FL_y | (cm) | | - | 0.50 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} | (W/cm ²) | 210 | | | | |
| | Power | (dB) | | 0 | - | 0 | - | - |
| | Beam angle | (deg) | | 0 | - | 0 | - | - |
| | Sample vol. position | (mm) | | 10 | - | 10 | - | - |
| | Sample vol. | (mm) | | 1.00 | - | 1.48 | - | - |
| | Scale | (m/sec) | | 0.40 | - | 0.40 | - | - |
| | Frequency | (MHz) | | 5.00 | - | 5.00 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: 9T-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|-----|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 0.63 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 26.53 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | # | - | 5.00 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.26 | - | # |
| | | Y (cm) | | - | 0.65 | - | # |
| Other Info | t_d | PD (μ sec) | # | | | | |
| | prr | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.14 | - | # |
| | | FL_y (cm) | | - | 0.50 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. position (mm) | | - | - | 300 | - | - |
| | Frequency (MHz) | | - | - | 5.00 | - | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: P2D-RS

Operating Mode: CW

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | | |
|------------------------------|---|----------------------|---|--------|------|---------------------|-------|-----------------|-----|-----|---|
| | | | | | scan | non-scan | | | | | |
| | | | | | | Maximum Index Value | | | | (a) | - |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (Mpa) | # | | | | | | |
| | P | | W_0 | (mW) | | - | 59.66 | | - | # | |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | | |
| | Z_s | | Z_1 | (cm) | | | | - | | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | | |
| | Z_b | | Z_{sp} | (cm) | | | | | - | | |
| | z at max $I_{\text{pi} \alpha}$ | | Z_{sp} | (cm) | # | | | | | | |
| | $d_{\text{eq}}(Z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - | | |
| | f_{awf} | | f_c | (MHz) | # | - | 2.00 | - | - | # | |
| | Dim of A_{aprt} | | | X (cm) | | - | 0.68 | - | - | # | |
| | | | | Y (cm) | | - | 1.15 | - | - | # | |
| Other Info | t_d | | PD (μsec) | # | | | | | | | |
| | p_{rr} | | PRF (Hz) | # | | | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | # | | | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | | | |
| | Focal Length | FL _x (cm) | | | - | 0.36 | - | | # | | |
| | | FL _y (cm) | | | - | 0.48 | - | | # | | |
| | $I_{\text{pi} \alpha}$ at max MI | | $I_{\text{PA},3} @$ MI _{max} (W/cm ²) | # | | | | | | | |
| Operator Control | Sample vol. position (mm) | | | - | - | 85 | - | - | - | | |
| | Frequency (MHz) | | | - | - | 2.00 | - | - | - | | |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: P6D-RS

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|-----|------|----------------------------------|----------------------------|-----|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | (a) | - | $A_{\text{aprt}} \leq 1$ 0.07 | $A_{\text{aprt}} > 1$ - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (Mpa) | # | | | | |
| | P | W_0 (mW) | | - | 2.84 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | # | - | 5.00 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.30 | - | # |
| | | Y (cm) | | - | 0.27 | - | # |
| Other Info | t_d | PD (μsec) | # | | | | |
| | p_{rr} | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (Mpa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.32 | - | # |
| | | FL_y (cm) | | - | 0.34 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm^2) | # | | | | |
| | Sample vol. position | (mm) | - | - | 20 | - | - |
| | Frequency | (MHz) | - | - | 5.00 | - | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav8

Operating Mode: 2D

| Index Label | | | | MI | TIS | | TIB non-scan | TIC | |
|------------------------------|---|-------------|---|-------|-------|----------|-----------------|-----|--------------------------|
| | | | | | scan | non-scan | | | |
| | | | | | | | | | $A_{\text{aprt}} \leq 1$ |
| Maximum Index Value | | | | 1.15 | 1.26 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 3.08 | | | | | |
| | P | | W_0 (mW) | | 36.62 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | Z_s | | Z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 0.55 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 7.13 | 7.25 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 0.70 | - | - | - | # |
| | | Y (cm) | | | 0.19 | - | - | - | # |
| Other Info | t_d | | PD (ms) | 0.25 | | | | | |
| | p_{rr} | | PRF (Hz) | 17 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 3.53 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.58 | - | - | | # |
| | | FL_y (cm) | | | 0.14 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 352 | | | | | |
| | | | | | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 1 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 8.00 | 8.00 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 60.00 | 75.00 | - | - | - | - |
| | Depth | | (mm) | 160 | 160 | - | - | - | - |
| | Focus | | (mm) | 99 | 99 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav8

Operating Mode: M-Mode

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 1.14 | - | 0.31 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 3.05 | | | | |
| | P | W_0 (mW) | | - | 8.73 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 0.55 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 7.13 | - | 7.38 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.70 | - | # |
| | | Y (cm) | | - | 0.19 | - | # |
| Other Info | t_d | PD (ms) | 0.25 | | | | |
| | p_{rr} | PRF (Hz) | 1000 | | | | |
| | p_r at max $I_{\text{pi},\alpha}$ | $P_r @ P_{\text{II,max}}$ (MPa) | 3.49 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II,max}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.58 | - | # |
| | | FL_y (cm) | | - | 0.16 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm ²) | 307 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | (deg) | 0 | - | 0 | - | - |
| | Frequency | (MHz) | 8.00 | - | 8.00 | - | - |
| | Depth | (mm) | 160 | - | 160 | - | - |
| | Focus | (mm) | 99 | - | 69 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: AcuNav8

Operating Mode: CMM

| Index Label | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|-------------|------|--|-----------------|-----|
| | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | 1.13 | - | 0.43 | - | (b) |
| Assoc. Acoustic Param. | IEC FDA Units | | | | | |
| | $P_{r,\alpha}$ $P_{r,3}$ (MPa) | 1.62 | | | | |
| | P W_0 (mW) | | - | 16.18 | - | # |
| | Min of $[P_\alpha(Z_s), I_{ta,\alpha}(Z_s)]$ Min of $[W_3(Z_1), I_{TA,3}(Z_1)]$ (mW) | | | - | | |
| | Z_s Z_1 (cm) | | | - | | |
| | Z_{bp} Z_{bp} (cm) | | | - | | |
| | Z_b Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ Z_{sp} (cm) | 0.55 | | | | |
| | $d_{eq}(Z_b)$ $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} f_c (MHz) | 4.95 | - | 4.95 | - | # |
| | Dim of A_{aprt} | X (cm) | - | 0.70 | - | # |
| | | Y (cm) | - | 0.19 | - | # |
| Other Info | t_d PD (ms) | 0.92 | | | | |
| | prr PRF (Hz) | 2000 | | | | |
| | p_r at max I_{pi} $P_r @ PII_{max}$ (MPa) | 1.78 | | | | |
| | d_{eq} at max I_{pi} $d_{eq} @ PII_{max}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | - | 0.58 | - | # |
| | | FL_y (cm) | - | 0.16 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI $I_{PA,3} @ MI_{max}$ (W/cm ²) | 94 | | | | |
| | Power (dB) | 0 | - | 0 | - | - |
| | PRF (Hz) | 2000 | - | 2000 | - | - |
| | ROI Span (mm) | 90 | - | 90 | - | - |
| | ROI Center (mm) | 20 | - | 20 | - | - |
| | Sample Volume (mm) | 0.73 | - | 0.73 | - | - |
| Operator Control | Frequency (MHz) | 5.00 | - | 5.00 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav8

Operating Mode: CFM

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|--|----------------------|-------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.17 | 0.98 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 0.96 | | | | | |
| | P | W_0 | (mW) | | 12.98 | - | | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | | |
| | z_s | z_1 | (cm) | | | | - | | |
| | z_{bp} | z_{bp} | (cm) | | | | - | | |
| | z_b | z_{sp} | (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | z_{sp} | (cm) | 0.55 | | | | | |
| | $d_{\text{eq}}(z_b)$ | $d_{\text{eq}}(z_{\text{sp}})$ | (cm) | | | | | - | |
| | f_{awf} | f_c | (MHz) | 4.95 | 5.00 | - | - | - | # |
| | Dim of A_{aprt} | X | (cm) | | 0.70 | - | - | - | # |
| | | Y | (cm) | | 0.19 | - | - | - | # |
| Other Info | t_d | PD | (ms) | 0.93 | | | | | |
| | p_{rr} | PRF | (Hz) | 496 | | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | 1.06 | | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - | |
| | Focal Length | FL_x | (cm) | | 0.58 | - | - | | # |
| | | FL_y | (cm) | | 0.14 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 30 | | | | | |
| Operator Control | Power | (dB) | 0 | 0 | - | - | - | - | |
| | Tilt | (deg) | 0 | 0 | - | - | - | - | |
| | Framerate | (index) | 0 | 0 | - | - | - | - | |
| | PRF | (Hz) | 500 | 500 | - | - | - | - | |
| | ROI Span | (mm) | 25 | 25 | - | - | - | - | |
| | ROI Center | (mm) | 40 | 40 | - | - | - | - | |
| | Sample Voume | (mm) | 0.73 | 0.73 | - | - | - | - | |
| | ROI Width | (deg or ratio to max width) | 15.00 | 15.00 | - | - | - | - | |
| | Frequency | (MHz) | 5.00 | 5.00 | - | - | - | - | |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav8

Operating Mode: PW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|------|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | 0.57 | - | 0.81 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 1.26 | | | | |
| | P | W_0 (mW) | | - | 25.50 | - | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | - | | |
| | Z_s | Z_1 (cm) | | | - | | |
| | Z_{bp} | Z_{bp} (cm) | | | - | | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{\text{pi},\alpha}$ | Z_{sp} (cm) | 0.51 | | | | |
| | $d_{\text{eq}}(Z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | 5.00 | - | 6.70 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.62 | - | # |
| | | Y (cm) | | - | 0.19 | - | # |
| Other Info | t_d | PD (ms) | 1.11 | | | | |
| | p_{rr} | PRF (Hz) | 2577 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 1.38 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.56 | - | # |
| | | FL_y (cm) | | - | 0.14 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI _{max} (W/cm ²) | 56 | | | | |
| | Power | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | (cm) | 0 | - | 0 | - | - |
| | Sample Volume Position | (mm) | 50 | - | 50 | - | - |
| | Sample Volume | (mm) | 1.02 | - | 2.52 | - | - |
| | Scale | (m/s) | 0.40 | - | 0.60 | - | - |
| | Frequency | (MHz) | 5.00 | - | 6.67 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav8

Operating Mode: CW

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|--|--|-----|------|--|-----------------|-----|
| | | | | scan | non-scan $A_{aprt} \leq 1$ $A_{aprt} > 1$ | | |
| Maximum Index Value | | | (a) | - | 0.10 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | # | | | | |
| | P | W_0 (mW) | | - | 4.19 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{TA,3}(Z_1)]$ (mW) | | | | - | |
| | Z_s | Z_1 (cm) | | | | - | |
| | Z_{bp} | Z_{bp} (cm) | | | | - | |
| | Z_b | Z_{sp} (cm) | | | | - | |
| | z at max $I_{pi\alpha}$ | Z_{sp} (cm) | # | | | | |
| | $d_{eq}(Z_b)$ | $d_{eq}(Z_{sp})$ (cm) | | | | - | |
| | f_{awf} | f_c (MHz) | # | - | 5.00 | - | # |
| | Dim of A_{aprt} | X (cm) | | - | 0.18 | - | # |
| | | Y (cm) | | - | 0.19 | - | # |
| Other Info | t_d | PD (ms) | # | | | | |
| | prr | PRF (Hz) | # | | | | |
| | p_r at max I_{pi} | $P_r @ PII_{max}$ (MPa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{eq} @$ PII_{max} (cm) | | | | - | |
| | Focal Length | FL_x (cm) | | - | 0.08 | - | # |
| | | FL_y (cm) | | - | 0.18 | - | # |
| Operator Control | $I_{pi\alpha}$ at max MI | $I_{PA,3} @$ MI_{max} (W/cm ²) | # | | | | |
| | Sample Volume Position (mm) | | - | - | 130 | - | - |
| | Frequency (MHz) | | - | - | 5.00 | - | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav™ 10/
SoundStar™ 3D 10FG / eco 10FG

Operating Mode: 2D

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC |
|------------------------------|---|-------------|---|-------|-------|----------|---|--------------------------|-----|
| | | | | | scan | non-scan | | | |
| | | | | | | | | $A_{\text{aprt}} \leq 1$ | |
| Maximum Index Value | | | | 1.03 | 1.07 | - | - | - | (b) |
| Assoc. Acoustic Param. | IEC | | FDA Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ (MPa) | 2.64 | | | | | |
| | P | | W_0 (mW) | | 34.69 | - | | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | | Min of $[W_{.3}(Z_1),$ $I_{\text{TA}.3}(Z_1)]$ (mW) | | | | - | | |
| | z_s | | z_1 (cm) | | | | - | | |
| | z_{bp} | | z_{bp} (cm) | | | | - | | |
| | z_b | | z_{sp} (cm) | | | | | - | |
| | z at max $I_{\text{pi } \alpha}$ | | z_{sp} (cm) | 0.72 | | | | | |
| | $d_{\text{eq}}(z_b)$ | | $d_{\text{eq}}(Z_{\text{sp}})$ (cm) | | | | | - | |
| | f_{awf} | | f_c (MHz) | 6.88 | 6.50 | - | - | - | # |
| | Dim of A_{aprt} | X (cm) | | | 0.70 | - | - | - | # |
| | | Y (cm) | | | 0.30 | - | - | - | # |
| Other Info | t_d | | PD (ms) | 0.22 | | | | | |
| | p_{rr} | | PRF (Hz) | 47 | | | | | |
| | p_r at max I_{pi} | | $P_r @ P_{\text{II} \text{max}}$ (MPa) | 3.13 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{\text{eq}} @$ $P_{\text{II} \text{max}}$ (cm) | | | | | - | |
| | Focal Length | FL_x (cm) | | | 0.62 | - | - | | # |
| | | FL_y (cm) | | | 0.22 | - | - | | # |
| | $I_{\text{pi } \alpha}$ at max MI | | $I_{\text{PA}.3} @$ MI_{max} (W/cm ²) | 305 | | | | | |
| Operator Control | Power | | (dB) | 0 | 0 | - | - | - | - |
| | Tilt | | (deg) | 0 | 0 | - | - | - | - |
| | Framerate | | (index) | 1 | 1 | - | - | - | - |
| | Frequency | | (MHz) | 8.00 | 8.00 | - | - | - | - |
| | Width | | (deg or ratio to max width) | 45.00 | 45.00 | - | - | - | - |
| | Depth | | (mm) | 160 | 160 | - | - | - | - |
| | Focus | | (mm) | 34 | 99 | - | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav™ 10/SoundStar™
3D 10FG / eco 10FG

Operating Mode: M-Mode

| Index Label | | | | MI | TIS | | | TIB non-scan | TIC | |
|------------------------------|--|--|---|----------------------|------|----------|------|-------------------|-----|----------------|
| | | | | | scan | non-scan | | | | |
| | | | | | | | | $A_{aprt} \leq 1$ | | $A_{aprt} > 1$ |
| Maximum Index Value | | | | 1.04 | - | 0.27 | - | - | (b) | |
| Assoc. Acoustic Param. | IEC | | FDA | Units | | | | | | |
| | $P_{r,\alpha}$ | | $P_{r,3}$ | (MPa) | 2.11 | | | | | |
| | P | | W_0 | (mW) | | - | 7.83 | | # | |
| | Min of $[P_\alpha(Z_s),$ $I_{ta,\alpha}(Z_s)]$ | | Min of $[W_{,3}(Z_1),$ $I_{TA,3}(Z_1)]$ | (mW) | | | | - | | |
| | Z_s | | Z_1 | (cm) | | | | - | | |
| | Z_{bp} | | Z_{bp} | (cm) | | | | - | | |
| | Z_b | | Z_{sp} | (cm) | | | | | - | |
| | z at max $I_{pi\alpha}$ | | Z_{sp} | (cm) | 0.72 | | | | | |
| | $d_{eq}(z_b)$ | | $d_{eq}(Z_{sp})$ | (cm) | | | | | - | |
| | f_{awf} | | f_c | (MHz) | 4.75 | - | 7.25 | - | - | # |
| | Dim of A_{aprt} | | X | (cm) | | - | 0.70 | - | - | # |
| | | | Y | (cm) | | - | 0.30 | - | - | # |
| Other Info | t_d | | PD | (ms) | 0.60 | | | | | |
| | p_{rr} | | PRF | (Hz) | 1000 | | | | | |
| | p_r at max I_{pi} | | $P_r@ P_{II_{max}}$ | (MPa) | 2.37 | | | | | |
| | d_{eq} at max I_{pi} | | $d_{eq}@$ $P_{II_{max}}$ | (cm) | | | | | - | |
| | Focal Length | | FL_x | (cm) | | - | 0.58 | - | | # |
| | | | FL_y | (cm) | | - | 0.24 | - | | # |
| | $I_{pi\alpha}$ at max MI | | $I_{PA,3}@$ MI _{max} | (W/cm ²) | 95 | | | | | |
| Operator Control | Power | | | (dB) | 0 | - | 0 | - | - | - |
| | Beam Angle | | | (deg) | 0 | - | 0 | - | - | - |
| | Frequency | | | (MHz) | 4.44 | - | 8.00 | - | - | - |
| | Depth | | | (mm) | 160 | - | 160 | - | - | - |
| | Focus | | | (mm) | 29 | - | 79 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: AcuNav™ 10/SoundStar™
3D 10FG / eco 10FG

Operating Mode: CMM

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|----------------------|------|------|--------------------------|-----------------------|-----|
| | | | | | scan | non-scan | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | | 1.05 | - | 0.28 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 1.12 | | | | |
| | P | W_0 | (mW) | | - | 6.24 | | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | |
| | z_s | z_1 | (cm) | | | | - | |
| | z_{bp} | z_{bp} | (cm) | | | | - | |
| | z_b | z_{sp} | (cm) | | | | | - |
| | z at max $I_{\text{pi},\alpha}$ | z_{sp} | (cm) | 1.77 | | | | |
| | $d_{\text{eq}}(z_b)$ | $d_{\text{eq}}(z_{\text{sp}})$ | (cm) | | | | | - |
| | f_{awf} | f_c | (MHz) | 4.15 | - | 6.65 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.70 | - | # |
| | | Y | (cm) | | - | 0.30 | - | # |
| Other Info | t_d | PD | (ms) | 1.17 | | | | |
| | p_{rr} | PRF | (Hz) | 6024 | | | | |
| | p_r at max $I_{\text{pi},\alpha}$ | $P_r @ P_{\text{II,max}}$ | (MPa) | 1.44 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | | - |
| | Focal Length | FL_x | (cm) | | - | 0.52 | - | # |
| | | FL_y | (cm) | | - | 0.26 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 46 | | | | |
| | Power | (dB) | | 0 | - | 0 | - | - |
| | PRF | (Hz) | | 6024 | - | 4000 | - | - |
| | ROI Span | (mm) | | 30 | - | 200 | - | - |
| | ROI Center | (mm) | | 20 | - | 20 | - | - |
| | Sample Volume | (mm) | | 0.74 | - | 0.77 | - | - |
| | Frequency | (MHz) | | 4.00 | - | 6.67 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav™ 10/SoundStar™
3D 10FG / eco 10FG

Operating Mode: CFM

| Index Label | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|-------|-------|--------------------------|-----------------------|-----|
| | | | | scan | non-scan | | |
| Maximum Index Value | | | 1.05 | 1.35 | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | (b) |
| Assoc. Acoustic Param. | IEC | FDA Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ (MPa) | 2.01 | | | | |
| | P | W_0 (mW) | | 37.23 | - | | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{.3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ (mW) | | | | - | |
| | z_s | z_1 (cm) | | | | - | |
| | z_{bp} | z_{bp} (cm) | | | | - | |
| | z_b | z_{sp} (cm) | | | | | - |
| | z at max $I_{\text{pi } \alpha}$ | z_{sp} (cm) | 0.72 | | | | |
| | $d_{\text{eq}}(z_b)$ | $d_{\text{eq}}(z_{\text{sp}})$ (cm) | | | | | - |
| | f_{awf} | f_c (MHz) | 6.60 | 6.60 | - | - | - |
| | Dim of A_{aprt} | X (cm) | | 0.70 | - | - | - |
| | | Y (cm) | | 0.30 | - | - | - |
| Other Info | t_d | PD (ms) | 0.45 | | | | |
| | p_{rr} | PRF (Hz) | 250 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II}_{\text{max}}}$ (MPa) | 2.36 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ (cm) | | | | | - |
| | Focal Length | FL_x (cm) | | 0.60 | - | - | # |
| | | FL_y (cm) | | 0.22 | - | - | # |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} (W/cm ²) | 144 | | | | |
| | Power | (dB) | 0 | 0 | - | - | - |
| | Tilt | (deg) | 0 | 0 | - | - | - |
| | Framerate | (index) | 1 | 1 | - | - | - |
| | PRF | (Hz) | 250 | 250 | - | - | - |
| | ROI Span | (mm) | 110 | 110 | - | - | - |
| | ROI Center | (mm) | 35 | 35 | - | - | - |
| | Sample Voume | (mm) | 0.42 | 0.42 | - | - | - |
| | ROI Width | (deg or ratio to max width) | 15.00 | 15.00 | - | - | - |
| | Frequency | (MHz) | 6.67 | 6.67 | - | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Acoustic information

Transducer Model: AcuNav™ 10/SoundStar™
3D 10FG / eco 10FG

Operating Mode: PW

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|--|----------------------|------|------|--------------------------|-----------------------|-----|
| | | | | | scan | non-scan | | |
| | | | | | | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ | |
| Maximum Index Value | | | | 0.63 | - | 1.01 | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | 1.38 | | | | |
| | P | W_0 | (mW) | | - | 31.89 | | # |
| | Min of $[P_{\alpha}(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_{,3}(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | |
| | z_s | z_1 | (cm) | | | | - | |
| | z_{bp} | z_{bp} | (cm) | | | | - | |
| | z_b | z_{sp} | (cm) | | | | | - |
| | z at max $I_{\text{pi},\alpha}$ | z_{sp} | (cm) | 0.67 | | | | |
| | $d_{\text{eq}}(z_b)$ | $d_{\text{eq}}(z_{\text{sp}})$ | (cm) | | | | | - |
| | f_{awf} | f_c | (MHz) | 4.13 | - | 6.65 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.62 | - | # |
| | | Y | (cm) | | - | 0.30 | - | # |
| Other Info | t_d | PD | (ms) | 1.06 | | | | |
| | p_{rr} | PRF | (Hz) | 2092 | | | | |
| | p_r at max I_{pi} | $P_r @ P_{\text{II,max}}$ | (MPa) | 1.52 | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II,max}}$ | (cm) | | | | | - |
| | Focal Length | FL_x | (cm) | | - | 0.52 | - | # |
| | | FL_y | (cm) | | - | 0.22 | - | # |
| Operator Control | $I_{\text{pi},\alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | 44 | | | | |
| | Power | | (dB) | 0 | - | 0 | - | - |
| | Beam Angle | | (cm) | 0 | - | 0 | - | - |
| | Sample Volume Position | | (mm) | 20 | - | 80 | - | - |
| | Sample Volume | | (mm) | 0.97 | - | 0.97 | - | - |
| | Scale | | (m/s) | 0.40 | - | 0.40 | - | - |
| | Frequency | | (MHz) | 4.00 | - | 6.67 | - | - |

Notes:

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Transducer Model: AcuNav™ 10/SoundStar™
3D 10FG / eco 10FG

Operating Mode: CW

| Index Label | | | | MI | TIS | | TIB non-scan | TIC |
|------------------------------|---|---|----------------------|-----|------|--|-----------------|-----|
| | | | | | scan | non-scan $A_{\text{aprt}} \leq 1$ $A_{\text{aprt}} > 1$ | | |
| Maximum Index Value | | | | (a) | - | 0.17 - | - | (b) |
| Assoc. Acoustic Param. | IEC | FDA | Units | | | | | |
| | $P_{r,\alpha}$ | $P_{r,3}$ | (MPa) | # | | | | |
| | P | W_0 | (mW) | | - | 7.14 | - | # |
| | Min of $[P_\alpha(Z_s),$ $I_{\text{ta},\alpha}(Z_s)]$ | Min of $[W_3(Z_1),$ $I_{\text{TA},3}(Z_1)]$ | (mW) | | | | - | |
| | Z_s | Z_1 | (cm) | | | | - | |
| | Z_{bp} | Z_{bp} | (cm) | | | | - | |
| | Z_b | Z_{sp} | (cm) | | | | | - |
| | z at max $I_{\text{pi } \alpha}$ | Z_{sp} | (cm) | # | | | | |
| | $d_{\text{eq}}(z_b)$ | $d_{\text{eq}}(Z_{\text{sp}})$ | (cm) | | | | | - |
| | f_{awf} | f_c | (MHz) | # | - | 5.00 | - | # |
| | Dim of A_{aprt} | X | (cm) | | - | 0.18 | - | # |
| | | Y | (cm) | | - | 0.30 | - | # |
| Other Info | t_d | PD | (ms) | # | | | | |
| | p_{rr} | PRF | (Hz) | # | | | | |
| | p_r at max $I_{\text{pi } \alpha}$ | $P_r @ P_{\text{II}_{\text{max}}}$ | (MPa) | # | | | | |
| | d_{eq} at max I_{pi} | $d_{\text{eq}} @$ $P_{\text{II}_{\text{max}}}$ | (cm) | | | | | - |
| | Focal Length | FL_x | (cm) | | - | 0.08 | - | # |
| | | FL_y | (cm) | | - | 0.22 | - | # |
| Operator Control | $I_{\text{pi } \alpha}$ at max MI | $I_{\text{PA},3} @$ MI_{max} | (W/cm ²) | # | | | | |
| | Sample Volume Position | | (mm) | - | - | 160 | - | - |
| | Frequency | | (MHz) | - | - | 5.00 | - | - |

Notes:

(a) This index is not required for this operating mode; see section 4.1.3.1 of the "Output Display Standard" (NEMA UD-3).

(b) This probe is not intended for transcranial or neonatal cephalic use.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed.

Chapter 4

Electromagnetic Compatibility

Vivid q N is intended for use in the electromagnetic environment specified in the tables below.


The user of Vivid q N should assure that the device is used in such an environment.

Guidance and manufacturer's declaration

| Guidance and manufacturer's declaration – electromagnetic emissions | | |
|---|------------|--|
| Emissions test | Compliance | Electromagnetic environment - guidance |
| RF emission CISPR 11 EN55011 | Group 1 | Vivid q N uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. Vivid q N system is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emission CISPR 11 EN55011 | Class A | |
| Harmonic emission EN/IEC 61000-3-2 | Class A | |
| Voltage fluctuations/ flicker emissions EN/IEC 61000-3-3 | Complies | |

| Guidance and manufacturer's declaration – electromagnetic immunity | | | |
|---|---|---|--|
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
| Electrostatic discharge (ESD) EC/IEC 61000-4-2 | ± 6 kV contact ± 8 kV air | ± 6 kV ± 8 kV | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrostatic transient / burst EC/IEC 61000-4-4 | ± 2 kV for power-supply lines ± 1 kV for input/output lines | ± 2 kV ± 1 kV | Mains power quality should be that of a typical commercial or hospital environment. |
| Surge EC/IEC 61000-4-5 | ± 1 kV differential mode ± 2 kV common mode | ± 1 kV ± 2 kV | Mains power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations on power supply input lines EC/IEC 61000-4-11 | < 5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec | Compliance for all test levels. Controlled shutdown with return to pre-disturbance condition after operator's intervention. (Power-on switch) | Mains power quality should be that of a typical commercial or hospital environment. If the user of Vivid q N requires continued operation during power mains interruptions, it is recommended that Vivid q N is powered from an uninterruptible power supply or a battery. |

| | | | |
|---|-------|---------------------|---|
| Power frequency (50/60 Hz) magnetic field EC/IEC 61000-4-8 | 3 A/m | 3A/m 50 and 60Hz | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment. |
| NOTE U_T is the a. c. mains voltage prior to application of the test level. | | | |

| Guidance and manufacturer's declaration – electromagnetic immunity. | | | |
|--|---------------------------------|--------------------------|--|
| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment – guidance ^c |
| Conducted RF EC/IEC 61000-4-6 | 3 Vrms 150 kHz to 80 MHz | 3 Vrms [V ₁] | <p>Portable and mobile RF communications equipment should be used no closer to any part of Vivid q N, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E_1} \right] \sqrt{P} \quad \begin{matrix} 80 \text{ MHz to} \\ 800 \text{ MHz} \end{matrix}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to 2,5 GHz}$ |
| Conducted RF IEC 61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3V/m [E ₁] | <p>where <i>p</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <i>d</i> is the recommended separation distance in metres (m).^b</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
| <p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.</p> | | | |

- ^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Vivid q N is used exceeds the applicable RF compliance level above, Vivid q N should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating Vivid q N.
- ^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.
- ^c See examples of calculated separation distances in next table.

| Recommended separation distances between portable and mobile RF communications equipment and Vivid q N | | | |
|---|--|---|--|
| Vivid q N is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Vivid q N can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Vivid q N as recommended below, according to the maximum output power of the communications equipment | | | |
| Rated maximum output of transmitter W | Separation distance according to frequency of transmitter m | | |
| | 150 kHz to 80 MHz $d = [\frac{3,5}{V_1}] \sqrt{P}$ | 80 MHz to 800 MHz $d = [\frac{3,5}{E_1}] \sqrt{P}$ | 800 MHz to 2,5 GHz $d = [\frac{7}{E_1}] \sqrt{P}$ |
| 0.01 | 0.12 | 0.12 | 0.23 |
| 0.1 | 0.38 | 0.38 | 0.73 |
| 1 | 1.2 | 1.2 | 2.3 |
| 10 | 3.8 | 3.8 | 7.3 |
| 100 | 12 | 12 | 23 |
| <p>For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p>NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p> | | | |



GE Medical Systems

*GE Medical Systems-Americas: Fax 414.544.3384
P.O. Box 414; Milwaukee, Wisconsin 53201-0414, U.S.A
GE Medical Systems-Europe: Fax 33.1.40.93.33.33
Paris, France
GE Medical Systems-Asia: Fax 65.291.7006
Singapore
GE Vingmed Ultrasound: Fax: +47 3302 1350
Horten, Norway*